








Bloodline/Tubing systems for blood purification

GENERAL INFORMATION

Refer to product or packaging label for:

- Type / Code of Bloodline/Tubing system
- Manufacturer and model of the machine for which the product is intended
- Applicable symbols

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Intended use: Bloodline/Tubing systems are intended for single use only for extracorporeal blood purification.

Contraindications: Special contraindications are unknown. General contraindications for extracorporeal blood purification are applicable.

Materials: Bloodline/Tubing systems: medical grade soft-PVC; connectors and other components: PP, PC, PVC, ABS, PE, PA, Silicone, PBT, PTFE, PHE, TPE, PS, PETA, PETG.

One or more components contain in a concentration > 0.1 mass% according to Article 33 and 59 (1, 10) of Regulation (EC) No. 1907/2006 („REACH“):

- 2-ethylhexyl 10-ethyl-4,4-dioctyl-7-oxo-8-oxa-3,5-dithia-4-stannatetradecanoate (DOTE)

DIRECTIONS FOR USE

The Bloodline/Tubing systems are designed for machines for extracorporeal blood purification multiFiltrate or multi (Version: multiFiltratePRO) (refer to label) and shall be used only after appropriate instruction or training. Use only Bloodline/Tubing systems approved for the selected therapy.

Refer to the Instructions for Use of the machine during set-up, treatment and termination. Also the Instructions for Use of other disposables such as the haemofilter, plasmapher, or adsorbents must be taken into account.

Colour codes should be followed and used in line with the corresponding markings on the machine. Do not use if the sterile package is damaged or the protective or closure caps are not in place. Do not use after use-by date (refer to label).

Unpack and connect the Bloodline/Tubing systems aseptically. For hygienic and functional reasons, it is recommended that tubing systems are inserted immediately before preparation only. Keep preparation and circulation times before the treatment as short as possible and in compliance with applicable guidelines.

Make sure all tubes and chambers are correctly inserted into the respective holders, avoiding twisting, kinking or occlusion of the tubes.

Ensure that all caps and connections are secure. To ensure correct functionality of the bubble catchers, fill to about 1 cm below the top. This avoids fluid contacting the hydrophobic filter as well as air entering the extracorporeal circuit. In the event that fluid reaches and comes into contact with the hydrophobic filter within the Transducer Protector (TP) the Bloodline/Tubing systems must be replaced. With bloodlines that have the possibility of an additional connection intended to be connected to the pressure measurement system, this can be used to connect a new pressure measurement line (available as an accessory in the Fresenius Medical Care product range). Do not push back the fluid in the pressure measurement line by using a syringe. This might damage the membrane in the TP and thus lead to contamination. In the event that fluid might have passed through the TP, check the machine for contamination after treatment is completed. If contamination has occurred, the machine must be taken out of service and disinfected according to the manufacturer's recommendations before further use.

Fill and rinse the Bloodline/Tubing systems in accordance with the Instructions for Use of the machine or supplemental Instructions for Use for additional options, and training, if applicable.

Use only needles with a beveled tip and a diameter equal or smaller than 20 gauge (outer needle diameter to be 0.9 mm or less) to puncture injection sites (if present).

The infusion administration lines should be clamped, except when needed. Handle pumps and pump clips in accordance to the Instructions for Use of the machine.

To ensure a secure connection between patient access and bloodline, hold and screw the coloured (blue, red) coupling nut on the bloodline only. Do not apply the screwing torque to the inner part of the connector. After the connection, check that the components are firmly screwed together.

CAUTION The Bloodline/Tubing systems are intended for single use only. Re-use may be hazardous to both patient and operator. Cleansing solutions and disinfectants may damage materials employed for the Bloodline/Tubing systems.

Fresenius Medical Care Bloodlines/Tubing systems are designed to withstand the maximum and minimum manufacturer's recommended pressures and flow rates generated in use with the respective multiFiltrate or multiFiltratePRO.

The Bloodline/Tubing systems must be replaced after the maximum usage time or respectively after a specific blood volume pumped as indicated on the primary packaging.

Excessive negative pressure may cause partial collapse of the pump segment resulting in an actual blood flow substantially less than indicated on the machine.

Minimum temperature of use: 18°C (64°F). Careful attention must be given to the installation and blood level of the bubble catchers (see above). Make sure that due to pressure fluctuation neither the blood contacts the hydrophobic filter nor air enters the extracorporeal circuit. Check the compatibility and safety of the connections (e.g. Luer-Lock) to the patient, to all connected disposables (e.g. haemofilter) and to the solution bags. Make sure that the components/tubes in direct interface with the machine are fitting properly.

Inspect the Bloodline/Tubing system for kinking and leaks during the priming phase and treatment, taking corrective measures (e.g. tightening Luer-Lock connection) or exchanging it as necessary.

The plastics used can be incompatible with drugs or disinfectants (e.g. connectors made of polycarbonate can develop cracks when in contact with aqueous solutions with pH > 10).

If nutritional solutions are administered into the bloodline, wetting of the luer-lock connection to the feeding system with lipidic fluids can weaken the properties of the plastic material used. Make sure that during the connection of the feeding system, the infusion line close to the luer-lock connection site remains completely free of nutrition solution.

Change the Bloodline/Tubing system if a visible clot formation („white strip“) is seen in the bloodline downstream of the return line bubble catcher, for instance at the point where the calcium infusion line merges into the bloodline.

WASTE UTILIZATION Used products shall be subject to disposal and destruction in accordance with the current legislation of the country in which the Bloodline/Tubing system is used. If required, the waste should be disposed by an accredited organization in accordance with the applicable law.

WARRANTY The manufacturer shall not be liable for any misuse, improper handling, non compliance with instructions for use and cautionary notes, and for any damage incurred subsequent to the manufacturer's delivery of the Bloodline/Tubing systems.

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