



Shanghai United Imaging Healthcare Co., Ltd  
% Shumei Wang  
QM & RA VP  
No. 2258 Chengbei Rd., Jiading Industrial District  
Shanghai, 201807  
CHINA

January 26, 2020

Re: K193241  
Trade/Device Name: uMI 550  
Regulation Number: 21 CFR 892.1200  
Regulation Name: Emission computed tomography system  
Regulatory Class: Class II  
Product Code: KPS, JAK  
Dated: November 22, 2019  
Received: November 25, 2019

Dear Shumei Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.  
Director  
Division of Radiological Health  
OHT7: Office of In Vitro Diagnostics  
and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K193241

Device Name

uMI 550

Indications for Use (Describe)

The uMI 550 PET/CT is a diagnostic imaging system that combines two existing imaging modalities - PET and CT. The quantitative distribution information of PET radiopharmaceuticals within the patient body measured by PET can assist healthcare providers in assessing the metabolic and physiological functions. CT provides diagnostic tomographic anatomical information as well as photon attenuation information for the scanned region. The accurate registration and fusion of PET and CT images provides anatomical reference for the findings in the PET images.

This system is intended to be operated by qualified healthcare professionals to assist in the detection, localization, diagnosis, staging, restaging, treatment planning and treatment response evaluation for diseases and disorders in, but not limit to, oncology, cardiology and neurology.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) SUMMARY

**1. Date of Preparation:** January 22, 2020

**2. Sponsor Identification**

**Shanghai United Imaging Healthcare Co.,Ltd.**

No.2258 Chengbei Rd. Jiading District, 201807, Shanghai, China

Contact Person: Shumei Wang

Position: QM&RA VP

Tel: +86-021-67076888-6776

Fax: +86-021-67076889

Email: [shumei.wang@united-imaging.com](mailto:shumei.wang@united-imaging.com)

**3. Identification of Proposed Device**

**Trade Name:** uMI 550

**Common Name:** Emission Computed Tomography System

**Model(s):** uMI 550

**Regulatory Information**

**Regulation Number:** 21 CFR 892.1200

**Regulation Name:** Emission Computed Tomography System

**Regulatory Class:** II

**Product Code:** KPS, JAK,

**Review Panel:** Radiology

**4. Identification of Predicate Device(s)**

**Predicate Device**

**510(k) Number:** K182237

**Device Name:** uMI 550

**Model(s):** uMI 550

**Regulatory Information**

**Regulation Number:** 21 CFR 892.1200

**Regulation Name:** Emission Computed Tomography System

**Regulatory Class:** II

**Product Code:** KPS, JAK,

**Review Panel:** Radiology

**5. Device Description**

The uMI 550 PET/CT system is a combined multi-slice X-Ray Computed Tomography and Positron Emission Tomography scanner. This system is intended

to be operated by qualified healthcare professionals for performing diagnostic imaging examinations. The spatial alignment and precise image registration between PET and CT ensure the PET and CT images of the same region can be fused accurately for reading. PET measures the distribution of PET radiopharmaceuticals inside the human body quantitatively. CT produces the anatomical information of the same scanned region, and provides accurate localization for the findings in the PET images. The attenuation information contained in the CT images can be utilized in the PET image reconstruction to ensure quantitation accuracy.

The uMI 550 PET/CT system also includes a patient table, a workstation with associated software installed. The software is used for patient management, data management, scan control, image reconstruction and image reading. All patient images produced by the system conform to the DICOM 3.0 standard.

The uMI 550 PET/CT has been previously cleared by FDA via K182237. The modifications performed on the uMI 550 PET/CT (K182237) in this submission are due to the addition of HYPER Iterative and Auto-Planbox function. Meanwhile the sensitivity specification has been updated. HYPER Iterative allows more iterations while remains the image noise at an acceptable level by incorporating a noise control term into the objective function. HYPER Iterative can achieve high image contrast and quantification accuracy. Auto-Planbox plan the scan range by recognizing body parts on CT scout image. It locates the different body parts based on anatomy characteristic. The scan range is generated to cover the whole body parts according to protocol selection. This function will simplify scanning process, which will be convenient for user operation.

The modifications are listed as follow:

- Introduce HYPER Iterative function.
- Introduce Auto-Planbox function.
- Update sensitivity specification.

|                 | uMI 550 PET/CT previous submission(K182237) | uMI 550 PET/CT this submission        | Remark                                      |
|-----------------|---|---------------------------------------|---|
| Sensitivity     | >/=10cps/kBq @0cm<br>>/=10cps/kBq @10cm     | >/=9cps/kBq @0cm<br>>/=9cps/kBq @10cm | Sensitivity specification has been updated. |
| HYPER Iterative | Not available                               | Available                             | New function has been added                 |
| Auto-Planbox    | Not available                               | Available                             | New function has been added                 |

## **6. Indications for Use**

The uMI 550 PET/CT is a diagnostic imaging system that combines two existing imaging modalities - PET and CT. The quantitative distribution information of PET radiopharmaceuticals within the patient body measured by PET can assist healthcare providers in assessing the metabolic and physiological functions. CT provides diagnostic tomographic anatomical information as well as photon attenuation information for the scanned region. The accurate registration and fusion of PET and CT images provides anatomical reference for the findings in the PET images.

This system is intended to be operated by qualified healthcare professionals to assist in the detection, localization, diagnosis, staging, restaging, treatment planning and treatment response evaluation for diseases and disorders in, but not limit to, oncology, cardiology and neurology.

## **7. Comparison of Technological Characteristics with the Predicate Devices**

The technology characteristics of the modified uMI 550, reflected in this 510(k) submission, do not alter the scientific technology of the devices and are substantially equivalent to those of the predicate devices.

Table 1 below provides a comparison of the technological characteristics of the proposed device in comparison to the predicate device.

Table 1 Comparison of Technological Characteristics

| ITEM           | Proposed Device<br>uMI 550  | Predicated Device<br>uMI 550(K182237)   | Remark |
|----------------|---|---|--------|
| <b>General</b> |   |   |        |
| Product Code   | KPS, JAK  | KPS, JAK  | Same   |
| Regulation No. | 21 CFR 892.1200   | 21 CFR 892.1200   | Same   |
| Class          | Class II  | Class II  | Same   |
| Intended Use   | <p>The uMI 550 PET/CT is a diagnostic imaging system that combines two existing imaging modalities - PET and CT. The quantitative distribution information of PET radiopharmaceuticals within the patient body measured by PET can assist healthcare providers in assessing the metabolic and physiological functions. CT provides diagnostic tomographic anatomical information as well as photon attenuation information for the scanned region. The accurate registration and fusion of PET and CT images provides anatomical reference for the findings in the PET images.</p> <p>This system is intended to be operated by qualified healthcare professionals to assist in the detection, localization, diagnosis, staging, restaging, treatment planning and treatment response evaluation for diseases and disorders in, but not limit to, oncology, cardiology and neurology.</p> | <p>The uMI 550 PET/CT is a diagnostic imaging system that combines two existing imaging modalities - PET and CT. The quantitative distribution information of PET radiopharmaceuticals within the patient body measured by PET can assist healthcare providers in assessing the metabolic and physiological functions. CT provides diagnostic tomographic anatomical information as well as photon attenuation information for the scanned region. The accurate registration and fusion of PET and CT images provides anatomical reference for the findings in the PET images.</p> <p>This system is intended to be operated by qualified healthcare professionals to assist in the detection, localization, diagnosis, staging, restaging, treatment planning and treatment response evaluation for diseases and disorders in, but not limit to, oncology, cardiology and neurology.</p> | Same   |



| <b>Specifications</b>    |  |                                |  |                                |                   |                        |      |
|--------------------------|--|--------------------------------|--|--------------------------------|-------------------|------------------------|------|
| PET Specifications       |  |                                |  |                                |                   |                        |      |
| Scintillator material    | LYSO   |                                | LYSO   |                                | Same              |                        |      |
| Scintillator dimensions  | 2.76mm×2.76mm×16.3mm                                 |                                | 2.76mm×2.76mm×16.3mm                                 |                                | Same              |                        |      |
| Number of detector rings | 84   |                                | 84   |                                | Same              |                        |      |
| Number of image planes   | 167  |                                | 167  |                                | Same              |                        |      |
| Axial field of view      | 24cm   |                                | 24cm   |                                | Same              |                        |      |
| Image matrix sizes       | 128×128, 150×150, 192×192, 256×256, 512×512, 600×600 |                                | 128×128, 150×150, 192×192, 256×256, 512×512, 600×600 |                                | Same              |                        |      |
| Coincidence window       | 4.0ns  |                                | 4.0ns  |                                | Same              |                        |      |
| Spatial Resolution       | Axial FWHM@1cm                                       | <math>\leq 3.5\text{mm}</math> |  | <math>\leq 3.5\text{mm}</math> |                   | Same                   |      |
|                          | Transaxial FWHM@1cm                                  | <math>\leq 3.5\text{mm}</math> |  | <math>\leq 3.5\text{mm}</math> |                   |                        |      |
|                          | Axial FWHM@10cm                                      | <math>\leq 4.0\text{mm}</math> |  | <math>\leq 4.0\text{mm}</math> |                   |                        |      |
|                          | Transaxial FWHM@10cm                                 | <math>\leq 4.0\text{mm}</math> |  | <math>\leq 4.0\text{mm}</math> |                   |                        |      |
| Sensitivity              | >math>\geq 9\text{cps/kBq}</math>                    |                                | >math>\geq 10\text{cps/kBq}</math>                   |                                | Note No.1         |                        |      |
| NECR Peak Value          | >math>\geq 90\text{kcps}@13\text{kBq/cc}</math>      |                                | >math>\geq 90\text{kcps}@13\text{kBq/cc}</math>      |                                | Same              |                        |      |
| Peak True Count Rate     | >math>\geq 300\text{kcps}@27\text{kBq/cc}</math>     |                                | >math>\geq 300\text{kcps}@27\text{kBq/cc}</math>     |                                | Same              |                        |      |
| PET Scatter Fraction     | <math>\leq 0.44</math>                               |                                | <math>\leq 0.44</math>                               |                                | Same              |                        |      |
| Count Rate Bias          | <math>\leq \pm 5\%</math>                            |                                | <math>\leq \pm 5\%</math>                            |                                | Same              |                        |      |
| Image Quality            | Sphere   | Acceptance Value               |  | Sphere                         | Acceptance Value  |                        | Same |
|                          |  | Contrast Recovery              | Background variability                               |                                | Contrast Recovery | Background variability |      |
|                          | 10mm   | $\geq 45.0\%$                  | $< 7.5\%$  | 10mm                           | $\geq 45.0\%$     | $< 7.5\%$              |      |
|                          | 13mm   | $\geq 55.0\%$                  | $< 7.0\%$  | 13mm                           | $\geq 55.0\%$     | $< 7.0\%$              |      |
|                          | 17mm   | $\geq 65.0\%$                  | $< 7.0\%$  | 17mm                           | $\geq 65.0\%$     | $< 7.0\%$              |      |

|  |                                    |        |                                    |               |        |       |  |
|--|------------------------------------|--------|------------------------------------|---------------|--------|-------|--|
|  | 22 mm                              | ≥72.0% | <7.0%                              | 22 mm         | ≥72.0% | <7.0% |  |
|  | 28 mm                              | ≥65.0% | <7.0%                              | 28 mm         | ≥65.0% | <7.0% |  |
|  | 37 mm                              | ≥70.0% | <7.0%                              | 37 mm         | ≥70.0% | <7.0% |  |
|  | Long Residual                      | ≤16.0% | <7.0%                              | Long Residual | ≤16.0% | <7.0% |  |
| <b>CT Specifications</b>                                       |                                    |        |                                    |               |        |       |  |
| Scan Regime  | Continuous Rotation                |        | Continuous Rotation                |               | Same   |       |  |
| Scan Modes   | Topo<br>Axial Scan<br>Helical Scan |        | Topo<br>Axial Scan<br>Helical Scan |               | Same   |       |  |
| Detector Material  | Solid-state GOS                    |        | Solid-state GOS                    |               | Same   |       |  |
| Z-plane coverage   | 22mm                               |        | 22mm                               |               | Same   |       |  |
| Size of detector element in Z-plane                            | 0.55mm                             |        | 0.55mm                             |               | Same   |       |  |
| Number of element per row                                      | 864                                |        | 864                                |               | Same   |       |  |
| Number of detector row   | 40                                 |        | 40                                 |               | Same   |       |  |
| Maximum slices generated per rotation (multi-slice capability) | 80                                 |        | 80                                 |               | Same   |       |  |
| Minimum slice thickness  | 0.55mm                             |        | 0.55mm                             |               | Same   |       |  |
| Maximum sampling rate  | Up to 4800 views per 360°          |        | Up to 4800 views per 360°          |               | Same   |       |  |
| Tube anode storage capacity                                    | 5.3MHU                             |        | 5.3MHU                             |               | Same   |       |  |
| Maximum cooling rate   | 815 kHU/min                        |        | 815 kHU/min                        |               | Same   |       |  |
| Focal spot size  | 0.5x1.0mm<br>1.0x1.0mm             |        | 0.5x1.0mm<br>1.0x1.0mm             |               | Same   |       |  |
| Power  | 50 kW                              |        | 50 kW                              |               | Same   |       |  |
| mA Range   | 10-420mA                           |        | 10-420mA                           |               | Same   |       |  |
| kV Settings  | 70, 80, 100, 120, 140              |        | 70, 80, 100, 120, 140              |               | Same   |       |  |
| Aperture   | 700mm                              |        | 700mm                              |               | Same   |       |  |
| Rotation speed   | Up to 0.5 sec per 360° rotation    |        | Up to 0.5 sec per 360° rotation    |               | Same   |       |  |
| Image Spatial Resolution                                       | High mode:<br>>20 lp/cm @ MTF 0%   |        | High mode:<br>>20 lp/cm @ MTF 0%   |               | Same   |       |  |

|  |  |  |      |
|--|--|--|------|
|  | 16.5±1.7 lp/cm @<br>MTF10%<br>11.5±1.2 lp/cm @<br>MTF50%   | 16.5±1.7 lp/cm @<br>MTF10%<br>11.5±1.2 lp/cm @<br>MTF50%   |      |
| Image Noise                              | 3.0±0.5 HU at 120 kV, 5<br>mm slice thickness,<br>CTDIvol 28.9 mGy   | 3.0±0.5 HU at 120 kV, 5<br>mm slice thickness,<br>CTDIvol 28.9 mGy   | Same |
| CT Number<br>Display Range               | -1024 ~+8191 HU  | -1024 ~+8191 HU  | Same |
| Scan Field of<br>View                    | Up to 500 mm   | Up to 500 mm   | Same |
| Reconstruction<br>Field of View          | 40mm-500mm<br>40mm-600mm with<br>extend FOV  | 40mm-500mm<br>40mm-600mm with<br>extend FOV  | Same |
| Image Matrix                             | Up to 1024 x 1024  | Up to 1024 x 1024  | Same |
| Reconstructed<br>slice thickness         | 0.55mm,1.1mm,2.2mm,<br>5.5mm,11mm (axial)<br>0.55-10mm(helical)  | 0.55mm,1.1mm,2.2mm,<br>5.5mm,11mm (axial)<br>0.55-10mm(helical)  | Same |
| Pitch                                    | 0.1~2.0  | 0.1~2.0  | Same |
| Maximum<br>continuous<br>exposure time   | Up to 100seconds   | Up to 100seconds   | Same |
| <b>Others Specifications</b>             |  |  |      |
| Display                                  | 24inch, 1920 x 1200  | 24inch, 1920 x 1200  | Same |
| Horizontal motion<br>range               | 2160 mm  | 2160 mm  | Same |
| Vertical motion<br>range                 | 550 mm-945 mm from<br>the floor  | 550 mm-945 mm from<br>the floor  | Same |
| Maximum<br>horizontal<br>scannable range | Standard 1700mm with<br>up to full 2000mm.   | Standard 1700mm with<br>up to full 2000mm.   | Same |
| Table Horizontal<br>Speed                | Up to 200mm/sec  | Up to 200mm/sec  | Same |
| Vertical speed                           | Up to 40 mm/sec  | Up to 40 mm/sec  | Same |
| Table Horizontal<br>Position accuracy    | ±0.25mm  | ±0.25mm  | Same |
| Table Maximum<br>table load              | 250kg  | 250kg  | Same |
| <b>Safety</b>                            |  |  |      |
| Electrical Safety                        | ANSI AAMI ES60601-<br>1:2005/(R)2012 and<br>A1:2012,<br>C1:2009/(R)2012 and<br>A2:2010/(R)2012                         | ANSI AAMI ES60601-<br>1:2005/(R)2012 and<br>A1:2012,<br>C1:2009/(R)2012 and<br>A2:2010/(R)2012                         | Same |
| EMC                                      | IEC 60601-1-2 Medical<br>electrical equipment -<br>Part 1-2: General<br>requirements for basic<br>safety and essential | IEC 60601-1-2 Medical<br>electrical equipment -<br>Part 1-2: General<br>requirements for basic<br>safety and essential | Same |

|                  |   |  |      |
|------------------|---|--|------|
|                  | performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests  | performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests   |      |
| Biocompatibility | Patient Contact Materials were tested and demonstrated no cytotoxicity (ISO 10993-5), no evidence for irritation and sensitization (ISO 10993-10).  | Patient Contact Materials were tested and demonstrated no cytotoxicity (ISO 10993-5), no evidence for irritation and sensitization (ISO 10993-10). | Same |
| Clinical         | Comparing with the predicate device, the difference of clinical image is due to the addition of HYPER Iterative function. The comparison sample images showed improved image contrast, resolution and increased SNR.  |  |      |
| Justification    |   |  |      |
| Note ID          | Justification   |  |      |
| Note 1           | Since activity measurement factor has been updated according to well-counter vendor's notice. The sensitivity specification, which depends on activity measurement, shall be updated accordingly. Provide the lower system sensitivity for PET imaging which would be compensated by a little longer scan time. The difference of the system sensitivity from the two devices does not affect the safety and effectiveness. |  |      |

## 8. Performance Data

The following performance data were provided in support of the substantial equivalence determination.

### Non-Clinical Testing

Non-clinical testing including dosimetry and image performance tests were conducted for the uMI 550 during the product development.

UNITED IMAGING HEALTHCARE claims conformance to the following standards and guidance:

### Electrical Safety and Electromagnetic Compatibility (EMC)

Electrical Safety and Electromagnetic Compatibility (EMC) testing were conducted on the uMI 550 in accordance with the following standards:

- ANSI AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012 (Consolidated Text) Medical electrical equipment - Part 1: General requirements for basic safety and essential performance (IEC 60601-1:2005, MOD)
- IEC 60601-1-2 Edition 4.0 2014-02 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests

- IEC 60601-2-44 Edition 3.2: 2016 Medical electrical equipment - Part 2-44: Particular requirements for the basic safety and essential performance of x-ray equipment for computed tomography
- IEC 60825-1 Edition 2.0 2007-03 Safety of laser products - Part 1: Equipment classification, and requirements [Including: Technical Corrigendum 1 (2008), Interpretation Sheet 1 (2007), Interpretation Sheet 2 (2007)]

#### **Product Particular Standards**

- IEC 60601-1-3 Edition 2.1 2013-04 Medical electrical equipment - Part 1-3: General requirements for basic safety and essential performance - Collateral Standard: Radiation protection in diagnostic X-ray equipment
- IEC 61223-3-5 First edition 2004-08 Evaluation and routine testing in medical imaging departments - Part 3-5: Acceptance tests -Imaging performance of computed tomography X-ray equipment [Including: Technical Corrigendum 1 (2006)]
- NEMA XR 25-2010: Computed Tomography Dose Check
- NEMA XR 28-2013 Supplemental Requirements for User Information and System Function Related to Dose in CT
- NEMA XR 29-2013: Standard Attributes on CT Equipment Related to Dose Optimization and Management

#### **Performance Verification**

- NEMA NU 2-2012 Performance Measurements of Positron Emission Tomographs
- Clinical Image Evaluation for sample clinical images
- AEC Test for AEC performance study
- Performance Evaluation for HYPER Iterative function
- Performance Evaluation for Auto-Planbox function

Bench testing was performed to support substantial equivalence and the product performance claims. For HYPER Iterative function, the evaluation and analysis used the identical raw datasets obtained on UIH's uMI 550 system and then applied both OSEM and HYPER Iterative. The resultant images were compared for:

- Quantification accuracy and signal to noise ratio (SNR) using the NEMA IQ phantom
- Lesion detectability using patient data
- Influence of regularization parameter selection using patient data
- Effectiveness for large weight patient
- Image contrast improvement on brain imaging

Bench testing showed that HYPER Iterative can be used in PET body and brain imaging, can provide the flexibility to achieve user-preferable image by adjusting the regularization parameter, and its images have demonstrated improved image quality, higher quantification accuracy, and improved image SNR than OSEM images.

For Auto-Planbox function, the evaluation and analysis used 16 group scout images and compared Auto-Planbox recognition with manual annotation body parts position. Bench testing results showed that Auto-Planbox can locate head center position, body

start position and body end position with proper accuracy and can assist operator to plan the scan range.

### **Software**

- NEMA PS 3.1-3.20(2011): Digital Imaging and Communications in Medicine (DICOM)
- IEC 62304: Medical Device Software - software life cycle process
- Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices
- Content of Premarket Submissions for Management of Cybersecurity in Medical Devices

### **Biocompatibility**

- ISO 10993-10 Third Edition 2010-08-01 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
- ISO 10993-5 Third edition 2009-06-01 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity

### **Other Standards and Guidances**

- ISO 14971: Medical Devices – Application of risk management to medical devices
- Code of Federal Regulations, Title 21, Part 820 - Quality System Regulation
- Code of Federal Regulations, Title 21, Subchapter J - Radiological Health
- Laser Products - Conformance with IEC 60825-1 and IEC 60601-2-22; Guidance for Industry and FDA Staff (Laser Notice No. 50)
- Provision for Alternate Measure of the Computed Tomography Dose Index (CTDI) to Assure Compliance with the Dose Information Requirements of the Federal Performance Standard for Computed Tomography

### **Software Verification and Validation**

Software documentation for a Moderate Level of Concern software per FDA' Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" is included as a part of this submission.

The risk analysis was completed and risk control was implemented to mitigate identified hazards. The testing results show that all the software specifications have met the acceptance criteria. Verification and validation testing of the proposed device was found acceptable to support the claim of substantial equivalence.

UNITED IMAGING HEALTHCARE conforms to the Cybersecurity requirements by implementing a process of preventing unauthorized access, modification, misuse or denial of use, or unauthorized use of information that is stored, accessed, or transferred from a medical device to an external recipient. Cybersecurity information in accordance with guidance document "Content of Premarket Submissions for Management of Cybersecurity in Medical Devices" is included in this submission.

### **Clinical Image Evaluation**

The reader study used a total of 20 retrospectively collected clinical cases. The raw data from each of these cases was reconstructed with both OSEM and HYPER Iterative. Each image was read by 3 American board-certified nuclear medicine

physicians who provided an assessment of image contrast according to 3-point scale (“+” represents better contrast image for small lesion detection; “-” represents worse contrast image for small lesion detection; “=” represents similar contrast image for small lesion detection), and image quality according to a 5-point scale (1= images are unreadable and cannot be used for diagnosis, 2=insufficient image quality and may affect diagnosis; 3=image quality is barely ok for diagnosis; 4=image quality is sufficient for diagnosis; 5=wonderful image quality for diagnosis). The results of the study indicated that HYPER Iterative has better image contrast than OSEM and the image quality is sufficient for clinical diagnosis.

### **Clinical Testing**

No Clinical Study is included in this submission.

### **9. Conclusions**

Based on the comparison and analysis above, the proposed device has same intended use, similar performance, equivalence safety and effectiveness as the predicate device. The differences above between the proposed device and predicate device do not affect the intended use, technology characteristics, safety and effectiveness. And no issues are raised regarding to safety and effectiveness. The proposed device is determined to be Substantially Equivalent (SE) to the predicate device.