



January 9, 2022

Careglove Global SDN. BHD.  
Siti Ahmad  
Official Correspondent  
Careglove Global SDN BHD  
Lot 17479, Lorong Senawang 2/3, Off Jalan Senawang 3,  
Senawang Industrial Estate  
Seremban, Negeri Sembilan Darul Khusus 70450  
Malaysia

Re: K211666

Trade/Device Name: Powder Free Nitrile Examination Gloves, Blue Chemotest  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-Powdered Patient Examination Glove  
Regulatory Class: Class I, reserved  
Product Code: LZA, LZC  
Dated: November 25, 2021  
Received: December 6, 2021

Dear Siti Ahmad:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray III, PhD  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K211666

Device Name

POWDER FREE NITRILE EXAMINATION GLOVES, BLUE, CHEMOTEST

Indications for Use (Describe)

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

These gloves were tested for use with Chemotherapy Drugs as per ASTM D6978-05 Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs:

Chemotherapy Drugs Permeation

The following chemicals have been tested with these gloves.

Chemotherapy Drugs Concentration Breakthrough Detection Time in Minutes

\*Carmustine (BCNU) 3.3 mg/ml 22.6

Cisplatin 1.0 mg/ml >240 min

Cyclophosphamide (Cytosan) 20.0 mg/ml >240 min

Dacarbazine 10.0 mg/ml >240 min

Doxorubicin HCL 2.0 mg/ml >240 min

Etoposide 20.0 mg/ml >240 min

Fluorouracil 50.0 mg/ml >240 min

Ifosfamide 50.0 mg/ml >240 min

Mitoxantrone 2 mg/ml >240 min

Paclitaxel 6.0 mg/ml >240 min

\*Thio Tepa 10.0 mg/ml >43.9

Vincristine Sulfate 1.0 mg/ml >240 min

\*WARNING: Not recommended for use with Carmustine and Thiotepa. Please note that following drugs have extremely low permeation times:

Carmustine (BCNU): 22.6 minutes and Thiotepa: 43.9 minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
*PRAStaff@fda.hhs.gov*

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

**510(K) SUMMARY**  
**K211666**

Applicant: **CAREGLOVE GLOBAL SDN BHD**

Location Lot 17479, Lorong Senawang 2/3  
Off Jalan Senawang 3,  
Senawang Industrial Estate,  
70450 Seremban,  
Negeri Sembilan Darul Khusus,  
Malaysia.

Phone No. (60) 6 6782377 Fax No. (60) 6 6785377

Contact Person: Lim Kwee Shyan

Summary Preparation Date: 5<sup>th</sup> January, 2021

**Device Information**

Trade Name: POWDER FREE NITRILE EXAMINATION GLOVES, BLUE, CHEMOTEST

Common Name: POWDER FREE NITRILE EXAMINATION GLOVES

Classification Name: Patient Examination Gloves

Product Code: LZA, LZC

Regulation: 21 CFR 880.6250

**Predicate Device Identification**

510(k) Number: K162858,

Device Name: Careplus Powder Free Nitrile Examination Glove, Blue, Chemotest

**Device Description**

It is the powder-free variation of the class I latex patient examination gloves made by on-line polymer-coating and mild on-line chlorination process. The process modifies the surface characteristics and causes it to remain tack-free without the use of any dusting or donning powder.

**Intended Use of Device**

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

**Comparison of technological characteristics between the predicate and subject devices.**

The Powder Free Nitrile Examination Gloves – Blue, Chemotest, are summarized with the following technological characteristic compared to ASTM D6319 or equivalent standards.

Characteristic	Standard	Specification	Subject Device Powder Free Nitrile Examination Gloves, Blue, Chemotest) K211666	Predicate Device (Careplus Powder Free Nitrile Examination Glove, Blue, Chemotest) K162858	Remarks
Product Code	-	-	LZA, LZC	LZA, LZC	Same
Intended Use	-	-	Intended for medical purposes that is worn on the examiner's hand to prevent contamination	Intended for medical purposes that is worn on the examiner's hand to prevent	Same

			between patient and examiner	contamination between patient and examiner	
Design	-	-	Powder Free, Non-Sterile, Ambidextrous, Beaded Cuff	Powder Free, Non-Sterile, Ambidextrous, Beaded Cuff	Same
Indications for Use	-	-	A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner	A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner	Same
Construction	-	-	Ambidextrous, Polymer Coated or Chlorinated, Powder Free Nitrile	Ambidextrous, Polymer Coated or Chlorinated, Powder Free Nitrile	Same
Color Description	-	-	Blue	Blue	Similar
Material	-	-	Nitrile	Nitrile	Same
Single Use	-	-	Yes	Yes	Same
Packaging	-	-	Packed in Dispenser Boxes	Packed in Dispenser Boxes	Same
Chemo Drugs Claim	-	-	Chemo Claim	Chemo Claim	Same

Summary of non-clinical performance test result

Characteristic	Standard	Specification	Subject Device Powder Free Nitrile Examination Gloves, Blue, Chemotest) K211666	Predicate Device (Careplus Powder Free Nitrile Examination Glove, Blue, Chemotest) K162858	Remarks
<u>Dimension</u> Length (size: XSmall), mm Length (size: Small), mm Length (size: Medium), mm Length (size: Large), mm Length (size: XLarge), mm  Thickness (palm), mm Thickness (finger), mm  Width (size: XSmall), mm Width (size: Small), mm Width (size: Medium), mm Width (size: Large), mm Width (size: XLarge), mm	ASTM D6319-19	220 min 220 min 230 min 230 min 230 min  0.05 min 0.05 min  70 ± 10 80 ± 10 95 ± 10 110 ± 10 120 ± 10	Meet 220mm min Meet 220mm min Meet 230mm min Meet 230mm min Meet 230mm min  Meet 0.05mm min Meet 0.05mm min  Meet 70 ± 10 mm Meet 80 ± 10 mm Meet 95 ± 10 mm Meet 110 ± 10 mm Meet 120 ± 10 mm	Meet 220mm min Meet 220mm min Meet 230mm min Meet 230mm min Meet 230mm min  Meet 0.05mm min Meet 0.05mm min  Meet 70 ± 10 mm Meet 80 ± 10 mm Meet 95 ± 10 mm Meet 110 ± 10 mm Meet 120 ± 10 mm	Same
<u>Physical Properties</u> (Before Ageing) i) Tensile Strength (MPa) ii) Ultimate Elongation (%)	ASTM D6319-19	Min. 14 Min. 500	Meet 14MPa min. Meet 500% min	Meet 14MPa min. Meet 500% min	Similar



			- >240 min Dacarbazine (10 mg/ml) - >240 min Doxorubicin HCL (2.0 mg/ml) - >240 min Etoposide (20 mg/ml) - >240 min Flurouracil (50 mg/ml) - >240 min Ifosfamide (50.0 mg/ml) - >240 min Mitoxantrone (2 mg/ml) - >240 min Paclitaxel (6.0 mg/ml) - >240 min ThioTepa (10.0 mg/ml) - 43.9 Vincristine Sulfate (1.0 mg/ml) - >240 min	(Cytoxan) (20mg/ml) - >240 min Dacarbazine (10.0 mg/ml) - >240 min Doxorubicin HCL (2.0 mg/ml) - >240 min Etoposide (20.0 mg/ml) - >240 min Flurouracil (50.0 mg/ml) - >240 min Ifosfamide (50 mg/ml) - >240 min Mitoxantrone (2.0 mg/ml) - >240 min Paclitaxel (6.0 mg/ml) - >240 min ThioTepa (10.0 mg/ml) - >40.5 min Vincristine Sulfate (1.0 mg/ml) - >240 min	
--	--	--	---	--	--

Following is a table showing the actual measured parameters of the gloves (e.g. length, thickness, width, physical properties, etc.) as compare to ASTM. All data meets the standard reference requirement.

Test	Standard	Acceptance Criteria	Result																								
Freedom From Holes	ASTM D6319-19 ASTM D5151-19	Meet requirement inspection level G-I, AQL 2.5	Pass																								
Dimension	ASTM D6319-19	<table border="1"> <thead> <tr> <th>Size</th> <th>XSmall</th> <th>Small</th> <th>Medium</th> <th>Large</th> <th>XLarge</th> </tr> </thead> <tbody> <tr> <td>Length, min. mm</td> <td colspan="2">220</td> <td colspan="3">230</td> </tr> <tr> <td>Thickness, min. mm</td> <td colspan="5">0.05</td> </tr> <tr> <td>Width, ± 10 mm</td> <td>70</td> <td>80</td> <td>95</td> <td>110</td> <td>120</td> </tr> </tbody> </table>	Size	XSmall	Small	Medium	Large	XLarge	Length, min. mm	220		230			Thickness, min. mm	0.05					Width, ± 10 mm	70	80	95	110	120	Pass
Size	XSmall	Small	Medium	Large	XLarge																						
Length, min. mm	220		230																								
Thickness, min. mm	0.05																										
Width, ± 10 mm	70	80	95	110	120																						
Physical Properties	ASTM D6319-19	<table border="1"> <thead> <tr> <th colspan="2">Before Aging</th> <th colspan="2">After Accelerated Aging</th> </tr> <tr> <th>Tensile Strength</th> <th>Ultimate Elongation</th> <th>Tensile Strength</th> <th>Ultimate Elongation</th> </tr> </thead> <tbody> <tr> <td>14 MPa min.</td> <td>500 % min.</td> <td>14 MPa min.</td> <td>400 % min.</td> </tr> </tbody> </table>	Before Aging		After Accelerated Aging		Tensile Strength	Ultimate Elongation	Tensile Strength	Ultimate Elongation	14 MPa min.	500 % min.	14 MPa min.	400 % min.	Pass												
Before Aging		After Accelerated Aging																									
Tensile Strength	Ultimate Elongation	Tensile Strength	Ultimate Elongation																								
14 MPa min.	500 % min.	14 MPa min.	400 % min.																								
Residual Powder Content	ASTM D6319-19 ASTM D6124-06	Not more than 2 mg per glove	Pass																								

Summary of Clinical Testing

Not applicable

Conclusions:

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission K211666, the Powder Free Nitrile Examination Glove Blue Chemotest is a safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K162858.