



April 26, 2021

Ningbo Green Textile Co., Ltd.  
% Boyle Wang  
Official Correspondent  
Shanghai Truthful Information Technology Co., Ltd.  
Room 608, No. 738 Shangcheng Rd., Pudong  
Shanghai, Shanghai 200120  
China

Re: K202647  
Trade/Device Name: Surgical Mask  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical Apparel  
Regulatory Class: Class II  
Product Code: FXX,  
Dated: March 19, 2021  
Received: March 23, 2021

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD  
Acting Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K202647

Device Name  
Surgical Mask

### Indications for Use (Describe)

The Surgical Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# 510(k) Summary

## K202647

This 510(k) summary is being submitted in accordance with requirements of 21 CFR 807.92.

### **1.0 Submitter's Information**

Name: NINGBO GREEN TEXTILE CO.,LTD.  
Address:NO.498,3RD XINGCI ROAD,HANGZHOU BAY NEW ZONE,  
NINGBO,ZHEJIANG,CHINA 315336  
Tel: 86-574-63541998  
Fax: 86- 574-63528444  
Contact: Carol  
Date of Preparation: Mar.19,2021

### **Designated Submission Correspondent**

Mr. Boyle Wang  
Shanghai Truthful Information Technology Co., Ltd.  
Room 608, No. 738 Shangcheng Rd., Pudong Shanghai, 200120 China  
Tel: +86-21-50313932  
Email: [Info@truthful.com.cn](mailto:Info@truthful.com.cn)

### **2.0 Device Information**

Trade name: Surgical Mask  
Common name: Surgical Face Mask  
Classification name: Surgical Face Mask  
Model: Green2020-32

### **3.0 Classification**

Production code: FXX  
Regulation number: 21CFR 878.4040  
Classification: Class II  
Panel: Surgical Apparel

### **4.0 Predicate Device Information**

Manufacturer: WUHAN DYMEX HEALTHCARE CO., LTD.  
Device: SURGICAL FACE MASK

510(k) number: K182515

## **5.0 Device Description**

The Surgical Masks are single use, three-layer, flat-pleated style with ear loops and nose piece. The Surgical Masks are manufactured with three layers, the inner and outer layers are made of spun-bond polypropylene, and the middle layer is made of melt blown polypropylene filter. The ear loops are held in place over the users' mouth and nose by two elastic ear loops welded to the facemask. The elastic ear loops are not made with natural rubber latex. The nose piece in the layers of facemask is to allow the user to fit the facemask around their nose, which is made of Polyvinyl chloride coated iron wire. The surgical masks will be provided in blue. The surgical masks are sold non-sterile and are intended to be single use, disposable devices.

## **6.0 Indication for Use Statement**

The Surgical Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

## **7.0 Comparison to the Predicate Device**

**Table 1 General Comparison**

<b>Item</b>	<b>Subject Device K202647</b>	<b>Predicate Device K182515</b>	<b>Remark</b>
Product Name	Surgical Mask	Surgical Face Mask	--
Product Code	FXX	FXX	Same
Regulation No.	21 CFR 878.4040	21 CFR 878.4040	Same
Class	II	II	Same

Intended Use& Indications for use		The Surgical Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood	The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control	Same
		and body fluids. This is a single use, disposable device(s), provided non-sterile.	practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.	
Design features		Ear Loops, 3 layers	Ear Loops, 3 layers	Same
Mask Styles		Flat pleated	Flat pleated	Same
Material	Outer facing layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Middle layer	Melt blown polypropylene filter	Melt blown polypropylene filter	Same
	Inner Facing layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Nose piece	Polyvinyl chloride coated iron wire	Malleable polyethylene wire	Different
	Ear loops	not made with natural rubber latex	not made with natural rubber latex	Same
Color		Blue	Blue	Same
Dimension (Length)		17.5 cm +/- 1cm	17.5 cm +/- 1cm	Same
Dimension (Width)		9.5 cm +/- 1cm	9.5 cm +/- 1cm	Same
OTC use		Yes	Yes	Same
Sterility		Non-Sterile	Non-Sterile	Same
Single Use		Yes	Yes	Same
ASTM F2100 Level		Level 2	Level 2	Same

Biocompatibility			
Cytotoxicity	Comply with ISO 10993-5; Under the conditions of the study, not cytotoxicity effect.	Comply with ISO 10993-5	Same
Sensitization	Comply with ISO 10993-10; Under the conditions of the study, not an sensitization	Comply with ISO 10993-10	Same
Irritation	Comply with ISO 10993-10; Under the conditions of the study, not an irritant.	Comply with ISO 10993-10	Same

## 8.0 Non-Clinical Test Conclusion

The proposed device was tested and conformed to the following standards and the requirements stated in the Guidance for Industry and FDA Staff: Surgical Masks – Premarket Notification [510(k)] Submission issued on March 5, 2004.

**Table 2 - Performance Testing**

Test Method	Subject Device	Acceptance Criteria	Result
Fluid Resistance Performance ASTM F1862	31 Out of 32 pass at 120 mmHg(16.0 kPa)	29 Out of 32 pass at 120 mmHg (16.0 kPa)	Pass
Particulate Filtration Efficiency ASTM F2299	Average 99.6%	≥ 98%	Pass
Bacterial Filtration Efficiency ASTM F2101	Average 99.75%	≥ 98%	Pass

Differential Pressure (Delta P) MILM-36954C	Average 2.2 mmH <sub>2</sub> O/cm <sup>2</sup>	< 6.0 mmH <sub>2</sub> O/cm <sup>2</sup>	Pass
Flammability 16 CFR 1610	Class 1	Class 1	Pass

**Table 3 - Biocompatibility Testing**

Test Method	Subject Device	Result
Cytotoxicity	Under the conditions of the study, the device is noncytotoxic.	Pass
Irritation	Under the conditions of the study, the device is nonirritating.	Pass
Sensitization	Under the conditions of the study, the device is nonsensitizing	Pass

**9.0 Clinical Test Conclusion**

No clinical study is included in this submission.

**10.0 Conclusion**

The conclusion drawn from the non-clinical tests demonstrates that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device in K182515.