



March 3, 2020

Dentsply Sirona Inc.  
Karl Nittinger  
Vice President Corporate Regulatory Affairs  
221 West Philadelphia Street, Suite 60W  
York, Pennsylvania 17401

Re: K193064  
Trade/Device Name: Atlantis® suprastructures  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: Class II  
Product Code: NHA  
Dated: December 3, 2019  
Received: December 4, 2019

Dear Karl Nittinger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Srinivas Nandkumar, Ph.D.  
Director  
DHT1B: Division of Dental Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K193064

Device Name

Atlantis® suprastructures

Indications for Use (Describe)

Atlantis® suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function.

Atlantis® suprastructures are intended for attachment to a minimum of two (2) implants. Atlantis® suprastructures are compatible with following implants and abutments:

Manufacturer/	Implant system
BioHorizons	<p><b>Implant Level</b> Internal/ Tapered 3.5, 4.5, 5.7</p> <p><b>Abutment Level</b> Multi-unit abutment</p>
Biomet 3i	<p><b>Implant Level</b> Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 Certain™ XP 4/5, 5/6</p> <p><b>Abutment Level</b> Low Profile Abutment</p>
Camlog	<p><b>Implant Level</b> Screw-line Implant 3.3 Screw-line/Root-line Implant 3.8, 4.3 Screw-line/Root-line Implant 5.0, 6.0</p>
Dentsply Sirona	<p><b>Implant Level</b></p> <p><b>Astra Tech Implant System®</b> Astra Tech Implant EV 3.0 – Green Astra Tech Implant EV 3.6 – Purple Astra Tech Implant EV 4.2 – Yellow Astra Tech Implant EV 4.8 – Blue Astra Tech Implant EV 5.4 – Brown Astra Tech Implant EV Profile 4.2 Astra Tech Implant EV Profile 4.8 OsseoSpeed® TX 3.0 – Yellow OsseoSpeed® TX 3.5, 4.0 – Aqua OsseoSpeed® TX 4.5, 5.0 – Lilac OsseoSpeed® TX Profile 4.5/5.0</p> <p><b>Xive®</b> Xive® S 3.0 Xive® S 3.4, 3.8, 4.5, 5.5</p> <p><b>Abutment Level</b></p> <p><b>Astra Tech Implant System®</b> UniAbutment EV Angled Abutment EV Multibase Abutment EV 20° UniAbutment 45° UniAbutment Angled Abutment 20°</p> <p><b>Ankylos®</b> Ankylos® Balance Base Narrow D4.2 Ankylos® Balance Base D5.5</p>

	<b>Xive®</b> Xive® TG 3.4, 3.8, 4.5 Xive® MP 3.4, 3.8, 4.5 Xive® MP 5.5
Keystone	<b>Implant Level</b> PrimaConnex SD 3.3, 3.5 PrimaConnex RD 4.0, 4.1 PrimaConnex WD 5.0 Genesis 3.8, 4.5, 5.5, 6.5
MIS	<b>Implant Level</b> Internal Hex Narrow SEVEN – M4 Internal Hex Standard SEVEN – M4 Internal Hex Wide SEVEN – M4 Conical Connection Narrow C1 – V3 Conical Connection Standard C1 – V3 Conical Connection Wide C1 <b>Abutment Level</b> Multi-Unit Abutment
Nobel Biocare	<b>Implant Level</b> NobelActive NP 3.5 NobelActive RP 4.3, 5.0 NobelActive WP 5.5 NobelReplace NP 3.5 NobelReplace RP 4.3 NobelReplace WP 5.0 NobelReplace 6.0 Branemark System NP 3.3 Branemark System RP 3.75, 4.0 <b>Abutment Level</b> Multi-Unit Abutment RP
Straumann	<b>Implant Level</b> Bone Level 3.3 NC Bone Level 4.1, 4.8 RC Standard/Standard Plus 4.8 RN (Regular Neck) Standard/Standard Plus 4.8 WN (Wide Neck) <b>Abutment Level</b> RN Abutment Level WN Abutment Level Bone Level Multi-Base Abutment D3.5, 4.5 Bone Level Multi-Base Angled Abutment Screw-Retained Abutment 3.5 Screw-Retained Abutment 4.6
Zimmer Dental	<b>Implant Level</b> Tapered Screw-Vent™ 3.5, 4.5, 5.7 Screw Vent™ 3.3, 3.7, 4.5 <b>Abutment Level</b> Tapered Abutment

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

**510(k) SUMMARY**  
**For K193064**  
**Atlantis<sup>®</sup> suprastructures**

**1. Submitter Information:**

Dentsply Sirona Inc.  
221 West Philadelphia Street  
Suite 60W  
York, PA 17401

Contact Person: Karl Nittinger  
Telephone Number: 717-849-4424  
Fax Number: 717-849-4343  
Email: [karl.nittinger@dentsplysirona.com](mailto:karl.nittinger@dentsplysirona.com)

Date Prepared: 03 March 2020

**2. Device Name:**

- Proprietary Name: Atlantis<sup>®</sup> suprastructures
- Classification Name: Endosseous dental implant abutment
- CFR Number: 21 CFR 872.3630
- Device Class: Class II
- Product Code: NHA

**3. Predicate Device:**

Predicate Device for subject Atlantis suprastructures made by milling:

Predicate Device Name	510(k)	Company Name
Multibase Abutments EV and ATLANTIS <sup>™</sup> Suprastructures	K163350	Dentsply Sirona

Predicate Device for subject Atlantis suprastructures made by additive manufacturing:

Predicate Device Name	510(k)	Company Name
ATLANTIS <sup>™</sup> Suprastructures	K163398	Dentsply Sirona

Reference Devices:

<b>Reference Device Name</b>	<b>510(k)</b>	<b>Company Name</b>
ATLANTIST™ ISUS Implant Suprastructures	K160207	Dentsply Sirona
BioHorizons Tapered Internal Implant System	K071638	BioHorizons Implant Systems Inc.
BioHorizons Simple Solutions with Laser-Lok®	K100985	BioHorizons Implant Systems, Inc.
MIS Internal Hex Dental Implant System	K180282	MIS Implants Technologies Ltd. (Dentsply Sirona)
MIS C1 Narrow Platform Conical Connection Implant System, MIS C1 Wide Platform Conical Connection Abutments	K172505	MIS Implants Technologies Ltd. (Dentsply Sirona)
MIS V3 Conical Connection Dental Implant System	K163349	MIS Implants Technologies Ltd. (Dentsply Sirona)
Conical Connection Implants	K112162	MIS Implants Technologies Ltd. (Dentsply Sirona)
NobelActive Wide Platform (WP)	K133731	Nobel Biocare AB
Various Branemark System Dental Implant Products	K022562	Nobel Biocare AB

4. Description of Device:

Atlantis suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function. The subject Atlantis suprastructures include new compatible interfaces of the currently marketed Atlantis suprastructures made by milling (K163350) or additive manufacturing techniques (K163398) for the following abutments and implants:

- BioHorizons Multi-unit abutment
- MIS Internal Hex Narrow SEVEN – M4
- MIS Internal Hex Standard SEVEN – M4
- MIS Internal Hex Wide SEVEN – M4
- MIS Conical Connection Narrow C1 – V3
- MIS Conical Connection Standard C1 – V3
- MIS Conical Connection Wide C1
- MIS Multi-Unit Abutment
- Nobel Biocare NobelActive WP 5.5
- Nobel Biocare Branemark System NP 3.3
- Nobel Biocare Branemark System RP 3.75, 4.0

The design of the subject device is derived from patient dental models and completed by Dentsply Sirona technicians using computer-assisted design (CAD) according to the clinician’s prescription. The final CAD design of the Atlantis suprastructures are fabricated using additive manufacturing (AM) to produce a customized, patient-specific device. The Atlantis suprastructures subject of this bundled premarket notification are fabricated by milling or by additive manufacturing techniques.

The milled Atlantis suprastructures are composed of commercially pure titanium (CPTi) or cobalt chrome alloy and are available in following design types: Bar, Bridge, Hybrid and 2 in 1. Additional design limitations regarding minimum required segment cross-section, maximum span between implants and maximum cantilever extension, have been introduced for Atlantis suprastructures made by milling.

Atlantis suprastructures made by additive manufacturing are provided as Bridge and Hybrid types with optional mechanical retention (pin or cell retention) on the surface. The manufacturing of the Bridges and Hybrids of the Atlantis suprastructures by an additive manufacturing technique start from a titanium alloy and a cobalt-chrome alloy in powder form.

Milled variants of the subject Atlantis® suprastructures are offered in versions composed of unalloyed titanium conforming to ASTM F67 (*Standard Specification for Unalloyed Titanium for Surgical Implant Applications*) and cobalt chromium alloy (CoCr) conforming to ISO 22674 (*Dentistry – Metallic materials for fixed and removable restorations and appliances*). Variants of the subject Atlantis® suprastructures which are fabricated utilizing additive manufacturing are manufactured using titanium alloy powder or cobalt chromium alloy powder conforming to ISO 22674 (*Dentistry – Metallic materials for fixed and removable restorations and appliances*).

Fixation screws for use with the subject Atlantis® suprastructures are composed of titanium alloy confirming to ASTM F136 (*Standard Specification for Wrought Titanium-6 Aluminum-4 Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications*).

Screws are available for all compatible implants and abutments to screw the Atlantis suprastructures into the implant or onto the abutment.

Labeling is modified by providing a separate compatibility chart which lists all implants and abutments compatible with the Atlantis suprastructures.

##### 5. Indications for Use:

Atlantis® suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function.

Atlantis® suprastructures are intended for attachment to a minimum of two (2) implants. Atlantis® suprastructures are compatible with the implants and abutments listed in the Atlantis® Suprastructures Compatibility Chart:



<b>Manufacturer/</b>	<b>Implant system</b>
BioHorizons	<b>Implant Level</b>
	Internal/ Tapered 3.5, 4.5, 5.7
	<b>Abutment Level</b>
	Multi-unit abutment
Biomet 3i	<b>Implant Level</b>
	Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0
	Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5
	Certain™ XP 4/5, 5/6
	<b>Abutment Level</b>
	Low Profile Abutment
Camlog	<b>Implant Level</b>
	Screw-line Implant 3.3
	Screw-line/Root-line Implant 3.8, 4.3
	Screw-line/Root-line Implant 5.0, 6.0
Dentsply Sirona	<b>Implant Level</b>
	<b>Astra Tech Implant System®</b>
	Astra Tech Implant EV 3.0 – Green
	Astra Tech Implant EV 3.6 – Purple
	Astra Tech Implant EV 4.2 –
	Yellow Astra Tech Implant EV 4.8
	– Blue Astra Tech Implant EV 5.4 –
	Brown Astra Tech Implant EV
	Profile 4.2 Astra Tech Implant EV
	Profile 4.8 OsseoSpeed® TX 3.0 –
	Yellow OsseoSpeed® TX 3.5, 4.0 –
	Aqua
	OsseoSpeed® TX 4.5, 5.0 – Lilac
	OsseoSpeed® TX Profile 4.5/5.0
	<b>Xive®</b>
	Xive® S 3.0
	<b>Abutment Level</b>
<b>Astra Tech Implant System®</b>	
UniAbutment EV	
Angled Abutment EV	
Multibase Abutment EV	
20° UniAbutment	
45° UniAbutment	
Angled Abutment 20°	
<b>Ankylos®</b>	
Ankylos® Balance Base Narrow D4.2	
Ankylos® Balance Base D5.5	
	<b>Xive®</b>
Xive® TG 3.4, 3.8, 4.5	
Xive® MP 3.4, 3.8, 4.5	
Xive® MP 5.5	

Keystone	<b>Implant Level</b>
	PrimaConnex SD 3.3, 3.5 PrimaConnex RD 4.0, 4.1 PrimaConnex WD 5.0 Genesis 3.8, 4.5, 5.5, 6.5
	<b>Abutment Level</b>
	Multi-Unit Abutment
MIS	<b>Implant Level</b>
	Internal Hex Narrow SEVEN – M4 Internal Hex Standard SEVEN – M4 Internal Hex Wide SEVEN – M4 Conical Connection Narrow C1 – V3 Conical Connection Standard C1 – V3 Conical Connection Wide C1
	<b>Abutment Level</b>
	Multi-Unit Abutment
Nobel Biocare	<b>Implant Level</b>
	NobelActive NP 3.5 NobelActive RP 4.3, 5.0 NobelActive WP 5.5 NobelReplace NP 3.5 NobelReplace RP 4.3 NobelReplace WP 5.0 NobelReplace 6.0
	<b>Abutment Level</b>
	Multi-Unit Abutment RP
Straumann	<b>Implant Level</b>
	Bone Level 3.3 NC Bone Level 4.1, 4.8 RC Standard/Standard Plus 4.8 RN (Regular Neck) Standard/Standard Plus 4.8 WN (Wide Neck)
	<b>Abutment Level</b>
	RN Abutment Level WN Abutment Level Bone Level Multi-Base Abutment D3.5, 4.5 Bone Level Multi-Base Angled Abutment Screw-Retained Abutment 3.5 Screw-Retained Abutment 4.6
Zimmer Dental	<b>Implant Level</b>
	Tapered Screw-Vent™ 3.5, 4.5, 5.7 Screw Vent™ 3.3, 3.7, 4.5
	<b>Abutment Level</b>
	Tapered Abutment

6. Substantial Equivalence:

An overview of the similarities and differences between the subject and predicate devices is given in Table 5.1 and 5.2 below.

The indications for use of the subject and predicate Atlantis suprastructures (K163350, K163398) are identical with respect to the clinical indications. The indications for use of the subject and predicate devices only differ regarding the compatible implants and abutments which are added by this 510(k) premarket notification. Design types, material, manufacturing techniques (milling, additive manufacturing), screw types (standard screw, Angulated Screw Access (ASA) as well as mechanical retention methods on the suprastructure surface (cell, pin) of the subject Atlantis suprastructures are the same as cleared for the predicate Atlantis suprastructures made by milling

(K163350) or additive manufacturing (K163398). The additional design limitations for some design types of milled Atlantis suprastructures (Non-Standard bar, 2 in1 primary structure, Hybrid, Bridge) are within the design parameters cleared for the predicate device (K163350). The design parameters cleared for the predicate Atlantis suprastructures made by additive manufacturing (K163398) remain unchanged for the subject device.

#### 7. Non-Clinical Performance Data:

Following non-clinical test data and analyses are included to support substantial equivalence:

- Cross-sectional material analysis of the subject Atlantis suprastructure interfaces with the compatible implants and abutments, comparison to existing worst case interface geometry and reference to results of fatigue testing (K160207) conducted according to test methods of ISO 14801: *Dental-implants Dynamic Fatigue Test for Endosseous Dental Implants. ISO 14801*.
- Geometric measurement data and statistical compatibility analysis of OEM implant bodies, OEM abutments, and OEM fixation screws, to support the compatibility of the subject Atlantis suprastructure interfaces with the compatible implants and abutments (K071638, K100985, K133731, K022562).
- Reference to process validation testing (K163398) for Atlantis suprastructures made by additive manufacturing.
- Reference by equivalence to validation of moist heat sterilization parameters conducted according to ISO 17665-1 *Sterilization of health care products -- Moist heat -- Part 1: Requirements for the development, validation, and routine control of a sterilization process for medical devices* and ISO17665-2 *Sterilization of health care products - Moist heat - Part 2: Guidance on the application of ISO 17665-1*.
- Reference by equivalence to biocompatibility data of predicate devices (K163350, K163398).

#### 8. Clinical Performance Data:

No human clinical data was included in this premarket notification to support the substantial equivalence of the subject Atlantis suprastructures.

#### 9. Conclusion Regarding Substantial Equivalence:

The information included in this bundled 510(k) submission supports the substantial equivalence of the subject Atlantis suprastructures. The subject Atlantis suprastructures have the same intended use as the legally marketed predicate devices cleared under premarket notifications K163350 and K163398. The subject device also has similar indications for use and incorporates the same fundamental technology as the predicate devices (K163350, K163398). Results of technical analyses and compatibility studies to verify the performance of the subject Atlantis suprastructures support a determination of substantial equivalence.

**Table 5.1: Indications for Use for the subject and the predicate devices**

<b>Subject Device</b>	<b>Indications for Use</b>								
<p>Dentsply Sirona</p> <p>Atlantis® suprastructures</p>	<p>Atlantis® suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function. Atlantis® suprastructures are intended for attachment to a minimum of two (2) implants.</p> <p>The Atlantis® suprastructures are compatible with the implants and abutments listed in the Atlantis® Suprastructures Compatibility Chart:</p> <table border="1" data-bbox="375 426 1466 877"> <thead> <tr> <th data-bbox="375 426 565 478"><b>Manufacturer/ Company</b></th> <th data-bbox="565 426 1466 478"><b>Implant system</b></th> </tr> </thead> <tbody> <tr> <td data-bbox="375 478 565 594">BioHorizons</td> <td data-bbox="565 478 1466 594"> <p><b>Implant Level</b> Internal/ Tapered 3.5, 4.5, 5.7</p> <p><b>Abutment Level</b> Multi-unit abutment</p> </td> </tr> <tr> <td data-bbox="375 594 565 762">Biomet 3i</td> <td data-bbox="565 594 1466 762"> <p><b>Implant Level</b> Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 Certain™ XP 4/5, 5/6</p> <p><b>Abutment Level</b> Low Profile Abutment</p> </td> </tr> <tr> <td data-bbox="375 762 565 877">Camlog</td> <td data-bbox="565 762 1466 877"> <p><b>Implant Level</b> Screw-line Implant 3.3 Screw-line/Root-line Implant 3.8, 4.3 Screw-line/Root-line Implant 5.0, 6.0</p> </td> </tr> </tbody> </table>	<b>Manufacturer/ Company</b>	<b>Implant system</b>	BioHorizons	<p><b>Implant Level</b> Internal/ Tapered 3.5, 4.5, 5.7</p> <p><b>Abutment Level</b> Multi-unit abutment</p>	Biomet 3i	<p><b>Implant Level</b> Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 Certain™ XP 4/5, 5/6</p> <p><b>Abutment Level</b> Low Profile Abutment</p>	Camlog	<p><b>Implant Level</b> Screw-line Implant 3.3 Screw-line/Root-line Implant 3.8, 4.3 Screw-line/Root-line Implant 5.0, 6.0</p>
<b>Manufacturer/ Company</b>	<b>Implant system</b>								
BioHorizons	<p><b>Implant Level</b> Internal/ Tapered 3.5, 4.5, 5.7</p> <p><b>Abutment Level</b> Multi-unit abutment</p>								
Biomet 3i	<p><b>Implant Level</b> Certain™ 3.25, 4.0, 4/3, 5.0, 5/4, 6.0 Certain™ Prevail™ 3/4/3 (3.4), 4/5/4, 5/6/5, 4/3, 5/4, 6/5 Certain™ XP 4/5, 5/6</p> <p><b>Abutment Level</b> Low Profile Abutment</p>								
Camlog	<p><b>Implant Level</b> Screw-line Implant 3.3 Screw-line/Root-line Implant 3.8, 4.3 Screw-line/Root-line Implant 5.0, 6.0</p>								
	<p>Dentsply Sirona</p> <p><b>Implant Level</b></p> <p><b>Astra Tech Implant System®</b> Astra Tech Implant EV 3.0 – Green Astra Tech Implant EV 3.6 – Purple Astra Tech Implant EV 4.2 – Yellow Astra Tech Implant EV 4.8 – Blue Astra Tech Implant EV 5.4 – Brown Astra Tech Implant EV Profile 4.2 Astra Tech Implant EV Profile 4.8 OsseoSpeed® TX 3.0 – Yellow OsseoSpeed® TX 3.5, 4.0 – Aqua OsseoSpeed® TX 4.5, 5.0 – Lilac OsseoSpeed® TX Profile 4.5/5.0</p> <p><b>Xive®</b> Xive® S 3.0 Xive® S 3.4, 3.8, 4.5, 5.5</p> <p><b>Abutment Level</b></p> <p><b>Astra Tech Implant System®</b> UniAbutment EV Angled Abutment EV Multibase Abutment EV 20° UniAbutment 45° UniAbutment Angled Abutment 20°</p> <p><b>Ankylos®</b> Ankylos® Balance Base Narrow D4.2 Ankylos® Balance Base D5.5</p> <p><b>Xive®</b> Xive® TG 3.4, 3.8, 4.5 Xive® MP 3.4, 3.8, 4.5 Xive® MP 5.5</p>								

**Table 5.1: Indications for Use for the subject and the predicate devices (cont.)**

<b>Subject Device</b>	<b><u>Indications for Use</u></b>	
	Keystone Dental	<b>Implant Level</b> PrimaConnex SD 3.3, 3.5 PrimaConnex RD 4.0, 4.1 PrimaConnex WD 5.0 Genesis 3.8, 4.5, 5.5, 6.5
	MIS	<b>Implant Level</b> Internal Hex Narrow SEVEN – M4 Internal Hex Standard SEVEN – M4 Internal Hex Wide SEVEN – M4 Conical Connection Narrow C1 – V3 Conical Connection Standard C1 – V3 Conical Connection Wide C1 <b>Abutment Level</b> Multi-Unit Abutment
	Nobel Biocare	<b>Implant Level</b> NobelActive NP 3.5 NobelActive RP 4.3, 5.0 NobelActive WP 5.5 NobelReplace NP 3.5 NobelReplace RP 4.3 NobelReplace WP 5.0 NobelReplace 6.0 Branemark System NP 3.3 Branemark System RP 3.75, 4.0 <b>Abutment Level</b> Multi-Unit Abutment RP
	Straumann	<b>Implant Level</b> Bone Level 3.3 NC Bone Level 4.1, 4.8 RC Standard/Standard Plus 4.8 RN (Regular Neck) Standard/Standard Plus 4.8 WN (Wide Neck) <b>Abutment Level</b> RN Abutment Level WN Abutment Level Bone Level Multi-Base Abutment D3.5, 4.5 Bone Level Multi-Base Angled Abutment Screw-Retained Abutment 3.5 Screw-Retained Abutment 4.6
	Zimmer Dental	<b>Implant Level</b> Tapered Screw-Vent™ 3.5, 4.5, 5.7 Screw Vent™ 3.3, 3.7, 4.5 <b>Abutment Level</b> Tapered Abutment

Table 5.1: Indications for Use for the subject and the predicate devices (cont.)

Predicate Devices	Indications for Use		
Dentsply Sirona  Multibase Abutments EV and ATLANTIS™ Suprastructures  K163350	ATLANTIS™ Suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function. ATLANTIS™ Suprastructures are intended for attachment to a minimum of two (2) implants and are indicated for compatibility with the following implant and abutment systems:		
	<b>Implants:</b>		
	Manufacturer	Name of Implant	Size
	Biomet 3i	Certain	3.25, 4/3 – Prevail ¾/3, 4/3
		Certain	4.0, 5/4 – Prevail 4/5/4, 5/4
		Certain	5.0, XP 4/5 – Prevail 5/6/5, 6/5
		Certain	6.0, XP 5/6
	BioHorizons	Internal/Tapered	3.5, 4.5, 5.7
	Camlog	Screw-line Implant	3.3
		Screw-line / Root-line Implant	3.8, 4.3, 5.0, 6.0
	DENTSPLY Implants	XiVE	S 3.0, S 3.4, S 3.8, S 4.5, S 5.5
		OsseoSpeed™ TX	3.0, 3.5/4.0, 4.5/5.0
		Osseospeed™ Profile TX	4.5/5.0
		Osseospeed™ EV	3.0, 3.6, 4.2, 4.8, 5.4
	Keystone Dental	PrimaConnex	SD 3.3/3.5
		PrimaConnex	RD 4.0/4.1
		PrimaConnex	WD 5.0
		Genesis	3.8, 4.5, 5.5/6.5
	Nobel Biocare	NobelActive	NP 3.5 – RP 4.3, 5.0
		NobelReplace	NP-3.5 – RP 4.3 – WP 5.0 – 6.0
	Straumann	Bone Level	3.3 NC – 4.1, 4.8 NC
		Standard Plus	3.5 NN
		Standard / Standard Plus	4.8 RN – 4.8 WN
	Zimmer Dental	Tapered Screw Vent / Screw Vent*	Tapered S-V 3.5/S-V 3.3, 3.7 / Tapered S-V 4.5/ S-V 4.5
		Tapered Screw Vent	5.7
	<b>Abutments:</b>		
	Manufacturer	Name of Abutment	
	Biomet 3i	Low Profile Abutment	
	DENTSPLY	ATIS Uni Abutment EV	
		ATIS UniAbutment 20°, ATIS Un	
		ATIS Angled Abutment EV	
		ATIS Angled Abutment 20°	
		ANKYLOS Balance Base Narrow D4.2, Balance Base D5.5	
ATIS Multibase Abutment EV			
XiVE MP 3.4, MP 3.8, MP 4.5, MP 5.5			
Nobel Biocare	Multi-Unit Abutment RP		
Straumann	Bone Level Multi-Base Angled Abutment		
	Bone Level Multi-Base Abutment D3.5, D4.5		
	RN Abutment Level, WN Abutment Level		
	Screw-Retained Abutment 3.5, 4.6		
Zimmer Dental	Tapered Abutment		

Table 5.1: Indications for Use for the subject and the predicate devices (cont.)

Predicate Devices	Indications for Use
<p>Dentsply Sirona</p> <p>ATLANTIS™ Suprastructures</p> <p>K163398</p>	<p>ATLANTIS suprastructures are indicated for attachment to dental implants or abutments in the treatment of partially or totally edentulous jaws for the purpose of restoring chewing function.</p> <p>ATLANTIS suprastructures are intended for attachment to a minimum of two (2) implants and are indicated for compatibility with the following implant and abutment systems:</p> <p><b>Implants:</b></p> <ul style="list-style-type: none"> <li>• Biomet 3i Certain 3.25, 4/3 - Prevail 3/4/3, 4/3</li> <li>• Biomet 3i Certain 4.0, 5/4 - Prevail 4/5/4, 5/4</li> <li>• Biomet 3i Certain 5.0,XP4/5 - Prevail 5/6/5, 6/5</li> <li>• Biomet 3i Certain 6.0, XP 5/6</li> <li>• BioHorizons Internal/Tapered 3.5, 4.5, 5.7</li> <li>• Camlog Screw-line Implant 3.3</li> <li>• Camlog Screw-line / Root-line Implant 3.8, 4.3, 5.0, 6.0</li> <li>• DENTSPLY Implants XiVE S 3.0, S 3.4, S 3.8, S 4.5, S 5.5</li> <li>• DENTSPLY Implants OsseoSpeed TX 3.0, 3.5/4.0, 4.5/5.0</li> <li>• DENTSPLY Implants OsseoSpeed Profile TX 4.5, 5.0</li> <li>• DENTSPLY Implants OsseoSpeed EV 3.0, 3.6, 4.2, 4.8, 5.4</li> <li>• DENTSPLY Implants OsseoSpeed Profile EV 4.2, 4.8</li> <li>• Keystone Dental PrimaConnex SD 3.3, 3.5</li> <li>• Keystone Dental PrimaConnex RD 4.0, 4.1</li> <li>• Keystone Dental PrimaConnex WD 5.0</li> <li>• Keystone Dental Genesis 3.8, 4.5, 5.5/6.5</li> <li>• Nobel Biocare NobelActive NP 3.5 - RP 4.3, 5.0</li> <li>• Nobel Biocare NobelReplace NP 3.5 - RP 4.3 - WP 5.0 - 6.0</li> <li>• Straumann Bone Level 3.3 NC - 4.1, 4.8 RC</li> <li>• Straumann Standard Plus 3.5 NN</li> <li>• Straumann Standard/Standard Plus 4.8 RN - 4.8 WN</li> <li>• Zimmer Dental Tapered S-V 3.5/ S-V 3.3, 3.7 / Tapered S-V 4.5/ S-V 4.5</li> <li>• Zimmer Dental Tapered Screw-Vent 5.7</li> </ul> <p><b>Abutments:</b></p> <ul style="list-style-type: none"> <li>• Biomet 3i Low Profile Abutment</li> <li>• DENTSPLY Implants ATIS Uni Abutment EV</li> <li>• DENTSPLY Implants ATIS UniAbutment 20°, ATIS UniAbutment 45°</li> <li>• DENTSPLY Implants ATIS Angled Abutment EV, ATIS Angled Abutment 20°</li> <li>• DENTSPLY Implants Ankylos Balance Base Narrow D4.2, Balance Base D5.5</li> <li>• DENTSPLY Implants XiVE MP 3.4, MP 3.8, MP 4.5, MP 5.5</li> <li>• DENTSPLY Implants XiVE TG 3.4, TG 3.8, TG 4.5</li> <li>• Nobel Biocare Multi-Unit Abutment RP</li> <li>• Straumann Bone Level Multi-Base Angled Abutment</li> <li>• Straumann Bone Level Multi-Base Abutment D3.5, D4.5</li> <li>• Straumann RN Abutment Level, WN Abutment Level</li> <li>• Straumann Screw-Retained Abutment 3.5</li> <li>• Straumann Screw-Retained Abutment 4.6</li> <li>• Zimmer Dental Tapered Abutment</li> </ul>

Table 5.2: Similarities and Differences between the subject and the predicate devices

	Subject device	Predicate Devices	
	Dentsply Sirona Atlantis® suprastructures	Dentsply Sirona Multibase Abutments EV and ATLANTIS™ Suprastructures K163350	Dentsply Sirona ATLANTIS™ Suprastructures K163398
<b>Design</b>			
Restoration	multi-unit	multi-unit	multi-unit
Design type	Bar, Bridge, Hybrid, 2 in 1	Bar, Bridge, Hybrid, 2 in 1	Bridge, Hybrid
Prosthesis attachment	screw-retained, friction fit (2in1)	screw-retained, friction fit (2in1)	screw-retained, friction fit (2in1)
Connection suprastructure – implant/abutment	internal, external	internal, external	internal, external
Design Technology	CAD/CAM (Milling, Additive Manufacturing (AM))	CAD/CAM (Milling)	CAD/CAM (Additive Manufacturing) (AM)
<u>Design parameters Standard Bar Micro/ Standard Bar Macro</u>	20 mm/ 40 mm 6 mm/ 8 mm 2,5 mm <sup>2</sup> / 4.98 mm <sup>2</sup>	20 mm/ 40 mm 6 mm/ 8 mm 2,5 mm <sup>2</sup> / 4.98 mm <sup>2</sup>	N/A
<u>Design parameters Non-Standard bar, 2 in1 primary structure, Hybrid (milled), Bridge (milled):</u>	40 mm 15 mm 2.5 mm <sup>2</sup> / 4.98 mm <sup>2</sup> / 4.0 mm height, 3.0 mm width*	40 mm 15 mm 4.0 mm height, 3.0 mm width	N/A
<u>Design parameters Bridge and Hybrid (AM)</u>	40 mm 15 mm 4.0 mm height, 3.0 mm width	N/A	40 mm 15 mm 4.0 mm height, 3.0 mm width

\*2.5 mm<sup>2</sup>: if span between implants 0 to 20 mm or cantilever 0 to 6 mm), 4.98 mm<sup>2</sup>: if span between implants >20 to 40 mm or cantilever >6 to 8 mm), 4.0 mm height, 3.0 mm width: if cantilever >8 to 15 mm



Table 5.2: Similarities and Differences between the subject and the predicate devices (cont.)

	Subject device	Predicate Devices	
	Dentsply Sirona Atlantis® suprastructures	Dentsply Sirona Multibase Abutments and ATLANTIS™ Suprastructures K163350	Dentsply Sirona ATLANTIS™ Suprastructures K163398
Mechanical retention	no retention Cell or pin retention (optional for AM)	no retention	Cell or pin retention (optional)
Screws	Standard screws Angulated Screw Access (ASA) screws	Standard screws Angulated Screw Access (ASA) screws	Standard screws Angulated Screw Access (ASA) screws
<b><u>Material</u></b>			
Abutment/ suprastructure	Commercially pure titanium (block) Titanium alloy (powder) Cobalt-chrome alloy (block, powder)	Commercially pure titanium (block) Titanium alloy (powder) Cobalt-chrome alloy (block, powder)	Commercially pure titanium (block) Titanium alloy (powder) Cobalt-chrome alloy (block, powder)
Prosthetic screw	Titanium alloy block	Titanium alloy block	Titanium alloy block