



April 14, 2020

GA Health Company Limited
% Rafael Aguila
Responsible Third-Party Official
Accelerated Device Approval Services, LLC
6800 S.W. 40th Street, Ste. 444
Ludlum, FL 33155-3708

Re: K200388
Trade/Device Name: Vacutore Air/Water Bottle Tubing
Vacutore CO2 Source Tubing with Luer Connector
Regulation Number: 21 CFR 876.1500
Regulation Name: Endoscope and accessories
Regulatory Class: II
Product Code: OCX, FCX
Dated: April 1, 2020
Received: March 30, 2020

Dear Rafael Aguila:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Shanil P. Haugen, Ph.D.
Assistant Director
DHT3A: Division of Renal, Gastrointestinal,
Obesity and Transplant Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K200388

Device Name

Vacutore Air/Water Bottle Tubing, Vacutore CO2 Source Tubing with Luer Connector

Indications for Use (Describe)

The Vacutore Air/Water Bottle Tubing is to connect an air/CO2 source, a sterile water source (water bottle), and an endoscope to supply air/CO2 and water during gastrointestinal endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.

The Vacutore CO2 Source Tubing with Luer Connector is intended to be used with a carbon dioxide (CO2) source with the purpose of supplying CO2 to the endoscope during endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

The 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of the Safe Medical Devices Act of 1990 and 21 CFR 807.92.

1. Submission Sponsor

Submitter's Name: GA Health Company Limited
Submitter's Address: Unit 18, 21/F, Metropole Square
2 On Yiu Street ,Shatin, N.T,
Hong Kong, CHINA

2. Sponsor Contact

Contact Person: Cindy Ye
Chief Executive Officer
Telephone: +852 2833 9010
Email: cindy@andorate.com

3. Date Prepared

13th April 2020

4. Device Identification

Trade Device Name: Vacutore[®] Air/Water Bottle Tubing
Common Device Name: Air/Water Tubing for Endoscope
Classification Name: OCX - Endoscopic Irrigation/Suction System
Regulation Number: 876.1500
Classification: 2
Regulation Name: Endoscope and accessories.

Trade Device Name: Vacutore[®] CO₂ Source Tubing with Luer Connector
Common Device Name: CO₂ Source Tubing
Classification Name: FCX - Insufflator, Automatic Carbon-Dioxide For Endoscope
Regulation Number: 876.1500
Classification: 2
Regulation Name: Endoscope and accessories.

5. Predicate Device Identification

Predicate Device 510(k) No.: K093665
Predicate Device Trade Name: Endo SmartCap[™]
Predicate Device Product Code: FAJ

Predicate Device 510(k) No.: K172916
Predicate Device Trade Name: FUJIFILM Water Tank
Predicate Device Product Code: FDF, FDS

There are two predicate devices as mentioned above. The primary predicate device is Endo SmartCap™ (K093665). The secondary predicate device is FUJIFILM Water Tank (K172916).

6. Device Description:

The Vacutore® Air/Water Bottle Tubing and CO₂ Source Tubing with Luer Connector (Hereafter called CO₂ Source Tubing) are intended for 24-hour multi-patient use. Air/water bottle tubing is supplied in sterile and CO₂ Source Tubing is supplied in non-sterile. Table 1 shows the components included in the application.

Table 1 – Components included in the application

Components	Qty	Classification Name	Regulation Number	Classification
Vacutore® Air/Water Bottle Tubing (GAR088US)	1	OCX - Endoscopic Irrigation/Suction System	876.1500	2
Vacutore® Air/Water Bottle Tubing (GAR090US)	1	OCX - Endoscopic Irrigation/Suction System	876.1500	2
Vacutore® CO ₂ Source Tubing with Luer Connector (GAR047)	1	FCX - Insufflator, Automatic Carbon-Dioxide For Endoscope	876.1500	2

The Vacutore® Air/Water Bottle Tubing is manufactured for use in conjunction with sterile water bottle, and together with Fujifilm 500/600 and 700 series endoscopes. The air/water bottle tubing is individually packed in sealed package, sold as a sterile device. The air/water bottle tubing is designed to be attached to the air/water port of the endoscopes to provide irrigation through the air/water channel to the distal end of endoscope.

The Vacutore® CO₂ Source Tubing with Luer Connector is manufactured for use in conjunction with air/water bottle tubing and CO₂ insufflator. The CO₂ Source tubing is individually packed in sealed package, sold as a non-sterile device. The CO₂ Source Tubing is designed to be attached to the air/water bottle tubing and the outlet of the CO₂ insufflator to provide irrigation and CO₂ insufflation through the air/water channel to the distal end of endoscope.

According to the Medical Device Recalls database in FDA website, no product recall was found for the predicate device. According to the FDA MAUDE database, safety issues on fluid leakage and backflow were identified. In the performance test, the leakage test are conducted for verification.

There were no prior submissions for the Vacutore® Air/Water Bottle Tubing and Vacutore® CO₂ Source Tubing with Luer Connector.

7. Intended Use:

The Vacutore® Air/Water Bottle Tubing is to connect an air/CO₂ source, a sterile water source (water bottle), and an endoscope to supply air/CO₂ and water during gastrointestinal endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.

The Vacutore[®] CO₂ Source Tubing with Luer Connector is intended to be used with a carbon dioxide (CO₂) source with the purpose of supplying CO₂ to the endoscope during endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.

8. Technological Characteristics

Table 2, 3 and 4 summarize the Vacutore[®] Air/Water Bottle Tubing and Vacutore[®] CO₂ Source Tubing with Luer Connector technological characteristics as compared to the predicate device from Endo SmartCap[™] and FUJIFILM Water Tank.

Table 2 Summary of design, features and principles of operation between the Vacutore[®] Air/Water Bottle Tubing (GAR088US) technological characteristics as compared to the predicate devices.

Specification	Predicate Device	Proposed Device	Substantial Equivalence
Device name	Endo SmartCap [™]	Vacutore [®] Air/Water Bottle Tubing	N/A
K number	K093665	---	N/A
Manufacturer	Medivators, Inc.	GA Health Company Limited	N/A
Product code	FAJ	OCX	N/A
Classification	2	2	Identical
Regulation No	876.1500	876.1500	Identical
Regulation Name	Endoscope and accessories	Endoscope and accessories	Identical
Supplied Sterile	Yes	Yes	Identical
Compatibility	Fujinon [®] Endoscopes	Fujifilm 500/600 series endoscope	Substantial Equivalent
Indications for Use	ENDO SMARTCAP [™] Tubing is intended to be used with an air or CO ₂ source and/ or pump along with a sterile water source to supply air or CO ₂ and sterile water to an gastrointestinal endoscope during endoscopic procedures	The Vacutore [®] Air/Water Bottle Tubing is to connect an air/CO ₂ source, a sterile water source (water bottle), and an endoscope to supply air/CO ₂ and water during gastrointestinal endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device	Substantial Equivalent
Environment of Use	Hospital and or clinics	Hospital and or clinics	Identical
Material	Methyl methacrylate-acrylonitrile-butadiene-styrene copolymer, Polycarbonate, Polyethylene, Polyvinyl Chloride, Thermoplastic Elastomer	Polyvinyl Chloride, Silicone, Polyoxymethylene, Polycarbonate, Acrylonitrile Butadiene Styrene	Substantial Equivalent
Packaging	Each model packed separately in a seal pouched	Each model packed separately in a seal pouched	Identical

510(k) Premarket Notification: Traditional
Air/Water Bottle Tubing and CO₂ Source Tubing with Luer Connector
GA Health Company Limited
CONFIDENTIAL

Manufacturing method	Injection moulding and extrusion	Injection moulding and extrusion	Substantial Equivalent
Sterilization	EO gas	EO gas	Identical
Shelf Life	Three years	One year	Substantial Equivalent

Table 3 Summary of design, features and principles of operation between the Vacutore® Air/Water Bottle Tubing (GAR090US) technological characteristics as compared to the predicate devices.

Specification	Predicate Device	Proposed Device	Substantial Equivalence
Device name	FUJIFILM Water Tank	Vacutore® Air/Water Bottle Tubing	N/A
K number	K172916	---	N/A
Manufacturer	Fujifilm Corporation	GA Health Company Limited	N/A
Product code	FDS, FDF	OCX	N/A
Classification	2	2	Identical
Regulation No	876.1500	876.1500	Identical
Regulation Name	Endoscope and accessories	Endoscope and accessories	Identical
Supplied Sterile	No	Yes	Substantial Equivalent*
Compatibility	FUJIFILM 700 system scopes	Fujifilm 700 series endoscope	Substantial Equivalent
Indications for Use	The Water Tank Model WT-603 is intended for use in combination with FUJIFILM 700 system scopes to deliver air and water through the endoscope under the management of a physician in medical facilities. Do not use this product for any other purpose.	The Vacutore® Air/Water Bottle Tubing is to connect an air/CO ₂ source, a sterile water source (water bottle), and an endoscope to supply air/CO ₂ and water during gastrointestinal endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.	Substantial Equivalent
Environment of Use	Hospital and or clinics	Hospital and or clinics	Identical
Material	Silicone Rubber, Poly (4-methyl-1-pentene)	Polyvinyl Chloride, Silicone, Polyoxymethylene, Polycarbonate, Acrylonitrile Butadiene Styrene	Substantial Equivalent
Packaging	Each model packed in a PE bag with a carton box	Each model packed separately in a seal pouched	Substantial Equivalent
Manufacturing method	Injection moulding	Injection moulding and extrusion	Substantial Equivalent
Sterilization	Supplied non-sterile Sterilized in a daily basis	EO gas	Substantial Equivalent*
Shelf Life	One year	One year	Identical

*Remark: The predicate device is supplied as non-sterile and the end-user should be sterilized prior using, while the subject device is sold as sterile in which no sterilization is needed by the end-user prior using.

Table 4 Summary of design, features and principles of operation between the Vacutore® CO₂ Source Tubing with Luer Connector (GAR047) technological characteristics as compared to the predicate devices.

Specification	Predicate Device	Proposed Device	Substantial Equivalence
Device name	Endo SmartCap™ CO₂ Source Tubing	Vacutore® CO₂ Source Tubing with Luer Connector	N/A
K number	K093665	---	N/A
Manufacturer	Medivators, Inc.	GA Health Company Limited	N/A
Product code	FAJ	FCX	N/A
Classification	2	2	Identical
Regulation No	876.1500	876.1500	Identical
Regulation Name	Endoscope and accessories	Endoscope and accessories	Identical
Supplied Sterile	No	No	Identical
Indications for Use	ENDO SMARTCAP™ CO ₂ Source Tubing with Luer Input is intended to be used with a CO ₂ insufflator along with ENDO SMARTCAP™ Irrigation Tubing or ENDOGATOR™ Hybrid Irrigation Tubing to supply CO ₂ to a GI endoscope during GI endoscopic procedures	The Vacutore® CO ₂ Source Tubing with Luer Connector is intended to be used with a carbon dioxide (CO ₂) source with the purpose of supplying CO ₂ to the endoscope during endoscopic procedures. This device is intended to be used for not more than 24 hours. It is a 24-hour multi-patient use device.	Substantial Equivalent
Environment of Use	Hospital and or clinics	Hospital and or clinics	Identical
Material	Polycarbonate, Polyvinyl Chloride, high impact polystyrene, Polyethylene terephthalate, polytetrafluoroethylene	Polyvinyl Chloride, polypropylene, Polycarbonate, polytetrafluoroethylene,	Substantial Equivalent
Packaging	Each model packed separately in a PE bag	Each model packed separately in a seal pouched	Identical
Manufacturing method	Injection moulding	Injection moulding	Substantial Equivalent
Shelf Life	Three years	Three years	Substantial Equivalent

9. Non-Clinical Performance Data

The bench testing was performed to support substantial equivalence. The following testing were performed on GA Health samples from initial production lots, including sterilization.

9.1 Performance Test

9.1.1 Vacutore[®] Air/Water Bottle Tubing

- 9.1.1.1 Assembling Integrity
- 9.1.1.2 Endoscope Compatibility
- 9.1.1.3 Compatibility with Bottle
- 9.1.1.4 Two-Way Valve Integrity
- 9.1.1.5 Flow Clamp Test
- 9.1.1.6 Air Flow Test
- 9.1.1.7 Water Flow Test
- 9.1.1.8 24-hour use test

9.1.2 Vacutore[®] CO₂ Source Tubing with Luer Connector

- 9.1.2.1 Assembling Integrity
- 9.1.2.2 Compatibility with CO₂ Insufflator
- 9.1.2.3 Compatibility with air/water bottle tubing
- 9.1.2.4 CO₂ delivery test
- 9.1.2.5 Water delivery test
- 9.1.2.6 Air leak test
- 9.1.2.7 Water backflow test

9.2 Sterilization

Vacutore[®] Air/Water Bottle Tubing is sold in sterile package. The subject device has been sterilized in a validated EO sterilization cycle. The EO sterilization cycle has a Sterility Assurance Level (SAL) of 10⁻⁶. EO residuals on the components are below the maximum levels defined in ANSI/AAMI/ISO 10993-7:2008 Biological Evaluation of Medical Devices - Part 7: Ethylene Oxide sterilization residuals. The Vacutore[®] Air/Water Bottle Tubing, and the predicate devices, are not labeled as pyrogen-free because they do not have any blood or cerebrospinal fluid contact.

Vacutore[®] CO₂ Source Tubing with Luer Connector is supplied as non-sterile.

9.3 Shelf Life

The Vacutore[®] Air/Water Bottle Tubing has a one (1) year shelf life and the Vacutore[®] CO₂ Source Tubing with Luer Connector has a three (3) year shelf life, based on the design and existing sterile barrier data from existing packaging. The subject devices are packaged in a paper/film pouch respectively like other sterile products currently manufactured. These pouches have been tested by Sanitation Environment Technology Institute, Soochow University, which is an CNAS accredited laboratory. The tests conducted including accelerated aging, seal strength, dye penetration, microbial barrier properties, vacuum leak test and the sterility test. The performance test was conducted after the accelerated aging process. The test result can imply that the air/water bottle tubing can provide and maintain a sterile barrier and its intended performance for at least one (1) year. It also implies that the CO₂ source tubing maintains its intended performance for at least three (3) years.

9.4 Biocompatibility

The biocompatibility of the Vacutore[®] Air/Water Bottle Tubing and Vacutore[®] CO₂ Source Tubing with Luer Connector were conducted in accordance with the FDA guideline "Use of International Standard ISO 10993-1, "Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process"". It included the following tests.

- Acute Systemic Toxicity Test
- In Vitro Cytotoxicity
- Skin sensitization
- Irritation

The air/water bottle tubing and CO₂ source tubing are classified as surface device with mucosal membrane contact for a limited duration (not more than 24 hours). The test result shows that the air/water bottle tubing is biocompatible.

9.5 Simulated Use Testing

Simulated use testing of the subject devices was conducted. The test results can imply that the subject devices provide sufficient backflow prevention with the intended devices and verified 24-hour multi-patient use.

10. Clinical Testing

Similar devices have been on the market for many years with proven safety and efficacy for the use of the device. These devices have no direct patient contact. Based on this history and the use of the device, clinical testing was not necessary to support substantial equivalence data. The non-clinical testing performed supports safety and efficacy of the devices and provides data to show substantial equivalence to the predicate device.

11. Conclusion

Vacutore[®] Air/Water Bottle Tubing and Vacutore[®] CO₂ Source Tubing with Luer Connector has the same intended use as the predicate devices.

Based on the technological characteristics and overall performance of the devices in bench testing, It is believed that no significant differences exist between the subject devices and the predicate devices.

The Vacutore[®] Air/Water Bottle Tubing and Vacutore[®] CO₂ Source Tubing with Luer Connector do not raise any new issues of safety and effectiveness.

From a clinical perspective and comparing design specifications, the Vacutore[®] Air/Water Bottle Tubing and Vacutore[®] CO₂ Source Tubing with Luer Connector, and the predicate device are substantially equivalent.