



August 13, 2021

Weony (Shenzhen) Technology Co., Ltd.
Lucy Yan
Consultant
3rd Floor B, Building 19, HeYi BeiFangYongFa
Science & Technology Park, HeYi Community, ShaJing Street,
Shenzhen, GuangDong 518104
China

Re: K210671

Trade/Device Name: Digital Blood Pressure Monitor WBP Series: WBP101, WBP102, WBP103,
WBP104, WBP105, WBP106, WBP107

Regulation Number: 21 CFR 870.1130

Regulation Name: Noninvasive Blood Pressure Measurement System

Regulatory Class: Class II

Product Code: DXN

Dated: July 13, 2021

Received: July 23, 2021

Dear Lucy Yan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

LCDR Stephen Browning
Assistant Director
Division of Cardiac Electrophysiology, Diagnostics
and Monitoring Devices
Office of Cardiovascular Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K210671

Device Name

Digital Blood Pressure Monitor WBP Series: WBP101, WBP102, WBP103, WBP104, WBP105, WBP106, WBP107

Indications for Use (Describe)

The device intended to measure the diastolic, systolic blood pressures and pulse rate for adult population in home and hospital facilities by using a non-invasive oscillometric technique with a single upper arm cuff (22-36 cm).

The device detects the appearance of irregular heart beats during measurement and gives a warning signal with readings.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

This summary of 510(k) information is submitted as required by requirements of SMDA and 21 CFR §807.92.

1 Administrative Information

Submission Date	Jan. 8, 2021
Manufacturer information	<p>Submitter's Name: WEONY (SHENZHEN) TECHNOLOGY CO., LTD.</p> <p>Address: Address: 3rd Floor B, Building 19, HeYi BeiFangYongFa Science & Technology Park, HeYi Community, ShaJing Street, BaoAn District, ShenZhen, 518104, P.R. China</p> <p>Contact person: Autumn Liu TEL: 86-755-86057437 E-Mail: autumn.liu@weony-sz.com</p>
Submission Correspondent	<p>Contact person: Ms Lucy.Yan E-Mail: Lucy.yan@aivikon.com 3rd Floor B, Building 19, HeYi BeiFangYongFa Science & Technology Park, HeYi Community, ShaJing Street, BaoAn District, ShenZhen, 518104, P.R. China</p>
Establishment registration number	NA

2 Device Information

Common name of the device	System, Measurement, Blood-Pressure, Non-Invasive
Trade name of the device	Digital Blood Pressure Monitor WBP Series
Type/Model of the device	WBP101, WBP102, WBP103, WBP104, WBP105, WBP106, WBP107
Classification information	<p>Classification panel: Cardiovascular</p> <p>Classification name: System, Measurement, Blood-Pressure, Non-Invasive</p> <p>Regulation Number: 870.1130</p> <p>Device Class: II</p> <p>Product Code: DXN</p>
type of submission	510(k) Traditional

3 Predicate Device Information

Sponsor:	Truly Instrument Limited
Device:	Truly Automatic Arm Blood Pressure Monitor
510(K) Number:	K091434

4 Device Descriptions

Weony Digital Blood Pressure Monitor WBP Series are designed to measure the systolic and diastolic blood pressure and pulse rate of an individual (at least 18 or above) by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm. Our method to define systolic and diastolic pressure is similar to the auscultatory method but uses an electronic pressure sensor rather than a stethoscope and mercury manometer. The sensor converts tiny alterations in cuff pressure to electrical signals, by analyzing those signals to define the systolic and diastolic blood pressure and calculating pulse rate, which is a well-known technique in the market called the "oscillometric method".

The main components of the Blood Pressure Monitor are the main unit and cuff unit. ABS is used to outer housing of the main unit. The preformed cuff unit, which is applicable to arm circumference approximately between 220 and 360 mm, includes the inflatable bladder and nylon shell. All models of the arm blood pressure monitor use a single size of cuff. The device consists of the microprocessor, the pressure sensor, the operation keys, the pump, the electromagnetic deflation control valve, and the LCD. The subject devices are powered by four AA alkaline batteries or adapter.

The device has irregular heart beat (IHB) indicator which compares the longest and the shortest time intervals of detected pulse waves to mean time interval and displays a warning signal with the reading to indicate the detection of irregular pulse rhythm when the difference of the time intervals is over a specified range.

5 Intended Use/ Indications for Use

The subject device intended to measure the diastolic, systolic blood pressures and pulse rate of an adult individual in hospitals and home environments by using a non-invasive oscillometric technique with a single upper arm cuff (22-36 cm).

The device detects the appearance of irregular heart beats during measurement and gives a warning signal with readings.

6 SE Comparison

Table 1. Substantial Equivalence Comparison

Characteristics	Subject device	Predicate device (K091434)	Remark
Device Name	Digital Blood Pressure Monitor	Truly Automatic Arm Blood Pressure Monitor DB series	NA
Device Model	WBP101, WBP102, WBP103, WBP104, WBP105, WBP106, WBP107.	DB21, DB22, DB23, DB31, DB32, DB61M, DB62M, DB63M, DB71M	NA
Manufacturer	WEONY (SHENZHEN) TECHNOLOGY CO., LTD.	Truly Instrument Limited	NA
Intended Use/ Indication for Use	The subject device intended to measure the diastolic, systolic blood pressures and pulse rate of an adult individual in hospitals and home environments by using a non-invasive oscillometric technique with a single upper arm cuff (22-36 cm). The device detects the appearance of irregular heart beats during measurement and gives a warning signal with readings.	Truly Automatic Arm Blood Pressure Monitor DB series, Models DB21, DB22, DB23, DB31, DB32, DB61M, DB62M, DB63M, DB71M are a series device intended to measure the systolic and diastolic blood pressure and pulse rate of an adult individual, by using a non-invasive technique in which an inflatable cuff is wrapped around the upper arm. The devices features include irregular pulse rhythm detection during measurement, and display a warning signal with the reading once the irregular heartbeat is detected.	SE
Intended Population	adult	adult	same
Intended Anatomical site	upper arm	upper arm	same
Prescription & OTC	OTC	OTC	same
Working Principle	Oscillometric method	Oscillometric method	same
Pressurization Source	Automatic internal pump	Automatic internal pump	same
Internal Power supply	4pcs "AA" alkaline Batteries	4- size "AA" alkaline Batteries	same
Memory Function	2 × 90 memories (SYS, DIA, Pulse)	DB21: 2×60; DB22: 2×50; DB23: 4×99; DB31: 2×60; DB32: 1×99; DB61M: 4×99; DB62M: 4×99; DB63M: 4×99; DB71M: 4×99	SE
Cuff Size	220mm~360mm	220mm~340mm	Similar Note01
Measuring range	Pressure: 0 to 299 mmHg (in 1 mmHg increment);	Pressure: (20mmHg~280mmHg)	Similar Note02
	Pulse: 40 to 180 beat/minute	Pulse rate range (40-195) beats/minute	

Measuring resolution	1 mmHg	1 mmHg	same
Accuracy	Pressure: ± 3 mmHg; Pulse: $\pm 5\%$	Pressure: ± 3 mmHg; Pulse $\pm 5\%$.	same
Display Type	LCD digital display	LCD digital display	same
Cuff attachment method	By plastic host connected to monitor	By plastic host connected to monitor	same
Irregular Heartbeat Detection	The subject devices have the IHB function.	DB22, DB23, DB61M, DB62M, DB63M, DB71M have the IHB feature.	same
Operating Environment	10~40°C,	10~40°C,	same
	15%~90%RH	15%~90%RH	
Materials	Patient contact materials of the cuff have been tested in accordance with ISO 10993 tested in accordance with ISO 10993 and FDA guidance	Patient contact materials of the cuff have been tested in accordance with ISO 10993 tested in accordance with ISO 10993 and FDA guidance	Same

Note01: The subject devices have the larger arm circumference than predicate device, but the subject devices have been tested by ISO81060-2.

Note02: The subject device has a different measuring range of pressure and pulse from the predicate device, but the subject devices have been validated all the full claimed range.

The subject device is as same as predicate device in Working Principle, Intended patient population, intended application site, measuring accuracy. Only their Cuff size, measuring ranges are a little bit different (refer to Note01 to Note 02) which had been validated. Als, the differences would not raise any safety or effectiveness issue based on tests in this submission.

Thus, the subject device is Substantially Equivalent (SE) to the predicate device which is legally marketed in US.

7 Brief discussions of the non-clinical tests

The subject device conforms to the following guidances and standards:

- ✧ Non-Invasive Blood Pressure (NIBP) Monitor Guidance
- ✧ IEC 60601-1:2005+A1:2012: Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance;
- ✧ IEC 60601-1-2:2014 Medical Electrical Equipment - Part 1-2: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirements and Tests.
- ✧ IEC 60601-1-11: 2010 Medical Electrical Equipment - Part 1-11: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Requirements for Medical Electrical Equipment and Medical Electrical Systems

Used in The Home Healthcare Environment;

- ✧ ISO 10993-5: 2009 /(R)2014 Biological evaluation of medical devices – Part 5: Tests for In Vitro cytotoxicity;
- ✧ ISO 10993-10: 2010 Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization;
- ✧ IEC 80601-2-30: 2013 Medical electrical equipment - Part 2-30: Particular Requirements for the Basic Safety and Essential Performance of Automated Non-invasive Sphygmomanometers

8 Brief discussions of clinical tests

- ✧ ISO 81060-2:2018 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type;

In this clinical investigation, 85 patients (47 males and 38 females) participated in the clinical study. Same arm sequential method was adopted during the clinical testing. The manual Mercury Sphygmomanometer was used as a reference device. All the subjects were volunteer to take part in the clinical study, all the subjects completed the clinical study without any AE or side-effect.

The results showed the accuracy of the blood pressure monitor is within acceptable scope specified in ISO 81060-2.

9 Other information (such as required by FDA guidance)

No other information.

11 Conclusions

The subject device:

Digital Blood Pressure monitor manufactured by WEONY (SHENZHEN) TECHNOLOGY CO., LTD. is respectively substantially equivalent to the predicate device Arm Blood Pressure Monitor manufactured by Truly Instrument Limited (K091434).