



March 9, 2022

Guangdong Xingcan Brothers Medical Technology Co.,Ltd  
% Cassie Lee  
Manager  
Share Info (Guangzhou) Medical Consultant Ltd.  
No. 1919-1920, Building D3, Minjie Plaza, Shuixi Road  
Huangpu District  
Guangzhou, Guangdong  
China

Re: K213851

Trade/Device Name: Disposable Medical Nitrile gloves  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LZA  
Dated: December 6, 2021  
Received: December 10, 2021

Dear Cassie Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray, III, Ph.D.  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K213851

Device Name  
Disposable Medical Nitrile gloves

### Indications for Use (Describe)

The Disposable Medical Nitrile gloves intended to be worn on the hands of examiners to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## **K213851 510(k) Summary**

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirement of 21 CFR 807.92.

### **1. Date of the summary prepared: December 6, 2021**

### **2. Submitter's Information**

Sponsor Name: Guangdong Xingcan Brothers Medical Technology Co.,Ltd

Address: Room B10, 4th Floor, No. 137 (plant A1), Pacific Industrial Zone, Xintang Town, Zengcheng District, Guangzhou

Establishment Registration Number: Applying

Post Code: 511300

Contact name: Suhai Yin

Phone: +86 139 2426 4846

Fax: +86 755-8654 6919

E-mail: 1273722763@qq.com

### **Application Correspondent:**

Contact Person: Ms. Cassie Lee

Company: Share Info (Guangzhou) Medical Consultant Ltd.

Address: No. 1919-1920, Building D3, Minjie Plaza, Shuixi Road, Huangpu District, Guangzhou, China

Tel: +86 20 8200 6973

Email: [regulatory@share-info.com](mailto:regulatory@share-info.com)

### **3. Subject Device Information**

Type of 510(k): Traditional

Common Name: Polymer Patient Examination Glove

Classification Name: Non-powdered patient examination glove

Trade Name: Disposable Medical Nitrile gloves

Specifications: Small (S), Medium (M), Large (L), X-large (XL)

Review Panel: General Hospital

Product Code: LZA

Regulation Number: 21 CFR 880.6250

Regulatory Class: Class I

### **4. Predicate Device Information**

Sponsor: Guang Dong Kingfa SCI. & TECH.CO., LTD.  
 Common Name: Polymer Patient Examination Glove  
 Classification Name: Non-Powdered Patient Examination Glove  
 Trade Name: Patient Examination Gloves  
 510(k) Number: K203593  
 Review Panel: General Hospital  
 Product Code: LZA  
 Regulation Number: 21 CFR 880.6250  
 Regulatory Class: Class I

**5. Device Description**

The proposed devices are powder-free nitrile examination gloves, provided as non-sterile and disposable devices. The proposed devices are mainly made from nitrile and there are four sizes, including small (S), medium (M), large (L), X-large (XL) for optional. The gloves are composed of Butadiene-Acrylonitrile Copolymers (CAS No.9003-18-3) and blue colorant Pigment Blue (CAS No.147-14-8). The examination glove is a smooth surface and a rolled rim at the cuff edge.

The Medical nitrile examination gloves meet the specifications in ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application.

**6. Intended Use / Indications for Use**

The Disposable Medical Nitrile gloves intended to be worn on the hands of examiners to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.

**7. Technological Characteristics Comparison**

Elements of Comparison	Subject Device	Predicate Device	Result
Company	Guangdong Xingcan Brothers Medical Technology Co.,Ltd	Guang Dong Kingfa SCI. & TECH.CO., LTD.	--
510 (k) Number	K213851	K203593	--
Trade Name	Disposable Medical Nitrile gloves	Patient Examination Gloves	--
Product Code	LZA	LZA	Same
Classification Name	Non-Powdered Patient Examination Glove	Non-Powdered Patient Examination Glove	Same
Classification	Class I	Class I	Same
Regulation No.	21 CFR 880.6250	21 CFR 880.6250	Same

Elements of Comparison	Subject Device	Predicate Device	Result
Indications For Use	The Disposable Medical Nitrile gloves intended to be worn on the hands of examiners to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.	The nitrile examination glove is intended to be worn on the hands of examiners to prevent contamination between patient and examiner. This is a single-use, powder-free, non-sterile device.	Same
Material of Use	Nitrile	Nitrile	Same
Color	Blue	Blue	Same
Size (ASTM D6319-19)	Small, Medium, Large, X Large	Small, Medium, Large, X Large	Same
Sterilization	Non-sterile	Non-sterile	Same
Usage	Single usage	Single usage	Same
Dimensions (ASTM D6319-19)	<b>Length:</b> For S: $\geq 220$ mm For M/L/XL: $\geq 230$ mm	Length: S (220mm min) M (230mm min) L (230mm min) XL (230mm min)	Same
	<b>Width:</b> For S: $80 \pm 10$ mm For M: $95 \pm 10$ mm For L: $110 \pm 10$ mm For XL: $120 \pm 10$ mm	Palm width Small ( $80 \pm 10$ mm) Medium ( $95 \pm 10$ mm) Large ( $110 \pm 10$ mm) X large ( $120 \pm 10$ mm)	
Physical Properties (ASTM D6319-19)	Meets requirements of the ASTM D6319-19 <b>Before Aging:</b> Tensile Strength: $\geq 14$ Mpa Ultimate Elongation: $\geq 500\%$	Meets requirements of the ASTM D6319-19 <b>Before Aging:</b> Tensile Strength: Min 14 Mpa Ultimate Elongation: Min 500%	Same
	Meets requirements of the ASTM D6319-19 <b>After Aging:</b>	Meets requirements of the ASTM D6319-19 <b>After Aging:</b>	

<b>Elements of Comparison</b>	<b>Subject Device</b>	<b>Predicate Device</b>	<b>Result</b>
	Tensile Strength: $\geq 14$ Mpa Ultimate Elongation: $\geq 400\%$	Tensile Strength: Min 14Mpa Ultimate Elongation: Min 400%	
Thickness (ASTM D6319-19)	Palm min. 0.05 mm Finger min. 0.05 mm	Palm min. 0.05 mm Finger min. 0.05 mm	Same
Powder Free (ASTM D6319-19)	$\leq 2$ mg/glove	$\leq 2$ mg/glove	Same
Freedom from Holes (Water Tight -1000 ml)- ASTM D6319-19 (Cross Reference D5151)	Meets requirements of the ASTM D6319-19	Meets requirements of the ASTM D6319-19	Same
Biocompatibility – Skin Sensitization (ISO 10993-10:2010)	Under the conditions of the study not a sensitizer	Under the conditions of the study not a sensitizer	Same
Biocompatibility – Skin Irritation (ISO 10993-10:2010)	Under the conditions of study not an irritant	Under the conditions of study not an irritant	Same
Biocompatibility – Acute Systemic Toxicity (ISO 10993- 11: 2017)	Cytotoxicity is assessed via rationale. Under the condition of acute systemic toxicity test, the test article did not show acute systemic toxicity in vivo.	Cytotoxicity is assessed via rationale. Under the condition of acute systemic toxicity test, the test article did not show acute systemic toxicity in vivo.	Same

## 8. Test Summary

### 2.1 Summary of Non-Clinical Performance Testing

#### 1) Performance Testing Summary:

Test Method	Test Purpose	Acceptance Criteria	Test Results	Conclusion
ASMT D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application – Physical Dimensions Test	To determine the width, length, and thickness of the gloves	<b>Width:</b> For S: 80±10mm For M: 95±10mm For L: 110±10mm For XL: 120±10mm  <b>Length:</b> For S: ≥220mm For M: ≥230mm For L: ≥230mm For XL: ≥230mm	<b>For Lot 20210825B:</b> <b>Width:</b> For S: 83~86mm For M: 91~92mm For L: 104-106mm For XL: 112-114mm  <b>Length:</b> For S: 234~236mm For M: 235~237mm For L: 237~239mm For XL: 237~239mm  <b>For Lot 20210827A:</b> <b>Width:</b> For S: 83~85mm For M: 91~93mm For L: 104-106mm For XL: 111-113mm  <b>Length:</b> For S: 234~236mm For M: 235~237mm For L: 237~239mm For XL: 237~239mm  <b>For Lot 20210830A:</b> <b>Width:</b> For S: 83~85mm For M: 90~92mm For L: 104-106mm For XL: 112-115mm  <b>Length:</b> For S: 234~236mm	Passed



			<p>For M: 235~237mm  For L: 237~240mm  For XL: 237~239mm</p>	
		<p><b>Thickness:</b>  Finger: <math>\geq 0.05\text{mm}</math>  Palm: <math>\geq 0.05\text{mm}</math></p>	<p><b>For Lot 20210825B:</b>  <b>Finger:</b>  For S/M: Pass at 0.11mm  For L/XL: Pass at 0.12mm  <b>Palm:</b>  For S/M/L/XL: Pass at 0.08mm</p> <p><b>For Lot 20210827A:</b>  <b>Finger:</b>  For S: Pass at 0.11mm  For M/L/XL: Pass at 0.12mm  <b>Palm:</b>  For S/M/L/XL: Pass at 0.08mm</p> <p><b>For Lot 20210830A:</b>  For S/M/L: Pass at 0.11mm  For XL: Pass at 0.12mm  <b>Palm:</b>  For S/M/L/XL: Pass at 0.08mm</p>	

ASMT D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application – Physical Dimensions Test	To determine the tensile strength and ultimate elongation before and after acceleration aging	<b>Before Aging:</b> Tensile Strength: ≥14Mpa Ultimate Elongation: ≥500% <b>After Aging:</b> Tensile Strength: ≥14Mpa Ultimate Elongation: ≥400%	<b>For all three lots:</b> <b>Before Aging:</b> Tensile Strength: ≥14Mpa Ultimate Elongation: ≥500% <b>After Aging:</b> Tensile Strength: ≥14Mpa Ultimate Elongation: ≥400%	Passed
ASTM D6319-19 (ASTM D5151-11) Standard Test Method for Detection of Holes in Medical Gloves	To determine the holes in the gloves	AQL 2.5	<b>For all three lots:</b> Pass at AQL 2.5	Passed
ASMT D6319-19 (ASTM D6124-11) Standard Test Method for Residual Powder on Medical Gloves	To determine the residual powder in the gloves	< 2.0 mg/glove	<b>For Lot 20210825B:</b> For S/M/L/XL: Pass at 0.22 mg/glove  <b>For Lot 20210827A:</b> For S/M/L/XL: Pass at 0.19 mg/glove  <b>For Lot 20210830A:</b> For S/M/L/XL: Pass at 0.20 mg/glove	Passed

## 2) Biocompatibility Testing Summary:

Test Method	Test Purpose	Acceptance Criteria	Test Results	Conclusion
ISO 10993-10 Biological evaluation of medical devices – Part 10: Tests for	To evaluate the potential intracutaneous reactivity caused by	Under the conditions of study not an irritation	Under the conditions of study not an irritation	Passed

skin irritation and skin sensitization	intracutaneously inject the extract to rabbits			
ISO 10993-10 Biological evaluation of medical devices – Part 10: Tests for skin irritation and skin sensitization	To determine the skin sensitization potential in guinea pigs.	Under the conditions of the study not a sensitization	Under the conditions of the study not a sensitization	Passed
ISO 10993-11:2017 Biological evaluation of medical devices – Part 11: Tests for acute systemic toxicity	The test item was evaluated for acute systemic toxicity in ICR mouse	Under the conditions of the study no systemic toxicity	Under the condition of acute systemic toxicity test, the test article did not show acute systemic toxicity in vivo.	Passed

**2.1 Summary of Clinical Performance Test**

No clinical study is included in this submission.

**9. Final Conclusion:**

The subject device is a safe, as effective, and perform as well as or better than the legally marketed predicated K203191.