



July 27, 2022

Shandong Maida Medical Technology Co.,Ltd.  
% Boyle Wang  
General Manager  
Shanghai Truthful Information Technology Co., Ltd.  
Room608,No.738,Shangcheng Rd.,Pudong  
Shanghai, Shanghai 200120  
China

Re: K221271

Trade/Device Name: Disposable Nitrile Powder-Free Examination Gloves (Tested for Use with  
Chemotherapy Drugs)

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA, LZC

Dated: April 18, 2022

Received: May 2, 2022

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmnmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, Ph.D.  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K221271

Device Name

Disposable Nitrile Powder-Free Examination Gloves (Tested for Use with Chemotherapy Drugs)

Indications for Use (Describe)

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs.

Chemotherapy Drug	Concentration	Breakthrough Detection Time in Minutes
Carmustine (BCNU)	3.3 mg/ml (3,300 ppm)	22.8(23.7, 22.8, 23.1) Minutes
Cisplatin	1.0 mg/ml(1,000 ppm)	> 240 Minutes
Cyclophosphamide (Cytosan)	20.0 mg/ml(20,000 ppm)	> 240 Minutes
Dacarbazine	10 mg/ml (10,000 ppm)	> 240 Minutes
Doxorubicin	2.0 mg/ml(2,000 ppm)	> 240 Minutes
Etoposide	20.0 mg/ml(20,000 ppm)	> 240 Minutes
Fluorouracil	50.0 mg/ml(50,000 ppm)	>240 Minutes
Paclitaxel	6.0 mg/ml(6,000 ppm)	>240 Minutes
Thio Tepa	10.0 mg/ml(10,000 ppm)	46.8(48.2, 48.6, 46.8) Minutes

Please note that the following drugs have low permeation times: Carmustine (BCNU) 3.3 mg/ml 22.8 Minutes (min.); Thio Tepa 10.0 mg/ml 46.8 Minutes (min.). Warning: Please do not use with Carmustine (BCNU) and Thiotepea.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# **K221271 510(k) Summary**

This summary of 510(k) is being submitted in accordance with 21 CFR 807.92.

## **1.0 Submitter's Information**

Name: Shandong Maida Medical Technology Co.,Ltd.

Address: Room 102, Eastern building, No.166, South 1st Road, Development zone, Dongying, Shandong, China

Contact: Room 102, Eastern building, No.166, South 1st Road, Development zone, Dongying, Shandong, China

Date of Preparation: 2022.04.25

## **Designated Submission Correspondent**

Mr. Boyle Wang

Shanghai Truthful Information Technology Co., Ltd.

Room 1801, No. 161 Lujiazui East Rd., Pudong Shanghai, 200120 China

Tel: +86-21-50313932

Email: Info@truthful.com.cn

## **2.0 Device Information**

Trade name: Disposable Nitrile Powder-Free Examination Gloves (Tested for Use with Chemotherapy Drugs)

Common name: Patient Examination Gloves

Classification name: Non-powdered patient examination glove

Model(s): S, M, L, XL

## **3.0 Classification**

Production code: LZA, LZC

Regulation number: 21CFR880.6250

Classification: Class I

Panel: General Hospital

## **4.0 Predicate Device Information**

Manufacturer: Ever Growth (Vietnam) Co., Ltd.

Device: Disposable Powder Free Nitrile Examination Glove, Tested For Use With Chemotherapy Drugs, Disposable Powder Free Nitrile Examination Glove, Tested For Use With Chemotherapy Drugs, Orange Color

510(k) number: K190860

## **5.0 Device Description**

The subject device is single use, disposable gloves intended for medical

purposes to be worn on the examiner's hands to prevent contamination between patient and examiner. The gloves are powder-free, ambidextrous with beaded cuff, blue colored, nitrile, and tested for use with chemotherapy drugs. The gloves are offered in four sizes: extra-small, small, medium, large, and extra-large.

## 6.0 Indication for Use

A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 *Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs*

Chemotherapy Drug	Concentration	Breakthrough Detection Time in Minutes (minutes)
Carmustine (BCNU)	3.3 mg/ml (3,300 ppm)	22.8(23.7, 22.8, 23.1)
Cisplatin	1.0 mg/ml(1,000 ppm)	> 240
Cyclophosphamide (Cytoxan)	20.0 mg/ml(20,000 ppm)	> 240
Dacarbazine	10 mg/ml (10,000 ppm)	> 240
Doxorubicin	2.0 mg/ml(2,000 ppm)	> 240
Etoposide	20.0 mg/ml(20,000 ppm)	> 240
Fluorouracil	50.0 mg/ml(50,000 ppm)	> 240
Paclitaxel	6.0 mg/ml(6,000 ppm)	> 240
Thio Tepa	10.0 mg/ml(10,000 ppm)	46.8(48.2, 48.6, 46.8)

Please note that the following drugs have low permeation times:

Carmustine (BCNU) 3.3 mg/ml 22.8 Minutes (min.);

Thio Tepa 10.0 mg/ml 46.8 Minutes (min.).

Warning: Please do not use with Carmustine (BCNU) and Thiotepa.

## 7.0 Technological Characteristic Comparison Table

**Table1-General Comparison**

Item	Subject Device Pending	Predicate Device (K190860)	Comparison
Product Code	LZA,LZC	LZA,LZC	Same
Regulation No.	21CFR880.6250	21CFR880.6250	Same
Class	I	I	Same

Intended Use	A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs.	A patient examination gloves is a disposable device intended for medical purpose that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. In addition, these gloves were tested for use with chemotherapy drugs in accordance with ASTM D6978-05 Standard Practice for Assessment of Medical gloves to Permeation by Chemotherapy Drugs.	Same
Powdered or Powered free	Powdered free	Powdered free	Same
Design Feature	Ambidextrous	Ambidextrous	Same
Sterility	Non-Sterile	Non-Sterile	Same
Labeling Information	Single-use indication, powder free, device color, device name, glove size and quantity, Non-Sterile, a statement of standard ASTM D6978-05 compliance and a summary of the testing results.	Single-use indication, powder free, device color, device name, glove size and quantity, Non-Sterile, a statement of standard ASTM D6978-05 compliance and a summary of the testing results.	Similar

**Table2 Device Dimensions Comparison**

Predicate Device(K190860)	Designation	Size					Tolerance
		XS	S	M	L	XL	
	Length, mm	230	230	230	230	230	min
	Width, mm	70	80	95	110	120	±10
	Thickness, mm:						
	Finger	0.05					min
	Palm	0.05					min
Subject Device	Designation	Size					Tolerance
		S	M	L	XL		
	Length, mm	220	230	230	230	min	
	Width, mm	80	95	110	120	±10	

	Thickness, mm:		
	Finger	0.05	min
	Palm	0.05	min
Remark	Different		

Analysis: The physical dimensions are different with that of the predicate, but they all meet the requirements of ASTM D6319-19.

**Table3 Performance Comparison**

Item		Subject device	Predicate device (K190860)	Comparison	
Colorant		Blue	White, Orange	Different 1	
Physical Properties	Before Aging	Tensile Strength	14MPa, min	14MPa, min	Same
		Ultimate Elongation	500% min	500% min	Same
	After Aging	Tensile Strength	14MPa, min	14MPa, min	Same
		Ultimate Elongation	400%min	400%min	Same
	Comply with ASTM D6319		Comply with ASTM D6319		Same
Freedom from Holes	Be free from holes when tested in accordance with ASTMD5151 AQL=2.5		Be free from holes when tested in accordance with ASTMD5151 AQL=2.5	Same	
Powder Content	0.15-0.19 mg per glove, Meet the requirements of ASTM D6124		Meet the requirements of ASTM D6124	Same	
Chemotherapy Drugs Tested with Minimum Breakthrough Detection Time as Tested per ASTM D 6978	Carmustine (BCNU) 3.3 mg/ml: 22.8 Minutes (min.)		Carmustine (BCNU) 3.3 mg/ml: White:11.8 Minutes; Orange:31.6Minutes	Similar	
	Cisplatin 1.0 mg/ml: > 240 Minutes		Cisplatin 1.0 mg/ml: > 240 Minutes	Same	
	Cyclophosphamide (Cytoxan) 20.0 mg/ml: > 240 Minutes		Cyclophosphamide (Cytoxan) 20.0 mg/ml: > 240 Minutes	Same	
	Dacarbazine 10 mg/ml:> 240 Minutes		Dacarbazine (DTIC) 10.0 mg/ml:	Same	
			>240 Minutes		

Doxorubicin 2.0 mg/ml: > 240 Minutes	Doxorubicin Hydrochloride 2.0 mg/ml: >240 Minutes	Similar
Etoposide 20.0 mg/ml: > 240 Minutes	Etoposide (Toposar) 20.0 mg/ml: >240 Minutes	Similar
Fluorouracil 50.0 mg/ml: > 240 Minutes	Fluorouracil 50.0 mg/ml: >240 Minutes	Same
Paclitaxel 6.0 mg/ml: >240 Minutes	Paclitaxel (Taxol) 6.0 mg/ml: >240 Minutes	Same
Thio Tapa 10.0 mg/ml: 46.8Minutes (min.)	Thio-Tepa 10.0 mg/ml: White:16.9 Minutes; Orange: 72.5 Minutes	Similar

**Analysis:**

Different 1: The color of the subject device is different of that of the predicate. Biocompatibility testing was successfully completed for the subject device.

**Table4 Safety Comparison**

Item		Proposed device	Predicated device	Comparison
Material		Nitrile	Nitrile	SAME
Biocompatibility	Irritation	Under the conditions of the study, not an irritant	Comply with ISO10993-10	SAME
	Sensitization	Under conditions of the study, not a sensitizer.		
	Cytotoxicity	Under the conditions of the study, the device is potentially cytotoxic	Comply with ISO10993-5	Analysis3
	Systemic toxicity	Under the conditions of the study, the device does not elicit a systemic toxicity response in the model animal.	Complies with ISO 10993-11 Third edition 2017-09	
Label and Labeling		Meet FDA's Requirement	Meet FDA's Requirement	SAME

Analysis3: The proposed device is potentially cytotoxic, but all proposed devices are conducted the systemic toxicity test.



## **8.0 Summary of Non-Clinical Testing**

### **Biocompatibility Testing**

The biocompatibility evaluation for Disposable Nitrile Powder-Free Examination Gloves (Tested for Use with Chemotherapy Drugs) was conducted in accordance with the following standards:

ISO 10993-10:2010 *Biological Evaluation of Medical Devices - Part 10: Tests for Irritation And Skin Sensitization.*

ISO 10993-5:2009 *Biological Evaluation of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity*

ISO 10993-11 Third edition 2017-09 *Biological evaluation of medical devices - Part 11: Tests for systemic toxicity*

### **Performance Testing (Bench)**

Physical performance qualities of the proposed device were evaluated per ASTM D6319-10, *Standard Specification for Disposable Nitrile Powder-Free Examination Gloves for Medical Application.*

Permeation testing was conducted to support the addition of the labeling claim: *Tested for use with chemotherapy drugs.* In addition, the proposed device was tested according to ASTM D6978-05 (Reapproved 2019), *Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs*, in which minimum breakthrough times were determined for a wide range of chemotherapy drugs.

In summary, the performance testing of the subject device was conducted to adequately demonstrate the effectiveness of the device in accordance with the relevant test methods cited below:

- ASTM D6124-06 (Reapproved 2017), *Standard Test Method for Residual Powder on Medical Gloves*
- ASTM D5151-19, *Standard Test Method for Detection of Holes in Medical Gloves.*
- ASTM D6319-19, *Standard Specification for Disposable Nitrile Powder-Free Examination Gloves for Medical Application.*
- ASTM D6978-05 (Reapproved 2019) ,*Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.*

No.	Name of the Test Methodology / Standard	Purpose	Acceptance Criteria	Results
1	ISO 10993-10:2010 Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.	This part of ISO 10993 assesses possible contact hazards from chemicals released from medical devices, which may produce skin and mucosal irritation, eye irritation or skin sensitization.	Skin Sensitization Test: provided grades less than 1, otherwise sensitization.	All grades are 0.  All animals were survived and no abnormal signs were observed during the study.
2			Skin Irritation Test: If the primary irritation index is 0-0,4, the response category is Negligible. 0,5-1,9 means slight 2-4,9 means moderate 5-8 means severe	The primary irritation index is 0.  The response of the proposed device was categorized as negligible under the test condition
3	ISO 10993-5:2009 Biological Evaluation Of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity	This part of ISO 10993 describes test methods to assess the in vitro cytotoxicity of medical devices.	The viab.% of the 100% extract of the test article is the final result, and if viability is reduced to <70% of the blank, it has cytotoxic potential.	Viab.% of 100% test article extract is 17.1%  It means the proposed device have potential toxicity to L-929 in the MTT method
4	ISO 10993-11: 2017 Biological evaluation of medical devices — Part 11: Tests for systemic toxicity	To evaluate the potential for medical device materials to cause adverse systemic reactions.	Within the monitoring period (72 h), if the toxicosis response of testing group is not greater than that of control group, the testing sample is regarded as acceptable.	There was no evidence of systemic toxicity from the extract.
5	ASTM D6124-06 (Reapproved 2017), Standard Test Method for Residual Powder on Medical Gloves	This standard is designed to determine the amount of residual powder (or filter-retained mass) found on medical gloves	powder residue limit of 2.0 mg	0.15-0.19 mg /glove

6	ASTM D5151-06(Reapproved 2015), Standard Test Method for Detection of Holes in Medical Gloves.	This test method covers the detection of holes in medical gloves.	Samples number: 125 gloves AQL: 2.5 (ISO 2859) Criterion $\leq 7$ gloves for water leakage	no glove water leakage found
7	ASTM D6319-10(Reapproved 2015), Standard Specification For Nitrile Examination Gloves For Medical Application.	This specification covers certain requirements for nitrile rubber gloves used in conducting medical examinations and diagnostic and therapeutic procedures.	<p>Sterility: no need</p> <p>Freedom from holes: pl. Refer to No. 5 in table 5</p> <p>Dimensions: S: width <math>80 \pm 10</math> mm Length <math>\geq 220</math> mm M: width <math>95 \pm 10</math> mm Length <math>\geq 230</math> mm L: width <math>110 \pm 10</math> mm Length <math>\geq 230</math> mm XL: width <math>120 \pm 10</math> mm Length <math>\geq 230</math> mm</p> <p>Thickness: Finger <math>\geq 0.05</math> mm Palm <math>\geq 0.05</math> mm</p> <p>Physical properties: Before aging Tensile strength <math>\geq 14</math> MPa Ultimate Elongation <math>\geq 500\%</math> After Accelerated Aging Tensile strength <math>\geq 14</math> MPa Ultimate Elongation <math>\geq 400\%</math></p> <p>Powder-free Residue: pl. Refer to No. 4 in table 5</p>	<p>N.A.</p> <p>Please refer to No. 5 in table 5</p> <p><b>Lot no.:</b>210515</p> <p>Dimensions: S: width: 85-87 mm Length 247-253 mm M: width 88-96 mm Length 242-257 mm L: width 90-99 mm Length 240-254 mm XL: width 110-115 mm Length 245-253 mm</p> <p>Thickness: Finger 0.09-0.21 mm Palm 0.06-0.15 mm</p> <p>Physical properties: Before aging Tensile strength 14.1-22.5 MPa Ultimate Elongation 503.274% - 670.613% After Accelerated Aging Tensile strength 14.6-20.8 MPa Ultimate Elongation 411.403% - 592.683%</p> <p>Powder-free Residue: pl. Refer to No. 4 in table 5</p> <p><b>Lot no.:</b>210518</p> <p>Dimensions: S: width: 84-86 mm Length 248-256 mm M: width 95-96 mm Length 237-266 mm L: width 105-108 mm</p>

				<p>Length 257-262 mm  XL: width 114-117 mm  Length 252-262 mm  Thickness:  Finger 0.10-0.12mm  Palm 0.07-0.08mm</p> <p>Physical properties:  Before aging  Tensile strength 15.2-29.8 MPa  Ultimate Elongation 500.492% - 593.853%  After Accelerated Aging  Tensile strength 14.1-23.8MPa  Ultimate Elongation 451.751% - 597.368%</p> <p>Powder-free Residue:  pl. Refer to No. 4 in table 5</p> <p><b>Lot no.:</b>210520  Dimensions:  S: width: 84-87 mm  Length 244-257 mm  M: width 93-98 mm  Length 245-260 mm  L: width 104-110mm  Length 250-263 mm  XL: width 114-119 mm  Length 252-260 mm  Thickness:  Finger 0.10-0.12 mm  Palm 0.06-0.08 mm</p> <p>Physical properties:  Before aging  Tensile strength 14.4-23.9MPa  Ultimate Elongation 501.484% - 547.660%  After Accelerated Aging  Tensile strength 14.2-23.9 MPa  Ultimate Elongation 492.901% -</p>
--	--	--	--	--

				599.996%  Powder-free Residue: pl. Refer to No. 4 in table 5
8	ASTM D 6978	Chemotherapy Drugs Tested with Minimum Breakthrough Detection Time		Carmustine (BCNU) 3.3 mg/ml: 22.8 Minutes (min.)  Cisplatin 1.0 mg/ml: > 240 Minutes  Cyclophosphamide (Cytosan) 20.0 mg/ml: > 240 Minutes  Dacarbazine 10 mg/ml: > 240 Minutes  Doxorubicin 2.0 mg/ml: > 240 Minutes  Etoposide 20.0 mg/ml: > 240 Minutes  Fluorouracil 50.0 mg/ml: >240 Minutes  Paclitaxel 6.0 mg/ml: >240 Minutes  Thio Tega 10.0 mg/ml: 46.8Minutes (min.)

**9.0 Summary of Clinical Testing**

Clinical testing is not needed for this device.

## **10.0 Conclusion**

The conclusions drawn from the nonclinical tests demonstrate that the subject device, Disposable Nitrile Powder-Free Examination Gloves (Tested for Use with Chemotherapy Drugs) is as safe, as effective, and performs as well as or better than the legally marketed predicated device under K190860.