

216*150*130



95% PARTICULATE CUP RESPIRATOR P2
AND MEDICAL MASK

PERFORMANCE TESTED TO :
AS NZS 1716:2012 P2
BS EN 149:2001 FFP2 NR



95% PARTICULATE CUP RESPIRATOR P2
AND MEDICAL MASK

PERFORMANCE TESTED TO :
AS NZS 1716:2012 P2
BS EN 149:2001 FFP2 NR



95% PARTICULATE CUP RESPIRATOR P2
AND MEDICAL MASK

PERFORMANCE TESTED TO :
AS NZS 1716:2012 P2
BS EN 149:2001 FFP2 NR

INSTRUCTIONS FOR USE:

- 1) Wash your hands thoroughly before and after putting on a respirator.
- 2) To work as intended, a mask needs to fit your face properly. To ensure the correct fit:
 - a. Adjust the mask to cover your mouth, nose and chin.
 - b. Position the elastic head loops over head to secure the mask firmly in place.
 - c. The mask is designed to mould to the shape of the bridge of your nose and face to allow you to breathe normally.
- 3) Avoid touching the mask frequently as this may reduce its level of protection. If you must touch your face then you should wash your hands first and after doing so.
- 4) When removing the mask, unhook the Head loops. **DO NOT REMOVE THE MASK BY TOUCHING THE OUTER FRONT PANEL**, as it may be contaminated.
- 5) After removing, dispose of masks in accordance with local policies and procedures.
- 6) Disposable particulate cup respirators are for single use only, and **SHOULD NOT** be used for more than four hours under any circumstances. If the respirator is damaged or soiled, replace it immediately.



WARNINGS:
DO NOT touch the inner surface of the respirator. DO NOT use if the packaging is damaged. DO NOT use in highly toxic atmosphere. DO NOT use when oxygen is less than 19.5% and in explosive atmospheres. If breathing becomes difficult or dizziness occurs, return to an area where you can safely remove the respirator. Facial hair may reduce the effectiveness of the respirator, as it may not form an adequate seal on the face. Failure to follow all instructions, could reduce the effectiveness of this respirator and may result in sickness or death.

STORAGE CONDITIONS:

Store in clean, dry conditions between +5°C to +40°C and at < 60% relative humidity.
Use original packaging during storage or transport. Do not use if packaging is damaged.

Single use only



MFD 12/2020

EXP 12/2022

20 PCS

LOT MH2020-SRC46-12-01

CODE SM-RC201/20

ARTG NUMBER: 347410

This product is NOT made with natural rubber latex.



LA-009 V1

SINGLE USE
20 PCS
HEADLOOP



SINGLE USE
20 PCS
HEADLOOP



SINGLE USE
20 PCS
HEADLOOP

