



June 11, 2021

Zhuhai Gaoge Medical Technology Co., Ltd.  
% Olivia Meng  
Regulatory Affairs Manager  
Guangzhou Osmunda Medical Device Technical Service Co., Ltd.  
8-9th Floor, R&D Building, No. 26 Qinglan Street,  
Panyu District  
Guangzhou, Guangdong 510006  
China

Re: K203732

Trade/Device Name: Surgical Face Mask  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical apparel  
Regulatory Class: Class II  
Product Code: FXX  
Dated: April 25, 2021  
Received: April 30, 2021

Dear Olivia Meng:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD  
Acting Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K203732

Device Name  
Surgical Face Mask

### Indications for Use (Describe)

The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary K203732

In accordance with 21 CFR 807.92 the following summary of information is provided:

### 1. SUBMITTER

Zhuhai Gaoge Medical Technology Co., Ltd.

3rd floor, flat 5, No. 6 JinYuan one street, TangJiaWan town, Zhuhai National Hi-tech industrial development zone, Zhuhai city, Guangdong province, China

Phone: +86-0756-8803221

Fax: +86-0756-8307281

Primary Contact Person:	Olivia Meng Regulatory Affairs Manager Guangzhou Osmunda Medical Device Technical Service Co., Ltd. Tel: (+86)-20-6231 6262 Fax: (+86)-20-8633 0253
-------------------------	---

Secondary Contact Person:	Xiangliang Cheng Quality Manager Zhuhai Gaoge Medical Technology Co., Ltd. Phone: +86-0756-8803221 Fax: +86-0756-8307281
---------------------------	--

Date prepared	Apr 25 <sup>th</sup> , 2021
---------------	-----------------------------

### 2. DEVICE

Device Name:	Surgical Face Mask
--------------	--------------------

Common name:	Mask, Surgical
--------------	----------------

Model:	GK103
--------	-------

Regulation number	21 CFR 878.4040
-------------------	-----------------

Regulation Class:	II
-------------------	----

Product Code:	FXX
---------------	-----

### 3. PREDICATE DEVICE

K153496, Disposable Surgical Face Mask

This predicate has not been subject to a design-related recall.

#### 4. DEVICE DESCRIPTION

The Surgical Face Mask is designed and manufactured by Zhuhai Gaoge Medical Technology Co., Ltd. It is non-sterile and for single use.

The Surgical Face Mask has one model, GK103, that is the ear-loop style. It is made of three-layer nonwovens, ear loops and nose piece. Inner and outer layers are made of spun-bond polypropylene, and the middle layer is made of melt blown polypropylene filter. The ear loops are held in place over the users' mouth and nose by two elastic ear loops welded to the mask. The ear loops are made of polyester. The nose piece in the layers of mask is to allow the user to fit the mask around their nose, which is made of iron wire wrapped with plastic.

It is a self-inhalation filter mask, which works by filtering the air containing harmful substances through the filter material of the mask before being inhaled or exhaled.

The product is level 2 according to ASTM F2100-19. The main parameters of the product are listed as followed:

- Bacterial filtration efficiency (BFE)  $\geq 98\%$
- Sub-micron particle filtration efficiency  $\geq 98\%$
- Different pressure:  $< 6.0 \text{ mm H}_2\text{O}/\text{cm}^2$
- Flammability: class 1
- Resistance to penetration by synthetic blood: 120 mmHg

#### 5. INDICATIONS FOR USE

The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

6. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

Item	Proposed device	Predicate device	Comparison result
Manufacturer	Zhuhai Gaoge Medical Technology Co., Ltd.	Xiantao Rayxin Medical Products Co., Ltd.	NA
510K Number	K203732	K153496	NA
Product Common Name	Surgical Face Mask	Disposable Surgical Face Mask	NA
Intended Use	The Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.	The Disposable Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.	Same
Mask style	Flat pleated	Flat pleated	Same
Design feature	Ear loop, 3 layers	Ear Loops, Tie-On, 3 layers	Similar
Material of outer facing layer	Spun-bond polypropylene	Spun-bond polypropylene	Same
Material of middle layer	Melt blown polypropylene filter	Melt blown polypropylene filter	Same
Material of inner facing layer	Spun-bond Polypropylene	Spun-bond Polypropylene	Same
Nose piece	Galvanized iron wire wrapped with polyethylene(PE)	Malleable aluminum wire	Similar
Attachment	Ear loops: Polyester	Ear loops: Polyester	Similar
		Tieon: Spun-bond Polypropylene	
Color	Blue	Blue	Same
Dimension (Length × Width)	17.5 cm × 9.5 cm	17.5 cm × 9.5 cm	Same
OTC use	Yes	Yes	Same
Sterility	Non-sterile	Non-sterile	Same
Single use	Yes	Yes	Same
ASTM F 2100 level	Level 2	Level 2	Same
Fluid Resistance	For 3 non-consecutive lots,	32 out of 32 pass	Same

Performance ASTM F1862	32 out of 32 pass at 120 mmHg	at 120 mmHg	
Particulate Filtration Efficiency ASTM F2299	For 3 non-consecutive lots, 32 out of 32 pass, average at 98.76%	98.46%	Similar
Bacterial Filtration Efficiency ASTM F2101	For 3 non-consecutive lots, 32 out of 32 pass, average at 99.9%	98.7%	Similar
Differential Pressure (Delta P) ASTM F2100	For 3 non-consecutive lots, 32 out of 32 pass, average at 3.7 mmH <sub>2</sub> O/cm <sup>2</sup>	4.2 mmH <sub>2</sub> O/cm <sup>2</sup>	Similar
Flammability 16 CFR 1610	For 3 non-consecutive lots, 32 out of 32 pass, Class 1 Non-Flammable	Class 1 Non-Flammable	Same
<b>Biocompatibility</b>			
Cytotoxicity	Under the conditions of the study, not cytotoxicity effect	Under the conditions of the study, not cytotoxicity effect	Same
Irritation	Under the conditions of the study, not an irritant	Under the conditions of the study, not an irritant	Same
Sensitization	Under conditions of the study, not a sensitizer	Under conditions of the study, not a sensitizer	Same

The proposed Surgical Face Mask and the predicate device is identical in the intended use, mask style, design feature, color, ASTM F2100 level and biocompatibility, and similar only in ear loops' material. So the proposed device is identical to the predicate device.

## 7. PERFORMANCE DATA

The following performance data were provided to demonstrate that the subject device met the acceptance criteria of the test method or standard.

### **Biocompatibility testing**

The biocompatibility evaluation for the Surgical Face Mask was conducted in accordance with the International Standard ISO 10993-1:2018, "Biological evaluation of medical

devices - Part 1: Evaluation and testing within a risk management process” as recognized by FDA. The biocompatible testing included the following tests:

- Cytotoxicity - (ISO 10993-5: 2009)
- Sensitization - (ISO 10993-10:2010)
- Skin Irritation - (ISO 10993-10:2010)

### Performance testing

Performance testing was conducted on the Surgical Face Mask. All of the tested parameters met the predefined acceptance criteria.

Test item & Test Methods	Test purpose	Acceptance criteria	Results
<b>Flammability</b> ASTM F2100-19 16 CFR Part 1610-2008	Testing the characteristics of a material that pertain to its relative ease of ignition and relative ability to sustain combustion.	Class 1 ASTM F2100	Pass
<b>Bacterial Filtration Efficiency</b> ASTM F2100-2019 9.1 ASTM F2101-2019	Testing the effectiveness of medical face mask material in preventing the passage of aerosolized bacteria.	Level 2: ≥ 98% ASTM F2100	Pass
<b>Different Pressure, mm H<sub>2</sub>O/cm<sup>2</sup></b> ASTM F2100-19 9.2	Measuring the pressure of dropping across a medical face mask material.	Level 2: < 6.0 ASTM F2100	Pass
<b>Sub-Micron Particle Filtration Efficiency</b> ASTM F2100-2019 9.3 ASTM F2299/F2299M-2017	Testing the efficiency of the filter material in capturing aerosolized particles smaller than one micron.	Level 2: ≥ 98% ASTM F2100	Pass
<b>Resistance to Penetration by synthetic blood</b> ASTM F2100-2019 9.4 ASTM F1862/F1862-2017	Testing the efficiency of resistance to penetration by synthetic blood.	Level 2: pass at 120 mmHg ASTM F2100	Pass

### 8. CONCLUSION

The conclusions drawn from the nonclinical tests demonstrate that the device is as safe, as effective, and performs as well as or better than the legally marketed device.