



August 4, 2020

Leontyne Banks
Regulatory Affairs Specialist
Three Lakes Drive
Northfield, Illinois 60093

Re: K193666

Trade/Device Name: Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy and Fentanyl) - Regular Cuff; Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA, OPJ

Dated: June 10, 2020

Received: June 10, 2020

Dear Leontyne Banks:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

CAPT Elizabeth Claverie, M.S.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K193666

Device Name

Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff

Indications for Use (Describe)

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

| Chemo Drugs Tested: | Breakthrough Time (Minutes) |
|--|-----------------------------|
| Arsenic Trioxide 1 mg/ml (1,000 ppm) | >240 |
| Azacitidine 25 mg/ml (25,000 ppm) | >240 |
| Bendamustine 5 mg/ml (5,000 ppm) | >240 |
| Bleomycin 15 mg/ml (15,000 ppm) | >240 |
| Bortezomib 1 mg/ml (1,000 ppm) | >240 |
| Busulfan 6 mg/ml (6,000 ppm) | >240 |
| Carboplatin 10.0 mg/ml (10,000 ppm) | >240 |
| Carfilzomib 2.0 mg/ml (2,000) | >240 |
| Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) | 33.1 min |
| Cetuximab 2.0 mg/ml (2,000) | >240 |
| Cisplatin 1.0 mg/ml (1,000 ppm) | >240 |
| Cyclophosphamide (Cytosan) 20.0 mg/ml (20,000 ppm) | >240 |
| Cytarabine 100 mg/ml (100,000 ppm) | >240 |
| Cytovene 10 mg/ml (10,000 ppm) | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm) | >240 |
| Daunorubicin 5 mg/ml (5,000 ppm) | >240 |
| Decitabine 5 mg/ml (5,000 ppm) | >240 |
| Docetaxel 10.0 mg/ml (10,000 ppm) | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) | >240 |
| Epirubicin (Ellence) 2.0 mg/ml (2,000 ppm) | >240 |
| Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) | >240 |
| Fludarabine 25 mg/ml (25,000 ppm) | >240 |
| Fluorouracil 50.0 mg/ml (50,000 ppm) | >240 |
| Fulvestrant 50.0 mg/ml (50,000 ppm) | >240 |
| Gemcitabine (Gemzar) 38 mg/ml (38,000 ppm) | >240 |
| Idarubicin 1 mg/ml (1,000 ppm) | >240 |
| Ifosfamide 50.0 mg/ml (50,000 ppm) | >240 |
| Irinotecan 20.0 mg/ml (20,000 ppm) | >240 |
| Mechlorethamine HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Melphalan 5 mg/ml (5,000 ppm) | >240 |
| Methotrexate 25 mg/ml (25,000 ppm) | >240 |
| Mitomycin C 0.5 mg/ml (500 ppm) | >240 |
| Mitoxantrone 2.0 mg/ml (2,000 ppm) | >240 |
| Oxaliplatin 2.0 mg/ml (2,000 ppm) | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) | >240 |
| Paraplatin 10 mg/ml (10,000 ppm) | >240 |
| Pemetrexed Disodium 25.0 mg/ml (25,000 ppm) | >240 |
| Pertuzumab 30 mg/ml (30,000 ppm) | >240 |
| Raltitrexed 0.5 mg/ml (500 ppm) | >240 |

| | |
|---|----------|
| Retrovir 10.0 mg/ml (10,000 ppm) | >240 |
| Rituximab 10 mg/ml (10,000 ppm) | >240 |
| Temsirolimus 25.0 mg/ml (25,000 ppm) | >240 |
| Thiotepa 10.0 mg/ml (10,000 ppm) | 69.2 min |
| Topotecan HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Trastuzumab 21.0 mg/ml (21,000 ppm) | >240 |
| Trisonex 1 mg/ml (1,000 ppm) | >240 |
| Vinblastine 1 mg/ml (1,000 ppm) | >240 |
| Vincristine Sulfate 1.0 mg/ml (1,000 ppm) | >240 |
| Vinorelbine 10 mg/ml (10,000 ppm) | >240 |
| Zoledronic Acid 0.8 mg/ml (800 ppm) | >240 |
| Fentanyl Citrate, 100mcg/2mL | >240 |

Please note that the following drug has low permeation time:

Carmustine (BCNU) (3.3 mg/ml) 33.1 minutes

Thiotepa (THT) (10.0mg/ml) 69.2 minutes

CAUTION: Testing showed an average breakthrough time of 33.1 minutes with Carmustine and 69.2 minutes with Thiotepa

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
 Food and Drug Administration
 Office of Chief Information Officer
 Paperwork Reduction Act (PRA) Staff
 PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Indications for Use

510(k) Number (if known)

K193666

Device Name

Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff

Indications for Use (Describe)

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

| Chemotherapy Drug Tested | Breakthrough Time (in minutes) |
|---------------------------------------|--------------------------------|
| Arsenic Trioxide 1 mg/ml | >240 |
| Azacitidine 25 mg/ml | >240 |
| Bendamustine 5 mg/ml | >240 |
| Bleomycin 15 mg/m | >240 |
| Bortezomib 1 mg/ml | >240 |
| Busulfan 6 mg/ml | >240 |
| Carfilzomib 2 mg/ml | >240 |
| Carmustine (BCNU) 3.3 mg/ml | 59.4 min |
| Cetuximab | >240 |
| Cisplatin 1.0 mg/ml | >240 |
| Cyclophosphamide (Cytoxan) 20.0 mg/ml | >240 |
| Cytarabine Hydrochloride 100 mg/ml | >240 |
| Cytovene 10 mg/ml | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml | >240 |
| Daunorubicin 5 mg/ml | >240 |
| Decitabine 5 mg/ml | >240 |
| Docetaxel 10.0 mg/ml | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml | >240 |
| Epirubicin (Ellence) 2.0 mg/ml | >240 |
| Etoposide (Toposar) 20.0 mg/ml | >240 |
| Fludarabine 25 mg/ml | >240 |
| Fluorouracil 50.0 mg/ml | >240 |
| Fulvestrant 50.0 mg/ml | >240 |
| Gemcitabine (Gemzar) 38 mg/ml | >240 |
| Idarubicin 1 mg/ml | >240 |
| Ifosfamide 50.0 mg/ml | >240 |
| Irinotecan 20.0 mg/ml | >240 |
| Mechlorethamine HCl 1.0 mg/ml | >240 |
| Melphalan 5 mg/ml | >240 |
| Methotrexate 25.0 mg/ml | >240 |
| Mitomycin C 0.5 mg/ml | >240 |
| Mitoxantrone 2mg/ml | >240 |
| Oxaliplatin 5.0 mg/ml | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml | >240 |
| Paraplatin (Carboplatin) 10 mg/ml | >240 |
| Pemetrexed Disodium 25.0 mg/ml | >240 |
| Pertuzumab 30 mg/ml | >240 |
| Raltitrexed 0.5 mg/ml | >240 |
| Retrovir 10.0 mg/ml | >240 |

| | |
|-------------------------------|----------|
| Rituximab 10 mg/ml | >240 |
| Temsirolimus 25.0 mg/ml | >240 |
| Thiotepa 10.0 mg/ml | 68.2 min |
| Topotecan HCl 1.0 mg/ml | >240 |
| Trastuzumab 21.0 mg/ml | >240 |
| Trisonex 1 mg/ml | >240 |
| Vinblastine 1 mg/ml | >240 |
| Vincristine Sulfate 1.0 mg/ml | >240 |
| Vinorelbine 10 mg/ml | >240 |
| Zoledronic Acid 0.8 mg/ml | >240 |
| Fentanyl Citrate, 100mcg/2mL | >240 |

Please note that the following drug has low permeation time:

Carmustine (BCNU) (3.3 mg/ml) 59.4 minutes

Thiotepa (THT) (10.0 mg/ml) 68.2 minutes

CAUTION: Testing showed an average breakthrough time of 59.4 minutes with Carmustine and 68.2 minutes with Thiotepa

Type of Use (*Select one or both, as applicable*)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
 Food and Drug Administration
 Office of Chief Information Officer
 Paperwork Reduction Act (PRA) Staff
 PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

510(k) SUMMARY [AS REQUIRED BY 21CFR807.92(c)]

Submitter / 510(k) Sponsor

Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093
Registration Number: 1417592
Applicant Contact: Leontyne Banks
Regulatory Affairs Specialist

Summary Preparation Date

July 31, 2020

Type of 510(k) Submission

Traditional

Device Name / Classification

Name of Device: Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff; Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff

Proprietary Name: Medline Powder-Free Blue Nitrile Patient Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff; Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff

Common Name: Patient Examination Glove, Specialty

Classification Name: Polymer Patient Exam Glove, Medical Gloves with Chemotherapy/Fentanyl Labeling Claims

Product Code: LZA, OPJ, QDO

Classification Panel: General Hospital

Regulation #: 21 CFR 880.6250

Predicate Device

YTY Industry (Manjung): Non-Sterile, Powder-Free Nitrile Examination Gloves (Cobalt Blue) Tested for use with Chemotherapy Drugs – K111248



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

Device Description

The Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff are non-sterile, single use, disposable gloves intended for medical purposes to be worn on the hands of examiners to prevent contamination between a patient and an examiner. The gloves are nitrile, powder-free, ambidextrous, and blue-colored with a beaded cuff. The proposed device is offered in two versions: Regular cuff (VS311) and Extended cuff (VS711).

TABLE 1: The Medline Powder-Free Blue Nitrile Examination Gloves – Regular Cuff Sizes

| Model Number | Size |
|--------------|-------------|
| VS311XS | Extra-Small |
| VS311S | Small |
| VS311M | Medium |
| VS311L | Large |
| VS311XL | Extra-Large |

TABLE 2: The Medline Powder-Free Blue Nitrile Examination Gloves – Extended Cuff Sizes

| Model Number | Size |
|--------------|-------------|
| VS711XS | Extra-Small |
| VS711S | Small |
| VS711M | Medium |
| VS711L | Large |
| VS711XL | Extra-Large |

The gloves are designed and manufactured in accordance with the ASTM D6319-10 Standard Specification for Nitrile Examination Gloves for Medical Application and are tested for use with chemotherapy drugs as well as for use with Fentanyl per ASTM D6978-05 (Reapproved 2019).

Indications for Use – Regular Cuff (VS311XS, VS311S, VS311M, VS311L and VS311XL)

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner’s hand to prevent contamination between patient and examiner. These gloves were tested



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

for use with chemotherapy drugs per ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs

| Chemo Drugs Tested: | Breakthrough Time (Minutes) |
|--|-----------------------------|
| Arsenic Trioxide 1 mg/ml (1,000 ppm) | >240 |
| Azacitidine 25 mg/ml (25,000 ppm) | >240 |
| Bendamustine 5 mg/ml (5,000 ppm) | >240 |
| Bleomycin 15 mg/ml (15,000 ppm) | >240 |
| Bortezomib 1 mg/ml (1,000 ppm) | >240 |
| Busulfan 6 mg/ml (6,000 ppm) | >240 |
| Carboplatin 10.0 mg/ml (10,000 ppm) | >240 |
| Carfilzomib 2.0 mg/ml (2,000) | >240 |
| Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) | 33.1 min |
| Cetuximab 2.0 mg/ml (2,000) | >240 |
| Cisplatin 1.0 mg/ml (1,000 ppm) | >240 |
| Cyclophosphamide (Cytosan) 20.0 mg/ml (20,000 ppm) | >240 |
| Cytarabine 100 mg/ml (100,000 ppm) | >240 |
| Cytovene 10 mg/ml (10,000 ppm) | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm) | >240 |
| Daunorubicin 5 mg/ml (5,000 ppm) | >240 |
| Decitabine 5 mg/ml (5,000 ppm) | >240 |
| Docetaxel 10.0 mg/ml (10,000 ppm) | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) | >240 |
| Epirubicin (Ellence) 2.0 mg/ml (2,000 ppm) | >240 |
| Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) | >240 |
| Fludarabine 25 mg/ml (25,000 ppm) | >240 |
| Fluorouracil 50.0 mg/ml (50,000 ppm) | >240 |
| Fulvestrant 50.0 mg/ml (50,000 ppm) | >240 |
| Gemcitabine (Gemzar) 38 mg/ml (38,000 ppm) | >240 |
| Idarubicin 1 mg/ml (1,000 ppm) | >240 |
| Ifosfamide 50.0 mg/ml (50,000 ppm) | >240 |
| Irinotecan 20.0 mg/ml (20,000 ppm) | >240 |
| Mechlorethamine HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Melphalan 5 mg/ml (5,000 ppm) | >240 |
| Methotrexate 25 mg/ml (25,000 ppm) | >240 |
| Mitomycin C 0.5 mg/ml (500 ppm) | >240 |
| Mitoxantrone 2.0 mg/ml (2,000 ppm) | >240 |
| Oxaliplatin 2.0 mg/ml (2,000 ppm) | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) | >240 |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | |
|---|----------|
| Paraplatin 10 mg/ml (10,000 ppm) | >240 |
| Pemetrexed Disodium 25.0 mg/ml (25,000 ppm) | >240 |
| Pertuzumab 30 mg/ml (30,000 ppm) | >240 |
| Raltitrexed 0.5 mg/ml (500 ppm) | >240 |
| Retrovir 10.0 mg/ml (10,000 ppm) | >240 |
| Rituximab 10 mg/ml (10,000 ppm) | >240 |
| Temsirolimus 25.0 mg/ml (25,000 ppm) | >240 |
| Thiotepa 10.0 mg/ml (10,000 ppm) | 69.2 min |
| Topotecan HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Trastuzumab 21.0 mg/ml (21,000 ppm) | >240 |
| Trisonex 1 mg/ml (1,000 ppm) | >240 |
| Vinblastine 1 mg/ml (1,000 ppm) | >240 |
| Vincristine Sulfate 1.0 mg/ml (1,000 ppm) | >240 |
| Vinorelbine 10 mg/ml (10,000 ppm) | >240 |
| Zoledronic Acid 0.8 mg/ml (800 ppm) | >240 |
| Fentanyl Citrate, 100mcg/2mL | >240 |

Please note that the following drugs have low permeation times:

- Carmustine (BCNU) (3.3 mg/ml) 33.1 minutes
- Thiotepa (10.0 mg/ml) 69.2 minutes

CAUTION: Testing showed an average breakthrough time of 33.1 minutes with Carmustine and 69.2 minutes with Thiotepa.

Indications for Use – Extended Cuff (VS711XS, VS711S, VS711M, VS711L and VS711XL)

A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner’s hand to prevent contamination between patient and examiner.

These gloves were tested for use with chemotherapy drugs per ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs.

| Chemo Drugs Tested: | Breakthrough Time (Minutes) |
|--------------------------|-----------------------------|
| Arsenic Trioxide 1 mg/ml | >240 |
| Azacitidine 25 mg/ml | >240 |
| Bendamustine 5 mg/ml | >240 |
| Bleomycin 15 mg/m | >240 |
| Bortezomib 1 mg/ml | >240 |
| Busulfan 6 mg/ml | >240 |



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

| | |
|---------------------------------------|----------|
| Carfilzomib 2 mg/ml | >240 |
| Carmustine (BCNU) 3.3 mg/ml | 59.4 min |
| Cetuximab | >240 |
| Cisplatin 1.0 mg/ml | >240 |
| Cyclophosphamide (Cytosan) 20.0 mg/ml | >240 |
| Cytarabine Hydrochloride 100 mg/ml | >240 |
| Cytovene 10 mg/ml | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml | >240 |
| Daunorubicin 5 mg/ml | >240 |
| Decitabine 5 mg/ml | >240 |
| Docetaxel 10.0 mg/ml | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml | >240 |
| Epirubicin (Ellence) 2.0 mg/ml | >240 |
| Etoposide (Toposar) 20.0 mg/ml | >240 |
| Fludarabine 25 mg/ml | >240 |
| Fluorouracil 50.0 mg/ml | >240 |
| Fulvestrant 50.0 mg/ml | >240 |
| Gemcitabine (Gemzar) 38 mg/ml | >240 |
| Idarubicin 1 mg/ml | >240 |
| Irinotecan 20.0 mg/ml | >240 |
| Mechlorethamine HCl 1.0 mg/ml | >240 |
| Melphalan 5 mg/ml | >240 |
| Methotrexate 25.0 mg/ml | >240 |
| Mitomycin C 0.5 mg/ml | >240 |
| Oxaliplatin 5.0 mg/ml | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml | >240 |
| Paraplatin (Carboplatin) 10 mg/ml | >240 |
| Pemetrexed Disodium 25.0 mg/ml | >240 |
| Pertuzumab 30 mg/ml | >240 |
| Raltitrexed 0.5 mg/ml | >240 |
| Retrovir 10.0 mg/ml | >240 |
| Rituximab 10 mg/ml | >240 |
| Temsirolimus 25.0 mg/ml | >240 |
| Thiotepa 10.0 mg/ml | 68.2 min |
| Topotecan HCl 1.0 mg/ml | >240 |
| Trastuzumab 21.0 mg/ml | >240 |
| Trisonex 1 mg/ml | >240 |
| Vinblastine 1 mg/ml | >240 |
| Vincristine Sulfate 1.0 mg/ml | >240 |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | |
|------------------------------|------|
| Vinorelbine 10 mg/ml | >240 |
| Zoledronic Acid 0.8 mg/ml | >240 |
| Fentanyl Citrate, 100mcg/2mL | >240 |

Please note that the following drugs have low permeation times:
 Carmustine (BCNU) (3.3 mg/ml) 59.4 minutes
 Thiotepa (10.0 mg/ml) 68.2 minutes

CAUTION: Testing showed an average breakthrough time of 59.4 minutes with Carmustine and 68.2 minutes with Thiotepa.

Comparison of Technological Characteristics

TABLE 3: Comparison of Proposed – VS311 and Predicate Device

| Device Characteristic | Proposed Device VS311 | Predicate Device | Comparison Analysis |
|-------------------------|--|---|---------------------|
| Product Name | Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) Regular Cuff | Non-Sterile, Powder-Free Nitrile Examination Gloves (Cobalt) Tested for use with chemotherapy drugs | N/A |
| 510(k) Reference | K193666 | K111248 | N/A |
| Product Owner | Medline Industries, Inc. | YTY Industry (Manjung) SDN, BHD | N/A |
| Product Code | LZA, OPJ, QDO | LZA, LZC | Similar |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|-------------------------------|---|---|----------------|
| Intended Use | A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner. These gloves were tested for use with Chemotherapy drugs and Fentanyl Citrate as per ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Medical Gloves of Permeation by Chemotherapy Drugs | A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or fingers to prevent contamination between patient and examiner. These glove were tested for use with chemotherapy drugs. | Same |
| Regulation Number | 21 CFR 880.6250 | 21 CFR 880.6250 | Same |
| Materials | Powder-Free Nitrile | Powder-Free Nitrile | Same |
| Color | Blue | Blue | Same |
| Sizes | x- small, small, medium, large, x-large | Medium, Large, Extra-large | Similar |
| Dimensions – Length | Complies with: <u>ASTM D6319-10</u> 220mm min | Complies with: <u>ASTM D6319-10</u> 220mm min. | Same |
| Dimensions – Width | Complies with: <u>ASTM D6319-10</u> 70mm min | Complies with: <u>ASTM D6319-10</u> 70mm min | Same |
| Dimensions – Thickness | Complies with: <u>ASTM D6319-10</u> Palm – 0.05mm min. Finger – 0.05mm min. | Complies with: <u>ASTM D6319-10</u> Palm – 0.05mm min. Finger – 0.05mm min. | Same |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|----------------------------------|---|---|------------------|
| Physical Properties | Complies with: <u>ASTM D6319-10</u> minimum: Tensile Strength: Before Aging ≥14 MPa, min. After Aging ≥14 MPa, min. | Complies with: <u>ASTM D6319-10</u> minimum: Tensile Strength: Before Aging ≥14 MPa, min. After Aging ≥14 MPa, min. | Same |
| | Elongation: Before Aging 500%, min. After Aging 400%, min. | Elongation: Before Aging 500% min. After Aging 400% min. | |
| Freedom from Holes | Complies with: ASTM D6319-10 and ASTM D5151-06 G-1, AQL 1.5 | Complies with: ASTM D6319-10 and ASTM D5151-06 G-1, AQL 2.5 | Same |
| Powder or Powder-Free | Powder-Free | Powder-Free | Same |
| Residual Powder | Complies with: <u>ASTM D6319-10</u> <2mg per glove | Complies with: <u>ASTM D6319-05</u> <2mg per glove | Same |
| Contact Durations | Limited ≤24 hours | Limited ≤24 hours | Same |
| Biocompatibility | AAMI/ANSI/ISO 10993-10: Not a skin irritant Not a skin sensitizer | AAMI/ANSI/ISO 10993-10: Not a skin irritant Not a skin sensitizer | Same |
| Sterility | Non-sterile | Non-sterile | Same |
| Rx Only or OTC | Over the Counter | Over the Counter | Same |
| Tested Chemotherapy Drugs | Aresenic Trioxide 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| | Azacididine (Vidaza) 25.0 mg/ml (25,000 ppm) >240 min. | | Different |
| | Bendamustine 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| | Bleomycin 15.0 mg/ml (15,000 ppm) >240 min. | | Different |
| | Bortezomib 1.0 mg/ml (1,000 ppm) >240 min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|---|---|------------------|
| Busulfan 6.0 mg/ml (6,000 ppm) >240 min. | | Different |
| Carboplatin 10mg/ml 10,000 ppm >240 min | Carboplatin 10mg/ml 10,000 ppm >240 min | Same |
| Carfilzomib 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) 33.1 min. | Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) 1.82 min | Similar |
| Cetuximab 2.0 mg/ml (2,000 ppm) >240 min | | Different |
| Cisplatin 1.0 mg.ml (1,000 ppm) >240 min. | Cisplatin 1.0 mg/ml (1,000 ppm) >240 min. | Same |
| Cyclophosphamide (Cytosan) 20.0 mg/ml (20,000 ppm) >240 min. | Cyclophosphamide (Cytosan) 20.0 mg/ml (20,000 ppm) >240 min. | Same |
| Cytarabine 100.0 mg/ml (100,000 ppm) >240 min. | | Different |
| Cytovene 10.0 mg/ml (10,000 ppm) >240 min. | | Different |
| Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm) >240 min. | Dacarbazine (DTIC) 10 mg/ml (10,000 ppm) >240 min. | Same |
| Daunorubicin 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Decitabine 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Docetaxel 10.0 mg/ml (10,000 ppm) >240 min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|---|--|------------------|
| Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) >240 min. * | Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) >240 min | Same |
| Epirubicin (Ellence) 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) >240 min. * | Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) >240 min. | Same |
| Fludarabine 25.0 mg/ml (25,000 ppm) >240 min. | | Different |
| Fluorouracil 50.0 mg/ml (50,000 ppm) >240 min. | Fluorouracil 50.0 mg/ml (50,000 ppm) >240 min. | Same |
| Fulvestrant 50.0 mg/ml (50,000 ppm) >240 min. | | Different |
| Gemcitabine (Gemzar) 38.0 mg/ml (38,000 ppm) >240 min. | | Different |
| Idarubicin 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| Ifosfamide 50.0 mg/ml (50,000 ppm) >240 min. | Ifosfamide 50.0 mg/ml (50,000 ppm) >240 min. | Same |
| Irinotecan 20.0 mg/ml (20,000 ppm) >240 min. | | Different |
| Mechlorethamine HCl 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| Melphalan 5.0 mg/ml (5,000 ppm) >240 min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|--|--|------------------|
| Methotrexate 25.0 mg/ml (25,000 ppm) >240 min. * | Methotrexate 25.0 mg/ml (25,000 ppm) >240 min. | Same |
| Mitomycin C 0.5 mg/ml (500 ppm) >240 min. | Mitomycin C 0.5 mg/ml (2,000 ppm) >240 min. | Same |
| Mitoxantrone 2.0 mg/ml (2,000 ppm) >240 min. | Mitoxantrone 2.0 mg/ml (2,000 ppm) >240 min. | Same |
| Oxaliplatin 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) >240 min. | Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) >240 min. | Same |
| Paraplatin 10.0 mg/ml (10,000 ppm) >240 min. | | Different |
| Pemetrexed Disodium 25.0 mg/ml (25,000 ppm) >240 min. | | Different |
| Pertuzumab 30.0 mg/ml (30,000 ppm) >240 min. | | Different |
| Raltitrexed 0.5 mg/ml (500 ppm) 240> min. | | Different |
| Retrovir 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| Rituximab 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| Temsirolimus 25.0 mg/ml (25,000 ppm) 240> min. | | Different |
| Thiotepa 10.0 mg/ml (10,000 ppm) 69.2 min. | Thiotepa 10.0 mg/ml (10,000 ppm) 0.93 min | Similar |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|-------------------------|---|---|-----------|
| | Topotecan HCl 1.0 mg/ml (1,000 ppm) 240> min. | | Different |
| | Trastuzumab 21.0 mg/ml (21,000 ppm) 240> min. | | Different |
| | Trisonex 1.0 mg/ml (1,000 ppm) 240> min. | | Different |
| | Vinblastine 1 mg/ml (1,000 ppm) 240> min. | | Different |
| | Vincristine Sulfate 1.0 mg/ml (1,000 ppm) >240 min. | Vincristine Sulfate 1.0 mg/ml (1,000 ppm) >240 min. | Same |
| | Vinorelbine 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| | Zoledronic Acid 0.8 mg/ml (800 ppm) 240> min. | | Different |
| Fentanyl Testing | Fentanyl Citrate, 100.0 mcg/2mL >240 min. | | Different |

TABLE 4: Comparison of Proposed – VS711 and Predicate Device

| Device Characteristic | Proposed Device VS711 | Predicate Device | Comparison Analysis |
|-------------------------|---|---|---------------------|
| Product Name | Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) Extended Cuff | Non-Sterile, Powder-Free Nitrile Examination Gloves (Cobalt) Tested for use with chemotherapy drugs | N/A |
| 510(k) Reference | K193666 | K111248 | N/A |
| Product Owner | Medline Industries, Inc. | YTY Industry (Manjung) SDN, BHD | N/A |
| Product Code | LZA, OPJ, QDO | LZA, LZC | Similar |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|-----------------------------------|--|---|-----------------------|
| <p>Intended Use</p> | <p>A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner’s hand to prevent contamination between patient and examiner. These gloves were tested for use with Chemotherapy drugs and Fentanyl Citrate as per ASTM D6978-05 (Reapproved 2019) Standard Practice for Assessment of Medical Gloves of Permeation by Chemotherapy Drugs</p> | <p>A patient examination glove is a disposable device intended for medical purposes that is worn on the examiner’s hand or fingers to prevent contamination between patient and examiner.</p> <p>These glove were tested for use with chemotherapy drugs.</p> | <p>Same</p> |
| <p>Regulation Number</p> | <p>21 CFR 880.6250</p> | <p>21 CFR 880.6250</p> | <p>Same</p> |
| <p>Materials</p> | <p>Powder-Free Nitrile</p> | <p>Powder-Free Nitrile</p> | <p>Same</p> |
| <p>Color</p> | <p>Blue</p> | <p>Blue</p> | <p>Same</p> |
| <p>Sizes</p> | <p>x- small, small, medium, large, x-large</p> | <p>Medium, Large, Extra-large</p> | <p>Similar</p> |
| <p>Dimensions – Length</p> | <p>Complies with: <u>ASTM D6319-10</u> 230mm min</p> | <p>Complies with: <u>ASTM D6319-10</u> 230mm min.</p> | <p>Same</p> |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|----------------------------------|---|---|------------------|
| Dimensions – Width | Complies with: <u>ASTM D6319-10</u> 70mm min | Complies with: <u>ASTM D6319-10</u> 70mm min | Same |
| Dimensions – Thickness | Complies with: <u>ASTM D6319-10</u> Palm – 0.05mm min. Finger – 0.05mm min. | Complies with: <u>ASTM D6319-10</u> Palm – 0.05mm min. Finger – 0.05mm min. | Same |
| Physical Properties | Complies with: <u>ASTM D6319-10</u> minimum: Tensile Strength: Before Aging ≥14 MPa, min. After Aging ≥14 MPa, min. | Complies with: <u>ASTM D6319-10</u> minimum: Tensile Strength: Before Aging ≥14 MPa, min. After Aging ≥14 MPa, min. | Same |
| | Elongation: Before Aging 500%, min. After Aging 400%, min. | Elongation: Before Aging 500% min. After Aging 400% min. | |
| Freedom from Holes | Complies with: ASTM D6319-10 and ASTM D5151-06 G-1, AQL 1.5 | Complies with: ASTM D6319-10 and ASTM D5151-06 G-1, AQL 2.5 | Same |
| Powder or Powder-Free | Powder-Free | Powder-Free | Same |
| Residual Powder | Complies with: <u>ASTM D6319-10</u> <2mg per glove | Complies with: <u>ASTM D6319-05</u> <2mg per glove | Same |
| Contact Durations | Limited ≤24 hours | Limited ≤24 hours | Same |
| Biocompatibility | AAMI/ANSI/ISO 10993-10: Not a skin irritant Not a skin sensitizer | AAMI/ANSI/ISO 10993-10: Not a skin irritant Not a skin sensitizer | Same |
| Sterility | Non-sterile | Non-sterile | Same |
| Rx Only or OTC | Over the Counter | Over the Counter | Same |
| Tested Chemotherapy Drugs | Aresenic Trioxide 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| | Azacitidine (Vidaza) 25.0 mg/ml (25,000 ppm) >240 min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|---|---|------------------|
| Bendamustine 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Bleomycin 15.0 mg/ml (15,000 ppm) >240 min. | | Different |
| Bortezomib 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| Busulfan 6.0 mg/ml (6,000 ppm) >240 min. | | Different |
| Carboplatin 10mg/ml 10,000 ppm >240 min | Carboplatin 10mg/ml 10,000 ppm >240 min | Same |
| Carfilzomib 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) 59.4 min. | Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) 1.82 min | Similar |
| Cetuximab 2.0 mg/ml (2,000 ppm) >240 min | | Different |
| Cisplatin 1.0 mg/ml (1,000 ppm) >240 min. | Cisplatin 1.0 mg/ml (1,000 ppm) >240 min. | Same |
| Cyclophosphamide (Cytoxan) 20.0 mg/ml (20,000 ppm) >240 min. | Cyclophosphamide (Cytoxan) 20.0 mg/ml (20,000 ppm) >240 min. | Same |
| Cytarabine 100.0 mg/ml (100,000 ppm) >240 min. | | Different |
| Cytovene 10.0 mg/ml (10,000 ppm) >240 min. | | Different |
| Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm) >240 min. | Dacarbazine (DTIC) 10 mg/ml (10,000 ppm) >240 min. | Same |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|--|---|------------------|
| Daunorubicin 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Decitabine 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Docetaxel 10.0 mg/ml (10,000 ppm) >240 min. | | Different |
| Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) >240 min. * | Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) >240 min | Same |
| Epirubicin (Ellence) 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) >240 min. * | Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) >240 min. | Same |
| Fludarabine 25.0 mg/ml (25,000 ppm) >240 min. | | Different |
| Fluorouracil 50.0 mg/ml (50,000 ppm) >240 min. | Fluorouracil 50.0 mg/ml (50,000 ppm) >240 min. | Same |
| Fulvestrant 50.0 mg/ml (50,000 ppm) >240 min. | | Different |
| Gemcitabine (Gemzar) 38.0 mg/ml (38,000 ppm) >240 min. | | Different |
| Idarubicin 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| Irinotecan 20.0 mg/ml (20,000 ppm) >240 min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | |
|---|--|------------------|
| Mechlorethamine HCl 1.0 mg/ml (1,000 ppm) >240 min. | | Different |
| Melphalan 5.0 mg/ml (5,000 ppm) >240 min. | | Different |
| Methotrexate 25.0 mg/ml (25,000 ppm) >240 min. * | Methotrexate 25.0 mg/ml (25,000 ppm) >240 min. | Same |
| Mitomycin C 0.5 mg/ml (500 ppm) >240 min. | Mitomycin C 0.5 mg/ml (2,000 ppm) >240 min. | Same |
| Oxaliplatin 2.0 mg/ml (2,000 ppm) >240 min. | | Different |
| Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) >240 min. | Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) >240 min. | Same |
| Paraplatin 10.0 mg/ml (10,000 ppm) >240 min. | | Different |
| Pemetrexed Disodium 25.0 mg/ml (25,000 ppm) >240 min. | | Different |
| Pertuzumab 30.0 mg/ml (30,000 ppm) >240 min. | | Different |
| Raltitrexed 0.5 mg/ml (500 ppm) 240> min. | | Different |
| Retrovir 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| Rituximab 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| Temsirolimus 25.0 mg/ml (25,000 ppm) 240> min. | | Different |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|-------------------------|---|---|------------------|
| | Thiotepa (THT) 10.0 mg/ml (10,000 ppm) 68.2 min. | Thiotepa 10.0 mg/ml (10,000 ppm) 0.93 min | Similar |
| | Topotecan HCl 1.0 mg/ml (1,000 ppm) 240> min. | | Different |
| | Trastuzumab 21.0 mg/ml (21,000 ppm) 240> min. | | Different |
| | Trisonex 1.0 mg/ml (1,000 ppm) 240> min. | | Different |
| | Vinblastine 1 mg/ml (1,000 ppm) 240> min. | | Different |
| | Vincristine Sulfate 1.0 mg/ml (1,000 ppm) >240 min. | Vincristine Sulfate 1.0 mg/ml (1,000 ppm) >240 min. | Same |
| | Vinorelbine 10.0 mg/ml (10,000 ppm) 240> min. | | Different |
| | Zoledronic Acid 0.8 mg/ml (800 ppm) 240> min. | | Different |
| Fentanyl Testing | Fentanyl Citrate, 100.0 mcg/2mL >240 min. | | Different |

Summary of Non-Clinical Testing:

Biocompatibility Testing

Table 5 below summarize the biocompatibility results for the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and the Medline Powder-Free Blue Nitrile Examination Gloves-(Tested for use with Chemotherapy Drugs and Fentanyl) – Extended Cuff



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

Table 5: Biocompatibility Testing – VS311 & VS711

| Name of Test /Citation | Purpose | Acceptance Criteria | Results |
|--|---------------------------|---------------------|---|
| ISO 10993-10: 2010 Biological Evaluation of Medical Devices, Part 10-Tests for Irritation and Skin Sensitization | Irritation Testing | Pass/Fail | Pass Under the conditions of the study, the subject device is not a primary skin irritant. |
| ISO 10993-10:2010. Biological Evaluation of Medical Devices, Part 10: Tests for Irritation and Skin Sensitization. | Sensitization Testing | Pass/Fail | Pass Under the conditions of the study, the subject device is not a primary skin sensitizer. |
| ISO 10993-5:2009 Biological Evaluation of Medical Devices, Part 5: Tests for <i>In Vitro</i> Cytotoxicity | Cytotoxicity Testing | Pass/Fail | Failed Under the conditions of the study, the subject device is cytotoxic. |
| ISO 10993-11: 2017 Biological Evaluation of Medical Devices Part 11: Tests for Systemic Toxicity. | Systemic Toxicity Testing | Pass/Fail | Pass Under the conditions of the study, the subject device is not toxic. |

Table 6: Quality Assurance Testing For Finished Gloves

| Characteristics | Reference Test Method & Sampling | Inspection Level | Inspection Level | Acceptance Criteria |
|--|--|------------------|------------------|------------------------|
| Dimension (length, width, thickness) | ASTM D6319 – 10 (2015) | S-2 | 4.0 | Meet Acceptance Number |
| Tensile strength and Ultimate elongation | ASTM D 6319 – 10 (2015) | S-2 | 4.0 | Meet Acceptance Number |
| Water Leak | ASTM D 5151-06 (2015) and ASTM D 6319-10 | G-1 | 1.5 | Meet Acceptance Number |
| Powder Residue | ASTM D6124 - 06 (2011) | N/A | N/A | 2 mg/glove |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

Physical Dimensions

Thirteen (13) glove samples of each size (extra-small, small, medium, large and extra-large) were tested for physical dimensions according to ASTM D6319-10. All samples from each of the five size gloves met the specifications for physical dimensions.

Tables 7 through 16 below summarize the physical dimension results for the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff.

TABLE 7: ASTM D6319-10 Test Results – Extra-Small, VS311XS (Regular Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 220mm min. | 247mm | Yes |
| Palm Width | 70 ± 10mm | 78mm | Yes |
| Thickness Finger | 0.05mm min. | 0.16mm | Yes |
| Thickness Palm | 0.05mm min. | 0.11mm | Yes |

TABLE 8: ASTM D6319-10 Test Results – Small, VS311S (Regular Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 220mm min. | 248mm | Yes |
| Palm Width | 80 ± 10mm | 86mm | Yes |
| Thickness Finger | 0.05mm min | 0.16mm | Yes |
| Thickness Palm | 0.05mm min | 0.11mm | Yes |

TABLE 9: ASTM D6319-10 Test Results – Medium, VS311M (Regular Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 250mm | Yes |
| Palm Width | 95 ± 10mm | 97mm | Yes |
| Thickness Finger | 0.05mm | 0.15mm | Yes |
| Thickness Palm | 0.05mm | 0.11mm | Yes |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

TABLE 10: ASTM D6319-10 Test Results – Large, VS311L (Regular Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 250mm | Yes |
| Palm Width | 110 ± 10mm | 107mm | Yes |
| Thickness Finger | 0.05mm | 0.16mm | Yes |
| Thickness Palm | 0.05mm | 0.11mm | Yes |

TABLE 11: ASTM D6319-10 Test Results – Extra-Large, VS311XL (Regular Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 251mm | Yes |
| Palm Width | 120 ± 10mm | 119mm | Yes |
| Thickness Finger | 0.05mm | 0.15mm | Yes |
| Thickness Palm | 0.05mm | 0.11mm | Yes |

TABLE 12: ASTM D6319-10 Test Results – Extra-Small, VS711XS (Extended Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 220mm min. | 302mm | Yes |
| Palm Width | 70 ± 10mm | 78mm | Yes |
| Thickness Finger | 0.05mm min. | 0.17mm | Yes |
| Thickness Palm | 0.05mm min. | 0.12mm | Yes |

TABLE 13: ASTM D6319-10 Test Results – Small, VS711S (Extended Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 220mm | 300mm | Yes |
| Palm Width | 80 ± 10mm | 89mm | Yes |
| Thickness Finger | 0.05mm min. | 0.18mm | Yes |
| Thickness Palm | 0.05mm min. | 0.12mm | Yes |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

TABLE 14: ASTM D6319-10 Test Results – Medium, VS711M (Extended Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 301mm | Yes |
| Palm Width | 95 ± 10mm | 99mm | Yes |
| Thickness Finger | 0.05mm | 0.18mm | Yes |
| Thickness Palm | 0.05mm | 0.13mm | Yes |

TABLE 15: ASTM D6319-10 Test Results – Large, VS711L (Extended Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 299mm | Yes |
| Palm Width | 110 ± 10mm | 108mm | Yes |
| Thickness Finger | 0.05mm | 0.16mm | Yes |
| Thickness Palm | 0.05mm | 0.12mm | Yes |

TABLE 16: ASTM D6319-10 Test Results – Extra-Large, VS711XL (Extended Cuff)

| Characteristic | Device Specification | Average Device Performance (n=13) | Complies with Standard (Y/N) |
|------------------|----------------------|-----------------------------------|------------------------------|
| Glove Length | 230mm min. | 298mm | Yes |
| Palm Width | 120 ± 10mm | 115mm | Yes |
| Thickness Finger | 0.05mm | 0.18mm | Yes |
| Thickness Palm | 0.05mm | 0.23mm | Yes |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

Freedom from Holes

Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff were tested for freedom from holes as per ASTM D5151-06, Standard Test Method for Detection of Holes in Medical Gloves and the FDA 1000ml Water Leak Test per 21 CFR 800.20. The results for all sizes (XS/S/M/L/XL) are summarized in Tables below.

TABLE 17: ASTM D5151 Freedom from Holes Test Results: PF Blue Nitrile Glove – Regular Cuff (VS311)

| Characteristic | Specification | Device Performance | Complies with Standard (Y/N) |
|--------------------|---|--|------------------------------|
| Freedom from Holes | ASTM D 5151-06 (2015) and ASTM D 6319-10, G-1, AQL 1.5 | XS/Small/Medium – 1/200 Large – 4/400 X-Large -2/200 | Y |

TABLE 18: ASTM D5151 Freedom from Holes Test Results: PF Blue Nitrile Glove-Extended Cuff (VS711)

| Characteristic | Specification | Device Performance | Complies with Standard (Y/N) |
|--------------------|---|---|------------------------------|
| Freedom from Holes | ASTM D 5151-06 (2015) and ASTM D 6319-10, G-1, AQL 1.5 | XS/Small/Medium, Large, X-Large -0/200 | Y |

Powder Content

Residual powder content for the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff was confirmed in accordance with ASTM D6124-06, Standard Test Method for Residual Powder on Medical Gloves. The average results for all sizes (XS/S/M/L/XL) are summarized below.

TABLE 19: ASTM D6124 Residual Powder Test Results: PF Blue Nitrile Glove – Regular Cuff (VS311)

| Characteristic | Specification | Device Performance | Complies with Standard (Y/N) |
|--------------------------|---------------------|--|------------------------------|
| Powder Content per Glove | Less Than 2mg/glove | X-Small – 0.16 mg/glove Small – 0.20 mg/glove Medium – 0.16 mg/glove Large – 0.16 mg/glove X-Large – 0.18 mg/glove | Y |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | | | |
|--|--|--|--|
| | | | |
|--|--|--|--|

TABLE 20: ASTM D6124 Residual Powder Test Results: PF Blue Nitrile Glove-Extended Cuff (VS711)

| Characteristic | Specification | Device Performance | Complies with Standard (Y/N) |
|--------------------------|---------------------|--|------------------------------|
| Powder Content per Glove | Less Than 2mg/glove | X-Small- 0.16 Small- 0.20 Medium- 0.16 Large- 0.18 X-Large- 0.18 | Y |

Permeation Testing

Permeation testing was conducted to support the addition of the labeling claim: Tested for use with chemotherapy drugs and Fentanyl. The gloves were tested according to ASTM D6978, Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. Minimum breakthrough times were determined for a wide range of chemotherapy drugs and Fentanyl. A summary of the minimum breakthrough times is provided in Table 21 and Table 22 (below).

Please note that the following drug, Carmustine and Thiotepa have a low permeation time. Labeling will include caution statement in the presence of these chemotherapy drugs.

TABLE 21: Chemotherapy Drugs/Fentanyl Tested – Regular Cuff (VS311XS, VS311S, VS311M, VS311L and VS311XL)

Summarized in Table 21 below are the drugs the proposed devices have been tested for use with per ASTM D6978-05 (Reapproved 2019) *Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs*.

| Chemotherapy Drug Tested | Breakthrough Time(minutes) |
|---------------------------------------|----------------------------|
| Aresenic Trioxide 1 mg/ml (1,000 ppm) | >240 |
| Azacitidine 25 mg/ml (25,000 ppm) | >240 |
| Bendamustine 5 mg/ml (5,000 ppm) | >240 |
| Bleomycin 15 mg/ml (15,000 ppm) | >240 |
| Bortezomib 1 mg/ml (1,000 ppm) | >240 |
| Busulfan 6 mg/ml (6,000 ppm) | >240 |
| Carboplatin 10.0 mg/ml (10,000 ppm) | >240 |



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

| | |
|--|----------|
| Carfilzomib 2.0 mg/ml (2,000) | >240 |
| Carmustine (BCNU) 3.3 mg/ml (3,300 ppm) | 33.1 min |
| Cetuximab 2.0 mg/ml (2,000) | >240 |
| Cisplatin 1.0 mg/ml (1,000 ppm) | >240 |
| Cyclophosphamide (Cytosan) 20.0 mg/ml (20,000 ppm) | >240 |
| Cytarabine 100 mg/ml (100,000 ppm) | >240 |
| Cytovene 10 mg/ml (10,000 ppm) | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml (10,000 ppm) | >240 |
| Daunorubicin 5 mg/ml (5,000 ppm) | >240 |
| Decitabine 5 mg/ml (5,000 ppm) | >240 |
| Docetaxel 10.0 mg/ml (10,000 ppm) | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml (2,000 ppm) | >240 |
| Epirubicin (Ellence) 2.0 mg/ml (2,000 ppm) | >240 |
| Etoposide (Toposar) 20.0 mg/ml (20,000 ppm) | >240 |
| Fludarabine 25 mg/ml (25,000 ppm) | >240 |
| Fluorouracil 50.0 mg/ml (50,000 ppm) | >240 |
| Fulvestrant 50.0 mg/ml (50,000 ppm) | >240 |
| Gemcitabine (Gemzar) 38 mg/ml (38,000 ppm) | >240 |
| Idarubicin 1 mg/ml (1,000 ppm) | >240 |
| Ifosfamide 50.0 mg/ml (50,000 ppm) | >240 |
| Irinotecan 20.0 mg/ml (20,000 ppm) | >240 |
| Mechlorethamine HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Melphalan 5 mg/ml (5,000 ppm) | >240 |
| Methotrexate 25 mg/ml (25,000 ppm) | >240 |
| Mitomycin C 0.5 mg/ml (500 ppm) | >240 |
| Mitoxantrone 2.0 mg/ml (2,000 ppm) | >240 |
| Oxaliplatin 2.0 mg/ml (2,000 ppm) | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml (6,000 ppm) | >240 |
| Paraplatin 10 mg/ml (10,000 ppm) | >240 |
| Pemetrexed Disodium 25.0 mg/ml (25,000 ppm) | >240 |
| Pertuzumab 30 mg/ml (30,000 ppm) | >240 |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | |
|---|---|
| Raltitrexed 0.5 mg/ml (500 ppm) | >240 |
| Retrovir 10.0 mg/ml (10,000 ppm) | >240 |
| Rituximab 10 mg/ml (10,000 ppm) | >240 |
| Temsirolimus 25.0 mg/ml (25,000 ppm) | >240 |
| Thiotepa 10.0 mg/ml (10,000 ppm) | 69.2 min |
| Topotecan HCl 1.0 mg/ml (1,000 ppm) | >240 |
| Trastuzumab 21.0 mg/ml (21,000 ppm) | >240 |
| Trisonex 1 mg/ml (1,000 ppm) | >240 |
| Vinblastine 1 mg/ml (1,000 ppm) | >240 |
| Vincristine Sulfate 1.0 mg/ml (1,000 ppm) | >240 |
| Vinorelbine 10 mg/ml (10,000 ppm) | >240 |
| Zoledronic Acid 0.8 mg/ml (800 ppm) | >240 |
| Fentanyl Tested | Breakthrough Time (in minutes) |
| Fentanyl Citrate, 100mcg/2mL | >240 |

TABLE 22: Chemotherapy/Fentanyl Drugs Tested – Extended Cuff (VS311XS, VS311S, VS311M, VS311L and VS311XL)

Summarized in 18 below, are the drugs the proposed devices have been tested for use with per ASTM D6978-05 (Reapproved 2019) *Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs*.

| Chemotherapy Drug Tested | Breakthrough Time (in minutes) |
|---------------------------------------|---|
| Arsenic Trioxide 1 mg/ml | >240 |
| Azacitidine 25 mg/ml | >240 |
| Bendamustine 5 mg/ml | >240 |
| Bleomycin 15 mg/m | >240 |
| Bortezomib 1 mg/ml | >240 |
| Busulfan 6 mg/ml | >240 |
| Carfilzomib 2 mg/ml | >240 |
| Carmustine (BCNU) 3.3 mg/ml | 59.4 min |
| Cetuximab | >240 |
| Cisplatin 1.0 mg/ml | >240 |
| Cyclophosphamide (Cytoxan) 20.0 mg/ml | >240 |
| Cytarabine Hydrochloride 100 mg/ml | >240 |
| Cytovene 10 mg/ml | >240 |
| Dacarbazine (DTIC) 10.0 mg/ml | >240 |



Medline Industries, Inc.
 Three Lakes Drive
 Northfield, IL 60093

| | |
|-------------------------------------|---|
| Daunorubicin 5 mg/ml | >240 |
| Decitabine 5 mg/ml | >240 |
| Docetaxel 10.0 mg/ml | >240 |
| Doxorubicin Hydrochloride 2.0 mg/ml | >240 |
| Epirubicin (Ellence) 2.0 mg/ml | >240 |
| Etoposide (Toposar) 20.0 mg/ml | >240 |
| Fludarabine 25 mg/ml | >240 |
| Fluorouracil 50.0 mg/ml | >240 |
| Fulvestrant 50.0 mg/ml | >240 |
| Gemcitabine (Gemzar) 38 mg/ml | >240 |
| Idarubicin 1 mg/ml | >240 |
| Irinotecan 20.0 mg/ml | >240 |
| Mechlorethamine HCl 1.0 mg/ml | >240 |
| Melphalan 5 mg/ml | >240 |
| Methotrexate 25.0 mg/ml | >240 |
| Mitomycin C 0.5 mg/ml | >240 |
| Oxaliplatin 5.0 mg/ml | >240 |
| Paclitaxel (Taxol) 6.0 mg/ml | >240 |
| Paraplatin (Carboplatin) 10 mg/ml | >240 |
| Pemetrexed Disodium 25.0 mg/ml | >240 |
| Pertuzumab 30 mg/ml | >240 |
| Raltitrexed 0.5 mg/ml | >240 |
| Retrovir 10.0 mg/ml | >240 |
| Rituximab 10 mg/ml | >240 |
| Temsirolimus 25.0 mg/ml | >240 |
| Thiotepa 10.0 mg/ml | 68.2 min |
| Topotecan HCl 1.0 mg/ml | >240 |
| Trastuzumab 21.0 mg/ml | >240 |
| Trisonex 1 mg/ml | >240 |
| Vinblastine 1 mg/ml | >240 |
| Vincristine Sulfate 1.0 mg/ml | >240 |
| Vinorelbine 10 mg/ml | >240 |
| Zoledronic Acid 0.8 mg/ml | >240 |
| Fentanyl Tested | Breakthrough Time (in minutes) |
| Fentanyl Citrate, 100mcg/2mL | >240 |



Medline Industries, Inc.
Three Lakes Drive
Northfield, IL 60093

Conclusion

In accordance with 21 CFR Part 807, and based on the information provided in this premarket notification, Medline Industries, Inc. concludes that The Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) – Regular Cuff and the Medline Powder-Free Blue Nitrile Examination Gloves (Tested for use with Chemotherapy Drugs and Fentanyl) - Extended Cuff are as safe, as effective, and performs as well as or better than the legally marketed predicate device, the YTY Industry (Manjung) Non-Sterile, Powder-Free Nitrile Examination Gloves (Cobalt Blue) Tested for use with Chemotherapy Drugs – K111248.