



June 9, 2020

Beijing Superlaser Technology Co., Ltd.  
% Ray Wang  
General Manager  
Beijing Believe-Med Technology Service Co., Ltd  
Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd.,  
FangShan District  
Beijing, 102401 CN

Re: K200042

Trade/Device Name: CO2 Laser System

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In  
Dermatology

Regulatory Class: Class II

Product Code: GEX

Dated: May 8, 2020

Received: May 12, 2020

Dear Ray Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Colin Kejing Chen  
Acting Assistant Director  
DHT4A: Division of General Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K200042

Device Name

CO2 Laser System

Indications for Use (Describe)

The CO2 Laser System is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

This 510(k) Summary is submitted in accordance with the requirements of 21 CFR Section 807.92.

1. Date of Preparation  
05/07/2020
2. Applicant Name and Address  
Beijing Superlaser Technology Co., Ltd.  
No.2, Zhongfu Street, Economic and Technological Industrial zone,  
Xihongmen Tower, Daxing District, Beijing, China.
3. Contact Person Information  
Shi Shuang  
Registration Specialist  
Tel: 86-10-81284899 to 806  
Fax: 86-10-81284899  
Email: 672257488@qq.com
4. Submission Correspondent  
Mr. Ray Wang  
Beijing Believe-Med Technology Service Co., Ltd.  
Rm.912, Building #15, XiYueHui, No.5, YiHe North Rd.,  
FangShan District, Beijing City, China, 102401  
Tel: +86-18910677558  
Fax: +86-10-56335780  
[Email: ray.wang@believe-med.com](mailto:ray.wang@believe-med.com)
5. Identification of Proposed Device  
Trade Name: CO2 Laser System  
Common Name: Powered Laser Surgical Instrument  
Model(s): SL-LC01  
  
Classification Name: Powered Laser Surgical Instrument  
Class: II  
Product Code: GEX  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology  
Review Panel: General & Plastic Surgery
6. Identification of the Primary Predicate  
510(k) Number: K192528  
Device Name: CO2 Laser Therapy System  
Manufacturer: Shangdong Huamei Technology Co., Ltd.
7. Device Description

The device emits CO2 laser at the wavelength of 10.6µm, which is the spectral line in the far infrared range. The device is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.

8. Indications for Use

The CO2 Laser System is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.

9. Substantially Equivalent (SE) Comparison

Tab 1 General Comparison

ITEM	Proposed Device	Predicate Device K192528	Remark
Product Code	GEX	GEX	SE
Regulation No.	21 CFR 878.4810	21 CFR 878.4810	SE
Class	2	2	SE
Where used	hospital	hospital	SE
Intended Use	The CO2 Laser System is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.	The CO2 Laser Therapy Machine is used for human tissue vaporization, coagulation in dermatology and plastic surgery, general surgery, gynecology, podiatry, dental and otorhinolaryngology.	SE

Tab 2 Performance Comparison

ITEM	Proposed Device		Predicate Device K192528		Remark
Maximum Power	30W		30W		SAME
Work Mode	Surgery (Single Pulse, Continuous, Pulse)		Surgery (Single Pulse, Continuous, Muti-Pulse)		SAME
Wavelength	10.6 um		10.6 um		SAME
Beam Delivery	7 joint light guide arm		7 knuckle arm key joints light arm		SAME
Aiming Beam	Red indicator light (650nm, ≤5 mW)		630-650nm red diode laser (≤5 mW)		SAME
Spot Size	0.5 mm		0.5 mm		SAME
Output Power	pulse,	1-30W	Single Pulse	1-30W	SAME
	single	1-30W	Muti-Pulse	1-30W	
	continuous	1-30W	Continuous	1-30W	
Pulse Duration	Single Pulse Mode: 1-1000 ms Pulse mode: 1-1000 ms		Single Pulse Mode: 0.1-1000 ms Muti-Pulse Mode: 0.1-1000 ms		SIMILAR
Control System	Touch screen, footswitch		Touch screen, footswitch		SAME
Laser Operation	Footswitch		Footswitch		SAME
Laser medium/energy source	CO2		CO2		SAME
Cooling System	Closed inner circulating water cooling		Air cooling		SIMILAR

<b>Cleaning Method</b>	70% isopropyl alcohol	70% medical alcohol	SAME
<b>Dimension</b>	37.5 cm x 29 cm x 113 cm	66*42*125cm(without light arm)	Analysis
<b>Weight</b>	40kg	80 kg	Analysis
<b>Power input</b>	110V 60Hz or 230V 50Hz	AC 110V/60Hz ;	SIMILAR

Tab 3 Safety Comparison

<b>Item</b>	<b>Proposed Device</b>	<b>Predicate Device K161925</b>	<b>Remark</b>
EMC, Electrical and Laser Safety			
Electrical Safety	Comply with IEC 60601-1, IEC 60601-2-22	Comply with IEC 60601-1, IEC 60601-2-22	SE
EMC	Comply with IEC 60601-1-2	Comply with IEC 60601-1-2	SE
Laser Safety	Comply with IEC 60601-2-22, IEC 60825	Comply with IEC 60601-2-22, IEC 60825	SE

### Analysis

The difference between proposed device and predicate device lies in the appearance (dimension, weight) and pulse duration. The difference will not affect the safety and effectiveness of proposed device in comparison to the predicate.

### 10. Non-Clinical Testing

Non clinical tests were conducted to verify that the proposed device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The test results demonstrated that the proposed device complies with the following standards:

- IEC 60601-1:2005/A1:2012 Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance;
- IEC 60601-2-22:2012, Medical Electrical Equipment - Part 2-22: Particular Requirements For Basic Safety And Essential Performance Of Surgical, Cosmetic, Therapeutic And Diagnostic Laser Equipment;
- IEC 60825-1: 2014, Safety of laser products - Part 1: Equipment classification and requirements.
- IEC 60601-1-2:2014, Medical electrical equipment- Part 1-2: General requirements for basic safety and essential performance- Collateral standard: Electromagnetic compatibility- Requirements and tests.

In Addition, the following non-clinical tests were performed to make sure that the device performs as intended:

- Software Validation & Verification Test
- Bench Performance Tests

### 11. Clinical Testing

No clinical study is performed to support substantial equivalence.

### 12. Conclusion

Based on the comparison and analysis above, the proposed device is determined to be Substantially Equivalent (SE) to the predicate device.