



May 19, 2021

Curexo, Inc.  
% Do Hyun Kim  
CEO  
BT Solutions, Inc.  
Unit 904, Eonju-ro 86-gil 5, Gangnam-gu  
Seoul, 06210  
Republic of Korea

Re: K201569  
Trade/Device Name: CUVIS-spine  
Regulation Number: 21 CFR 882.4560  
Regulation Name: Stereotaxic Instrument  
Regulatory Class: Class II  
Product Code: OLO  
Dated: April 16, 2021  
Received: April 19, 2021

Dear Do Hyun Kim:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For: Shumaya Ali, M.P.H.  
Assistant Director  
DHT6C: Division of Restorative, Repair, and Trauma  
Devices  
OHT6: Office of Orthopedic Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K201569

Device Name

CUVIS-spine

Indications for Use (Describe)

CUVIS-spine is intended for use as an aid for precisely locating anatomical structures and for the spatial positioning and orientation of guide bush to be used by surgeons for navigating and/or guiding compatible surgical instruments in open or percutaneous pedicle screw placement provided that the required markers and rigid patient anatomy can be identified on O-arm or C-arm.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

### General Information

Submitter/Applicant: CUREXO, INC.  
Address: 577, Gangnam-Daero, Seocho-Gu  
Seoul, 06530, Republic of Korea  
Tel: +82-31-788-7935  
Jungeun PARK  
Quality Management Representative  
Email: [jepark@curexo.com](mailto:jepark@curexo.com)

Applicant correspondent: Do Hyun Kim, BT Solutions, Inc.  
Address: Unit 904, Eonju-ro 86gil 5, Gangnam-gu  
Seoul 06210, Korea.  
Tel: +82-2-538-9140  
Email: [ceo@btsolutions.co.kr](mailto:ceo@btsolutions.co.kr)

Preparation Date: May 20, 2021

### Device Name and Code

Device Trade Name: CUVIS-spine  
Common Name: Pedicle Screw Guide System  
Classification Name: Stereotaxic instrument  
Product Code: OLO  
Regulation Number: 21 CFR 882.4560  
Classification: II  
Review Panel: Orthopedic

### Predicate Device

CUVIS-spine is substantially equivalent to the following legally marketed predicate devices

Table 1 Primary Predicate device

Applicant	Device Name	510(k) Number
Globus Medical Inc.	Excelsius GPS™	K171651

### **Device Description**

The CUVIS-spine is a mobile system mainly comprising the robotic arm, the main console and the staff console as an option. The robotic arm is positioned on the floor near the side of the surgical table. The location of the main console or the staff console is appropriately determined considering the user preference and the environments.

The CUVIS-spine is a pedicle screw guide system which consists of Robotic Arm, Main Console, Staff Console, Guide bush, Source Calibrator, Registration Tool, Registration Tool Adapter, Robotic Arm Drape, Tool Drape, Detector Drape, Marker Ball, Patient Marker, Marker Driver, Detector Calibrator, Dilator, Serration-tip Dilator, Drill Bit, Tapper, Stylet Tapper, Screwdriver, Instrument Container.

### **Indications / Intended Use**

CUVIS-spine is intended for use as an aid for precisely locating anatomical structures and for the spatial positioning and orientation of guide bush to be used by surgeons for navigating and/or guiding compatible surgical instruments in open or percutaneous pedicle screw placement provided that the required markers and rigid patient anatomy can be identified on O-arm or C-arm.

### **Technical Characteristics in Comparison to Predicate Devices**

CUVIS-spine is substantially equivalent to the following legally marketed predicate devices

	<b>Subject Device</b>	<b>Predicate Device</b>
	<b>CUVIS-spine(K201569)</b>	<b>Excelsius GPS™ (K171651)</b>
<b>Manufacturer</b>	CUREXO, INC.	Globus Medical Inc.
<b>Classification</b>	Class II	Class II
<b>Product Code</b>	OLO	OLO
<b>Regulation No.</b>	882.4560	882.4560

CUVIS-spine  
510(k) Summary

	<b>Subject Device</b>	<b>Predicate Device</b>
	<b>CUVIS-spine(K201569)</b>	<b>Excelsius GPS™ (K171651)</b>
<b>General Device description</b>	The CUVIS-spine is a mobile system mainly comprising the robotic arm, the main console and the staff console as an option. The robotic arm is positioned on the floor near the side of the surgical table. The location of the main console or the staff console is appropriately determined considering the user preference and the environments.	The EXCELSIUS GPS™ is a Robotic Positioning System that includes a computer controlled robotic arm, hardware, and software that enables real time surgical navigation and robotic guidance using radiological patient images(preoperative CT, intraoperative CT and fluoroscopy), using a dynamic reference base and positioning camera.
<b>Indications for use</b>	CUVIS-spine is intended for use as an aid for precisely locating anatomical structures and for the spatial positioning and orientation of guide bush to be used by surgeons for navigating and/or guiding compatible surgical instruments in open or percutaneous pedicle screw placement provided that the required markers and rigid patient anatomy can be identified on O-arm or C-arm.	The EXCELSIUS GPS™ is intended for use as an aid for precisely locating anatomical structures and for the spatial positioning and orientation of an instrument holder or guide tube to be used by surgeons for navigating and/or guiding compatible surgical instruments in open or percutaneous procedures provided that the required fiducial markers and rigid patient anatomy can be identified on CT scans or fluoroscopy. The system is indicated for the placement of spinal and orthopaedic bone screws.
<b>Technical Characteristics</b>		
<b>Principle of operation</b>	<ul style="list-style-type: none"> <li>- Intraoperative images</li> <li>- Patient registration</li> <li>- Surgical planning</li> <li>- Real-time tracking of navigated instruments</li> <li>- Guidance of instruments</li> </ul>	<ul style="list-style-type: none"> <li>- Intraoperative/preoperative images</li> <li>- Patient registration</li> <li>- Surgical planning</li> <li>- Real-time tracking of navigated instruments</li> <li>- Guidance of instruments</li> </ul>
<b>Image</b>	O-arm, C-arm	O-arm, C-arm, CT

CUVIS-spine  
510(k) Summary

	<b>Subject Device</b>	<b>Predicate Device</b>
	<b>CUVIS-spine(K201569)</b>	<b>Excelsius GPS™ (K171651)</b>
<b>Input Images</b>	- 3D intra-operative exam - 2D intra-operative exam	- 3D pre-operative exam - 3D intra-operative exam - 2D intra-operative exam
<b>Tracker</b>	Optical Tracking System	Optical Tracking System
<b>Guide</b>	Dilator and Tapper	Dilator and Tapper
<b>Target Tracking</b>	YES	YES
<b>Integrated Software</b>	- SRC (Control Software) - SPN (Planner Software)	Excelsius GPS Planning and Navigation Application Software
<b>Save/load Planning</b>	YES	YES
<b>Merge images functionality</b>	YES	YES
<b>Trajectory planning parameters</b>	- Entry point - Target point - Instrument length/diameter	- Entry point - Target point - Instrument length/diameter
<b>Localization means</b>	Optical system (infrared camera)	Optical system (infrared camera)
<b>Image-guided</b>	YES	YES
<b>Controller</b>	Forced-controlled movement allowing robot arm positioning (called hand guide function)	Force-controlled movement allowing robotic arm positioning
<b>Patient registration method</b>	- Intra-op CT: Registration tool - Fluoroscopy: Source calibrator	- Pre-op CT: Fluoroscopic to pre-op CT merge - Intra-op CT: Registration fixture - Fluoroscopy: Registration fixture
<b>Real time display of instrument position</b>	YES	YES
<b>Accessories</b>	- Registration instruments - Patient reference instruments - Surgical instruments	- Registration instruments - Patient reference instruments - Surgical instruments - End effector
<b>Performance data</b>		
<b>Electrical Safety and Electromagenti</b>	- IEC60601-1 - IEC60601-1-2	- IEC60601-1 - IEC60601-1-2

CUVIS-spine  
510(k) Summary

	<b>Subject Device</b>	<b>Predicate Device</b>
	<b>CUVIS-spine(K201569)</b>	<b>Excelsius GPS™ (K171651)</b>
<b>c compatibility</b>		
<b>Accuracy verification on anatomical landmarks</b>	Yes	Yes
<b>Biocompatibility</b>	The biocompatibility evaluation for the patient marker has been conducted in accordance with EN ISO 10993 standards.	The biocompatibility evaluation for EXCELSIUS GPS™ has been conducted in accordance with ISO 10993 standards.
<b>Parts of Contact with Patient</b>	Patient marker	Patient reference instruments (Patient marker)
<b>Nature of Body Contact</b>	Implant device / bone	Implant device / bone
<b>Patient Fixation</b>	Reference is fixed to patient's bony structure for tracking system	Reference is fixed to patient's bony structure for tracking system

**Performance Data**

Non-clinical tests: Pose accuracy and Repeatability of the CUVIS-spine were tested and validated. A cadaveric study of the proposed device had been performed for the robotic-assisted pedicle screw placement.

Biocompatibility were tested using following consensus standards:

- Tests for in vitro cytotoxicity were tested and evaluated according to the FDA-recognized consensus standard, ISO 10993-5.
- Tests for irritation, skin sensitization and intracutaneous reactivity were tested and evaluated according to the FDA-recognized consensus standard, ISO 10993-10.

Electromagnetic compatibility and electrical safety, etc, were tested using following consensus standards:

- Basic safety and essential performance of the CUVIS-spine is tested and evaluated according to the FDA-recognized consensus standard, ES 60601-1.
- Effect to the device by electromagnetic disturbances were tested and evaluated according to the FDA-recognized consensus standard IEC 60601-1-2.
- Risk management was recorded by referring to ISO 14971.



- Usability was documented by referring to IEC 60601-1-6.

### **Substantial Equivalence**

CUVIS-spine is not based on a new technology. The technological characteristics of the subject device is comparable to the predicate device for comparable indications for use. Thus, subject device CUVIS-spine is concluded to be substantially equivalent to the predicate device.

### **Conclusions**

On the basis of the information provided in this Summary, CUREXO, INC. believes that CUVIS-spine is substantially equivalent to legally commercialized predicate devices for the purposes of this 510 (k) submission.