



November 22, 2022

Beiang Air Tech LTD.  
% Tyra Chiu  
Regulatory Specialist  
Medical Wizdom, LLC  
12F.-4, No. 81, Sec. 2, Chang'an E. Rd., Zhongshan Dist.  
Taipei, 10491  
Taiwan

Re: K211507

Trade/Device Name: Airdog X5 Recirculating Air Cleaner (model KJ300F-X5)  
Regulation Number: 21 CFR 880.5045  
Regulation Name: Medical recirculating air cleaner  
Regulatory Class: Class II  
Product Code: FRF  
Dated: October 21, 2022  
Received: October 21, 2022

Dear Tyra Chiu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
**Clarence W. Murray III -S**

Clarence W. Murray, III, PhD  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K211507

Device Name  
Airdog X5 Recirculating Air Cleaner (model KJ300F-X5)

### Indications for Use (Describe)

The Airdog is a mobile air cleaner intended to be used to remove particles from the air for medical purposes. The device is intended for indoor use only.

The Airdog has been demonstrated to effectively inactivate H3N2 and reduce Staphylococcus albus by 4 log with L4 speed.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary

**Date** November 22, 2022

**Manufacturer/** BeiAng Air Tech Ltd.

**510(k) Owner** 175#, SONGBEI ROAD, SUZHOU INDUSTRIAL PARK,  
JIANGSU PROVINCE, 215000, CHINA

**Contact Person** Yan Zhang

Phone: +86-19951234257

E-mail: yan@beiantech.com

**Device Trade Name** Airdog X5 Recirculating Air Cleaner (model KJ300F-X5)

**Common Name** Medical Recirculating Air Cleaner

**Classification Name** Cleaner, Air, Medical Recirculating

**Device Class** II

**Review Panel** General Hospital

**Regelation Number** 880.5045

**Product Code** FRF

**Device Description  
and Technology  
Characteristics**

The Airdog X5 Recirculating Air Cleaner is a floor-standing air cleaner with an air quality sensor. The air from the room enters the X5 and passes a pre-filter, ionization frame, ionization field, collecting plates and carbon layer that captures particular matter and virus from the air. The X5 contains multiple PCB modules, air quality sensor, device status indicators, power switch and a display panel, it allows user to select between different operation parameter, including auto-mode, sleep mode and wind speed. The Ionization Wireframe and the Collecting Plate can be cleaned routinely and are reusable. X5 is powered from an AC wall outlet. It is intended to be used indoors only.

**Models** KJ300F-X5

**Indications for Use** The Airdog is a mobile air cleaner intended to be used to remove particles from the air for medical purposes. The device is intended for indoor use only.

The Airdog has been demonstrated to effectively inactivate *H3N2* and reduce *Staphylococcus albus* by 4 log with L4 speed.

**Predicate Device(s)** K200321  
Novaerus NV1050/ Novaerus US Inc

### Summary of Comparison and Technological Characteristics

	Proposed Device	Predicate Device	Comparison
<b>Device Name</b>	Airdog X5 Recirculating Air Cleaner(model KJ300F-X5)	Novaerus NV1050	-
510(k) #	K211507	K200321	-
Applicant	BeiAng Air Tech Ltd.	Novaerus US Inc	-
Product code	FRF (21 CFR 880.5045)	FRF (21 CFR 880.5045)	Same
Classification	II	II	Same
OTC use	YES	YES	Same
Intended Use	<p>The Airdog is a mobile air cleaner intended to be used to remove particles from the air for medical purposes. The device is intended for indoor use only.</p> <p>The Airdog has been demonstrated to effectively inactivate <i>H3N2</i> and reduce <i>Staphylococcus albus</i> by 4 log with L4 speed.</p>	<p>The Novaerus NV 1050 is intended as a room recirculating air cleaner. The system is used for filtering out and inactivating airborne particles form the air for medical purpose</p>	<p>Same.</p> <p>Both devices are room recirculating air cleaners for medical purpose.</p>
Use Environment	Indoor	Indoor	Same

Technology	Air from the room is passed through a plasma ion field to neutralize airborne micro-organisms. A negative charged dust collecting electrode traps the resulting debris and a third carbon	Air from the room is passed through a plasma field to inactivate airborne micro-organisms. A HEPA filter traps the resulting debris and an activated carbon filter absorbs any ozone	Similar. Both have the plasma field and carbon filter. The predicate device use HEPA filter to traps the resulting debris while the
	stage absorbs any ozone generated as a byproduct. (ETL test results also show that even if the carbon net is removed, the ozone is lower than the regulatory requirements)	generated as a byproduct of the plasma field	subject device uses a negative charged dust collecting electrode.
Power Source	100-240V~50-60Hz	110V AC	Different
Weight(kg)	10.7	51	Different
Dimension	25.59 inch (H)* 12.04 inch (W) *12.44 inch (L)	36.5 inch (H) x 19.0 inch (W) x 19.1inch (L)	Different
Reduction of biological agents	H3N2 Influenza virus inactivating by 99.99% with operation with L4 speed Staphylococcus albus reduced by 99.99% in 1 hour	Bacillus Globigii endospores and MS2 phage reduced by 3 log reduction in 10 minutes and 4 log reduction in 15 minutes when operating at full fan speed	Similar. Both are capable of reducing specified biological agents by 4 log.
Filtration of particles	produces a 4 log reduction in PM <sub>2.5</sub> particles in 120 minutes in a 30 m <sup>3</sup> chamber	4 log reduction in 0.5 to 2.0 $\mu$ m sized particles in 10 minutes in a 580ft <sup>3</sup> (16.4m <sup>3</sup> ) sealed room	Similar. Both are capable of reducing particles by 4 log
Ozone emitted	Meet the requirements of ul867 (<50ppb) for ozone and USA CRB	below 10 ppb	Similar. Both <50ppb
Standards used	IEC 60601-1: 2005/A1:2012 IEC 60601-1-2: 2014 IEC 62304 ed.1.1 ISO 14971	IEC 60601-1: 2005/A1:2012	Different.

**Discussion on Performance Data Non-Clinical Tests** Bench tests show that the X5 produces 4-log reduction of particles, 4-log elimination ratio of total bacteria counts, a 4-log inactivation in H3N2 Influenza virus and produce ozone emission which meets the requirements of ul867 for ozone and USA CRB.

Software verification and validation was conducted, and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices."

**Discussion on Clinical**

**Test Performed** Not applicable

**Conclusion** The conclusions drawn from the nonclinical tests demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device.