



January 20, 2022

KARL STORZ Endoscopy-America, Inc.
Mario Trujillo
Regulatory Affairs Specialist
2151 E. Grand Avenue
El Segundo, CA 90245

Re: K212458
Trade/Device Name: SSU Flexible Video-Uretero-Renoscope System
Regulation Number: 21 CFR§ 876.1500
Regulation Name: Endoscope and Accessories
Regulatory Class: II
Product Code: FGB
Dated: December 20, 2021
Received: December 21, 2021

Dear Mario Trujillo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Mark J. Antonino, M.S.
Assistant Director
DHT3B: Division of Reproductive,
Gynecology and Urology Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K212458

Device Name

SSU Flexible Video-Uretero-Renoscope System

Indications for Use (Describe)

The SSU Flexible Video-Uretero-Renoscope System is indicated for endoscopic examination in the urinary tract and can be used to examine the interior of the kidney, and using additional accessories, to perform various diagnostic and therapeutic procedures.

E-Box: the product serves as an adaptor for operating the flexible single-use videoscope on the compatible CCU.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

This 510(k) Summary is being submitted in accordance with the requirements of the Safe Medical Devices Act (SMDA) of 1990 and 21 CFR 807.92 and the FDA guidance document titled “The 510(k) Program: Evaluating Substantial Equivalence in Premarket Notifications [510(k)]” issued on July 28, 2014. All data included in this document is accurate and complete to the best of KARL STORZ SE & Co. KG knowledge.

Submitter:	KARL STORZ SE & Co. KG Dr.-Karl-Storz-Straße 34 78532 Tuttlingen, Germany
Contact:	Mario Trujillo Associate Regulatory Affairs Specialist Tel.: (424) 218-8481 Email: Mario.Trujillo@karlstorz.com
Date of Preparation:	September 7, 2021
Type of 510(k) Submission:	Traditional
Device Identification:	Trade Name: SSU Flexible Video-Uretero-Renoscope System Classification Name: Endoscope and Accessories (21 CFR Part 876.1500);
Regulatory Class:	II
Product Code:	FGB
Guidance Document:	Not Applicable
Predicate Device:	KARL STORZ Flexible Video-Uretero-Renoscope System (K141250).
Device Description:	The videoscopes in the modified Flexible Video-Uretero-Renoscope SSU System are sterile single-use, flexible video-endoscopes. The distal tip houses the CMOS (Complementary Metal Oxide Semiconductor) imaging sensor and the LED light source. The raw data captured at the distal tip CMOS imaging sensor is transferred to the E-Box adaptor, where it is converted to a standard NTSC (National Television System Committee) video signal by the PCB (Printed Circuit Board), which is then driven into one of the CCUs (C-MAC, C-HUB II, or X-LINK + IMAGE 1S Connect) for further processing and video formatting for output to a display monitor. The videoscopes and E-Box are powered by the CCUs through the connecting cords.
Intended Use:	The SSU Flexible Video-Uretero-Renoscope is intended for visualization purposes during urological procedures.

<p>Indications For Use:</p>	<p>The SSU Flexible Video-Uretero-Renoscope is indicated for endoscopic examination in the urinary tract and can be used to examine the interior of the kidney, and using additional accessories, to perform various diagnostic and therapeutic procedures.</p> <p>E-Box: the product serves as an adaptor for operating the flexible single-use videoscope on the compatible CCU.</p>																																																										
<p>Technological Characteristics:</p>	<p>Comparison Table: Subject vs. Predicate Devices</p> <table border="1" data-bbox="574 548 1190 1633"> <thead> <tr> <th data-bbox="574 548 792 695"></th> <th data-bbox="792 548 972 695">Subject Device Flexible Video-Uretero-Renoscope SSU System</th> <th data-bbox="972 548 1190 695">Predicate Device, K141250 Flexible Video-Uretero-Renoscope System</th> </tr> </thead> <tbody> <tr> <td data-bbox="574 695 792 751">Maximal Outer diameter Insertion Portion</td> <td data-bbox="792 695 972 751">3.2 mm</td> <td data-bbox="972 695 1190 751">Same as subject device</td> </tr> <tr> <td data-bbox="574 751 792 808">Outer diameter Insertion Tube</td> <td data-bbox="792 751 972 808">2.9 mm</td> <td data-bbox="972 751 1190 808">Same as subject device</td> </tr> <tr> <td data-bbox="574 808 792 846">Outer diameter Distal End</td> <td data-bbox="792 808 972 846">3.2 mm</td> <td data-bbox="972 808 1190 846">Same as subject device</td> </tr> <tr> <td data-bbox="574 846 792 884">Insertion portion length</td> <td data-bbox="792 846 972 884">700 mm</td> <td data-bbox="972 846 1190 884">675 mm</td> </tr> <tr> <td data-bbox="574 884 792 921">Working channel</td> <td data-bbox="792 884 972 921">Present</td> <td data-bbox="972 884 1190 921">Same as subject device</td> </tr> <tr> <td data-bbox="574 921 792 978">Inner diameter Working Channel</td> <td data-bbox="792 921 972 978">1.2 mm</td> <td data-bbox="972 921 1190 978">Same as subject device</td> </tr> <tr> <td data-bbox="574 978 792 1016">Tip deflection up/down</td> <td data-bbox="792 978 972 1016">270°/270°</td> <td data-bbox="972 978 1190 1016">Same as subject device</td> </tr> <tr> <td data-bbox="574 1016 792 1054">Field of view</td> <td data-bbox="792 1016 972 1054">110°</td> <td data-bbox="972 1016 1190 1054">90°</td> </tr> <tr> <td data-bbox="574 1054 792 1092">Direction of View</td> <td data-bbox="792 1054 972 1092">0°</td> <td data-bbox="972 1054 1190 1092">Same as subject device</td> </tr> <tr> <td data-bbox="574 1092 792 1129">Depth of Field</td> <td data-bbox="792 1092 972 1129">5 - 50 mm</td> <td data-bbox="972 1092 1190 1129">4 - 60mm</td> </tr> <tr> <td data-bbox="574 1129 792 1251">On-axis Resolution</td> <td data-bbox="792 1129 972 1251">12.5 Lp/mm at 5 mm 4.5 Lp/mm at 15 mm 1.25 Lp/mm at 50 mm</td> <td data-bbox="972 1129 1190 1251">11.0 Lp/mm at 4 mm 4.0 Lp/mm at 12 mm 1.0 Lp/mm at 60 mm</td> </tr> <tr> <td data-bbox="574 1251 792 1289">Chip type</td> <td data-bbox="792 1251 972 1289">CMOS</td> <td data-bbox="972 1251 1190 1289">Same as subject device</td> </tr> <tr> <td data-bbox="574 1289 792 1327">Chip location</td> <td data-bbox="792 1289 972 1327">Distal</td> <td data-bbox="972 1289 1190 1327">Same as subject device</td> </tr> <tr> <td data-bbox="574 1327 792 1365">Illumination source</td> <td data-bbox="792 1327 972 1365">LED</td> <td data-bbox="972 1327 1190 1365">Same as subject device</td> </tr> <tr> <td data-bbox="574 1365 792 1476">Compatible CCU</td> <td data-bbox="792 1365 972 1476">C-MAC C-HUB II Image 1S</td> <td data-bbox="972 1365 1190 1476">Image 1S</td> </tr> <tr> <td data-bbox="574 1476 792 1514">How device is provided</td> <td data-bbox="792 1476 972 1514">Sterile single-use</td> <td data-bbox="972 1476 1190 1514">Unsterile, reusable</td> </tr> <tr> <td data-bbox="574 1514 792 1593">EO Sterilization cycle</td> <td data-bbox="792 1514 972 1593">EO, Overpressure 2.7 bar absolute, 8.5 % ETO in 91.5 % CO2</td> <td data-bbox="972 1514 1190 1593">N/A</td> </tr> <tr> <td data-bbox="574 1593 792 1633">Sterilizing Agent</td> <td data-bbox="792 1593 972 1633">Ethylene Oxide (EO)</td> <td data-bbox="972 1593 1190 1633">N/A</td> </tr> </tbody> </table>			Subject Device Flexible Video-Uretero-Renoscope SSU System	Predicate Device, K141250 Flexible Video-Uretero-Renoscope System	Maximal Outer diameter Insertion Portion	3.2 mm	Same as subject device	Outer diameter Insertion Tube	2.9 mm	Same as subject device	Outer diameter Distal End	3.2 mm	Same as subject device	Insertion portion length	700 mm	675 mm	Working channel	Present	Same as subject device	Inner diameter Working Channel	1.2 mm	Same as subject device	Tip deflection up/down	270°/270°	Same as subject device	Field of view	110°	90°	Direction of View	0°	Same as subject device	Depth of Field	5 - 50 mm	4 - 60mm	On-axis Resolution	12.5 Lp/mm at 5 mm 4.5 Lp/mm at 15 mm 1.25 Lp/mm at 50 mm	11.0 Lp/mm at 4 mm 4.0 Lp/mm at 12 mm 1.0 Lp/mm at 60 mm	Chip type	CMOS	Same as subject device	Chip location	Distal	Same as subject device	Illumination source	LED	Same as subject device	Compatible CCU	C-MAC C-HUB II Image 1S	Image 1S	How device is provided	Sterile single-use	Unsterile, reusable	EO Sterilization cycle	EO, Overpressure 2.7 bar absolute, 8.5 % ETO in 91.5 % CO2	N/A	Sterilizing Agent	Ethylene Oxide (EO)	N/A
	Subject Device Flexible Video-Uretero-Renoscope SSU System	Predicate Device, K141250 Flexible Video-Uretero-Renoscope System																																																									
Maximal Outer diameter Insertion Portion	3.2 mm	Same as subject device																																																									
Outer diameter Insertion Tube	2.9 mm	Same as subject device																																																									
Outer diameter Distal End	3.2 mm	Same as subject device																																																									
Insertion portion length	700 mm	675 mm																																																									
Working channel	Present	Same as subject device																																																									
Inner diameter Working Channel	1.2 mm	Same as subject device																																																									
Tip deflection up/down	270°/270°	Same as subject device																																																									
Field of view	110°	90°																																																									
Direction of View	0°	Same as subject device																																																									
Depth of Field	5 - 50 mm	4 - 60mm																																																									
On-axis Resolution	12.5 Lp/mm at 5 mm 4.5 Lp/mm at 15 mm 1.25 Lp/mm at 50 mm	11.0 Lp/mm at 4 mm 4.0 Lp/mm at 12 mm 1.0 Lp/mm at 60 mm																																																									
Chip type	CMOS	Same as subject device																																																									
Chip location	Distal	Same as subject device																																																									
Illumination source	LED	Same as subject device																																																									
Compatible CCU	C-MAC C-HUB II Image 1S	Image 1S																																																									
How device is provided	Sterile single-use	Unsterile, reusable																																																									
EO Sterilization cycle	EO, Overpressure 2.7 bar absolute, 8.5 % ETO in 91.5 % CO2	N/A																																																									
Sterilizing Agent	Ethylene Oxide (EO)	N/A																																																									
<p>Non-Clinical Performance Data:</p>	<p>There are no performance standards or special controls developed under Section 514 of the FD&C Act for endoscopes. However, the subject device follows the FDA recognized consensus standards and is tested according to the following standards and FDA Guidance:</p> <ul style="list-style-type: none"> • Electrical Safety and EMC <ul style="list-style-type: none"> ○ IEC 60601-1 ○ IEC 60601-1-2 ○ IEC 60601-2-18 																																																										

	<ul style="list-style-type: none"> ○ IEC 62471 ○ ISO 10993 ○ ISO 8600 ● Performance Testing <ul style="list-style-type: none"> ○ Color Contrast Enhancement ○ Image intensity uniformity ○ Depth of field & Spatial Resolution ○ Distortion ○ Signal-to-Noise Ratio (SNR) & Sensitivity <p>Additional bench testing was performed to ensure the device met its design specifications. The bench testing performed verified and validated that the SSU Flexible Video-Uretero-Renoscope has met all its design specification and is substantially equivalent to its predicate devices.</p>
Clinical Performance Data:	Clinical testing was not required to demonstrate the substantial equivalence to the predicate devices. Non-clinical bench testing was sufficient to establish the substantial equivalence of the modifications.
Conclusion:	The conclusions drawn from the nonclinical tests demonstrate that the subject devices, the SSU Flexible Video-Uretero-Renoscope performs as well as or better than the predicate devices that are currently marketed for the same intended use.