



Dentium Co., Ltd (ICT Branch)
% Mr. Dave Kim
Medical Device Regulatory Affairs
Mtech Group
7707 Fannin St. Ste 200, V111
HOUSTON TX 77054

April 16, 2021

Re: K200270
Trade/Device Name: rainbow MCT
Regulation Number: 21 CFR 892.1750
Regulation Name: Computed tomography x-ray system
Regulatory Class: Class II
Product Code: OAS
Dated: March 10, 2021
Received: March 15, 2021

Dear Mr. Kim:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.
Director
Division of Radiological Health
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K200270

Device Name
rainbow MCT

Indications for Use (Describe)

rainbow MCT is a computed tomography x-ray system intended to produce 3D and panoramic diagnostic images of the maxillofacial areas for treatment planning for adult and pediatric patients. The device is operated and used by physicians, dentists, and x-ray technicians.

Rainbow 3D Image Viewer software features functions for acquiring, saving, searching, displaying, diagnosing and sending digital X-ray image data in dental practices and clinics.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary**K200270**

This 510(k) is being submitted in accordance with the requirements of 21 CFR §807.92.

1. Date Summary Prepared:

March 10,2021

2. Submitter's Identification:

Submitter's Name : Dentium Co., Ltd (ICT Branch)
 Submitter's Address: 76, Changnyong-daero 256beon-gil, Yeongtong-gu,
 Suwon-si, Gyeonggi-do, 16229
 Republic of Korea
 Submitter's Telephone: ++82-70-7098-6932

Contact person: Mr. Sang Woo Lee (swlee1@dentium.com)

Official Correspondent: Dave Kim (davekim@mtech-inc.net)
 (U.S. Designated agent)
 Address: 7707 Fannin St. Ste 200, Houston, TX 77054
 Telephone: +1- 713-467-2607

3. Device:

Trade /Proprietary Name: rainbow MCT
 Device: X-Ray, Tomography, Computed, Dental
 Regulation Description: Computed tomography x-ray system.
 Review Panel: Radiology
 Product Code: OAS
 Regulation Number 892.1750
 Device Class II

4. Predicate Device:

Legally Marketed Predicate Device Information:
 510(k) Number: K181432
 Trade /Proprietary Name: ProVecta 3D Prime with VistaSoft
 Device: X-Ray, Tomography, Computed, Dental
 Regulation Description: Computed tomography x-ray system.
 Review Panel Radiology
 Product Code OAS
 Regulation Number 892.1750
 Device Class II

5. Reference Device:

Legally Marketed Predicate Device Information:

510(k) Number:	K102196
Trade /Proprietary Name:	PaX-Zenith3D
Regulation Description:	Computed tomography x-ray system
Review Panel	Radiology
Product Code	OAS
Regulation Number	892.1750
Device Class	II

6. Device Description:

- rainbow MCT is a cone beam CT X-ray device for generating sectional images of dental images such as tooth, nasal cavity and temporomandibular joint. this is a medical diagnostic equipment designed to generate sectional images by placing X-ray source opposite to the imaging detector unit and rotating it around a patient. 2D images of the region of interest are reconstructed using a mathematical algorithm in 3 dimensional volumetric view and displayed on the computer monitor.

- The system is composed of X-ray generator, X-ray detector, X-ray collimator, main frame, rotation unit, PC and Monitor, etc. in compliance with US performance standard and regulatory requirement.

7. Indications for use:

rainbow MCT is a computed tomography x-ray system intended to produce 3D and panoramic diagnostic images of the maxillofacial areas for treatment planning for adult and pediatric patients. The device is operated and used by physicians, dentists, and x-ray technicians.

Rainbow 3D Image Viewer software features functions for acquiring, saving, searching, displaying, diagnosing and sending digital X-ray image data in dental practices and clinics.

8. Summary of the technological characteristics of the device compared to the predicate devices:

Summary of the Technological Characteristics

Descriptive Information		Rainbow MCT K200270	K181432, ProVecta 3D Prime with VistaSoft (DÜRR DENTAL)
Indications for Use		<p>- rainbow MCT is a computed tomography x-ray system intended to produce 3D and panoramic diagnostic images of the maxillofacial areas for treatment planning for adult and pediatric patients. The device is operated and used by physicians, dentists, and x-ray technicians.</p> <p>Rainbow 3D Image Viewer software features functions for acquiring, saving, searching, displaying, diagnosing and sending digital X-ray image data in dental practices and clinics.</p>	<p>ProVecta 3D Prime is computed tomography x-ray unit intended to generate 3D and panoramic X- ray images in dental radiography for adult and pediatric patients. It provides diagnostic details of the maxillofacial areas for a dental treatment. The device is operated and used by physicians, dentists, and x-ray technicians.</p> <p>Not intended for mammography use. The VistaSoft software features functions for recording, displaying, analyzing, diagnosing, managing and sending digital or digitized video and X-ray images in dental practices and specialist dental clinics. Not intended for mammography use</p>
Image Acquisition Modes		Panoramic and computed tomography	Panoramic and computed tomography
Input Voltage		AC 100-240 V, 50/60 Hz	AC 200-240V
Tube Voltage		60~100 kV	50-99 KV
Tube Current		4~12 mA	4~16mA
Focal Spot Size		0.5 x 0.5 mm	0,5 mm
Exposure Time		Max. 19 s	Max. 16.4s
Slice Width		0.1 mm min.	0.1 mm min.
Total Filtration		2.5 mm Al	2.8 mm Al
Chin Rest		Bite block, chin rest and headrest	Bite block, chin rest and headrest
Mechanical		Compact design	Compact design
Electrical		LDCP logic circuit (Low Dark Current Processing)	LDCP logic circuit (Low Dark Current Processing)
Software		Rainbow 3D ImageViewer, DICOM 3.0 Format compatible	VistaSoft, DICOM 3.0 compatible
Anatomical Sites		Maxillofacial	Maxillofacial
Image Receptor Note: CT and	CBCT	DTX3024	Xmaru1404CF
	Panoramic	DTX3024	Xmaru1404CF

panoramic image performance is identical because the sensors are identical.	MTF@ 1 lp/mm	50%	53%
	DQE @ 0.5 lp/mm	63%	64%
Size of Imaging Volume (cm)		DTX3024: Max. 10x8, 23x21	Xmaru1404CF : Max. 10x8.5
Pixel Resolution	CBCT	DTX3024: 5 lp/mm	2.5 lp/mm - 4x4 binning
	Panoramic	DTX3024: 5 lp/mm	2.5 lp/mm - 4x4 binning
Pixel Size	CBCT	DTX3024 :100 μm	Xmaru1404CF : 99 μm - 2x2 binning 198 μm - 4x4 binning
	Panoramic	DTX3024:100 μm	Xmaru1404CF : 99 μm - 2x2 binning 198 μm - 4x4 binning

9. Rainbow MCT detectors VS. PaX-Zenith 3D Detectors

		Rainbow MCT (K200270)	Pax-Zenith3D (K102196)		
Imaging detector for CT and panoramic mode	CBCT	DTX3024 (10x8, 23 x 21cm)	Xmaru2430CF (FOV 24 x 19cm)	Xmaru1524CF (FOV 15 x 16cm)	
	Panoramic	DTX3024			Xmaru1501CF
Imaging Performance	MTF@ 1 lp/mm	50%	53%	52%	50%
	DQE @ 0.5 lp/mm	63%	64%	45%	60%
Size of Imaging Volume (cm)					
Pixel Resolution	CBCT	5 lp/mm - 1x1	2.5 lp/mm - 4x4 binning		
	Panoramic	5 lp/mm - 1x1			5 lp/mm
Pixel Size	CBCT	100 μm	200 μm -	200 μm -	
	Panoramic	100 μm			100 μm

10. Discussion of Similarities and Differences:

Rainbow MCT dental computed tomography X-ray system described in this 510(k) is similar to the predicate device in its indications for use, performance, materials, and safety characteristics.

The differences include the digital X-ray imagers and image viewing software. The subject device offers FOV 10x8 & 23x21 cm whereas the predicate device and the reference device provide FOV 10.8.5 and 24x19 cm, respectively. Performance testing was conducted for the subject device to assess whether or not the parameter required for functionalities related to imaging properties of the dental X-ray device meets the designated acceptance criteria. The MTF, DQE and pixel resolution of the subject device performed similar or better than those of the predicate and reference device.

All test results were satisfactory.

11. Non-Clinical Data and Performance Testing

Electrical, mechanical, environmental safety and performance testing according to standard IEC 60601-1(2005, AMD 1: 2012), IEC 60601-1-3 (2008 + A1: 2013), IEC 60601-2-63 (2012,041: 2017))were performed, and EMC testing were conducted in accordance with standard IEC 60601-1-2 (2014).

rainbow MCT meets the provisions of NEMA PS 3.1-3.18, Digital Imaging and Communications in Medicine (DICOM) Set.

Non-clinical & Clinical considerations according to FDA Guidance “Guidance for the submissions of 510(k)’s for Solid State X-ray Imaging Devices” were performed. Acceptance test according to IEC 61223-3-4 and IEC 61223-3-5 was performed. All test results were satisfactory.

12. **Clinical Data**: Not required for a finding of substantial equivalence.

13. **Conclusion**:

In accordance with the Federal Food, Drug and Cosmetic Act, 21 CFR Part 807 and based on the similarity to the predicate device in terms of technology, performance and indications for use, Dentium Co., Ltd concludes that the rainbow MCT is substantially equivalent to ProVecta 3D Prime with VistaSoft (K181432), the predicate device, and Pax-Zenith3D (K102196), the reference device as described herein.

The differences between the new device and the predicate device shown in the comparison table above do not raise any new questions about safety and effectiveness and so we consider it substantially equivalent to the predicate device.