



April 21, 2021

Chongqing COE Display Technology Co., Ltd.
% James Tsai
Consultant
Shenzhen Joyantech Consulting Co., Ltd.
1713A, 17th Floor, Block A, Zhongguan Times Square,
Nanshan District
Shenzhen, Guangdong 518000
China

Re: K202759

Trade/Device Name: Disposable medical surgical mask
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical Apparel
Regulatory Class: Class II
Product Code: FXX
Dated: September 12, 2020
Received: September 21, 2020

Dear James Tsai:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray III, Ph.D.
Acting Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K202759

Device Name
Disposable medical surgical mask

Indications for Use (Describe)

The disposable medical surgical masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided as sterile and non-sterile.

Level 3, surgical mask model, sterile: CCSM-SC

Level 3, surgical mask model, non-sterile: CCSM-NC

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K202759

1. Administrative Information

Date of Summary prepared	April 17, 2021
Applicant information	Company: Chongqing COE Display Technology Co., Ltd. Company address: No. 3, Tonggui Avenue, Yufengshan Town, Yubei District, Chongqing, China Contact person: Liu Qingjun Phone: 86-023-67849398 Fax: 86-023-67843279 E-mail: liuqj@szcoe.com
Submission Correspondent	Shenzhen Joyantech Consulting Co., Ltd. Address: 1713A, 17th Floor, Block A, Zhongguan Times Square, Nanshan District, Shenzhen Contact person: James Tsai E-Mail: james_tsai@cefd.com; field@cefd.com

2. Device Information

Type of 510(k) submission:	Traditional
Trade Name:	Disposable medical surgical mask
Classification name:	Surgical Face Mask, Apparel
Review Panel:	General and plastic surgery devices
Product Code:	FXX
Common name	Surgical mask
Device Class:	II
Regulation Number:	878.4040

3. Predicate Device Information

Sponsor:	SAN-M PACKAGE CO., LTD.
Device trade name:	Surgical face mask (Ear loops and Tie-on)
Device Class:	II
510(K) Number:	K160269
Regulation name	Masks, Surgical

Production regulation:	21 CFR §878.4040
Product code:	FXX
Review Panel:	General and plastic surgery devices

4. Device Descriptions

The surgical mask consists of a mask body, a nose piece, and ear loops. The mask body is divided into three layers, the inner and outer layers are made of polypropylene fabrics; the middle layer is composed of Melt-blown cloth (polypropylene); the nose piece is made of polypropylene & iron (Fe), and the ear loops are made of nylon & spandex.

5. Indications for Use

The disposable medical surgical masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided as sterile and non-sterile.

6. Summary of Technological Characteristics

6.1 Comparisons for non-sterile mask

Comparison item	Proposed Device	Predicate Device	Remark
Manufacturer	Chongqing COE Display Technology Co., Ltd.	SAN-M PACKAGE CO., LTD.	Different
510k number	K202759	K160269	Different
Product name	Disposable medical surgical mask	Surgical face mask	Similar
Product model	CCSM-NC	EL 20000	Different
Product Code	FXX	FXX	Same
Classification	Class II (21 CFR 878.4040)	Class II (21 CFR 878.4040)	Same
Intended use & Indications for Use	The disposable medical surgical masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the	The surgical face masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood	Similar

	potential exposure to blood and body fluids. This is a single-use, disposable device, provided as sterile and non-sterile	and body fluids. This is a single-use, disposable device, provided non-sterile.	
Mask features	Ear Loops, flat Pleated	Ear Loops, flat Pleated	Same
Layers	3 layers	4 layers	Different
Outer layer	Polypropylene	Polypropylene	Same
Filter media	Melt-blown cloth (polypropylene)	Polypropylene spunbond Polypropylene meltblown	Similar
Inner layer	Polypropylene	Polypropylene	Same
Ear loops	Nylon & Spandex	Ear loops: Polyester, polyurethane	Similar
Nose piece	Polypropylene & Iron (Fe)	Polyethylene coated steel wire	Similar
Color	Blue	White or Blue	Same
Dimension (Length*Width)	Mask body: 17.5cm*9.5cm Nosepiece: length: 8cm-12cm; width: 3cm Ear loop: 18cm±0.9cm	17.5cm*9.5cm	Same
OTC use	Yes	Yes	Same
Sterility	Non-sterile	Non-sterile	Same
Use	Single-use, disposable	Single-use, disposable	Same
Performance level	Level 3	Level 3	Same
Fluid Resistance Performance ASTM F1862	32 Out of 32 Pass at 160mmHg	32 Out of 32 Pass at 160mmHg	Same
Particulate Filtration Efficiency ASTM F2299	Pass at $\geq 98.0\%$	Pass at 99.72%	Similar
Bacterial Filtration Efficiency ASTM F2101	Pass at $\geq 98\%$	Pass at 99%	Similar
Differential Pressure ASTM F2100-19	Pass at 4.6 mm H ₂ O/cm ²	Pass at 2.5 mmH ₂ O/cm ₂	Similar
Flammability 16 CFR 1610	Class I	Class I	Same

Cytotoxicity ISO 10993-5	Under the conditions of the study, the subject device was non-cytotoxic	Under the conditions of the study, the subject device was non-cytotoxic	Same
Sensitization ISO 10993-10	Under the conditions of the study, the subject device was non-sensitizing	Under the conditions of the study, the subject device was non-sensitizing	Same
Irritation ISO 10993-10	Under the conditions of the study, the subject device was non-irritating	Under the conditions of the study, the subject device was non-irritating	Same

The subject device only includes level 3 and ear loops type, it was covered by the scopes of predicate device, which includes level 1, level 2 and level 3, as well as types of ear loop.

Both proposed device and predicate device conform to ASTM F2100, the difference is the versions of the standard, the proposed device was tested by the latest version, and this difference will not raise any new issue of safety and effectiveness of the proposed device.

Although the materials, and layers of subject device are a little different from the predicate device, the performance and biocompatibility testing of the subject device meet all the requirements of standards of ASTM F2100 and ISO 10993-5 &-10. So, the differences between the predicate device and subject device will not raise any new issue of safety and effectiveness of the subject device.

6.2 Comparisons for sterile mask

Comparison item	Proposed Device	Predicate Device	Remark
Manufacturer	Chongqing COE Display Technology Co., Ltd.	SAN-M PACKAGE CO., LTD.	Different
510k number	K202759	K160269	Different
Product name	Disposable medical surgical mask	Surgical face mask	Similar
Product model	CCSM-SC	EL 20000	Different
Product Code	FXX	FXX	Same
Classification	Class II (21 CFR 878.4040)	Class II (21 CFR 878.4040)	Same
Intended use & Indications for Use	The disposable medical surgical masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control	The surgical face masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the	Similar

	practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided as sterile and non-sterile	potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.	
Mask features	Ear Loops, flat Pleated	Ear Loops, flat Pleated	Same
Layers	3 layers	4 layers	Different
Outer layer	Polypropylene	Polypropylene	Same
Filter media	Melt-blown cloth (polypropylene)	Polypropylene spunbond Polypropylene meltblown	Similar
Inner layer	Polypropylene	Polypropylene	Same
Ear loops	Nylon & Spandex	Ear loops: Polyester, polyurethane	Similar
Nose piece	Polypropylene & Iron (Fe)	Polyethylene coated steel wire	Similar
Color	Blue	White or Blue	Same
Dimension (Length*Width)	Mask body: 17.5cm*9.5cm Nosepiece: length: 8cm-12cm; width: 3cm Ear loop: 18cm±0.9cm	17.5cm*9.5cm	Same
OTC use	Yes	Yes	Same
Sterility	Sterile	Non-sterile	Different
Packaging material	Paper plastic bag	Not publicly available	Different
Sterilization method and S.A.L.	Sterilized by ethylene oxide gas, SAL=10 ⁻⁶	Not applied	Different
Residues	M _{EO} ≤4mg/d; M _{ECH} ≤9mg/d	Not applied	Different
Shelf life	No shelf life claim	No shelf life claim	Same
Use	Single-use, disposable	Single-use, disposable	Same
Performance level	Level 3	Level 3	Same
Fluid Resistance Performance ASTM F1862	32 Out of 32 Pass at 160mmHg	32 Out of 32 Pass at 160mmHg	Same
Particulate Filtration Efficiency ASTM F2299	Pass at ≥ 98.0%	Pass at 99.72%	Similar

Bacterial Filtration Efficiency ASTM F2101	Pass at $\geq 98\%$	Pass at 99%	Similar
Differential Pressure ASTM F2100-19	Pass at 4.8 mm H ₂ O/cm ²	Pass at 2.5 mmH ₂ O/cm ₂	Similar
Flammability 16 CFR 1610	Class I	Class I	Same
Cytotoxicity ISO 10993-5	Under the conditions of the study, the subject device was non-cytotoxic	Under the conditions of the study, the subject device was non-cytotoxic	Same
Sensitization ISO 10993-10	Under the conditions of the study, the subject device was non-sensitizing	Under the conditions of the study, the subject device was non-sensitizing	Same
Irritation ISO 10993-10	Under the conditions of the study, the subject device was non-irritating	Under the conditions of the study, the subject device was non-irritating	Same

The subject device only includes level 3 and ear loops type, it was covered by the scopes of predicate device, which includes level 1, level 2 and level 3, as well as types of ear loop.

Both proposed device and predicate device conform to ASTM F2100, the difference is the versions of the standard, the proposed device was tested by the latest version, and this difference will not raise any new safety and effectiveness of the proposed device.

Although the materials, layers and sterility status of subject device are a little different from the predicate device, the performance, biocompatibility and EO sterilization of the subject device meet all the requirements of standards of ASTM F2100, ISO 10993-5 &-10 and ISO 11135. So, the differences between the predicate device and subject device will not raise any new issue of safety and effectiveness of the subject device.

7. Brief discussion of Non-Clinical Testing

Surgical Face Mask conforms to the following standards:

ASTM F 2100-19, Standard Specification for Performance of Materials Use in Medical Face Masks.

ISO 10993-1:2018, Biological Evaluation Of Medical Devices - Part 1: Evaluation And Testing Within A Risk Management Process.

Bench testing

The bench testing of Surgical Face Mask includes the following tests:

- *Fluid Resistance Performance
- *Particulate Filtration Efficiency
- *Bacterial Filtration Efficiency
- *Differential Pressure

*Flammability

Biocompatibility testing

The biocompatibility evaluation for the proposed device was conducted in accordance with the International Standard ISO 10993-1: 2018 “Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing within a Risk Management Process,” as recognized by FDA. Biocompatibility evaluation included the following tests:

- Cytotoxicity
- Skin Sensitization
- Skin Irritation

Ethylene oxide Sterilization Validation

The proposed device is also provided for sterility, sterilization validation is performed, and the results show that the proposed device complies with the following standards:

- ISO 11135:2014 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices [Including: Amendment 1 (2018)]
- ISO 11737-1:2018 Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on product
- ISO 10993-7:2008 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals
- ASTM F88-2015 Standard Test Method for Seal Strength of Flexible Barrier Materials
- ASTM F1929-2015 Standard test method for detecting seal leaks in porous medical packaging by dye penetration.

Summary

Based on the non-clinical performance data as documented above in the device development, the proposed device has a safety and effectiveness profile that is similar to the predicate device, the testing results are summarized in the following table:

Name of Test Methodology (standard)	Purpose	Acceptance Criteria	Results	
			Sterile mask	Non-sterile mask
ASTM F 1862-17	Fluid Resistance Performance	29 out of 32 pass at 160 mmHg	Lot 1# pass at 160mmHg; Lot 2# pass at 160mmHg; Lot 3# pass at 160mmHg	Lot 1# pass at 160mmHg; Lot 2# pass at 160mmHg; Lot 3# pass at 160mmHg;
ASTM F2101-19	Bacterial Filtration Efficiency	≥ 98%	Lot 1# 99.6%-99.9%; Lot 2# 99.7%-99.9%; Lot 3# 99.6%-99.9%	Lot 1# 99.6%-99.9%; Lot 2# 99.6%-99.8%; Lot 3# 99.7%-99.9%

	Performance			
EN 14683: 2019	Differential Pressure (Delta-P)	< 6.0mm H ₂ O/cm ²	Lot 1# 4.0-5.2; Lot 2# 4.4-5.4; Lot 3# 3.6-5.2	Lot 1# 4.0-5.1; Lot 2# 4.1-5.0; Lot 3# 4.3-5.0
ASTM F2299- 2007	Particulate Filtration Efficiency Performance	≥ 98%	Lot 1# 98.0%-98.9%; Lot 2# 98.3%-98.9%; Lot 3# 98.3%-98.9%	Lot 1# 98.2%-98.9%; Lot 2# 98.5%-99.1%; Lot 3# 98.6%-99.0%
16 CFR Part 1610	Flammability	Class I	Lot 1# Class I; Lot 2# Class I; Lot 3# Class I	Lot 1# Class I; Lot 2# Class I; Lot 3# Class I
ISO 10993-5	Cytotoxicity	Non-cytotoxic	Under the conditions of the study, the subject device was non-cytotoxic	
ISO 10993-10	Sensitization	Non-sensitizing	Under the conditions of the study, the subject device was non-sensitizing	
ISO 10993-10	Irritation	Non-irritating	Under the conditions of the study, the subject device was non-irritating	

8. Brief discussion of clinical tests

No clinical tests were performed.

9. Conclusions

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission K202759, the Disposable medical surgical mask is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K160269.