



March 13, 2020

Trophy  
% Ms. Marie-Pierre Labat-Camy  
Global Regulatory Affairs Senior Manager  
4 Rue F. Pelloutier – Croissy-Beaubourg  
77435 Marne La Vallée, Cedex 2  
FRANCE

Re: K200183  
Trade/Device Name: CS 9600  
Regulation Number: 21 CFR 892.1750  
Regulation Name: Computed tomography x-ray system  
Regulatory Class: Class II  
Product Code: OAS  
Dated: January 22, 2020  
Received: January 24, 2020

Dear Ms. Marie-Pierre Labat-Camy:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting->

[combination-products](#)); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.  
Director  
Division of Radiological Health  
OHT7: Office of In Vitro Diagnostics  
and Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

**ATTACHMENT 7**

**INDICATIONS FOR USE STATEMENT**

## Indications for Use

510(k) Number (if known)

K200183

Device Name

CS 9600

Indications for Use (Describe)

The CS 9600 is extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.

The CS 9600 can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# Traditional 510(k) Submission for CS 9600

## 510(K) SUMMARY

K200183

### 1. Date 510(k) Summary prepared

March 11, 2020

### 2. Submitter information

#### **Applicant**

Trophy  
4 Rue F. Pelloutier – Croissy-Beaubourg  
77435 Marne La Vallée, Cedex 2  
France  
Phone: +33 1 64 80 85 26  
Contact person: Ms. Marie-Pierre LABAT-CAMY

#### **United States Sales Representative (U.S. Designated Agent):**

Carestream Dental LLC  
3625 Cumberland Boulevard, Suite 700,  
Atlanta, GA USA 30339  
Phone: +1 (470) 481-4619  
Contact: Mr. Sonny T. Nguyen

### 3. Device name and classification

Trade name: **CS 9600**  
Regulation number: 21 CFR 892.1750  
Regulation name (classification name): Computed tomography x-ray system  
Common name: Cone-beam computed tomography system  
Device Class: II  
Product Code: OAS

### 4. Predicate device and reference device

We consider the **CS 9600** to be similar in design, composition and function to the following primary predicate device and reference device introduced into commercial distribution after May 28, 1976:

Primary Predicate Device Name	510(k)	Company Name
CS 9600	K181136	Trophy

Regulation number: 21 CFR 892.1750  
Regulation name (classification name): Computed tomography x-ray system  
Common name: Cone-beam computed tomography system  
Device Class: II  
Product Code: OAS

Reference Device Name	510(k)	Company Name
CS 8100SC	K151087	Trophy

Regulation number: 21 CFR 872.1800  
Regulation name (classification name): Extraoral source x-ray system  
Common name: X-ray imaging device  
Device Class: II  
Product Code: MUH

# Traditional 510(k) Submission for CS 9600

## 5. Device Description

**CS 9600** is an extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.

The **CS 9600** can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.

**CS 9600** is a cone-beam computed tomography (CBCT) x-ray system. It means **CS 9600** rotates around the patient, capturing data using a cone-shaped x-ray beam. These data are used to reconstruct a two or a three-dimensional (3D) image of the following regions of the patient's anatomy: dental (teeth); oral and maxillofacial region (mouth, jaw and neck); ears, nose and throat region (ENT); cervical spine or wrist region.

Additional features such as low dose mode, scout image and metal artifact reduction are also provided by the **CS 9600**.

The **CS 9600** can also be upgraded with cephalometric modality. The cephalometric modality of the proposed device **CS 9600** is the same than the one available in the reference device K151087. The cephalometric mode works with a narrow beam linear scanning process called a "slot technique". The patient head is scanned in lines with a flat, fan-shaped x-ray beam.

## 6. Indication for use

**CS 9600** is an extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.

The **CS 9600** can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.

## 7. Substantial Equivalence

The subject device **CS 9600**, the primary predicate device K181136 and the reference device K151087 have the same intended use: extraoral x-ray systems that are intended to produce digital X-ray images at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.

Both the **CS 9600** and primary predicate device are intended as Cone-Beam Computed Tomography. The primary predicate device K181136 has the most similar indication for use and technological characteristics with the proposed device **CS 9600**.

## Traditional 510(k) Submission for CS 9600

The table below provides the summary of the technological characteristics of **CS 9600** compared to the primary predicate device.

	<b>Proposed device</b>	<b>Primary Predicate Device</b>
<b>COMPANY NAME MODEL NAME</b>	<b>Trophy CS 9600</b>	<b>Trophy CS 9600</b>
<b>510(K) Number</b>		<b>K181136</b>
<b>Indication for use</b>	<p><b>CS 9600</b> is an extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.</p> <p>The <b>CS 9600</b> can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.</p>	CS 9600 is an extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.
<b>Performance specification</b>	<ul style="list-style-type: none"> <li>• Panoramic modality</li> <li>• 3D modality</li> <li>• Cephalometric modality (optional)</li> </ul>	<ul style="list-style-type: none"> <li>• Panoramic modality</li> <li>• 3D modality</li> </ul>
<b>Rated line voltage</b>	100-240 Vac - 50/60 Hz	Same
<b>X-ray tube voltage</b>	60-90 kV 60-120 kV (in option)	Same
<b>X-ray tube current</b>	2-15 mA	Same
<b>X-ray tube</b>	DF-071G or OX/120-0307	Same
<b>Tube focal spot</b>	0.3 or 0.7 mm	Same
<b>Patient sizes</b>	4 patients sizes: child, small adult, medium adult, large adult	Same

The proposed **CS 9600** and the primary predicate device CS 9600 (K181136) have the same panoramic two-dimensional modality and three-dimensional 3D modality.

	<b>Proposed device Trophy CS 9600</b>	<b>Primary Predicate device Trophy CS 9600</b>
<b>510(K) Number</b>		<b>K181136</b>
<b>Two-dimensional modality: PANORAMIC</b>		
<b>Sensor technology</b>	CMOS	Same
<b>Image field</b>	6.4 x 140 mm (for adult patient size) 6.4 x 120 mm (for child patient size) 120 x 140 mm (for sinus one-shot exam)	Same
<b>Gray scale</b>	16384 – 14 bits	Same
<b>Magnification</b>	1.28	Same
<b>Radiological Exams</b>	<ul style="list-style-type: none"> <li>• Full panoramic</li> <li>• Segmented panoramic</li> <li>• Bitewing</li> <li>• Maxillary sinus</li> <li>• Lateral TMJ x2 and Lateral TMJx4</li> <li>• Sinus AP / PA / Lateral</li> <li>• Orthogonal panoramic</li> </ul>	Same
<b>Exposure time</b>	0.5 to 13 seconds	Same

## Traditional 510(k) Submission for CS 9600

<b>COMPANY NAME MODEL NAME</b>	<b>Proposed device Trophy CS 9600</b>	<b>Primary Predicate device Trophy CS 9600</b>
<b>510(K) Number</b>		<b>K181136</b>
<b>Three-dimensional modality: 3D</b>		
<b>Sensor technology</b>	CMOS	Same
<b>Field of View (cm) diameter x height</b>	<ul style="list-style-type: none"> <li>• 4 x 4</li> <li>• 5 x 5 (child 4 x 4)</li> <li>• 5 x 8</li> <li>• 6 x 6</li> <li>• 8 x 5</li> <li>• 8 x 8</li> <li>• 10 x 5 (child 8 x 5)</li> <li>• 10 x 10* (child 8 x 8)</li> </ul>	<ul style="list-style-type: none"> <li>• 12 x 5</li> <li>• 12 x 10*</li> <li>• 16 x 6</li> <li>• 16 x 10*</li> <li>• 16 x 12*</li> <li>• 16 x 17*</li> <li>*with tip of the volume</li> </ul>
<b>Radiological Exams</b>	<ul style="list-style-type: none"> <li>• Tooth/Teeth</li> <li>• Jaw (full, upper or lower)</li> <li>• Upper cervical spine</li> </ul>	<ul style="list-style-type: none"> <li>• TMJ</li> <li>• Face</li> <li>• ENT</li> <li>• Wrist</li> </ul>
<b>Gray scale</b>	16384 – 14 bits	Same
<b>Magnification</b>	1.4	Same
<b>Voxel size (µm)</b>	75, 150, 300 and 400	Same
<b>Exposure time</b>	3-20 seconds	Same
<b>OTHER INFORMATION</b>		
<b>Low dose mode</b>	Yes	Same
<b>Scout image</b>	Yes	Same
<b>3D Face Photo</b>	Yes in option (CS Face Scan)	Same
<b>Metal Artefact Reduction</b>	Yes in option (CS MAR)	Same
<b>Concerned anatomical sites</b>	<ul style="list-style-type: none"> <li>• Dento-maxillofacial area</li> <li>• ENT area</li> <li>• Cervical spine</li> <li>• Hand/wrist</li> </ul>	Same
<b>Sensor model</b>	CGF81	Same
<b>Sensor active area (mm)</b>	120 x 140	Same
<b>Pixel size (µm)</b>	100 x 100	Same
<b>Sensor resolution</b>	1200 x 1400 pixels	Same
<b>Limiting resolution</b>	5 lp/mm	Same
<b>MTF, X-ray (%) at 1 lp/mm, Typical</b>	60	Same
<b>DQE, X-ray (%) at 0 lp/mm, Typical</b>	60	Same

The proposed **CS 9600** and the primary predicate device CS 9600 (K181136) have the same panoramic two-dimensional modality and three-dimensional 3D modality.

The main difference between the proposed **CS 9600** and the primary predicate device K181136 is the difference in 2D imaging capabilities available with the addition in the proposed **CS 9600** of the optional two-dimensional cephalometric modality.

The proposed device **CS 9600** can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.



## Traditional 510(k) Submission for CS 9600

When the **CS 9600** is upgraded with cephalometric modality to produce cephalometric digital X-ray images, the cephalometric modality available is the same as the cephalometric modality available in previously cleared reference device K151087.

The following table provides the summary of the technological characteristics of the cephalometric modality of **CS 9600** compared to the reference device K151087.

	Proposed device	Reference Device
<b>COMPANY NAME</b>	<b>Trophy</b>	<b>Trophy</b>
<b>MODEL NAME</b>	<b>CS 9600</b>	<b>CS 8100SC</b>
<b>510(K) Number</b>		<b>K151087</b>
<b>Indication for use</b>	<p><b>CS 9600</b> is an extraoral system intended to produce two-dimensional and three-dimensional digital X-ray images of the dento-maxillofacial, ENT (Ear, Nose and Throat), cervical spine and wrist regions at the direction of healthcare professionals as diagnostic support for pediatric and adult patients.</p> <p>The <b>CS 9600</b> can be upgraded to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.</p>	<p>The CS 8100SC and the CS 8100SC Access are indicated to produce complete and segmented tomographic digital panoramic and <b>cephalometric digital X-ray images</b> to be used at the direction of healthcare professionals of the dento-maxillofacial region of the human anatomy as diagnostic support for pediatric and adult patients.</p> <p><b>This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.</b></p> <p>In addition, the CS 8100SC is provided with the panoramic radiological exam option which allows the display of one determined segment of the dental structures using the capability of the system to explore each slide during an exam.</p>
<b>Two-dimensional modality: CEPHALOMETRIC</b>		
<b>Sensor technology</b>	CMOS	CMOS
<b>Radiological Exams</b>	<ul style="list-style-type: none"> <li>• Lateral</li> <li>• Frontal AP or PA</li> <li>• Oblique</li> <li>• Submento-vertex</li> <li>• Carpus (optional)</li> </ul>	Same
<b>Field of View (cm) LxH</b>	18x18, 18x24 and 26x24	Same
<b>Gray scale</b>	16384 – 14 bits	Same
<b>Exposure time</b>	2.96 to 10 seconds	3 to 10 seconds
<b>Sensor model</b>	CGC95	Same
<b>Sensor active area (mm)</b>	131.6 x 6.4	Same
<b>Pixel size (µm)</b>	100 x 100	Same
<b>Sensor resolution</b>	1316 x 64 pixels	Same
<b>Limiting resolution</b>	5 lp/mm	Same
<b>MTF, X-ray (%) at 1 lp/mm, Typical</b>	60	Same
<b>DQE, X-ray (%) at 0 lp/mm, Typical</b>	60	Same

## Traditional 510(k) Submission for CS 9600

When the **CS 9600** is upgraded with cephalometric modality to produce cephalometric digital X-ray images, the proposed **CS 9600** and the reference device K151087 have the same cephalometric modality. Both systems are able to produce cephalometric digital X-ray images. This includes imaging the hand and wrist to obtain carpus image for growth and maturity assessment.

While the indication for use of primary predicate device K181136 and reference device K151087 are different, both devices have the same intended use, namely extraoral x-ray system.

The proposed device **CS 9600** combines the indication for use of both the primary predicate device K181136 with the part of indication for use related to the cephalometric modality of the reference device K151087.

### 8. Non-Clinical Performance Data

The Testing to verify the performance requirements of the subject device **CS 9600** was conducted and including in this premarket notification.

#### **Standards Conformance:**

EMC and Electrical Safety testing were performed respectively by LCIE and UL laboratory and found to meet all the requirements in standards IEC 60601-1: 2005 with A1 2012 (AAMI/ANSI ES60601-1:2005/R(2012)), IEC 60601-1-2: 2014, IEC 60601-1-3: 2008 with A1 2013, IEC 60601-2-63: 2017 and IEC 62304: 2006 with A1 2015.

**CS 9600** meets the provisions of NEMA PS 3.1-3.20, Digital Imaging and Communications in Medicine (DICOM) Set. The software contained in **CS 9600** has been validated according to FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices and FDA Guidance for the Content of Premarket Submissions for Management of Cybersecurity in Medical Devices.

The **CS 9600** provides design features and instructions for pediatric population according to the FDA Guidance on Pediatric Information for X-ray Imaging Device Premarket Notifications.

#### **Performance Testing:**

Since there is no change on the two-dimensional panoramic and three-dimensional 3D modalities of the proposed device **CS 9600** compared to the primary predicate device K181136, the performance testing was conducted on the new cephalometric modality added in the proposed device **CS 9600**. Following the FDA Guidance for the Submission of 510(k)'s for Solid State X-ray imaging Devices, for the new cephalometric imaging applications, clinical images representative of the range of the different cephalometric radiological exams were taken. The cephalometric images were reviewed by a qualified expert and were evaluated to be of acceptable clinical effectiveness for the proposed indications for use. The **CS 9600** set of images were deemed to be of a clinically usable diagnostic quality.

Non clinical and bench testing was conducted as part of design control to ensure the substantial equivalence of **CS 9600** with the primary predicate device K181136 and reference device K151087. **CS 9600** has been tested to ensure that the system as a whole operated in a safe and effective manner that is substantially equivalent to the primary predicate device and for the cephalometric modality to the predicate device. The results of the performance testing support substantial equivalence.

# Traditional 510(k) Submission for CS 9600

## 9. Conclusion

The comparison of characteristics supports substantial equivalence. **CS 9600** is as safe and effective as the primary predicate device K181136 and as the reference device K151087 for the cephalometric modality.

**CS 9600** is considered substantially equivalent to the primary predicate device K181136 because both have the most similar indication for use and technological characteristics. Both have the same panoramic two-dimensional modality and three-dimensional 3D modality. Cephalometric modality added in the proposed device **CS 9600** is identical to the cephalometric modality included in previously cleared reference device K151087. The compared technical features for cephalometric imaging technology, Field of View, and other basic characteristics are matching very closely. There is no difference in the new cephalometric modality available in the proposed device **CS 9600** and the cephalometric modality available on the reference device K151087. Therefore, the addition of the cephalometric modality to the **CS 9600** does not raise questions of safety and effectiveness and does not have any effect on performance in practice.