

August 12, 2022

Shandong Huge Dental Material Corporation Maggie Zheng Regulatory Affairs Manager No. 68 Shanhai Road, Donggang District Rizhao City, Shandong 276800 China

Re: K220680

Trade/Device Name: Denture Base Polymers Regulation Number: 21 CFR 872.3760

Regulation Name: Denture Relining, Repairing, Or Rebasing Resin

Regulatory Class: Class II

Product Code: EBI Dated: April 27, 2022 Received: May 2, 2022

Dear Maggie Zheng:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Michael E. Adjodha, M.ChE.
Assistant Director
DHT1B: Division of Dental and
ENT Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K220680					
Device Name					
Denture Base Polymers					
Indications for Use (Describe)					
Denture Base Polymers is used for fabrication and repair of partial and full denture base for patients with missing teeth.					
Type of Use (Select one or both, as applicable)					
CONTINUE ON A SEPARATE PAGE IF NEEDED.					

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



005_510 (k) Summary K220680

This summary of 510(k) substantial equivalence information is being submitted in accordance with the requirements of 21 C.F.R. 807.92.

1. Date Summary Prepared: July 20, 2022

2. Submitter Information:

Name Shandong Huge Dental Material Corporation

Address No. 68 Shanhai Road, Donggang District, Rizhao City, Shandong

Province, 276800, P.R. China

Telephone 086-633-2277268

Contact Person Ms. Maggie Zheng

Contact Title Regulatory Affairs Manager

E-mail zhengxy@hugedent.com

3. Device Name

Trade name: Denture Base Polymers

Common name: Dental base material

Classification name: resin, denture, relining, repairing, rebasing (21 C.F.R. 872.3760)

Regulatory Class: II

Product Code: EBI

4. Substantially Equivalent Device

Company Name	Device Name	510 (k) NO.	Substantially Equivalent (SESE) Decision Date	Product code	Remarks
VERTEX-DENTAL BV	Vertex Rapid Simplified	K102654	12/03/2010	EBI	Primary Predicate
VERTEX-DENTAL BV	Vertex Self Curing, Vertex Castavaria, Vertex Castapress	K102640	03/25/2011	EBI	Secondary Predicate

These predicate devices have not been subject to a design-related recall.

5. Description of Device

Denture Base Polymers contains powder and liquid, wherein the power is mainly composed of polymethyl methacrylate, and the main composition of the liquid is methyl methacrylate.



Denture Base Polymers is used for fabrication and repair of partial and full denture base for patients with missing teeth, and its extensive indications can meet the needs of the production of complete denture.

6. Indications for use

Denture Base Polymers is used for fabrication and repair of partial and full denture base for patients with missing teeth.

7. Summary of Physical and Chemical Properties Tests

• Chemical Composition:

The device has similar chemical composition as the predicate devices. All of them are mainly composed of Polymethyl methacrylate and Pigments (Powder) and Methyl methacrylate, Ethyleneglycol dimethacrylate (Liquid).

• Technological characteristics:

The device has the same technological characteristics as the predicate devices (Heat Curing and Self Curing). And the device is similar in specifications and shades as the predicate devices.

Properties:

The device has comparable physical and chemical properties as the predicate devices. (Meeting the requirements of ISO standards for the Dental base material, ISO 20795-1)

• Applications:

The device has similar indications for use as the sum of the predicate devices: fabrication and repair of partial and full denture base.

8. Technological Characteristics:

The new device, Denture Base Polymers, has the same design, main materials and chemical composition as the predicate device.

	New Device	Primary Predicate	Secondary Predicate	
Comparison Items		Vertex Rapid Simplified	Vertex Self Curing, Vertex	
	Denture Base Polymers	vertex Rapid Simplified	Castavaria, Vertex Castapress	
	K220680	K102654	K102640	
Regulatory Classifications	same	same	same	
2) Indications for use	This product is used for fabrication and repair of partial and full denture base for patients with missing teeth.	The Vertex Rapid Simplified is indicated for: 1. Fabrication of full dentures 2. Fabrication of partial dentures	Vertex cold-curing denture base materials are indicated for: 1. Manufacture of full and partial dentures 2. Repair of full and partial dentures 3. Rebasing of full and partial dentures 4. Relining of full and partial	



Comparison Items	Comparison Items		New Device		Primary Predicate		Secondary Predicate			
3) Contraindications NA							Castavaria, Vertex Castapress			
3) Contraindications NA NA PMMA, EDMA PMMA, MMA, EDMA PMMA,				K22068	0	K1026	54		640	
Poperties Physical parameters Physica	3)	Contraindica	tions	NΔ		NT A				
FINMA, MMA, EDMA Type 1: \$10 per 15 o. 100 pug/m³ Type 2: \$100 pm³ Type 2: \$100 pm³ Type 1: \$100 p					EDIG					
Physical parameters			_	PMMA, MMA	, EDMA	PMMA, MM	IA, EDMA		A, EDMA	
MPa Type 2 : ≥ 60 MPa Type 2 : ≥ 80 Type 2 : ≥ 80 Type 2 : ≥ 1500 MPa Type			Physi	cal			Solubility	methyl methacrylate		
Type 1 New Device Properties Prope					MPa Type 2: ≥ 60	< 32 ug/mm ³	$\begin{array}{c c} ug/mm^3 \\ Type \ 2: \le 8 \end{array}$.0 2.2% Type 2: ≤	2000 MPa Type 2: ≥	
Type 2 Type 2 Criteria per ISO 20795-1 SO 20795-1	5)	· •	New	Type 1	criteria per ISC	criteria per ISO	Meet the criteria per IS	criteria O per ISO	criteria per ISO	
Primary Predicate (K102654)			Device		criteria per ISC	criteria per ISO	criteria per IS	criteria O per ISO	criteria per ISO	
Secondary Predicate (K102640)				•	criteria per ISC	criteria per ISO	criteria per IS	criteria O per ISO	criteria per ISO	
7) Target Population dental patients dental patients dental patients 8) Anatomical Site on teeth on teeth 9) Where Used used in hospital, dental clinic and relevant places dental professional dental professional 10) Human Factors dental professional dental professional dental professional 11) Design same same same 12) Cautions similar similar similar 13) Standards Met same same same 14) Biocompatibility biocompatible biocompatible 15) Sterility Non-sterile Non-sterile					criteria per ISC	criteria per ISO	criteria per IS	criteria O per ISO	criteria per ISO	
8) Anatomical Site on teeth on teeth on teeth 9) Where Used used in hospital, dental clinic and relevant places dental professional dental professional 10) Human Factors dental professional same same 12) Cautions similar similar similar 13) Standards Met same same same 14) Biocompatibility biocompatible biocompatible 15) Sterility Non-sterile Non-sterile Non-sterile On teeth on teeth used in hospital, dental clinic and relevant places dental professional used in hospital, dental clinic and relevant places and relevant places and relevant places are used in hospital, dental clinic and relevant places are used in hospital, dental professional are used in hospital, dental clinic and relev	6)				similar		similar			
9) Where Used used in hospital, dental clinic and relevant places dental professional dental professional dental professional dental professional same same same 12) Cautions similar same same same 13) Standards Met same same same 14) Biocompatibility biocompatible biocompatible biocompatible Non-sterile Non-sterile Non-sterile	7)	Target Popul	lation	tion dental patients		dental patients		dental patients		
9) Where Used used in hospital, dental clinic and relevant places dental professional dental professional dental professional dental professional 11) Design same same same 12) Cautions similar similar same same 13) Standards Met same same same 14) Biocompatibility biocompatible biocompatible 15) Sterility Non-sterile Non-sterile Used in hospital, dental clinic and relevant places dental professional same same same same 14) Non-sterile	8)	Anatomical	Site	-		on teeth		on teeth		
10) Human Factorsdental professionaldental professionaldental professional11) Designsamesame12) Cautionssimilarsimilarsimilar13) Standards Metsamesamesame14) Biocompatibilitybiocompatiblebiocompatiblebiocompatible15) SterilityNon-sterileNon-sterileNon-sterile		used in hospital, dental clinic					used in hospital, dental clinic			
12) Cautions similar similar 13) Standards Met same same 14) Biocompatibility biocompatible biocompatible 15) Sterility Non-sterile Non-sterile Non-sterile	10)	_			•					
13) Standards Metsamesame14) Biocompatibilitybiocompatiblebiocompatible15) SterilityNon-sterileNon-sterile				-		<u>-</u>				
14) Biocompatibilitybiocompatiblebiocompatible15) SterilityNon-sterileNon-sterile	12) Cautions similar			similar		similar				
15) Sterility Non-sterile Non-sterile Non-sterile	13) Standards Met same			same		same				
15) Sterility Non-sterile Non-sterile Non-sterile	<u> </u>		biocompat			biocompatible		biocompatible		
	,		_				-			
	16) Chemical Safety			similar						

9. Summary of Biocompatibility

The new device, Denture Base Polymers, is substantially equivalent to the predicate devices that have been on the market for years and with no clinical adverse events. The formulation of new device does not contain any new or non-conventional chemicals compared to the legally marketed predicate device.



We selected our Denture Base Polymers (Model: Type 2 Class 1) as the representative in biocompatibility tests and those biocompatibility test reports can be used in the biological evaluation of Denture Base Polymers.

Biocompatibility tests were performed fully following the ISO 10993 standards. The test items include Cytotoxicity; Sensitization; Irritation or Intracutaneous Reactivity; Acute Systemic Toxicity; Material-Mediated Pyrogenicity; Subchronic Toxicity; Genotoxicity and Implantation.

10. Clinical Performance Data

Not applicable. Clinical performance testing has not been performed on the subject device.

11. Summary of Substantial Equivalence

As with the comparison shown in substantial equivalence discussion, these devices are same or similar in almost all aspects. The details of indications for use, composition of materials, physical properties, labeling, stability/shelf life and cautions are slightly different, but the minor differences between the new device and the predicate devices are out of significance and do not affect neither the general intended use nor substantial equivalence.

We conclude that our Denture Base Polymers is substantially equivalent to the predicate devices described herein.