



July 19, 2022

Siemens Medical Solutions USA, Inc.  
% Patricia Jones  
Regulatory Affairs Professional  
40 Liberty Boulevard  
MALVERN PA 19355

Re: K221516

Trade/Device Name: ARTIS icono (VE23) System  
Regulation Number: 21 CFR 892.1650  
Regulation Name: Image-intensified fluoroscopic x-ray system  
Regulatory Class: Class II  
Product Code: OWB, IZI, JAA, JAK  
Dated: May 23, 2022  
Received: May 25, 2022

Dear Patricia Jones:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Laurel Burk, Ph.D.  
Assistant Director  
Diagnostic X-Ray Systems Team  
DHT8B: Division of Radiological Imaging Devices  
and Electronic Products  
OHT8: Office of Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K221516

Device Name  
ARTIS icono (VE23) System

### Indications for Use (Describe)

ARTIS is a family of dedicated angiography systems developed for single and biplane diagnostic imaging and interventional procedures including, but not limited to, pediatric and obese patients.

Procedures that can be performed with the ARTIS family include cardiac angiography, neuro-angiography, general angiography, rotational angiography, multipurpose angiography and whole body radiographic/fluoroscopic procedures as well as procedures next to the table for i.e. patient extremities. This does not include projection radiography.

Additional procedures that can be performed include angiography in the operating room, image guided surgery by X-ray, by image fusion, and by navigation systems. The examination table as an integrated part of the system can be used for X-ray imaging, surgery and interventions.

ARTIS can also support the acquisition of position triggered imaging for spatial data synthesis.

The ARTIS family include also the software option DynaCT with following indications for use:

DynaCT is an X-ray imaging software option, which allows the reconstruction of two-dimensional images acquired with a standard angiographic C-arm device into a three-dimensional image format.

DynaCT is intended for imaging both hard and soft tissues as well as other internal body structures for diagnosis, surgical planning, interventional procedures and treatment follow-up.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**K221516**

**510(k) Summary: ARTIS icono (VE23)**

**Company:** Siemens Medical Solutions USA, Inc.  
40 Liberty Boulevard, 65  
Malvern, PA 19355

**Date Prepared:** July 7, 2022

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR § 807.92.

**1. General Information:**

**Importer / Distributor:**

Siemens Medical Systems USA, Inc.  
40 Liberty Boulevard, 65-1A  
Malvern, PA 19355

**Establishment Registration Number:** 2240869

**Manufacturing Site:**

Siemens Healthcare GmbH  
Siemensstr. 1  
91301 Forchheim, Germany

**Establishment Registration Number:** 3004977335

**2. Contact Person:**

Ms. Patricia D. Jones  
Regulatory Affairs Professional  
Siemens Medical Solutions USA, Inc.  
40 Liberty Boulevard,  
Malvern, PA 19355  
Phone: (678) 575-8832  
Email: [patricia.jones@siemens-healthineers.com](mailto:patricia.jones@siemens-healthineers.com)

**3. Device Name and Classification:**

<b>Trade Name:</b>	<b>ARTIS icono (VE23) System</b>
<b>Classification Name:</b>	Image-intensified fluoroscopic x-ray System
<b>Classification Panel:</b>	Radiology
<b>Regulation Number:</b>	21 CFR §892.1650
<b>Device Class:</b>	Class II
<b>Product Codes:</b>	OWB, IZI, JAA, JAK

**4. Legally Marketed Primary Predicate Device**

<b>Trade Name:</b>	<b>ARTIS icono (VE20) System</b>
<b>510(k) Clearance</b>	K190768
<b>Clearance Date</b>	September 12, 2019
<b>Classification Name:</b>	Image-intensified fluoroscopic x-ray System
<b>Classification Panel:</b>	Radiology
<b>Regulation Number:</b>	21 CFR §892.1650

**Device Class:** Class II  
**Product Code:** OWB  
**Subsequent Product Codes:** IZI, JAA, JAK  
**Total Product Life Cycle:** All product Recall incidents are considered during the Design Input phase of development to ensure the latest models will not be affected by any of the applicable issues.

**Secondary Predicate Device:**  
**Trade Name:** Artis Q/Q.zen System  
**510(k) Clearance:** K181407  
**Clearance Date:** August 15, 2018  
**Classification Name:** Image-intensified fluoroscopic x-ray System  
**Classification Panel:** Radiology  
**Regulation Number:** 21 CFR §892.1650  
**Device Class:** Class II  
**Product Codes:** OWB  
**Subsequent Product Codes :** IZI, JAA, JAK  
**Total Product Life Cycle:** All product Recall incidents are considered during the Design Input phase of development to ensure the latest models will not be affected by any of the applicable issues.

## 5. Device Description:

The ARTIS icono (VE23) system is a medical device that allows visualization of vessels within the human body. It is of the utmost importance to find the right projections so physician can navigate catheters and other devices safely. The ARTIS icono (VE23) system consist of a patient table and a multi-axis motorized c-arm that can be positioned around the patient and angulated in a double-oblique fashion iso-centering the region of interest between the x-ray tube and the flat panel detector. The x-ray generator is placed separately. The displays for visualizing the x-ray images are mounted at the ceiling with a movable display suspension system. System operation is executed via control modules table side so that the physician can move and position the table and c-arm adequately for best imaging while manipulate the catheters or other devices during x-ray. X-ray release is table side via a footswitch.

The ARTIS icono (VE23) ceiling mounted system covers the complete range of angiographic applications, cardiac angiography, neuro-angiography, general angiography, surgery and surgical angiography, multipurpose angiography, rotational angiography, and radiographic/fluoroscopic procedures.

The following components are configured to create the Ceiling configuration:

- (1) Ceiling Mounted C-arm, X-ray tube assembly and FD
- (2) Patient table
- (3) Display ceiling suspension with displays
- (4) Footswitch for releasing radiation

- (5) Control console for controlling the stand, patient table and imaging system
- (6) Ceiling Rails (two configurations)

Images and operating elements are displayed on screens. Depending on the ARTIS icono (VE23) system configuration, different display variants are used to visualize image and information content. Displays that visualize single images or large displays that are configurable to visualize multiple images and information content in various layouts are used.

Post processing can be done in the exam room or in the control room that offers monitors as well, with a footswitch location in the exam room or the control room. The ARTIS icono (VE23) ceiling System is capable of 2D and 3D imaging.

Other systems and software *syngo* Application Software, *syngo* X Workplace, Sensis, and or third-party systems may also be integrated into the ARTIS icono (VE23) screen configuration. Different screen configurations and layouts are possible in the examination room and in the control room.

The 510(k) submission, Subject device ARTIS icono with software version VE23 will support the following modifications made to the Subject Device in comparison to the Predicate Device:

#### **List of Modifications**

1. Updated system Software from VE20 to VE23
  - A. NOMSIE DSA (IQ) Customization of Contrast in DSA Images
  - B. New Generator (Polydoros ACX)
  - C. Peristepping / Perivision
2. Ceiling configuration, new ceiling mounted C-arm
3. Product Claims List
4. Update 510(k) Information

#### **6. Indications for Use:**

ARTIS is a family of dedicated angiography systems developed for single and biplane diagnostic imaging and interventional procedures including, but not limited to, pediatric and obese patients.

Procedures that can be performed with the ARTIS family include cardiac angiography, neuro-angiography, general angiography, rotational angiography, multipurpose angiography, and whole body radiographic/fluoroscopic procedures as well as procedures next to the table for i.e., patient extremities. This does not include projection radiography.

Additional procedures that can be performed include angiography in the operating room, image guided surgery by X-ray, by image fusion, and by navigation systems. The examination table as an integrated part of the system can be used for X-ray imaging, surgery, and interventions.

ARTIS can also support the acquisition of position triggered imaging for spatial data synthesis.

The ARTIS family include also the software option DynaCT with following indications for use:

DynaCT is an X-ray imaging software option, which allows the reconstruction of two-dimensional images acquired with a standard angiographic C-arm device into a three-dimensional image format.

DynaCT is intended for imaging both hard and soft tissues as well as other internal body structures for diagnosis, surgical planning, interventional procedures, and treatment follow-up.

**7. Substantial Equivalence:**

The ARTIS icono (VE23) System is substantially equivalent to the following predicate devices comparable properties provided in the below Table 1,

**Table 1: Predicate Device Comparable Properties for Subject Device Modifications:**

Predicate Device Name and Manufacturer	510(k) Number	Clearance Date	Comparable Properties
<b>Primary Predicate</b>	K190768	09/12/2019	<ul style="list-style-type: none"> <li>• Indications for use</li> <li>• Software Version VE20</li> <li>• Image Quality</li> <li>• Generator</li> </ul>
ARTIS icono (VE20) Siemens			
<b>Secondary Predicate</b>	K181407	08/15/2018	<ul style="list-style-type: none"> <li>• Ceiling configuration</li> <li>• Peri-Stepping &amp; Perivision</li> </ul>
Artis Q/Q.zen (Ceiling configuration) Siemens			

**8. Summary of Technological Characteristics of the Subject Device as Compared with the Predicate Device:**

The ARTIS icono (VE23) ceiling System is designed as a set of components (C-arm, X-ray tube and housing, flat panel detector, digital imaging system, collimator, generator etc.) to provide specialized angiography systems. Components used with ARTIS icono (VE23) System are either commercially available with current Siemens systems or include modifications to existing components. Technological differences between the Subject Device and the Predicate Device are provided in the **Table 2** below for all modifications.



**Table 2: Summary of Comparison of Technological Characteristics**

Modifications	Subject Device ARTIS icono (VE23)	Primary Predicate Device ARTIS icono (VE20) K190768	
<b>System Software</b>	1. Updated system Software from VE20 to VE23	System Software Version VE20	
	A. NOMSIE DSA (IQ) Customization of Contrast in DSA Images	Classical edge enhancement with filters and automatic windowing was used.	
	B. New Generator (Polydoros ACX)	Generator (Polydoros A100G)	
<b>Comparison Results</b>	<b>Comparable:</b> System software VE23 was updated to support NOMSIE DSA (IQ) Customization of Contrast in DSA Images and a new Generator. System Software modifications conforms to “Guidance for the content of Premarket submission for software in Medical Devices”. Provided in Volume 016 is all required software testing. Provided in Volume 020 are Bench Test Summaries and Test Reports. These software and hardware changes do not raise any new risks regarding the safety or effectiveness of the device. All test results met all acceptance criteria.		
Modifications	Subject Device ARTIS icono (VE23)	Predicate Device Artis Q/Q.zen (K181407)	Comparison Results
<b>System Software</b>	C. Added Peristepping / Perivision feature	Peristepping / Perivision	<b>Same:</b>
<b>Ceiling Configuration:</b>	2. Ceiling configuration, new ceiling mounted C-arm	Ceiling Configuration	<b>Comparable:</b> Nonclinical tests have been performed to assess the functionality of the subject device. Results of all conducted testing were found acceptable and do not raise any new issues of safety or effectiveness. All test results met all acceptance criteria.
<b>Claims:</b>	3. Product Claims List		
<b>Update 510(k) Information:</b>	4. Updated 510(k) information for predicate devices.		

**9. Nonclinical Performance Testing:**

Non-clinical tests were conducted for the ARTIS icono (VE23) ceiling system during product development.

The ARTIS icono (VE23) was certified by Siemens Healthcare GmbH Corporate Testing Laboratory to comply with the following standards for Electrical safety, performance, and Electromagnetic Compatibility:

- AAMI ANSI ES60601-1:2005/(R)2012
- IEC 60601-1-2:2020
- IEC 60601-1-3:2013
- IEC 60601-1-6:2020
- IEC 60825-1:2021 (Recognized 2007)
- TR 60878:2015



- IEC 62304:2015
- IEC 80001-1:2010
- IEC 60601-2-28:2017
- IEC 60601-2-43:2019
- IEC 60601-2-54:2018
- ISO 10993-1:2018
- ISO 14971:2019

The modifications described in this Premarket Notification are supported with verification and validation testing.

**Verification and Validation:**

Software Documentation for a Moderate Level of Concern software per FDA's Guidance Document "Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices" issued on May 11, 2005, and "Off-The-Shelf Software Use in Medical Devices" is also included as part of this submission. The performance data demonstrates continued conformance with special controls for medical devices containing software. Non-clinical tests were conducted on ARTIS icono System software (VE23 during product development.

The Risk analysis was completed, and risk control implemented to mitigate identified hazards. The testing results support that all the software specifications have met the acceptance criteria. Testing for verification and validation for the device was found acceptable to support the claims of substantial equivalence.

ARTIS icono System software (VE23 was tested and found to be safe and effective for intended users, uses, and use environments through the design control verification and validation process. The Human Factor Usability Validation showed that Human factors are addressed in the system test according to the operator's manual. Customer employees are adequately trained in the use of this equipment.

Siemens conforms to the cybersecurity requirements by implementing a process of preventing unauthorized access, modifications, misuse or denial of use, or the unauthorized use of information that is stored, accessed, or transferred from a medical device to an external recipient. Provided in this submission is a cybersecurity statement that considers IEC 80001-1:2010. The responsibility for compliance with IEC 80001-1-2010 is the hospital. Provided in the Software Section is the required cybersecurity information.

**Summary:**

Performance tests were conducted to test the functionality of ARTIS icono (VE23 System. These tests have been performed to assess the functionality of the subject device. Results of all conducted testing and clinical assessment were found acceptable and do not raise any new issues of safety or effectiveness.

**10. General Safety and Effectiveness Concerns:**

Instructions for use are included within the device labeling, and the information provided will enable the user to operate the device in a safe and effective manner.

Risk management is ensured via a hazard analysis, which is used to identify potential hazards. These potential hazards are controlled via software development, verification, and validation testing. To minimize electrical, mechanical and radiation hazards, Siemens adheres to recognized and established industry practice, and all equipment is subject to final performance testing. Furthermore, the operators are health care professionals familiar with and responsible for the evaluating and post processing of X-ray images.

**11. Conclusion as to Substantial Equivalence:**

The predicate devices were cleared based on non-clinical supportive information. Similar non-clinical test results demonstrate that the ARTIS icono (VE23) Ceiling System acceptance criteria are adequate for the intended use of the device. The comparison of technological characteristics, non-clinical performance data and software validation data demonstrates that the Subject Device is as safe and effective when compared to the Predicate Devices that is currently marketed for the same intended use.