

October 26, 2022

Humanwell Healthcare Group Medical Supplies Co., Ltd. % Ivy Wang
Technical Manager
Shanghai Sungo Management Consulting Company Limited
14th F, 1500# Century Avenue
Shanghai, 200122
China

Re: K222498

Trade/Device Name: Medical Examination Gloves (Nitrile) (XS, S, M, L, XL)

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA Dated: August 18, 2022 Received: August 18, 2022

#### Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

K222498 - Ivy Wang Page 2

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (<a href="DICE@fda.hhs.gov">DICE@fda.hhs.gov</a>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bifeng Qian, M.D., Ph.D.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

## DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023
See PRA Statement below.

Submission Number (if known)
K222498
Device Name
Medical Examination Gloves (Nitrile) (XS, S, M, L, XL)
Indications for Use (Describe)
The Medical Examination Glove (Nitrile) is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.
Type of Use (Select one or both, as applicable)
Prescription Use (Part 21 CFR 801 Subpart D)
CONTINUE ON A SEPARATE PAGE IS NEEDED

# This section applies only to requirements of the Paperwork Reduction Act of 1995. \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

#### 510(K) Summary

(As requirement by 21 CFR 807.92)

Date prepared: 2022-10-25

#### A. Applicant:

Name: Humanwell Healthcare Group Medical Supplies Co., Ltd.

Address: 169 Shendun Yi Lu, East Lake Hi-tech Development Zone, Wuhan, Hubei Province, People's Republic

of China

Contact: Ms. Aixue Yao
Title: Sales Manager
Tel: +86- 18627726855
Email: 674599374@qq.com

Submission Correspondent: Primary contact: Ms. Ivy Wang

Shanghai SUNGO Management Consulting Co., Ltd.

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: <a href="mailto:haiyu.wang@sungoglobal.com">haiyu.wang@sungoglobal.com</a> Secondary contact: Mr. Raymond Luo

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: fda.sungo@gmail.com

#### B. Device:

Trade Name: Medical Examination Gloves (Nitrile) (XS, S, M, L, XL) Common Name: Nitrile Patient Examination Gloves (Powder Free)

#### **Regulatory Information**

Classification Name: Polymer Patient Examination Glove

Classification: Class I Product code: LZA

Regulation Number: 21 CFR 880.6250

Review Panel: General Hospital

#### C. Predicate device:

K171422

Disposable Powder Free Nitrile Examination Glove, White/Blue/ Black/ Pink Color

Ever Global (Vietnam) Enterprise Corp

#### D. Indications for use of the device:

The Medical Examination Glove (Nitrile) is a disposable device intended for medical purposes that is worn on

the examiner's hand to prevent contamination between patient and examiner.

#### E. Device Description:

The Medical Examination Gloves (Nitrile) are non-sterile, single use only, disposable examination gloves intended for medical purposes to be worn by examiners to prevent contamination between the patient and the examiner. The gloves are blue color, powder free, nitrile ambidextrous gloves. The gloves are offered in sizes extra-small, small, medium, large, extra-large packed in a paper box.

The gloves are designed and manufactured in accordance with the ASTM D6319-19 standard.

#### F. Summary of Technological Characteristics

Table 1 General Comparison of Proposed and Predicate Devices

Device	Proposed Device	Predicate Device	Result
510K#	K222498	K171422	-
Product Name Medical Examination Gloves (Nitrile)		Disposable Powder Free Nitrile Examination Glove, White/Blue/Black/ Pink Color	-
Product Code	LZA	LZA	Same
Regulation Number	21 CFR 880.6250	21 CFR 880.6250	Same
Indications for use	The Medical Examination Glove (Nitrile) is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.	The Disposable Powder Free Nitrile Examination Glove, White/ Blue/ Black/ d Pink Color is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.	
Powder free	Yes	Yes	Same
Design feature	Ambidextrous	Ambidextrous	Same
Material	Nitrile	Nitrile	Same
OTC use	Yes	Yes	Same
Sterility	Non-sterile	Non-sterile	Same
Use	Singe use	Single use	
Label	Single-use, indication, powder free, device color, device name, glove size and quantity, Nitrile Glove Powder Free Blue, Non-Sterile	Single-use, indication, powder free, device color, device name, glove size and quantity, Disposable Powder Free Nitrile Examination Glove, Non-Sterile	Same

Table 2 Device dimension comparison

Duodiaata	Dosignation	Size				Tolerance	
Predicate	Designation	XS	S	М	L	XL	Tolerance
device (K171422)	Length, mm	230	230	230	230	230	min
	Width, mm	75	85	90	105	115	±5

	Thickness, mm						
	Finger		0.05				
	Palm		0.05				
	Designation	Size					Talawanaa
	Designation	XS	S	М	L	XL	Tolerance
Droposod	Length, mm	220	220	230	230	230	min
Proposed device	Width, mm	70	80	95	110	120	±10
device		Thickness, mm					
	Finger		0.05 min				
	Palm		0.05				
Result			Similar				

Analysis: The physical dimensions are different with that of the predicate, but they all meet the requirements of ASTM D6319-19, so the differences do not raise any new safety or performance questions.

Table 3 Performance comparison

Item			Proposed	Predicate device	Result
			device	(K171422)	
Colorant			Blue	White/ Blue/ Black/ Pink	Similar
Physical	Before aging	Tensile strength	14MPa, min	14MPa, min	Same
properties		Ultimate elongation	500%, min	500%, min	Same
	After aging	Tensile strength	14MPa, min	14MPa, min	Same
		Ultimate elongation	400%, min	400%, min	Same
	Comply with ASTM D6319			Comply with ASTM	Same
				D6319	
Freedom fro	m holes		Be free from holes	Be free from holes when	Same
			when tested in	tested in accordance	
			accordance with	with ASTM D5151	
			ASTM D5151	G-1, AQL 2.5	
			G-1, AQL 2.5		
Residual Powder			Meet the	Meet the requirements	Same
			requirements of	of ASTM D6124	
			ASTM D6124		

Analysis: The subject device is only available in a single color (blue,) however, the predicate is available in multiple colors (white, blue, black, pink). Biocompatibility testing was successfully completed for the subject device, demonstrating that any color differences do not affect the safety of the proposed device.

Table 4 Biocompatibility comparison

Item		Proposed device	Predicate device (K171422)	Result
Material		Nitrile	Nitrile	Same
Biocompatibility	Irritation	Under the	Comply with ISO	Same
	ISO 10993-10	conditions of the	10993-10	

	study, not an		
	irritant. Comply		
	with ISO 10993-10.		
Sensitization	Under the		
ISO 10993-10	conditions of the		
	study, not a		
	sensitizer. Comply		
	with ISO 10993-10.		
Acute systemic	Under the	/	Different
toxicity	conditions of the		
ISO 10993-11	study, the device is		
	non-toxic.		

Analysis: Toxicity information for the predicate device is not publicly available. This does not raise different safety or performance questions since the subject device has acceptable biocompatibility per the biocompatibility endpoint assessment.

#### **G.** Summary of Non-Clinical Testing

#### Biocompatibility

The following tests for the subject device were conducted to evaluated the biocompatibility of Medical Nitrile Examination Gloves:

- ISO 10993-10: Primary Skin Irritation
- ISO 10993-10: Dermal Sensitization
- ISO 10993-11: Acute Systemic Toxicity

Test	Purpose	Acceptance Criteria	Results
	<del> </del> '	'	
Irritation	The purpose of the test is to verify	Under the conditions of the	Pass
ISO 10993-10	that the proposed device extract is	study, not an irritant.	
	non-irritating.		
	g.		
Sensitization	The purpose of the test is to verify	Under the conditions of the	Pass
ISO 10993-10	that the proposed device extract is	study, not a sensitizer.	
	non-sensitizing		
Acute systemic	The purpose of the test is to verify	Under the conditions of the	Pass
toxicity	that the proposed device extract is	study, the device is non-toxic.	
ISO 10993-11	non-cytotoxic.		

#### Performance Testing

Physical performance testing of the proposed device was conducted as per ASTM D6319-19 *Standard Specification for Nitrile Examination Gloves for Medical Application*.

To summarize, the performance testing of the subject device was conducted to adequately demonstrate the effectiveness of the device in accordance with the relevant test methods cited below:

- ASTM D6319-19 Standard Specification for Nitrile Examination Gloves for Medical Application
- ASTM D6124-06 (Reapproved 2017) Standard Test Method for Residual Powder on Medical Gloves

### K222498

#### • ASTM D5151-19 Standard Test Method for Detection of Holes in Medical Gloves

Test Method	Purpose	Acceptance Criteria	Results
Dimensions (width)	The purpose of the	Width 70mm min	Pass
(thickness)	test is to evaluate	Length 220mm min	76mm min width
	the physical		234mm min length
	dimension of the	Palm – 0.05mm min.	Pass
	glove	Finger–0.05mm min.	Palm – 0.08mm min.
			Finger–0.09mm min
Physical properties	The purpose of the	Tensile Strength:	Pass
	test is to evaluate	Before Aging ≥14 MPa, min.	Tensile Strength:
	the tensile strength	After Aging ≥14 MPa, min.	Before Aging 24.44 MPa, min.
	and ultimate	Elongation:	After Aging 27.21 MPa, min.
	elongation before	Before Aging 500%, min.	Elongation:
	and after aging	After Aging 400%, min.	Before Aging 513%, min.
			After Aging 541%, min.
Freedom from holes	The purpose of the	No leakage at sampling level of	Pass
	test is to detect	G-1, AQL 2.5	No leakage, 80 of 80 passed of
	holes in the gloves		each size
Residual Powder	The purpose of the	<2mg per glove	Pass
	test is to detect the		average 0.84 mg per glove
	powder residue in		
	the glove		

#### **H.** Clinical Test Conclusion

No clinical study is included in this submission.

#### I. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission, the Medical Examination Gloves (Nitrile) are as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K171422.