



Bringing new life to organ transplantation™

## **OCS™ Lung System for the Preservation of Donor Lungs for Transplantation**

---

**May 17, 2017**

TransMedics, Inc.

Gastroenterology-Urology Devices Panel

## Introduction

---

**Waleed Hassanein, MD**

President and CEO






TransMedics, Inc.



# TransMedics is a Clinically Driven Organization


- Origin of OCS Technology: 1995-98 academic cardiothoracic surgery research project
- TransMedics founded in 2000: developed Organ Care System (OCS™) technology to maintain human solid organs in near-physiologic and functioning state to overcome limitations of cold ischemic storage – OCS Lung, Heart & Liver
- OCS platform approved outside of U.S. with ~800 successful human transplants performed globally to date in standard, extended, and DCD organ criteria

# Establishing High Level of Clinical Evidence – OCS Global Clinical Programs

	LUNG	HEART	LIVER
Standard Criteria	 <p>OCS™ Lung <b>INSPIRE</b> Trial</p>	 <p>OCS™ Heart <b>PROCEED II</b> Trial</p>	 <p>OCS™ Liver <b>PROTECT</b> Trial</p>
Extended Criteria	 <p>OCS™ Lung <b>EXPAND Lung</b> Trial</p>	 <p>OCS™ Heart <b>EXPAND Heart</b> Trial</p>	

# OCS System Designed to Address Limitations of Cold Ischemic Storage

**REDUCE ISCHEMIC INJURY**



**Warm Oxygenated  
Blood Perfusion**

**OPTIMIZE ORGAN CONDITION**



**Ventilation  
Recruitment**

**EX-VIVO FUNCTIONAL ASSESSMENT**



**Oxygenation, Vascular Resistance,  
& Airway Compliance**

# OCS Lung System: Integrated, Portable, Ex-vivo Lung Perfusion and Ventilation System



OCS Lung Console



OCS Lung Perfusion Module



OCS Lung Solution



# INSPIRE Demonstrated Assurance of Safety and Effectiveness of OCS Lung System

- Met primary effectiveness and safety endpoints
- Clinically significant reduction in PGD Grade 3 within 72 hours
- Similar safety profile of OCS Lung System to standard of care
- Other clinical benefits that will be further studied in post-market

## Proposed Indication for Use

The TransMedics<sup>®</sup> Organ Care System<sup>™</sup> (OCS) Lung System is a portable organ perfusion, ventilation, and monitoring medical device intended to preserve donor lungs in a near physiologic, ventilated, and perfused state for transplantation.

# Agenda

## Clinical Needs and Current Limitaitons of Cold Storage Preservation

### **Gabriel Loor, MD**

Associate Professor of Surgery; Director, Lung Transplantation  
Baylor College of Medicine

## Regulatory History

### **Waleed Hassanein, MD**

## INSPIRE Trial Design

### **Abbas Ardehali, MD**

Professor of Surgery and Medicine; Chief of Cardiothoracic Transplantation  
UCLA School of Medicine

## INSPIRE Trial Adjudication

### **John Wallwork, FRCS, FmedSCI**

Emeritus Professor, Cardiothoracic Surgery Papworth Hospital  
Past President, International Society for Heart and Lung Transplant (ISHLT)

## INSPIRE Trial Results

### **Gregor Warnecke, MD**

Prof. of Surgery; Director of Cardiothoracic Transplantation  
Hannover Medical School

## Training and Post-Market Studies

### **Waleed Hassanein, MD**

## Benefit-Risk Assessment

### **Dirk Van Raemdonck, MD, PhD**

Professor of Surgery, University Hospital Leuven Medical Center  
Co-chair of ISHLT PGD Working Group



# Additional Experts

## DSMB Chairman

### **Joshua Sonett, MD**

Professor of Surgery  
Chief, Thoracic Surgery  
Columbia University Medical Center

## Biostatistics

### **Christopher Mullin, MS**

Biostatistician  
3D Communications, LLC

## Device Design and Engineering

### **John Sullivan, MS**

Vice President, Engineering  
TransMedics, Inc.

## Clinical Operations

### **Tamer Khayal, MD**

Chief Medical Officer  
TransMedics, Inc.

## Regulatory

### **Miriam Provost, PhD**

**Christine Brauer, PhD**

**Robert Sheridan**

# Clinical Needs and Current Limitations of Cold Storage Preservation

---

## **Gabriel Loor, MD**

Associate Professor, Department of Surgery

Surgical Director of Lung Transplantation

Division of Cardiothoracic Transplantation and Circulatory Support

Michael E. DeBakey Department of Surgery

Baylor College of Medicine

# Transplantation is Gold Standard for Treating End-stage Lung Failure

- Without transplant, <50% patients alive in 1-2 years
- Lung transplantation provides:
  - Longer life expectancy
  - Improved functional status
  - Better quality of life

# Challenges in Lung Transplant Today

- Organ availability
- Older and sicker patients
- Preservation limitation and transplant logistics
- Primary Graft Dysfunction (PGD)
- Bronchiolitis Obliterans Syndrome (BOS)

## No Advancements in Organ Preservation for 30 Years

- In last 30 years, many advancements in lung transplantation:
  - Surgical techniques
  - Pre- and peri-operative care of recipients
  - Immunosuppressives
- No advancements in organ preservation beyond cold storage since dawn of organ transplantation


## Three Key Limitations of Cold Ischemic Storage



- Time-dependent ischemia / reperfusion injury
- No lung optimization capabilities
- No assessment of lung function

# Clinical Constraints of Ischemia in Lung Transplantation

Most U.S. Lung Transplants are <6 Hours



## Primary Graft Dysfunction (PGD) is Acute Lung Injury Associated with Reperfusion Injury

- PGD can occur within first 72 hours after transplant
  - Assessed at T0, 24, 48, and 72 hours post-transplant
- Short-term morbidity associated with PGD
  - Severe hypoxemia, lung edema, difficulty with ventilation, etc.
- ISHLT PGD Grading from 0-3 (0 = absent to 3 = severe)



# Reported PGD3 Incidence of 30.8% Within Initial 72 Hours

## Clinical Risk Factors for Lung Transplantation

Joshua M. Diamond<sup>1</sup>, James C. Lee<sup>1</sup>,  
Scarlett L. Bellamy<sup>2</sup>, David J. Lederer<sup>3</sup>,  
Sangeeta M. Bhorade<sup>8</sup>, Maria Crespo<sup>4</sup>,  
Jonathan Orens<sup>12</sup>, Ashish S. Shah<sup>13</sup>,  
David S. Wilkes<sup>15</sup>, Lorraine B. Ware<sup>14</sup>  
for the Lung Transplant Outcomes Group

<sup>1</sup>Pulmonary, Allergy, and Critical Care Division,  
Cardiovascular Surgery, and <sup>6</sup>Department of  
Philadelphia, Pennsylvania; <sup>4</sup>Division of Pulmonary  
of Physicians and Surgeons, New York, New York;  
Michigan; <sup>8</sup>Division of Pulmonary and Critical

ing results. Some explanations for these  
sample sizes; inconsistencies in PGD phenom-  
enon for multiple confounding variables; a  
prospective, single center, or administrative  
rigorous PGD definitions (5, 6).

In 2005, the International Society for  
plantation (ISHLT) standardized the PGD  
research on risk factors associated with the  
syndrome (7). Subsequent studies have  
construct validity of this definition with clinical  
logic markers of ALI severity (8, 9). In this  
identify donor, recipient, and perioperative  
using the ISHLT definition in a large, multicenter, prospective  
cohort study design.

## METHODS



### Study Design and Subject Selection

## RESULTS


There were 2,011 lung and heart-lung transplants performed at study centers during the study period. Of these, 1,255 patients were enrolled in the cohort study (Figure 1). There were no significant differences in sex or age, but there was more chronic obstructive pulmonary disease, less cystic fibrosis, and more SLT in the enrolled group (*see* Table E3). A total of 211 subjects (16.8%; 95% CI, 14.7–18.9) met criteria for grade 3 PGD, and 386 subjects (30.8%; 95% CI, 28.2–33.3) met the secondary PGD definition of grade 3 PGD at any time during the first 72 hours after transplantation.

were enrolled in the cohort study (Figure 1). There were no significant differences in sex or age, but there was more chronic obstructive pulmonary disease, less cystic fibrosis, and more SLT in the enrolled group (*see* Table E3). A total of 211 subjects (16.8%; 95% CI, 14.7–18.9) met criteria for grade 3 PGD, and 386 subjects (30.8%; 95% CI, 28.2–33.3) met the secondary PGD definition of grade 3 PGD at any time during the first 72 hours after transplantation.


# PGD3 Within First 48 Hours Correlates with Lower Long-term Survival, Higher BOS Rate



# PGD 3 at T0 Correlates with Long-Term BOS Rates



# PGD 2 or 3 at T0 Significant Risk Factor for BOS



# Clinical Need for Advancements in Lung Preservation for Transplantation

- Lung transplantation is gold standard for end-stage lung disease
- Lung preservation limited to cold ischemic storage for past 30+ years with inherent limitations
- PGD3 at any time point within 72 hours associated with poor patient outcomes
- Need for technology to improve lung preservation
  - Minimize ischemic injury
  - Optimize and assess lung during preservation

# Regulatory History and Protocol Design of OCS Lung INSPIRE Trial

---

**Waleed Hassanein, MD**


President and CEO

TransMedics, Inc.

# INSPIRE Trial Regulatory Background

- First RCT for lung preservation for transplantation
- Several challenges to be addressed in study design:
  - Endpoints, analysis population, timing of evaluations, etc.
  - Complex organ allocation and retrieval process

# INSPIRE Trial Regulatory Timeline Summary



## 24 Months of Complex Negotiations

- FDA outlined conditions for IDE approval
- TransMedics agreed to conditions to initiate INSPIRE

## Rationale for Protocol Amendment

- Published literature on early PGD
- Successful EXPAND Trial appeal to ODE on scientific merits of PGD3 within 72 hours<sup>1-4</sup>

1. Whitson BA et al., JHLT 2007; 26:1004-1011  
2. Daud et al., AMJ. Resp. & Crit. Care Med. 2007

3. Christie J et al., JHLT 2010; 29:2131-2137  
4. Huang H et al., AJT 2008: 245402462



# Protocol Design Topics for Clarification

**Protocol Non-Inferiority Design and Margin**


**Rationale for Protocol Amendment**

**Rationale for Administrative Extension Cohort**

# Non-Inferiority Trial Design

- Rationale for INSPIRE non-inferiority trial design:
  - Very common pivotal FDA trial design for device approval
  - Maintain current success rate of lung transplantation

# FDA Required 4% Non-Inferiority Margin



## Interpretation of Results with Conservative 4% NI Margin

- To our knowledge, 4% NI margin is narrowest used in pivotal device trial
- OCS had to perform at least 4-5% better than Control to meet NI margin

# Protocol Design Topics for Clarification

**Protocol Non-Inferiority Design and Margin**


**Rationale for Protocol Amendment**

**Rationale for Administrative Extension Cohort**

## INSPIRE Protocol Amendment

- Allow PGD component of the primary endpoint to comprehensively assess PGD3 within all time points rather than at single point post-transplantation (T72)
- Re-designate the Per-Protocol as primary analysis population for effectiveness

# PGD3 Within 72 Hours Clinically Appropriate Endpoint for Preservation Technology



- PGD3 within 72 hours is comprehensive and robust assessment of PGD3 post-lung transplantation<sup>1-4</sup>
- Captures early timepoints that may be impacted by preservation injury as compared at only at T72
- PGD assessed at every timepoint throughout INSPIRE Trial

## Per Protocol Analysis is Clinically Appropriate Population

- TransMedics consistently maintained that PP was appropriate primary analysis population:
  - FDA guidance on non-inferiority trial analysis<sup>1</sup>
  - PP assesses treatment effect when OCS and cold storage used as intended, eliminating confounding variables:
    - Not treated as randomized (e.g. OCS recipient transplanted using cold storage preserved lungs)
    - Major protocol violations (e.g. donor lung with pneumonia)

<sup>1</sup> US FDA - E9 Statistical Principles for Clinical Trials



# Protocol Design Topics for Clarification

**Protocol Non-Inferiority Design and Margin**

**Rationale for Protocol Amendment**

**Rationale for Administrative Extension Cohort**

## INSPIRE Trial Cohorts



- IDE approved 2 perfusion solutions in the OCS arm: OCS Lung Solution and LPD Solution
- Several investigators observed lung edema during preservation using LPD. TMDX notified FDA of this observation to seek advice
- Agreement to file an administrative extension (Admin Ext) to allow time to define plan
- OCS Solution subgroup identified as an important adjunct analysis

# Comprehensive Data Presentation

- **Cohorts**
  - **INSPIRE Cohort** (N=320, pre-specified sample size)
  - **Combined Cohort** (N=349, INSPIRE + administrative extension)
- **Effectiveness Analysis Populations**
  - Per-protocol (primary)
  - Modified ITT (supportive)
- **Results presented both Overall and for OCS Solution Subgroup**

## **INSPIRE Trial Design**

---

**Abbas Ardehali, MD**

Professor of Surgery and Medicine

William E. Connor Chair in Cardiothoracic Transplantation

Director, Heart and Lung Transplant Center

UCLA Medical Center

# Donor Eligibility Criteria Reflect Standard Lung Transplantation

Inclusion	Exclusion
<ul style="list-style-type: none"><li>▪ Age &lt;65 years old</li><li>▪ Normal gas exchange [<math>\text{PaO}_2 / \text{FiO}_2 \geq 300</math>] at time of final acceptance of donor lung</li><li>▪ No active lung disease</li><li>▪ Lung suitable for both OCS or cold storage</li></ul>	<ul style="list-style-type: none"><li>▪ Presence of moderate to severe traumatic lung injury</li><li>▪ Presence of confirmed active pneumonia</li><li>▪ Positive serology (Hep. B/C, HIV etc.)</li></ul>

# Recipient Eligibility Criteria Reflect Standard Lung Transplantation

Inclusion	Exclusion
<ul style="list-style-type: none"><li>▪ Registered double-lung transplant candidate</li><li>▪ Age <math>\geq</math> 18 years old</li></ul>	<ul style="list-style-type: none"><li>▪ Prior solid organ or bone marrow transplant</li><li>▪ Multi-organ transplant recipient</li><li>▪ Single lung recipient</li><li>▪ Chronic renal failure</li></ul>

# Primary Effectiveness Endpoint and Safety Endpoint

- Primary effectiveness composite endpoint (NI margin = 4%)
  - All-cause survival post transplant at day 30 and absence of PGD Grade 3 within first 72 hours
- Safety endpoint: mean # of lung-graft-related SAEs through 30 days post transplant (NI margin = 0.07 events)
  - Moderate to severe acute rejection
  - Respiratory failure
  - Bronchial anastomotic complications
  - Lung related infections
- 30-day window relevant to assessing preservation-related issues as compared to later timepoints which could be impacted by other variables

# Secondary and Other Clinical Endpoints

- Secondary endpoints:
  - PGD Grade 3 at 72 hours (NI margin = 5%)
  - PGD Grade 2 or 3 at 72 hours (NI margin = 7.5%)
  - Patient survival at day 30 (NI margin = 4%)
  
- Other endpoints:
  - Bronchiolitis Obliterans Syndrome (BOS)
  - ICU length of stay
  - Hospital length of stay
  - Ventilation time



# PGD Assessment According to ISHLT 2005 Consensus Statement<sup>1</sup>

Grade	PaO <sub>2</sub> /FiO <sub>2</sub>	Radiographic infiltrates consistent with pulmonary edema
0	>300	Absent
1	>300	Present
2	200-300	Present
3	<200	Present

- **Clinical Implementation of ISHLT Consensus Statement as follows:**
  - Intubated patients graded based on PaO<sub>2</sub>/FiO<sub>2</sub> ratio & chest x-ray
  - Extubated patients graded 0/1 based chest x-ray
  - ECMO graded Grade 3 except prophylactic ECMO for IPAH

## PGD Grading Discrepancy Examples


Recipient **intubated**, PF ratio <200 mmHg with clear chest X-ray reading :

- INSPIRE Grade = PGD 3
- FDA Grade = PGD 0


Recipient **extubated**, PF ratio <200 mmHg on nasal supplemental O<sub>2</sub>:

- INSPIRE Grade= PGD 0 or 1 based on Chest X-ray reading
- FDA Grade = PGD 3

# Complexities of Donor Lung Offer and Randomization Process



# CONSORT Diagram of INSPIRE Trial (Randomization to mITT)



## Categories of Screen Failures (n=58)

Screen Failure Type*	Definition	Control N=15	OCS N=43
<b>Donor Screen Failure, Transplanted Off Study</b>	Donor lungs <b>did not meet INSPIRE inclusion criteria</b> and were transplanted off study	<b>6</b>	<b>17</b>
<b>Donor Screen Failure, Remained on Waiting List at End of Study</b>	Initial donor lungs were <b>not accepted for transplantation</b> , patient remained randomized and waiting for a second offer at time of trial completion	<b>4</b>	<b>14</b>
<b>Logistics</b>	Logistical issues prevented use of randomized preservation method to be used	<b>1</b>	<b>10</b>
<b>Recipient</b>	Recipient found to be no longer eligible for inclusion in the trial on day of transplant	<b>4</b>	<b>2</b>

\*Adjudicated by independent medical monitor


## Extensive Analyses Performed To Understand Screen Failures

- **Multiple factors led to screen failures:**
  - Randomization prior to physical evaluation of donor lungs
  - Some lungs were not suitable for transplantation and randomized recipients remained on waiting list awaiting another donor offer
  - Transplant logistics
- No clear reason for imbalance between arms, however, this imbalance did not result in any measurable difference favoring the OCS arm in the patients analyzed in INSPIRE

## No Evidence of Donor Lung Characteristics Favoring OCS Arm in INSPIRE Trial

Donor Parameters	Control N=184	OCS N=165
Age (year) (mean $\pm$ SD)	40 $\pm$ 14	42 $\pm$ 14
Final PaO <sub>2</sub> /FiO <sub>2</sub> (mean $\pm$ SD)	432 $\pm$ 73	441 $\pm$ 79
Smoking >20 Pky in last 6 months	17%	18%
<b>Abnormal Findings on Donor Lung Visualization at Retrieval</b>	<b>26%</b>	<b>36%</b>
Lung contusions	1%	4%
Emphysematous blebs	1%	2%
Granulomas	0.5%	2%
Pneumonia	0%	1%
Major atelectasis	21%	24%
Excessive lung adhesions, or parenchymal tears	1%	6%

# CONSORT Diagram of INSPIRE Trial (mITT to PP)





## Patients with Pre-Specified Major Protocol Violations Included in mITT, Not PP Population

Major Protocol Violation*	Control n=4	OCS n=10
Donor lungs not eligible for inclusion (active pneumonia, severe COPD with large blebs, or no final donor PF ratio to confirm eligibility)	1	4
Failure to follow Instruction for Use (IFU)/Protocol	3	4
Patient transplanted with preservation method different than randomized due to user error	0	2

\*Adjudicated by independent medical monitor

## Recipient Characteristics Similar Between Arms

Recipient Characteristic	Control N=184	OCS N=165
Age (years), Mean $\pm$ SD	50 $\pm$ 14	50 $\pm$ 13
Female, %	36%	48%
BMI (kg/m <sup>2</sup> ), Mean $\pm$ SD	23 $\pm$ 4.1	23 $\pm$ 4.6
LAS Score, Mean $\pm$ SD	48 $\pm$ 18	51 $\pm$ 20
On ECMO on Transplant Day, %	5%	5%
Use of Intraoperative Cardiopulmonary Bypass	38%	40%
Secondary Pulmonary Hypertension, %	32%	40%
Primary Cause of Lung Failure, %		
COPD	28%	28%
IPF	34%	35%
Cystic Fibrosis	23%	21%
IPAH	4%	9%
Sarcoidosis	5%	3%

## INSPIRE Trial Methodology Summary

- INSPIRE RCT was successfully implemented in 21 international academic lung Tx. Centers in the complex field of lung transplantation
- PGD assessment followed the clinical implementation of the 2005 ISHLT Consensus Statement
- Screen failure imbalance did not result in any measurable difference favoring the OCS arm
- Largest body of prospective clinical evidence supporting use of EVLP in standard lung transplantation

## INSPIRE Trial Adjudication and Trial Oversight

---

**John Wallwork, FRCS, FmedSCI**

Emeritus Professor, Cardiothoracic Surgery

Papworth Hospital

Cambridge University, UK

President (1994-95), International Society for Heart and Lung  
Transplant (ISHLT)

## Medical Monitor Adjudication Process

- Adjudicated PGD scores according to the ISHLT 2005 consensus statement guidelines. This process was implemented in a blinded and consistent manner for both study groups.
- Adjudicated all Serious Adverse Events (SAEs) according to the protocol definitions, without changes to the protocol safety endpoint definition. This process was implemented in a blinded and consistent manner for both study groups
- There was no conflict between Medical Monitor role and my role on the DSMB

## **INSPIRE Trial Results**

---

**Gregor Warnecke, MD, PhD**

Vice Chairman

Director of Heart and Lung Transplantation

Hannover Medical School

# Outline of INSPIRE Trial Results

## Critical Transplant Times and OCS Perfusion Parameters

### Composite Primary Effectiveness Endpoint

### Components of Primary Composite Endpoint

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

### Adjunct Effectiveness Analysis

### Secondary Endpoints

### Safety

### Additional Clinical Endpoints

# INSPIRE Trial Definitions of Cross-Clamp and Ischemic Times

Definition	
<b>Cross-Clamp Time</b>	Time from aortic cross-clamp in donor to pulmonary artery cross-clamp removal in recipient
<b>Ischemic Time</b>	Time donor lung was not perfused with oxygenated blood




# OCS Lung Perfusion Impact on Ischemic Times During Transplantation

## Control - Cross Clamp/Ischemic Times Are Same



**Donor**

- From aortic cross-clamp in donor to pulmonary artery cross-clamp removal in recipient
- Time donor lung was not perfused with oxygenated blood



**Recipient**


## OCS - Ischemic Times Are Limited Due to OCS Perfusion

Ischemia




Oxygenated Perfusion

Ischemia

# OCS Significantly Reduced Ischemic Time on Donor Lungs – Combined Cohort



# Stable Perfusion Parameters & Lung Oxygenation on OCS Lung System – Combined Cohort



# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**

**Secondary Endpoints**


**Safety**

**Additional Clinical Endpoints**

# Primary Effectiveness Endpoint - PP: Composite of 30-Day Survival and Freedom from PGD3 Within 72 Hours


INSPIRE Cohort (N=320)

Combined Cohort (N=349)



\* met superiority test

# INSPIRE Trial Met Pre-specified Primary Effectiveness Endpoint



# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours


**Adjunct Effectiveness Analysis**

**Secondary Endpoints**

**Safety**


**Additional Clinical Endpoints**

# 30-Day Patient Survival – Combined Cohort







## 30-Day All Causes of Mortality – Combined Cohort



## 30-Day AND In-Hospital Causes of Mortality - Combined Cohort



## 30-Day AND In-Hospital Survival - Combined Cohort



# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**

**Secondary Endpoints**


**Safety**

**Additional Clinical Endpoints**

# OCS Resulted in Significant Reduction of PGD3 Within 72 Hours

## INSPIRE Cohort (N=320)

## Combined Cohort (N=349)



\* superiority test

# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**

**Secondary Endpoints**

**Safety**


**Additional Clinical Endpoints**

# Post-Hoc Adjunct Effectiveness Analysis – PP :

Composite of 30-Day and In-Hospital Survival & Freedom from PGD3 Within 72 Hours


INSPIRE Cohort (N=320)

Combined Cohort (N=349)



\* met superiority test

# Post-Hoc Adjunct Effectiveness Analysis Demonstrates Consistent Benefit of OCS





# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**

**Secondary Endpoints**

**Safety**

**Additional Clinical Endpoints**


# Rate of PGD Grade 3 At 72 Hours Similar Between Arms

		Point Estimate		Treatment Difference [Upper 95% CI]	P-value
OCS Overall		Control	OCS		
PP	INSPIRE Cohort	4.2%	2.1%	-2.1% ●	0.0002
	Combined Cohort	5.0%	3.9%	-1.1% ●	
mITT	INSPIRE Cohort	4.7%	5.3%	0.6% ●	0.037
	Combined Cohort	5.5%	6.7%	1.2% ●	0.072
<b>OCS Solution</b>					
PP	INSPIRE Cohort	4.2%	2.3%	-1.9% ●	0.001
	Combined Cohort	5.0%	5.1%	0.0% ●	
mITT	INSPIRE Cohort	4.7%	4.5%	-0.2% ●	0.028
	Combined Cohort	5.5%	6.8%	1.3% ●	0.110

-20%      -10%      0%      5%      10%      20%

← Supports Non-Inferiority      Does Not Support NI →

# Rate of PGD Grade 2 or 3 At 72 Hours Similar Between Arms



# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**

**Secondary Endpoints**

**Safety**


**Additional Clinical Endpoints**

# OCS Lung System Met Primary Safety Endpoint

INSPIRE Combined Cohort (n=349)	Control N=184	OCS N=164
Lung-graft related SAEs, n (%)	45 (24.5)	40 (24.4)
Mean ± SD	0.29 ± 0.54	0.26 ± 0.48
Non-Inferiority p-value		0.042
Type of Lung-graft related SAEs, n (%)		
Acute Rejection	4 (2)	2 (1)
Respiratory Failure*	16 (9)	23 (14)
Bronchial Anastomotic Complication	4 (2)	0
Major Pulmonary-Related Infection	29 (16)	18 (11)

\* Need for re-intubation, tracheostomy or the inability to discontinue ventilator support within 4 days post-transplant

# Survival Profile for Respiratory Failure Patients



## Overall Safety Profile (including Mortality) Similar Between OCS and Control

Patients	Control N=184	OCS N=164
Any Type of AE	83%	83%
Definitely Related	0%	0%
Probably Related	0%	1%
Possibly Related	3%	3%
Unlikely Related	31%	36%
Not Related	71%	69%
Any SAEs	63%	56%
Any Severe AEs	29%	31%
Death up to 24 months	16%	16%

# Outline of INSPIRE Trial Results

**Critical Transplant Times and OCS Perfusion Parameters**

**Composite Primary Effectiveness Endpoint**

**Components of Primary Composite Endpoint**

Short-Term Patient Survival

Freedom from PGD3 Within 72 hours

**Adjunct Effectiveness Analysis**


**Secondary Endpoints**

**Safety**


**Additional Clinical Endpoints**



# Ventilation Times, Lengths of ICU and Initial Hospital Stay Comparisons – Combined Cohort



# OCS Associated with Lower Incidence of BOS Through 24 Months – Combined Cohort PP



# at risk	OCS				
	Control	154	143	134	130
	180	162	153	143	56

# INSPIRE Trial Demonstrated Safety and Effectiveness of the OCS Lung System

- Met primary effectiveness endpoint and safety endpoint
- Significant reduction of PGD3 within 72 hours
- Significant reduction of ischemic time on donor lungs
- No additional safety risk associated with OCS compared to Control
- Favorable 2-year BOS results to be further evaluated in post-approval study

# Training Program and Post-Market Study Plan

---

**Waleed Hassanein, MD**

President and CEO

TransMedics, Inc.

# Clinical Training Infrastructure



- Dedicated 15,000 Sq. F. facility equipped with latest surgical and diagnostics equipment to replicate a retrieval environment
- 86 global academic and clinical institutions
- >400 health care professionals trained









## OCS Clinical Training and Support Program

- Initial Hands-On Clinical Training and Certification of Every New Clinical Center Starting an OCS Lung Program
- Dedicated OCS Lung iPad Training & Support Application
- 24 X 7 Phone and Text Messaging Hotline



# Post Approval Study Plan

---

## Two-Part Post-Approval Study Plan

➤ **Long-term Follow-up of INSPIRE Patients**

➤ **OCS Thoracic Organ Perfusion (TOP) Registry**

# Long-Term Follow-up of INSPIRE Patients

## Goal

- Assess impact of OCS Lung preservation on the incidence of BOS and survival for up to 5 years

## Data Collection

- Incidence of BOS at year 3, 4, and 5
- Survival at year 3, 4, and 5

## OCS Thoracic Organ Perfusion Registry

- **Goal:** Expand clinical evidence for OCS Lung System in standard criteria lung transplantation post market
- **Primary Clinical Objective:** 5-year survival compared to SRTR/OPTN data for historical controls in same time period of enrollment
- **Other Clinical Objectives:**
  - Incidence of PGD within initial 72 hours
  - Incidence of BOS-free survival up to 5 years

# Clinical Perspectives and Benefit-Risk Assessment

---

**Dirk Van Raemdonck, MD, PhD**

Director, Transplant Center


University Hospitals Leuven

Department of Thoracic Surgery

# OCS Lung System Provides Necessary Advance to Field of Lung Transplantation



**Standard of Care  
(1980s)**



**2004**






**2017**


# OCS Provides New Optimization and Monitoring Capabilities That are Not Possible with Cold Storage

CO-5

OCS System Designed to Address Limitations of Cold Ischemic Storage


REDUCE ISCHEMIC INJURY	OPTIMIZE ORGAN CONDITION	EX-VIVO FUNCTIONAL ASSESSMENT
		
Warm Oxygenated Blood Perfusion	Ventilation Recruitment	Oxygenation, Vascular Resistance, & Airway Compliance

- Improve quality of donor lung preservation
- Improve clinical decision making






# OCS Reduces Ischemic Time on Donor Lungs



- First preservation technology to reduce ischemia on donor lungs
- Ability to travel longer
- Better transplant procedure logistics



# OCS First Device To Demonstrate Significant Reduction of Most Severe Form of PGD



- First technology to reduce PGD
- Offers potential for better short- and long-term outcomes
- Encouraging BOS results to be further evaluated post approval

# INSPIRE Trial Demonstrates Reasonable Assurance of Safety for OCS Lung System

- Met safety endpoint
- Higher 30-day mortality in OCS was due to non-lung graft related causes
- Overall hospital mortality was similar between arms
- Favorable long-term safety profile with similar survival through 2 years

## Effectiveness of OCS Lung System Clearly Demonstrated in INSPIRE Trial

- OCS performed similar to or better than control on most effectiveness measures
- OCS overcomes many limitations of cold storage:
  - Reduces ischemic injury
  - Provides optimization and monitoring capabilities
- Positive benefit-risk profile

## OCS Lung Approval Would Enable Future Advancements in Lung Transplantation

- OCS Lung System is a paradigm shift in lung transplantation
- Important first step toward further advancements:
  - Improve long-term viability of donor lungs
  - Increase availability of donor lungs currently wasted due to limitations of cold storage
    - Would reduce mortality on the waiting list
- Providing OCS Lung to patients/physicians in US now critical to advancing field of lung transplantation



Bringing new life to organ transplantation™

## **OCS™ Lung System for the Preservation of Donor Lungs for Transplantation**

---

**May 17, 2017**


TransMedics, Inc.

Gastroenterology-Urology Devices Panel



# **ONSCREEN BACK-UP SLIDES**

---

## PGD3 Within 72 Hours Stratified by Ischemic Time Quartiles: Combined Cohort (N=349) – PP Population




# BOS-Free Status Through 24 Months by Control Ischemic Quartile Per Protocol Combined Cohort (N=349), Treatment Splits






## PGD 3 Within 72 Hours and BOS Free Survival Combined Cohort (N=349) – PP



# at risk	PGD 0-1	157	149	140	133	66
	PGD 3	81	68	65	62	33

## PGD 3 at T0 and BOS Free Survival Combined Cohort (N=349) – PP



# at risk	PGD 0-1 in T0-72	157	149	140	133	66
	PGD 3 at T0	58	50	47	46	24

# Logistics Screen Failure – Transplanted Off Study

## At 1<sup>st</sup> Donor Offer

### OCS

- 2 Device was not available for retrieval due to malfunction
- 1 No trained personnel to run OCS

### CONTROL

- 1 Randomization envelope not opened prior to retrieval

### OCS

- 1 pRBCs expired during simultaneous donor offers
- 1 Out of geographical zone donor offer in the UK
- 2 OCS solution used in first donor offer and out of stock
- 3 No trained personnel available to run OCS

## At 2<sup>nd</sup> or 3<sup>rd</sup> Donor Offer

# Screen Failures

## ■ mITT Definition

- All randomized patients for whom a matching lung has been harvested **and** determined to be eligible for preservation with either Control or OCS before any attempt has been made to preserve the lung with either Control or OCS.

## ■ Screen Failures not eligible for mITT

- Donor lungs turned down for transplant altogether (i.e. “dry runs)
- Donor lungs harvested but not eligible for INSPIRE according to inclusion/exclusion criteria
- Donor lungs harvested and eligible for INSPIRE but for which OCS or Control preservation could not be attempted (“Logistics Screen Failures”)
- Recipient screen failures were subjects who did not meet inclusion/exclusion criteria for INSPIRE at the time that a matching donor lung became available




**ITT: N=407**

## **INSPIRE Trial: Pre-Specified Primary Effectiveness Endpoint**

	Treatment Difference [Upper 95% CI]
ITT OCS Arm vs. Control	<b>-3.1% [4.7%]</b>


- Imputation Methodology
  - 9 US subjects had known outcomes off-study
  - Multiple imputation without adjustment was used to address unknown outcomes

# OCS Performed Better than Control Regardless of Cross-Clamp Time: Combined Cohort (N=349) – PP



# Pulmonary Function Test

## FEV1% Predicted from Discharge Through 24 Months – Combined Cohort (PP)




## Recipient Characteristics Similar Between Arms



Recipient Characteristic	Control N=184	OCS N=165
Age (years), Mean $\pm$ SD	50 $\pm$ 14	50 $\pm$ 13
Female, %	36%	48%
BMI (kg/m <sup>2</sup> ), Mean $\pm$ SD	23 $\pm$ 4.1	23 $\pm$ 4.6
LAS Score, Mean $\pm$ SD	48 $\pm$ 18	51 $\pm$ 20
On ECMO on Transplant Day, %	5%	5%
<b>Use of Intraoperative Cardiopulmonary Bypass</b>	<b>38%</b>	<b>40%</b>
Secondary Pulmonary Hypertension, %	32%	40%
Primary Cause of Lung Failure, %		
COPD	28%	28%
IPF	34%	35%
Cystic Fibrosis	23%	21%
IPAH	4%	9%
Sarcoidosis	5%	3%



## PGD 3 Within 72 Hours Associated with Significant Increase of Time on Ventilation and ICU Stay – Combined Cohort (N=349) – PP




# Incidence of Post-Transplant PGD Grade 3




**N=90 Subjects**

## PGD3 Within First 48 Hours Correlates with Lower Long- term Survival, Higher BOS Rate




Whitson BA et al., JHLT 2007; 26:1004-1011

## PGD 3 at T0 Correlates with Long-Term BOS Rates



Daud et al. AMJ. Resp. & Crit. Care Med. 2007

## PGD 2 or 3 at T0 Significant Risk Factor for BOS



DerHovanesian et al, Am J Transplant. 2016; 16(2): 640-649.



## Subject D: Protocol Violation (Ineligible donor lungs due to presence of active pneumonia)

- Evidence from Site-Entered CRF:
  - Eligibility Donor Form: Presence of active pulmonary disease
  - Donor Assessment Form: **Pulmonary edema possibly due to aspiration**; ongoing LLL consolidation collapse with focal areas of left upper lobe “pneumonic infiltrates”
  - Donor Lung Assessment Form: Mucoid/mucopurulent secretions; apical lung scarring
- FDA’s Control Subject Counter Example – Site-Entered CRF:
  - Donor Assessment Form: +++ Polymorphs; ++ RBCs; Yeast isolated; Coagulase negative staphylococcus isolated; Enterobacter cloacae complex isolate
  - **This counter example is not representative of ineligibility, because it only represents upper respiratory flora and not active pneumonia**

## Survival by PGD Category

	Survival Probability at 24 Months
PGD 0/1 at all timepoints	87.1%
PGD 3 at T0; PGD 0/1 thereafter	80.0%
PGD 3 at T24; PGD 0/1 thereafter	72.9%
PGD 3 at T48; PGD 0/1 thereafter	83.3%
PGD 3 at 72; PGD 0/1 thereafter	52.5%