



Medacta International S.A.
Chris Lussier
Senior Director, Quality and Regulatory
Medacta USA
3973 Delp Street
Memphis, Tennessee 38118

November 23, 2022

Re: K213459

Trade/Device Name: Glenoid Reconstruction System
Regulation Number: 21 CFR 888.3660
Regulation Name: Shoulder Joint Metal/Polymer Semi-Constrained Cemented Prosthesis
Regulatory Class: Class II
Product Code: PHX, KWS, HSD, MBF
Dated: October 26, 2022
Received: October 28, 2022

Dear Chris Lussier:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531 - 542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Victoria A. Lilling -S

Digitally signed by Victoria A.
Lilling -S
Date: 2022.11.23 19:03:24 -05'00'

Victoria Lilling, M.D.
Assistant Division Director
DHT6A: Division of Joint Arthroplasty Devices
OHT6: Office of Orthopedic Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K213459

Device Name

Glenoid Reconstruction System

Indications for Use (Describe)

The Reverse Shoulder Prosthesis is indicated for treatment of humeral fractures and for primary or revision total shoulder replacement in patients with a grossly deficient rotator cuff shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly deficient rotator cuff shoulder joint. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device. The glenoid baseplate is intended for cementless application with the addition of screws for fixation. The humeral stems are intended for cemented or cementless use.

The Glenoid Reconstruction System baseplate is intended for cementless application with the addition of polyaxial screws for primary stability. A Glenoid Reconstruction System central screw can be used to provide additional fixation.

The Reverse Shoulder Prosthesis – Short Humeral Diaphysis is indicated for primary total shoulder replacement in patients with grossly deficient rotator cuff shoulder joint with severe arthropathy.

The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device.

The glenoid baseplate is intended for cementless application with the addition of screws for fixation.

The humeral short stem is intended for cementless use.

The Glenoid Reconstruction System baseplate is intended for cementless application with the addition of polyaxial screws for primary stability. A Glenoid Reconstruction System central screw can be used to provide additional fixation.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

I. Submitter

Medacta International SA
Strada Regina
6874 Castel San Pietro (CH)
Switzerland
Phone (+41) 91 696 60 60
Fax (+41) 91 696 60 66

Contact Person: Stefano Baj, Regulatory and Compliance Director, Medacta International SA
Applicant Correspondent: Chris Lussier, Senior Director of Quality and Regulatory, Medacta USA
Date Prepared: October 25, 2021
Date Revised: November 23, 2022

II. Device

Device Proprietary Name:	Glenoid Reconstruction System
Common or Usual Name:	Shoulder Prosthesis, Reverse Configuration
Classification Name:	Shoulder joint metal/polymer semi-constrained cemented prosthesis
Primary Product Code	PHX
Secondary Product Code:	KWS, HSD, MBF
Regulation Number:	21 CFR 888.3660, 21 CFR 888.3690, 21 CFR 888.3670
Device Classification	II

III. Predicate Device

Substantial equivalence is claimed to the following devices:

Primary predicate device

- Aequalis Perform Reversed, K161742, Wright Medical

Additional predicate devices

- Comprehensive Reverse Shoulder, K120121, Biomet Manufacturing Corp.
- Reverse Shoulder Prosthesis, K100741, Encore Medical L.P.
- Encore Reverse Shoulder Prosthesis (RSP), K051075, Encore Medical L.P.

Reference device

- Medacta Shoulder System, K170452, Medacta International SA

IV. Device Description

The Glenoid Reconstruction System is a Medacta Shoulder System line extension to provide a larger product offering. It includes GRS baseplates and central screws, sterile implantable devices used to replace only the glenoid side of the gleno-humeral joint in a shoulder reverse configuration.

The GRS baseplate, intended for cementless application, is designed to be fixed into the glenoid bone by means of both a press-fit central post and Glenoid Polyaxial Screws. If desired, a GRS Central Screw can be used to provide additional stability. The GRS baseplate is designed to provide an interface for glenosphere coupling.

The subject baseplate is available in two taper diameters ($\text{Ø}24.5$ and $\text{Ø}27$) with two different lengths (20 and 30 mm) of the central post and 3 lateralization options (0 mm, +3 mm and +6 mm). It is made of Ti6Al7Nb according to ISO 5832-11 and double coated with Ti coating according to ASTM F1580-18 and HA coating according to ASTM F1185-03.

The GRS central screw is available in 4 different lengths, from 15 to 30 mm and it is made of Ti6Al4V according to ISO 5832-3.

V. Indications for Use

The Reverse Shoulder Prosthesis is indicated for treatment of humeral fractures and for primary or revision total shoulder replacement in patients with a grossly deficient rotator cuff shoulder joint with severe arthropathy or a previously failed joint replacement with a grossly deficient rotator cuff shoulder joint. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device. The glenoid baseplate is intended for cementless application with the addition of screws for fixation. The humeral stems are intended for cemented or cementless use.

The Glenoid Reconstruction System baseplate is intended for cementless application with the addition of polyaxial screws for primary stability. A Glenoid Reconstruction System central screw can be used to provide additional fixation.

The Reverse Shoulder Prosthesis- Short Humeral Diaphysis is indicated for primary total shoulder replacement in patients with grossly deficient rotator cuff shoulder joint with severe arthropathy. The patient's joint must be anatomically and structurally suited to receive the selected implant(s), and a functional deltoid muscle is necessary to use the device. The glenoid baseplate is intended for cementless application with the addition of screws for fixation. The humeral short stem is intended for cementless use.

The Glenoid Reconstruction System baseplate is intended for cementless application with the addition of polyaxial screws for primary stability. A Glenoid Reconstruction System central screw can be used to provide additional fixation.

VI. Comparison of Technological Characteristics

The subject and predicate devices (K161742 and K120121) are substantially equivalent with regards to the following characteristics:

- general design;
- lateralization (except the predicate K120121);
- fixation;

- biocompatibility;
- device usage;
- sterility;
- shelf-life; and
- packaging.

The subject implants differ from the predicate devices (K161742 and K120121) as follows:

- baseplate taper and central post diameters;
- baseplate material and finishing; and
- baseplate central post length and central screw lengths.

Discussion

Medacta International SA has not made any change to the indications for use, general design and shape, fixation, device usage, biocompatibility, sterility, shelf life, and packaging of the subject devices respect to the predicate devices.

Based on the comparison of technological characteristics and performance data provided within this submission, the data supports the substantial equivalence of the Glenoid Reconstruction System implants to the identified predicate devices.

VII. Performance Data

Based on the risk analysis, performance testing was conducted to written protocols. The following tests and rationales are provided in support of the substantial equivalence determination:

Non-Clinical Studies

- *DESIGN VALIDATION*
 - Glenoid Reconstruction System Design Validation
- *PERFORMANCE TESTING*
 - Fatigue testing on glenoid reconstruction system
 - Micromotion assessment on Glenoid reconstruction system according to ASTM F2028-17
 - GRS Central Screw Static test according to ASTM F543-17
 - Characterization Report “Y367” Titanium + “Osprovit” Hydroxyapatite double coating on GRS Glenoid baseplate component
 - Scanning electron microscopy pictures of the "Y367" Titanium + "Osprovit" HA implant surfaces of the GRS Glenoid baseplate
 - Cross sectioned area of “Y367” Titanium + “Osprovit” HA implant surfaces of the GRS Glenoid baseplate
 - Rationale comparison between features of the Hydroxyapatite Osprovit coating deposited on the Medacta GRS Glenoid baseplate and on planar samples made of Ti6Al7Nb, based on XRD analyses

- *PYROGENICITY* evaluation
- *BIOCOMBATIBILITY* evaluation
- *SHELF-LIFE* evaluation

Clinical Studies:

- No clinical studies were conducted.

VIII. Conclusion

The information provided above supports that the Glenoid Reconstruction Systems implants are substantially equivalent to the predicate devices.