



Qingdao Hisense Medical Equipment Co., Ltd.  
% Lu Zhonghao  
Quality Engineer  
No. 399 Songling Road, Laoshan District  
Qingdao, Shandong 266100  
CHINA

Re: K222132

November 8, 2022

Trade/Device Name: Hisense LCD monitor (HMD2G21S, HMD3G21S)  
Regulation Number: 21 CFR 892.2050  
Regulation Name: Medical image management and processing system  
Regulatory Class: Class II  
Product Code: PGY  
Dated: October 7, 2022  
Received: October 14, 2022

Dear Lu Zhonghao:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Jessica Lamb, Ph.D.  
Assistant Director  
Imaging Software Team  
DHT8B: Division of Radiological Imaging Devices  
and Electronic Products  
OHT8: Office of Radiological Health  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K222132

Device Name  
LCD Monitor( HMD2G21S,HMD3G21S)

### Indications for Use (Describe)

1) The 2MP Monochrome LCD Monitor HMD2G21S is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners.

The device does not support the display of mammography images for diagnosis.

2)The 3MP Monochrome LCD Monitor HMD3G21S is intended to be used in displaying and viewing digital images for diagnosis X-ray or MRI, etc. by trained medical practitioners.

The device does not support the display of mammography images for diagnosis.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

---

**Section5: 510(k) Summary**

**LCD monitor**

---

**Section5: 510(k) Summary****(K222132)****1. Applicant information**

|                 |   |
|-----------------|---|
| Date            | Oct. 7, 2022  |
| Submitter       | Qingdao Hisense Medical Equipment Co.,Ltd.<br>Address: No. 399 Songling Road, Laoshan District<br>266100, Qingdao, Shandong, P. R. China  |
| Contact Person: | <p><b>Primary Contact Person</b><br/>Lu zhonghao<br/>Quality Engineer<br/>Qingdao Hisense Medical Equipment Co.,Ltd.<br/>Mail:luzhonghao@hisense.com<br/>Tel: +86-532-83091111<br/>Fax:+86-532-83091111</p> <p><b>Secondary Contact Person</b><br/>Wu yalan<br/>RA Manager<br/>Qingdao Hisense Medical Equipment Co.,Ltd.<br/>Mail: Wu yalan@hisense.com<br/>Tel: +86-532-83091111<br/>Fax:+86-532-83091111</p> |

**2.Device information**

|                      |  |
|----------------------|--|
| Device Trade Name:   | Hisense LCD monitor HMD2G21S,HMD3G21S        |
| Common/Usual Name:   | 2M/3M Monochrome LCD Monitor                 |
| Classification       | II   |
| Classification Name: | Display, Diagnostic Radiology 21CFR 892.2050 |
| Product Code:        | PGY  |

**3. Predicate Device(s):**

|               |   |
|---------------|---|
| Trade Name    | JUSHA-M260G LCD Monitor                   |
| 510(k) Number | K183497                                   |
| Product Code  | PGY                                       |
| Manufacturer  | Nanjing Jusha Display Technology Co., Ltd |

|               |   |
|---------------|---|
| Trade Name    | JUSHA-M33C LCD Monitor                    |
| 510(k) Number | K141690                                   |
| Product Code  | PGY                                       |
| Manufacturer  | Nanjing Jusha Display Technology Co., Ltd |

**4. Device Description**

Hisense HMD2G21S, HMD3G21S LCD monitor complies with the DICOM Part 14 standard and is applicable to DSA, MRI, DR, CR, CT, and PET medical imaging. It is intended for trained medical practitioners and provides the image viewing and medical diagnostic functions.

**5. Intended Use/Indication for Use**

1) The 2MP Monochrome LCD Monitor HMD2G21S is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners.

The device does not support the display of mammography images for diagnosis.

2) The 3MP Monochrome LCD Monitor HMD3G21S is intended to be used in displaying and viewing digital images for diagnosis X-ray or MRI, etc. by trained medical practitioners.

The device does not support the display of mammography images for diagnosis.

**6. Substantial Equivalence Comparison**

**6.1 HMD2G21S and its predicate**

Table 01: General Comparison Table

| Elements of Comparison              | Proposed device   | Predicate Device  | Remarks |
|-------------------------------------|---|---|---------|
| 510(k) Number                       | K222132   | K183497   | -       |
| Manufacturer                        | Qingdao Hisense Medical Equipment Co., Ltd  | Nanjing Jusha Display Technology Co., Ltd   | -       |
| Device type/model                   | HMD2G21S  | JUSHA-M260G LCD Monitor   | -       |
| Intended use/<br>Indication for use | The 2MP Monochrome LCD Monitor HMD2G21S is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners.<br><br>The device does not support the display of mammography images for diagnosis. | JUSHA-M260G LCD Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners.<br><br>The device does not support the display of mammography images for diagnosis. | same    |
| Prescription or OTC                 | RX  | RX  | same    |

| Elements of Comparison | Proposed device   | Predicate Device  | Remarks |
|------------------------|---|---|---------|
| Electrical Safety      | Compliance with IEC 60601-1                             | Compliance with IEC 60601-1                             | same    |
| EMC                    | Compliance with IEC 60601-1-2                           | Compliance with IEC 60601-1-2                           | same    |
| Performance testing    | FDA guidance “Display devices for Diagnostic Radiology” | FDA guidance “Display devices for Diagnostic Radiology” | same    |

Table 02: Detailed Comparison Table

|   | Elements of Comparison    | Proposed device                            | Predicate Device                          | Remarks  |
|---|---------------------------|--|---|--|
| ID  | 510(k) Number             | K222132                                    | K183497                                   | -  |
|   | Manufacturer              | Qingdao Hisense Medical Equipment Co., Ltd | Nanjing Jusha Display Technology Co., Ltd | -  |
|   | Device type/model         | HMD2G21S                                   | JUSHA-M260G LCD Monitor                   | -  |
| <b>1.Display Performance/Specifications</b> |                           |  |   |  |
| 1.1   | Screen size               | 21,3 inches (540 mm)                       | 21.3"                                     | same   |
| 1.2   | Screen Technology         | Gray scale TFT LCD panel                   | Mono-TFT LCD Panel                        | same   |
| 1.3   | Viewing Angle             | Horizontal: 178°; vertical: 178° (CR ≥ 10) | Horizontal 178°, Vertical 178°            | same   |
| 1.4   | Pixel array               | 1600 pixels (H) x 1200 pixels (V)          | 1600 x 1200/1200 x 1600                   | same   |
| 1.5   | Display Area              | 432 (H) x 324 (V) (mm)                     | 432.0 (H) x 324.0(V) mm                   | same   |
| 1.6   | Pixel Pitch               | 0,270 (H) x 0,270 (V) (mm)                 | 0.27x0.27 mm                              | same   |
| 1.7   | Subpixel driving          | Not Applicable                             | Not Applicable                            | -  |
| 1.8   | Contrast Ratio            | 1800:1 (typ.)                              | 1400:1                                    | the specification of the proposed device is superior to that of the predicate device |
| 1.9   | Frame rate / refresh rate | 37.9~75kHz;60Hz                            | 37.9~75kHz;60Hz                           | same   |
| 1.10  | Maximum Brightness(typ)   | 1900 cd/m <sup>2</sup>                     | 1000cd/m <sup>2</sup>                     | the specification of the proposed  |

|                                      | Elements of Comparison       | Proposed device                             | Predicate Device                          | Remarks  |
|--------------------------------------|------------------------------|---|---|--|
| ID                                   | 510(k) Number                | K222132                                     | K183497                                   | -  |
|                                      | Manufacturer                 | Qingdao Hisense Medical Equipment Co., Ltd  | Nanjing Jusha Display Technology Co., Ltd | -  |
|                                      | Device type/model            | HMD2G21S                                    | JUSHA-M260G LCD Monitor                   | -  |
|                                      | )                            |   |   | device is superior to that of the predicate device |
| 1.11                                 | Recommended brightness for   | 400 cd/m <sup>2</sup>                       | 400cd/m <sup>2</sup>                      | same   |
| 1.12                                 | Backlight type               | LED   | LED                                       | same   |
| 1.13                                 | Ambient light sensing        | Built-in ambient light sensor               | Built in calibration sensor provided      | same   |
| 1.14                                 | Response Time                | 19 ms                                       | 16 ms                                     | The difference does not affect diagnosis.          |
| 1.15                                 | Aspect Ratio                 | 4:3   | 4:3                                       | same   |
| 1.16                                 | Luminance calibration        | Front-facing sensor Upper Computer Software | Built in calibration sensor provided      | Same, only difference in words.                    |
| 1.17                                 | Touch-screen                 | Not Applicable                              | Not Applicable                            | -  |
| <b>2.Video Signals</b>               |                              |   |   |  |
| 2.1                                  | Input Video Signal           | DVI-D x 1<br>DisplayPort x 1                | DVI-D x 1,<br>DisplayPort x 1             | same   |
| 2.2                                  | Output Signal                | DisplayPort x 1                             | DisplayPort x1                            | same   |
| 2.3                                  | Video bandwidth              | DVI: 215MHz<br>DisplayPort : 215MHz         | DVI: 215MHz<br>DisplayPort : 215MHz       | same   |
| <b>3.Power Related Specification</b> |                              |   |   |  |
| 3.1                                  | Power Requirements           | 24 V DC, 2,1 A                              | DC 12V                                    | this difference doesn't affect product's safety.   |
| 3.2                                  | Power Consumption/ save mode | 50.4W/Below 0.5 W                           | 50W/less than 0.5W                        | this difference doesn't affect product's safety.   |



|  | Elements of Comparison           | Proposed device  | Predicate Device   | Remarks  |
|--|----------------------------------|--|--|--|
| <b>ID</b>                                      | 510(k) Number                    | K222132  | K183497  | -  |
|  | Manufacturer                     | Qingdao Hisense Medical Equipment Co., Ltd                                     | Nanjing Jusha Display Technology Co., Ltd                            | -  |
|  | Device type/model                | HMD2G21S   | JUSHA-M260G LCD Monitor  | -  |
| 3.3  | Power Management                 | DVI DMPM<br>DisplayPort 1.2a   | DVI DMPM<br>DisplayPort 1.1a   | this difference doesn't affect product's safety.                   |
| <b>4.Miscellaneous Features/Specifications</b> |                                  |  |  |  |
| 4.1  | USB ports                        | 1 uplink port and 2 downlink ports/ Rev. 2.0                                   | 1 upstream (endpoint), 2 downstream/ Rev. 2.0                        | same   |
| 4.2  | Grayscale Tones(LUT)             | 14-bit:16384   | 14-bit:16384   | same   |
| 4.3  | User controls                    | Off the shelf  | Off the shelf  | same   |
| 4.4  | Software/Firm ware:              | Built-in embedded software   | Built-in embedded software   | same   |
| 4.5  | Dimensions w/o stand (W x H x D) | Without base:366 mm x 482 mm x 63 mm<br>With base:366 mm x 502-616 mm x 244 mm | Without stand: 382mm x490mm x77mm<br>With stand: 382mm x635mm x238mm | Different design scheme, the difference does not affect diagnosis. |
| 4.6  | Net weight                       | 5 kg (excluding the base)  | 7.5 kg(without stand)  | Different weight due to different components and parts             |
| 4.7  | VESA standard                    | 100 x 100 (mm)   | 100 x 100 (mm)   | same   |

6.2 HMD3G21S and its predicate

Table 03: General Comparison Table

| Elements of Comparison              | Proposed device  | Predicate Device   | Remarks |
|-------------------------------------|--|--|---------|
| 510(k) Number                       | K222132  | K141690  | /       |
| Manufacturer                        | Qingdao Hisense Medical Equipment Co., Ltd   | Nanjing Jusha Display Technology Co., Ltd  | /       |
| Device type/model                   | HMD3G21S   | JUSHA-M33C   | /       |
| Intended use/<br>Indication for use | The 3MP Monochrome LCD Monitor HMD3G21S is intended to be used in displaying and viewing digital images for diagnosis X-ray or MRI, etc. by trained medical practitioners.<br><br>The device does not support the display of mammography images for diagnosis. | JUSHA-M33C Monitor is intended to be used in displaying and viewing digital images for diagnosis of X-ray or MRI, etc. by trained medical practitioners.<br><br>The device does not support the display of mammography images for diagnosis. | same    |
| Prescription or OTC                 | RX   | RX   | same    |
| Electrical Safety                   | Compliance with IEC 60601-1  | Compliance with IEC 60601-1  | same    |
| EMC                                 | Compliance with IEC 60601-1-2  | Compliance with IEC 60601-1-2  | same    |
| Performance testing                 | FDA guidance “Display devices for Diagnostic Radiology”  | FDA guidance “Display devices for Diagnostic Radiology”  | same    |

Table 04: Detailed Comparison Table

| ID  | Elements of Comparison | Proposed device                            | Predicate Device                          | Remarks |
|---|------------------------|--|---|---------|
|   | 510(k) Number          | K222132                                    | K141690                                   | -       |
|   | Manufacturer           | Qingdao Hisense Medical Equipment Co., Ltd | Nanjing Jusha Display Technology Co., Ltd | -       |
|   | Device type/model      | HMD3G21S                                   | JUSHA-M33C                                | -       |
| <b>1.Display Performance/Specifications</b> |                        |  |   |         |
| 1.1   | Screen size            | 21.3 inches (540 mm)                       | 21.3inches                                | same    |
| 1.2   | Screen Technology      | Gray scale TFT LCD panel                   | Mono-TFT LCD Panel                        | same    |

|      | Elements of Comparison      | Proposed device                              | Predicate Device                          | Remarks  |
|------|-----------------------------|--|---|--|
| ID   | 510(k) Number               | K222132                                      | K141690                                   | -  |
|      | Manufacturer                | Qingdao Hisense Medical Equipment Co., Ltd   | Nanjing Jusha Display Technology Co., Ltd | -  |
|      | Device type/model           | HMD3G21S                                     | JUSHA-M33C                                | -  |
| 1.3  | Viewing Angle               | Horizontal: 178°; vertical: 178° (CR ≥ 10:1) | Horizontal 176°,Vertical 176°             | the specification of the proposed device is superior to that of the predicate device |
| 1.4  | Pixel array                 | 2048 pixels (H) x 1536 pixels (V)            | 2048 x1536/1536x 2048                     | same   |
| 1.5  | Display Area                | 433.15 (H) x 324.86 (V) (mm)                 | 433.152 (H) x 324.864 (V) (mm)            | same   |
| 1.6  | Pixel Pitch                 | 0.2115 (H) x 0.2115 (V) (mm)                 | 0.2115x0.2115 mm                          | same   |
| 1.7  | Subpixel driving            | Not Applicable                               | Not Applicable                            | -  |
| 1.8  | Contrast Ratio              | 1500:1 (typ.)                                | 1400:1                                    | the specification of the proposed device is superior to that of the predicate device |
| 1.9  | Frame Rate and Refresh Rate | 96.7kHz;60Hz                                 | 96.7kHz;60Hz                              | same   |
| 1.10 | Maximum Brightness(typ )    | 2000cd/m <sup>2</sup>                        | 1700cd/m <sup>2</sup>                     | the specification of the proposed device is superior to that of the predicate device |
| 1.11 | Recommended brightness      | 500d/cm <sup>2</sup>                         | 500cd/m <sup>2</sup>                      | same   |
| 1.12 | Backlight type              | LED  | LED                                       | same   |
| 1.13 | Ambient light sensing       | Built-in ambient light sensor                | Built in calibration sensor provided      | same   |

|  | Elements of Comparison       | Proposed device   | Predicate Device                              | Remarks  |
|--|------------------------------|---|---|--|
| ID   | 510(k) Number                | K222132   | K141690                                       | -  |
|  | Manufacturer                 | Qingdao Hisense Medical Equipment Co., Ltd                    | Nanjing Jusha Display Technology Co., Ltd     | -  |
|  | Device type/model            | HMD3G21S  | JUSHA-M33C                                    | -  |
| 1.14   | Response Time                | 28 ms (typ.)  | 40 ms   | the specification of the proposed device is superior to that of the predicate device |
| 1.15   | Aspect Ratio                 | 4:3   | 4:3   | same   |
| 1.16   | Luminance calibration        | Front-facing sensor<br>Upper Computer Software<br>Body sensor | Built in calibration sensor provided          | Same, only difference in words.  |
| 1.17   | Touch-screen                 | Not Applicable  | Not Applicable                                | -  |
| <b>2.Video Signals</b>                         |                              |   |   |  |
| 2.1  | Input Video Signal           | DVI-D x 1<br>DisplayPort x 1                                  | DVI-D x 1<br>DisplayPort x 1                  | same   |
| 2.2  | Output Signal                | DisplayPort x 1   | DisplayPort x1                                | same   |
| 2.3  | Video bandwidth              | DVI: 215MHz<br>DisplayPort : 215MHz                           | DVI: 215MHz<br>DisplayPort : 215MHz           | same   |
| <b>3.Power Related Specification</b>           |                              |   |   |  |
| 3.1  | Power Requirements           | 24 V DC, 1.7 A  | AC 100~240V 50~60Hz                           | this difference doesn't affect product's safety.                                     |
| 3.2  | Power Consumption/ save mode | 40.8 W/Below 0.5 W  | 45W/less than 3W                              | this difference doesn't affect product's safety.                                     |
| 3.3  | Power Management             | DVI DMPM<br>DisplayPort 1.2a                                  | DVI DMPM<br>DisplayPort 1.1a                  | this difference doesn't affect product's safety.                                     |
| <b>4.Miscellaneous Features/Specifications</b> |                              |   |   |  |
| 4.1  | USB ports                    | 1 uplink port and 2 downlink ports/ Rev. 2.0                  | 1 upstream (endpoint), 2 downstream/ Rev. 2.0 | same   |

|     | Elements of Comparison           | Proposed device   | Predicate Device   | Remarks  |
|-----|----------------------------------|---|--|--|
| ID  | 510(k) Number                    | K222132   | K141690  | -  |
|     | Manufacturer                     | Qingdao Hisense Medical Equipment Co., Ltd                                      | Nanjing Jusha Display Technology Co., Ltd                          | -  |
|     | Device type/model                | HMD3G21S  | JUSHA-M33C   | -  |
| 4.2 | Grayscale Tones(LUT)             | 16384   | 65536  | Different design scheme. The difference does not affect diagnosis. |
| 4.3 | User controls                    | Off the shelf   | Off the shelf  | same   |
| 4.4 | Software/Firmware:               | Built-in embedded software  | Built-in embedded software   | same   |
| 4.5 | Dimensions w/o stand (W x H x D) | Without base:366 mm x 482 mm x 63 mm;<br>With base:366 mm x 502-616 mm x 244 mm | Without stand:382mm x490mm x75mm<br>With stand:382mm x533mm x238mm | Different design scheme, the difference does not affect diagnosis. |
| 4.6 | Net weight                       | 6 kg(excluding the base)  | 7 kg(Without stand)  | Different weight due to different components and parts             |
| 4.7 | VESA standard                    | 100 x 100 (mm)  | 100 x 100 (mm)   | same   |

### Discussion of Differences

#### About HMD2G21S and JUSHA-M260G

- 1) The proposed device is superior to that of the predicate device in Contrast Ratio, Maximum Brightness. The proposed device has better displaying image quality. And the predicate device is better than proposed device in Response Time, this difference doesn't affect diagnosis.
- 2) Because of different design scheme, the Power, Dimension and Net weight is not same. But this difference doesn't affect the diagnosis.

**About HMD3G21S and JUSHA-M33C**

- 1) The proposed device is superior to that of the predicate device in Viewing Angle, Contrast Ratio, Maximum Brightness. The proposed device has better displaying image quality. And the predicate device is better than proposed device in Response Time and Grayscale Tones, but this difference doesn't affect diagnosis.
- 2) Because of different design scheme, the Power, Dimension and Net weight is not same. But This difference doesn't affect the diagnosis.

**7.Bench testing:**

The bench tests were performed on the proposed devices HMD2G21S, HMD3G21S in accordance with Assessment of Display Performance for Medical Imaging Systems by AAPM Task Group 18 (TG18 guideline). The detail test item as below.

| Measurements Guidance   | HMD2G21S  | HMD3G21S  |
|---|---|---|
| a. Spatial resolution   | By reporting modulation transfer function.  | By reporting modulation transfer function   |
| b. Pixel defects (maximum counts, allowed defect types, and locations)                                  | Maximum number allowed for each type.   | Maximum number allowed for each type.   |
| c. Artifacts  | Measure Artifacts with TG18   | Measure Artifacts with TG18   |
| d. Temporal response  | Measure the rise and fall time constants for 5– 95% and 40–60% luminance transitions. | Measure the rise and fall time constants for 5– 95% and 40–60% luminance transitions. |
| e. Luminance (maximum, minimum, achievable, and recommended)  | Measure the maximum, minimum, achievable, and recommended luminance.                  | Measure the maximum, minimum, achievable, and recommended luminance.                  |
| f. Conformance to a grayscale-to-luminance function (e.g., DICOM GSDF)                                  | Luminance Response by AAPM-TG18.  | Luminance Response by AAPM-TG18.  |
| g. Luminance at 30° and 45° in diagonal, horizontal, and vertical directions at center and four corners | NA  | NA  |
| h. Luminance uniformity or Mura test  | NA  | NA  |
| i. Stability of luminance and chromaticity response with temperature and time of operation or on-time   | NA  | NA  |
| j. Spatial noise  | NA  | NA  |
| k. Reflection coefficient   | NA  | NA  |
| l. Veiling glare or small-spot contrast   | NA  | NA  |

| Measurements Guidance                              | HMD2G21S  | HMD3G21S  |
|--|---|---|
| m. Color tracking (primary colors and color gamut) | NA, HMD2G21S is a monochrome LCD monitor, not color monitor | NA, HMD3G21S is a monochrome LCD monitor, not color monitor |
| n. Gray tracking (gray shades and white point)     | NA  | NA  |

**8.Summary of Non-Clinical Tests**

The Hisense LCD Monitor were evaluated for electrical, electromagnetic and performance, and have been found to comply with applicable standards as following:

- (1) EN 60601-1:2006+A1:2013+A12:2014 & IEC60601-1:2005+CORR.1:2006+CORR.2:2007+AM1:2012 Medical electrical equipment Part 1: General requirements for basic safety and essential performance
- (2) IEC 60601-1-2:2014 & EN 60601-1-2:2015

Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests.

The following quality assurance measures are applied to the development of the system:

- (1) Risk Management
  - (2) Requirement review and Design reviews
  - (3) Integration testing
  - (4) Performance testing
  - (5) Safety testing
- etc.

**9.Clinical Testing**

Not applicable.

The subject of this premarket submission did not require clinical studies to support substantial equivalence. The proposed device is Substantially Equivalent (SE) to the predicate device which is US legally market device. Therefore, the subject device is determined as safe and effectiveness.

**10.Conclusion**

Hisense considers the LCD Monitor HMD2G21S, HMD3G21S to be as safe, as effective, and performance is substantially equivalent to the predicate device.