



August 24, 2022

Neuronetics, Inc.
Amanda Pentecost, PhD
Regulatory Affairs Specialist
3222 Phoenixville Pike
Malvern, PA 19355

Re: K222230

Trade/Device Name: NeuroStar, NeuroStar TMS Therapy System, NeuroStar Advanced Therapy System, NeuroStar Advanced Therapy System for Mental Health

Regulation Number: 21 CFR 882.5805

Regulation Name: Repetitive Transcranial Magnetic Stimulation System

Regulatory Class: Class II

Product Code: OBP

Dated: July 21, 2022

Received: July 25, 2022

Dear Dr. Pentecost:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Pamela Scott
Assistant Director
DHT5B: Division of Neuromodulation
and Physical Medicine Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K222230

Device Name

NeuroStar D-Tect MT Accessory (81-00241-000)

Indications for Use (Describe)

The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Date Prepared: July 19, 2022

Applicant: Neuronetics, Inc.
3222 Phoenixville Pike
Malvern, PA 19355

Contact Person: Amanda Pentecost, PhD
Regulatory Affairs Specialist
Phone: 610-640-4202, ext. 1132
Email: amanda.pentecost@neurostar.com

Secondary Contact: Robin Fatzinger, RAC
Sr. Director, Regulatory Affairs
Phone: 610-640-4202, ext. 1064
Email: robin.fatzinger@neurostar.com

Device Trade Name: NeuroStar
NeuroStar TMS Therapy System
NeuroStar Advanced Therapy System
NeuroStar Advanced Therapy System for Mental Health

Device Common Name: Transcranial Magnetic Stimulator

Classification: 21 CFR 882.5805

Product Code: OBP

Predicate Device: NeuroStar Advanced Therapy System (K083538, K130233, K133408, K160703, K161519, K201158, K213543, and K220127)

Device Description / Technological Characteristics:

The NeuroStar Advanced Therapy System is a transcranial magnetic stimulation device. Specifically, it is a computerized, electromechanical medical device that produces and delivers non-invasive magnetic fields to induce electrical currents targeting specific regions of the cerebral cortex. Transcranial magnetic stimulation (TMS) is a non-invasive technique used to apply brief magnetic pulses to the brain. The pulses are administered by passing high currents through an electromagnetic coil placed adjacent to a patient's scalp. The pulses induce an electric field in the underlying brain tissue. When the induced field is above a certain threshold and is directed in an appropriate orientation relative the brain's neuronal pathway, localized axonal depolarizations are produced, thus activating neurons in the targeted brain region.

The NeuroStar Advanced Therapy System is an integrated system consisting of a combination of the following components:

- Mobile Console
- System Software
- Treatment Chair
- Ferromagnetic Treatment Coil
- Head Support System
- SenStar® Connect Treatment Link & SenStar® Treatment Link
- MT Cap
- TrakStar™ Patient Data Management System
- D-Tect™ MT Accessory

The proposed change to the NeuroStar Advanced Therapy System that is the subject of this 510(k) is the addition of the NeuroStar D-Tect™ MT Accessory, which is a non-sterile, multi-use device that provides an optional method to aid in the Motor Threshold (MT) hunt process. This device is used to provide indication and amplitude of hand movement during this process. The D-Tect™ MT Accessory includes a user display interface and a human hand interface with built-in sensors for measuring thumb and finger movements. The device is a standalone accessory that does not communicate or interact with the NeuroStar Advanced Therapy System, except for receiving pulse signals in order to sync data collection. Use of the device is optional, and provides an alternative, quantitative method to determine the MT compared to the standard of care, which includes a qualitative, visual assessment. The D-Tect™ MT Accessory is only used during the Motor Threshold (MT) determination process, which is only performed the first time a patient is seen by the doctor and is performed prior to the first treatment session.

Indications for Use:

The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.

Performance Standards:

The use of the D-Tect™ MT Accessory with the NeuroStar Advanced Therapy System has been tested and conforms to the following recognized consensus standards:

- ISO 10993-1
- ANSI AAMI ES60601-1
- IEC 60601-1-2
- IEC 60601-1-6

Non-clinical Testing:

The contents of this 510(k) complies with the FDA Guidance Document: "Class II Special Controls Guidance Document: Repetitive Transcranial Magnetic Stimulation (rTMS) Systems - Guidance for Industry and Food and Drug Administration Staff." Non-clinical performance testing was performed according to the standards listed above.

Additionally, usability testing was completed in accordance with IEC 60601-1-6 Edition 3.1 2013-10 and also following FDA Guidance Document: "Applying Human Factors and Usability Engineering to Medical Devices."

Clinical Testing:

There is no clinical testing required to support this submission.

Technological Comparison:

	NeuroStar Advanced Therapy System (Subject Device)	NeuroStar Advanced Therapy System (Predicate Device)	Explanation of Differences
Indications for Use	The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.	The NeuroStar Advanced Therapy System is indicated for the treatment of depressive episodes and for decreasing anxiety symptoms for those who may exhibit comorbid anxiety symptoms in adult patients suffering from Major Depressive Disorder (MDD) and who failed to achieve satisfactory improvement from previous antidepressant medication treatment in the current episode.	No Difference
Intended Use	Major Depressive Disorder (MDD) and comorbid anxiety symptoms	Major Depressive Disorder (MDD) and comorbid anxiety symptoms	No Difference
Anatomical Sites	Left dorsolateral prefrontal cortex	Left dorsolateral prefrontal cortex	No Difference
Target Population	Adult patients	Adult patients	No Difference
Clinical Setting	Inpatient and outpatient settings, including physician's offices and clinics, hospitals, and general medical/surgical hospitals	Inpatient and outpatient settings, including physician's offices and clinics, hospitals, and general medical/surgical hospitals	No Difference
Materials	Standard materials commonly used in the manufacture of electrical medical devices	Standard materials commonly used in the manufacture of electrical medical devices	No Difference
Biocompatibility	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	Patient-contacting device components use standard materials compliant with ISO 10993-1 that are commonly used in consumer products and medical device applications	No Difference

Design	Computerized, electromechanical medical devices that use integrated systems to deliver TMS therapy	Computerized, electromechanical medical devices that use integrated systems to deliver TMS therapy	No Difference
Energy Source	Power console with magnetic coil for delivery for magnetic energy	Power console with magnetic coil for delivery for magnetic energy	No Difference
Electrical Safety & EMC	IEC 60601-1 compliant IEC 60601-1-2 compliant	IEC 60601-1 compliant IEC 60601-1-2 compliant	No Difference
Sterility	No parts of the device, accessories or components are required to be sterilized	No parts of the device, accessories or components are required to be sterilized	No Difference
Coil Type	Ferromagnetic Iron core Internal cooling fan	Ferromagnetic Iron core Internal cooling fan	No Difference
Coil Positioning System	Integrated into Head Support System Laser-aided coil placement	Integrated into Head Support System Laser-aided coil placement	No Difference
Treatment Schedule	5 days per week for 4-6 weeks Total of 20-30 treatment sessions	5 days per week for 4-6 weeks Total of 20-30 treatment sessions	No Difference
Device Components	<ul style="list-style-type: none"> • Mobile Console • Ferromagnetic Coil for delivering treatment • Head Support System for coil positioning • MT Cap for coil positioning • D-Tect™ MT Accessory for MT location and level determination • Multi-use disposable for contact sensing and magnetic field quality control • TrakStar System for recording patient data 	<ul style="list-style-type: none"> • Mobile Console • Ferromagnetic Coil for delivering treatment • Head Support System for coil positioning • MT Cap for coil positioning • Multi-use disposable for contact sensing and magnetic field quality control • TrakStar System for recording patient data 	<p>Different</p> <p>The addition of the D-Tect™ MT Accessory is the subject of this special 510(k).</p>

%MT Range	25% to 140% MT	25% to 140% MT	No difference
Pulses per Second (PPS) Range	For treatment: 1 to 30 PPS For MT determination: 0.1-0.3 PPS	For treatment: 1 to 30 PPS For MT determination: 0.1-0.3 PPS	No difference
Induced Electric Field at 2 cm at 1.0 SMT	135 V/m (Nominal)	135 V/m (Nominal)	No difference
Pulse Type	Biphasic Sinusoid	Biphasic Sinusoid	No difference
Pulse Width	185 μ S (Nominal)	185 μ S (Nominal)	No difference
Treatment Protocols	<p>Standard Treatment: Level: 120% MT with allowable adjustments Repetition Rate: 10 PPS Stimulation Time: 4 s Inter-train Interval: As low as 11 s Session Duration: As low as 18.75 min Pulses per Session: 3000 Sessions per Week: 5</p>	<p>Standard Treatment: Level: 120% MT with allowable adjustments Repetition Rate: 10 PPS Stimulation Time: 4 s Inter-train Interval: As low as 11 s Session Duration: As low as 18.75 min Pulses per Session: 3000 Sessions per Week: 5</p>	No difference
	<p>NeuroBurst Treatment: Level: 80-120% MT with allowable adjustments Stimulation Time: 2 s Inter-train Interval: 8 s Pulses per Burst: 3 Interpulse Interval: 20 ms Session Duration: 3.3 min Pulses per Session: 600 Bursts per Second: 5 Amplitude: 0.22-2.08 SMT (\leq5% drop)</p>	<p>NeuroBurst Treatment: Level: 80-120% MT with allowable adjustments Stimulation Time: 2 s Inter-train Interval: 8 s Pulses per Burst: 3 Interpulse Interval: 20 ms Session Duration: 3.3 min Pulses per Session: 600 Bursts per Second: 5 Amplitude: 0.22-2.08 SMT (\leq5% drop)</p>	No difference
Treatment Level Range	<p>Standard Treatment: 0.22 to 2.08 SMT Calibrated linear output</p>	<p>Standard Treatment: 0.22 to 2.08 SMT Calibrated linear output</p>	No difference

	<p><u>NeuroBurst Treatment:</u> 0.22 to 1.9 SMT 80-120% MT ≤5% drop</p>	<p><u>NeuroBurst Treatment:</u> 0.22 to 1.9 SMT 80-120% MT ≤5% drop</p>	No difference
Stimulation Time Pulse Train Duration Range	<p><u>Standard Treatment:</u> 1 PPS: 1 to 600 s >1 PPS: 1 to 20 s</p>	<p><u>Standard Treatment:</u> 1 PPS: 1 to 600 s >1 PPS: 1 to 20 s</p>	No difference
	<p><u>NeuroBurst Treatment:</u> 1 to 10 s</p>	<p><u>NeuroBurst Treatment:</u> 1 to 10 s</p>	No difference
Inter-train Interval Range	<p><u>Standard Treatment:</u> 1 PPS: 0 to 600 s >1 PPS: 10 to 60 s</p>	<p><u>Standard Treatment:</u> 1 PPS: 0 to 600 s >1 PPS: 10 to 60 s</p>	No difference
	<p><u>NeuroBurst Treatment:</u> 1 to 60 s</p>	<p><u>NeuroBurst Treatment:</u> 1 to 60 s</p>	No difference
Pulses per Treatment Session	<p><u>Standard Treatment:</u> Nominal: 3000 Maximum: 5000</p>	<p><u>Standard Treatment:</u> Nominal: 3000 Maximum: 5000</p>	No difference
	<p><u>NeuroBurst Treatment:</u> Nominal: 600 Maximum: 2000</p>	<p><u>NeuroBurst Treatment:</u> Nominal: 600 Maximum: 2000</p>	No difference
Pulses per Burst (PPB)	<p><u>NeuroBurst Treatment:</u> 1 to 5</p>	<p><u>NeuroBurst Treatment:</u> 1 to 5</p>	No difference
Interpulse Interval	<p><u>NeuroBurst Treatment:</u> 20 to 2000 ms</p>	<p><u>NeuroBurst Treatment:</u> 20 to 2000 ms</p>	No difference
Bursts per Second (BPS)	<p><u>NeuroBurst Treatment:</u> 0.1 to 20.0 Hz</p>	<p><u>NeuroBurst Treatment:</u> 0.1 to 20.0 Hz</p>	No difference

Conclusion:

The NeuroStar Advanced Therapy System and the primary predicate device have the same indications for use and the same technological characteristics. The use of the optional D-Tect™ MT Accessory does not raise any new questions of safety or effectiveness.