



September 9, 2022

Steris

Anthony Piotrkowski  
Director Regulatory Affairs  
5960 Heisey Road  
Mentor, Ohio 44060

Re: K222440

Trade/Device Name: Vis-U-All Low Temperature Sterilization Pouches / Tubing, PRO-LITE  
Sterilization Trays

Regulation Number: 21 CFR 880.6850

Regulation Name: Sterilization Wrap

Regulatory Class: Class II

Product Code: FRG, KCT

Dated: August 11, 2022

Received: August 12, 2022

Dear Anthony Piotrkowski:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K222440

Device Name

Vis-U-All Low Temperature Sterilization Pouches/Tubing

Indications for Use (Describe)

The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:

- medical devices in a single or double pouch configuration
- trays\* containing medical devices in a single or double pouch configuration
- small items requiring surface sterilization in a single pouch configuration within a tray\*

NOTE: Trays must be legally marketed for use in the V PRO Low Temperature or STERRAD Sterilization Systems and contain a vent surface area to tray volume ratio  $\geq 0.135$  in-1 with the maximum number of instrument organizers installed.

to be sterilized in the:

- Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems
- Default Cycle of the STERRAD 100S Sterilizer
- Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers
- Express, Standard, Flex Scope and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers

\*STERRAD and ALLClear are trademarks of Advanced Sterilization Products

The pouches maintain the sterility of the enclosed devices until used.

When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.

Intended Sterilization Cycles and Intended Pouch Loads when Medical Devices are:

- Directly pouched
- Placed inside of a tray and the tray pouched

V-PRO 60 & s2 Lumen Cycle

- Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.
- Non-lumened devices including non-lumened rigid and semi-rigid endoscopes
- Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations:

o single or dual lumen devices

$\geq 0.77$  mm internal diameter (ID) and  $\leq 410$  mm in length

$\geq 1.8$  mm ID x  $\leq 542$  mm in length

o triple lumen devices

$\geq 1.2$  mm ID and  $\leq 275$  mm in length

$\geq 1.8$  mm ID and  $\leq 310$  mm in length

Or

$\geq 2.8$  mm ID and  $\leq 317$  mm in length

V-PRO 60 & s2 Non Lumen Cycle

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with

---

diffusion-restricted spaces such as the hinged portion of forceps and scissors.

#### V-PRO 60 & s2 Flexible Cycle

Load 1: One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:

- single or dual lumen device with lumens that are  $\geq 1$  mm ID and  $\leq 990$  mm in length

Load 2: Non-lumened devices including non-lumened rigid semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following configurations:

- $\geq 0.76$  mm ID and  $\leq 233$  mm in length
- $\geq 1.0$  mm ID and  $\leq 254$  mm in length
- $\geq 1.8$  mm ID and  $\leq 542$  mm in length

#### V-PRO s2 Fast Cycle

- Non-lumened devices including non-lumened rigid and semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

- Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations:

- o single or dual lumen devices
  - $\geq 0.77$  mm ID and  $\leq 410$  mm in length
  - $\geq 1.8$  mm ID  $\times \leq 542$  mm in length

- Triple channeled devices with stainless steel lumens that are either:

- $\geq 1.2$  mm ID and  $\leq 275$  mm in length
- $\geq 1.8$  mm ID and  $\leq 310$  mm in length

Or

- $\geq 2.8$  mm ID and  $\leq 317$  mm in length

#### V-PRO 1, 1 Plus, maX & maX 2 Lumen Cycle

- Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

- Non-lumened devices including non-lumened rigid and semi-rigid endoscopes

- Medical devices, including single, dual or triple channeled stainless steel lumens that are:

- $\geq 0.77$  mm ID and  $\leq 527$  mm in length
- $\geq 0.8$  mm ID and  $\leq 542$  mm in length
- $\geq 0.48$  mm ID and  $\leq 100$  mm in length

- Medical devices with Dead end lumens that are  $\geq 1.3$  mm ID and  $\leq 73$  mm in length

- Devices with rigid non-metallic lumens (such as those used in endoscope sheaths, take-apart forceps and trocars) that are:

- $\geq 3$  mm ID and  $\leq 298$  mm in length
- $\geq 4$  mm ID and  $\leq 424$  mm in length

#### V-PRO 1, 1 Plus, maX & maX2 Non Lumen Cycle

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

#### V-PRO maX and maX 2 Flexible Cycle

Load 1: Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.

The flexible endoscopes may contain either a single or dual channel lumen that is  $> 1$  mm ID and  $< 1050$  mm in length

Load 2:

- Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors

- Single, dual or triple channel stainless steel lumen that is  $\geq 0.48$  mm ID and  $\leq 100$  mm in length.

#### V-PRO maX 2 Fast Non Lumen Cycle

---

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

STERRAD 100S Default Cycle

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Metal and nonmetal lumened instruments with:

≥ 6 mm ID and ≤ 310 mm in length

Medical devices with a single stainless steel lumen with:

≥ 1 mm ID and ≤ 125 mm in length

≥ 2 mm ID and ≤ 250 mm in length

≥ 3 mm ID and ≤ 400 mm in length

STERRAD NX and NX with ALLClear Technology Standard Cycle

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical devices with a single stainless steel lumen with:

≥ 1 mm ID and ≤ 150 mm in length

≥ 2 mm ID and ≤ 400 mm in length

STERRAD NX and NX with ALLClear Technology Advanced Cycle

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical Devices, including most flexible endoscopes, with:

a single stainless steel lumen with:

≥ 1 mm ID and ≤ 500 mm in length

Single channel polyethylene and Teflon (polytetrafluoroethylene)

≥ 1 mm ID and ≤ 850 mm in length

STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical devices with a single stainless steel lumen with:

≥ 0.7 mm ID and ≤ 500 mm in length

STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical Devices, including most flexible endoscopes, with:

Single channel polyethylene and Teflon (polytetrafluoroethylene)

≥ 1 mm ID and ≤ 850 mm in length

STERRAD 100NX and 100NX with ALLClear Technology Express Cycle

Metal and nonmetal medical devices (surfaces sterilization only) and instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle

Medical devices including:

most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with ≥ 1 mm ID and ≤ 875 mm in length

accessory devices that are normally connected to a flexible endoscope during use

flexible endoscopes without lumens

---

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
*PRAStaff@fda.hhs.gov*

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## Indications for Use

510(k) Number (if known)

K222440

Device Name

PRO-LITE Sterilization Tray

Indications for Use (Describe)

The PRO-LITE Sterilization Trays are used to contain medical devices for sterilization in the following Cycles:

- Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO Low Temperature Sterilization Systems
- Default Cycle of the STERRAD®\* 100S Sterilizer
- Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers
- Standard, Flex Scope, Express and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers

\*STERRAD and ALLClear are trademarks of Advanced Sterilization Products

Prior to placing in the Sterilizer, the trays must either be:

- wrapped with a legally marketed sterilization wrap for use in the Sterilizers listed above
- or
- placed inside a legally marketed pouch for enclosing trays in the Sterilizers listed above.

The PRO-LITE Sterilization Trays are not intended to maintain sterility; they are intended to be used in conjunction with a validated, FDA-cleared sterilization wrap or pouch in order to maintain sterility of the enclosed medical instruments.

Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052

V-PRO 60 and s2 Lumen Cycle:

- Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors
  - Non-lumened devices including non-lumened rigid and semi-rigid endoscopes
  - Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes with the following configurations:
    - o Single or dual channeled devices with stainless steel lumens that are:
      - ≥ 0.77 mm ID and ≤ 410 mm in length
      - ≥ 1.8 mm ID x ≤ 542 mm in length
    - o Triple channeled devices with stainless steel lumens that are:
      - ≥ 1.2 mm ID and ≤ 275 mm in length
      - ≥ 1.8 mm ID and ≤ 310 mm in length
- or
- ≥ 2.8 mm ID and ≤ 317 mm in length

V-PRO 60 and s2 Non Lumen Cycle:

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

V-PRO 60 and s2 Flexible Cycle:

Load 1: One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:

- Single or dual lumen device with lumens that are ≥ 1 mm ID and ≤ 990 mm in length

Load 2: Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened

---

devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes with the following configurations:

- ≥ 0.76 mm ID and ≤ 233 mm in length
- ≥ 1.0 mm ID and ≤ 254 mm in length
- ≥ 1.8 mm ID and ≤ 542 mm in length

Intended Sterilization Cycle and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049

V-PRO s2 Fast Cycle:

- Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes, and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps and scissors.
- Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations:
  - o Single or dual channeled devices with stainless steel lumens
    - ≥ 0.77 mm ID and ≤ 410 mm in length
    - ≥ 1.8 mm ID and ≤ 542 mm in length
  - o Triple channeled devices with stainless steel lumens that are either
    - ≥ 1.2 mm ID and ≤ 275 mm in length
    - ≥ 1.8 mm ID and ≤ 310 mm in length
    - ≥ 2.8 mm ID and ≤ 317 mm in length

Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052, VP0053

V-PRO 1, 1 Plus, maX, and maX 2 Lumen Cycle:

- Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Medical devices, including single, dual or triple channeled stainless steel lumens that are:
  - ≥ 0.77 mm ID and ≤ 527 mm in length
  - ≥ 0.8 mm ID and ≤ 542 mm in length
  - ≥ 0.48 mm ID and ≤ 100 mm in length
- Medical devices with dead end stainless steel lumens that are ≥ 1.3 mm ID and ≤ 73 mm in length
- Devices with rigid non-metallic lumens (such as those used in endoscope sheaths, take-apart forceps and trocars) that are:
  - ≥ 3 mm ID and ≤ 298 mm in length
  - ≥ 4 mm ID and ≤ 424 mm in length

V-PRO 1 Plus, maX, and maX 2 Non Lumen Cycle:

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors

V-PRO maX, and maX 2 Flexible Cycle:

Load 1: Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain either a single or dual lumen that is ≥ 1 mm ID and ≤ 1050 mm in length.

Load 2:

- Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors
- Single, dual or triple channel stainless steel lumens that are ≥ 0.48 mm ID and ≤ 100 mm in length

V-PRO maX 2 Fast Non Lumen Cycle:

Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.

Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049,



**STERRAD 100S Default Cycle:**

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Metal and nonmetal lumened instruments with  
≥ 6 mm ID and ≤ 310 mm in length

Medical devices with a single stainless steel lumen with:

≥ 1 mm ID and ≤ 125 mm in length

≥ 2 mm ID and ≤ 250 mm in length

≥ 3 mm ID and ≤ 400 mm in length

Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049

**STERRAD NX and NX with ALLClear Technology Standard Cycle:**

Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical devices with a single stainless steel lumen with:

≥ 1 mm ID and ≤ 150 mm in length

≥ 2 mm ID and ≤ 400 mm in length

**STERRAD NX and NX with ALLClear Technology Advanced Cycle:**

Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors

Medical devices, including most flexible endoscopes, with:

o a single stainless steel lumen with:

≥ 1 mm ID and ≤ 500 mm in length

o single channel polyethylene and Teflon (polytetrafluoroethylene)

≥ 1mm ID and ≤ 850 mm in length

Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049, VP0051, VP0052, VP0053

**STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle:**

Metal and nonmetal medical devices including instruments with have diffusion-restricted spaces, such as the hinged portion of forceps and scissors

Medical devices with a single stainless steel lumen with:

≥ 0.7 mm ID and ≤ 500 mm in length

**STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle:**

Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

Medical devices, including most flexible endoscopes, with:

• Single channel polyethylene and Teflon (polytetrafluoroethylene)

≥ 1mm ID and ≤ 850 mm in length

**STERRAD 100NX and 100NX with ALLClear Technology Express Cycle:**

Metal and nonmetal devices surfaces and instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.

**STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle:**

Medical devices including:

• most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with ≥ 1 mm ID and ≤ 875 mm in length

• accessory devices that are normally connected to a flexible endoscope during use

---

• flexible endoscopes without lumens

---

Type of Use (*Select one or both, as applicable*)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

---

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



**510(k) Summary**  
**for**  
**Vis-U-All Low Temperature Sterilization Pouches/Tubing**

**Sponsor Facility**

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060  
Phone: (440) 354-2600  
Fax No: (440) 357-9198

Contact: Anthony Piotrkowski  
Director, Regulatory Affairs

Telephone: (440) 392-7437  
Fax No: (440) 357-9198  
e-mail: [tony\\_piotrkowski@steris.com](mailto:tony_piotrkowski@steris.com)

Submission Date: September 7, 2022

Premarket Notification Number: K222440

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

**1. Device Name**

Trade Name: Vis-U-All Low Temperature Sterilization Pouches/Tubing  
 Device Classification: Class II  
 Common/Usual Name: Sterilization pouch  
 Classification Name: Sterilization wrap  
 Classification Number: 21 CFR 880.6850  
 Product Code: FRG

**2. Predicate Device**

Vis-U-All Low Temperature Sterilization Pouches/Tubing, K183401

**3. Description of Device**

The proposed Vis-U-All Low Temperature Sterilization Pouches/Tubing is identical to the predicate and is a Tyvek/plastic film sterilization containment pouch designed for devices to be sterilized by the health care provider in V-PRO or STERRAD Low Temperature Sterilization Systems. As is the predicate device, the proposed device is available as a self-seal pouch, a heat-seal pouch, or heat-seal tubing. Available sizes and configurations are shown in **Table 1**.

**Table 1.** Sizes and Configurations of Vis-U-All Low Temperature Sterilization Pouches/Tubing

Type	Size *	Type	Size*	Type	Size*
Heat Seal Pouch	3 x 7	Self Seal Pouch	3 x 7	Tubing	3" x 100'
	4 x 9		4 x 9		4" x 100'
	4 x 12		4 x 12		6" x 100'
	4 x 22		4 x 22		9" x 100'
	6 x 10		6 x 10		14" x 100'
	8 x 12		8 x 12		
	10 x 15		10 x 15		
	12 x 18		12 x 18		
	8 x 21				
	8 x 27				
	9 x 27				
	11 x 22				
	12 x 27				

*\*Sizes are in inches unless specified otherwise*

The purpose of this submission is to qualify use of the Vis-U-All Low Temperature Sterilization Pouches and Tubing for extended claims in the V-PRO 60, s2, maX and maX 2 Sterilizer Cycles.

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

**4. Intended Use/ Indications for Use**

The Vis-U-All Low Temperature Sterilization Pouches/Tubing Intended use remains the same. The indications for use are being modified to expand lumen claims. A detailed comparison of the cleared and proposed indications for use are provided in **Table 2**.

**5. Comparison of the Modified Device to the Predicate**

The proposed and predicate devices are single use sterilization pouches for use in V-PRO Sterilizers. **Table 2** summarizes the difference between the proposed device and predicate device cleared under K183401.

**Table 2.** Technical Comparison to the Predicate.

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch ( <b>proposed</b> ) K222440	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
Intended Use / Indications for Use	<p>The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:</p> <ul style="list-style-type: none"> <li>• medical devices in a single or double pouch configuration</li> <li>• trays* containing medical devices in a single or double pouch configuration</li> <li>• small items requiring surface sterilization in a single pouch configuration within a tray*</li> </ul> <p>NOTE: Trays must be legally marketed for use in the V-PRO Low Temperature or STERRAD Sterilization Systems and contain a vent surface area to tray volume ratio <math>\geq 0.135 \text{ in}^{-1}</math> with the maximum number of instrument organizers installed.</p> <p>to be sterilized in the:</p> <ul style="list-style-type: none"> <li>• Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems</li> <li>• Default Cycle of the STERRAD 100S Sterilizer</li> <li>• Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers</li> <li>• Express, Standard, Flex Scope and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers</li> </ul>	<p>The Vis-U-All Low Temperature Sterilization Pouches/Tubing are sterilization containment pouches for use by health care providers to enclose:</p> <ul style="list-style-type: none"> <li>• medical devices in a single or double pouch configuration</li> <li>• trays containing medical devices in a single or double pouch configuration</li> <li>• small items requiring surface sterilization in a single pouch configuration within a tray</li> </ul> <p>NOTE: Trays must be legally marketed for use in the V-PRO Low Temperature or STERRAD Sterilization Systems and contain a vent surface area to tray volume ratio <math>\geq 0.135 \text{ in}^{-1}</math> with the maximum number of instrument organizers installed.</p> <p>to be sterilized in the:</p> <ul style="list-style-type: none"> <li>• Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO® Low Temperature Sterilization Systems</li> <li>• Default Cycle of the STERRAD 100S Sterilizer</li> <li>• Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers</li> <li>• Express, Standard, Flex Scope and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers</li> </ul>	<p>Identical except for:</p> <ul style="list-style-type: none"> <li>• Clarification of lumen claims for V-PRO lumen cycles</li> <li>• Clarification of endoscope claims for V-PRO flexible cycles</li> <li>• Removal of “stainless steel or titanium” from V-PRO Non Lumen cycles.</li> </ul>

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch ( <b>proposed</b> ) <b>K222440</b>	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
	<p>*STERRAD and ALLClear are trademarks of Advanced Sterilization Products</p> <p>The pouches maintain the sterility of the enclosed devices until used.</p> <p>When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.</p> <p>Intended Sterilization Cycles and Intended Pouch Loads when Medical Devices are:</p> <ul style="list-style-type: none"> <li>• Directly pouched</li> <li>• Placed inside of a tray and the tray pouched</li> </ul> <p><u>V-PRO 60 &amp; s2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm internal diameter (ID) and <math>\leq 410</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID x <math>\leq 542</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> </li> </ul> <p align="center">Or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul>	<p>*STERRAD and ALLClear are trademarks of Advanced Sterilization Products</p> <p>The pouches maintain the sterility of the enclosed devices until used.</p> <p>When used to enclose medical devices, the pouches are intended to contain the devices in such a manner as to leave a minimum of one inch between the devices and seal on all sides. When used to enclose a tray, the tray must fit loosely within the pouch.</p> <p>Intended Sterilization Cycles and Intended Pouch Loads when Medical Devices are:</p> <ul style="list-style-type: none"> <li>• Directly pouched</li> <li>• Placed inside of a tray and the tray pouched</li> </ul> <p><u>V-PRO 60 &amp; s2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm internal diameter (ID) and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> </li> </ul> <p align="center">or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch ( <b>proposed</b> ) <b>K222440</b>	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
	<p><u>V-PRO 60 &amp; s2 Non Lumen Cycle</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 &amp; s2 Flexible Cycle</u> <u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p> <ul style="list-style-type: none"> <li>• single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following configurations:</p> <ul style="list-style-type: none"> <li>• <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>• <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> <li>• <math>\geq 1.8</math> mm ID and <math>\leq 542</math> mm in length</li> </ul> <p><u>V-PRO s2 Fast Cycle</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid and semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations: <ul style="list-style-type: none"> <li>○ single or dual lumen devices <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID x <math>\leq 542</math> mm in length</li> </ul> </li> </ul> </li> <li>• Triple channeled devices with stainless steel lumens that are either:</li> </ul>	<p><u>V-PRO 60 &amp; s2 Non Lumen Cycle</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 &amp; s2 Flexible Cycle</u> <u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p> <ul style="list-style-type: none"> <li>• single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes, with the following dimensions:</p> <ul style="list-style-type: none"> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> <li>• <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>• <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> </ul> <p><u>V-PRO s2 Fast Cycle</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid and semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations: <ul style="list-style-type: none"> <li>▪ single or dual channeled devices with stainless steel lumens that are <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>• Triple channeled devices with stainless steel lumens that are either:</li> </ul>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch (proposed) K222440	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
	<ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p align="center">Or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX 2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual or triple channeled stainless steel lumens that are: <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 527</math> mm in length</li> <li>▪ <math>\geq 0.8</math> mm ID and <math>\leq 542</math> mm in length</li> <li>▪ <math>\geq 0.48</math> mm ID and <math>\leq 100</math> mm in length</li> </ul> </li> <li>• Medical devices with Dead end lumens that are <math>\geq 1.3</math> mm ID and <math>\leq 73</math> mm in length</li> <li>• Devices with rigid non-metallic lumens (such as those used in endoscope sheaths, take-apart forceps and trocars) that are: <ul style="list-style-type: none"> <li>▪ <math>\geq 3</math> mm ID and <math>\leq 298</math> mm in length</li> <li>▪ <math>\geq 4</math> mm ID and <math>\leq 424</math> mm in length</li> </ul> </li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX2 Non Lumen Cycle</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX and maX 2 Flexible Cycle</u></p>	<ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p align="center">Or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX 2 Lumen Cycle</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ <u>single or dual lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ <u>triple lumen devices</u> <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p align="center">or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 1, 1 Plus, maX &amp; maX2 Non Lumen Cycle</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX and maX 2 Flexible Cycle</u></p>	



**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch ( <b>proposed</b> ) <b>K222440</b>	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
	<p>Load 1: Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain either a single or dual channel lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</p> <p>Load 2:</p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Single, dual or triple channel stainless steel lumen that is <math>\geq 0.48</math> mm ID and <math>\leq 100</math> mm in length.</li> </ul> <p><u>V-PRO maX 2 Fast Non Lumen Cycle</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>STERRAD 100S Default Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Metal and nonmetal lumened instruments with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 6</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 125</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 250</math> mm in length</li> <li>• <math>\geq 3</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAD NX and NX with ALLClear Technology Standard Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p>	<p>Load 1: Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load. The flexible endoscopes may contain either:</p> <ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>• one lumen that is <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> <li>• and the other lumen that is <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p>Load 2: Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.</p> <p><u>V-PRO maX 2 Fast Non Lumen Cycle</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>STERRAD 100S Default Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Metal and nonmetal lumened instruments with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 6</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 125</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 250</math> mm in length</li> <li>• <math>\geq 3</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAD NX and NX with ALLClear Technology Standard Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	Modified Vis-U-All Low Temperature Sterilization Pouch (proposed) K222440	Vis-U-All Low Temperature Sterilization Pouch (K183401)	Comparison
	<p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 150</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAD NX and NX with ALLClear Technology Advanced Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p> <p>Medical Devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>• a single stainless steel lumen with: <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> </li> <li>• Single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 0.7</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical Devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>• Single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Express Cycle</u> Metal and nonmetal medical devices (surfaces sterilization only) and instruments which have diffusion-</p>	<p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 150</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAD NX and NX with ALLClear Technology Advanced Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p> <p>Medical Devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>• a single stainless steel lumen with: <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> </li> <li>• Single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 0.7</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical Devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>• Single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Express Cycle</u> Metal and nonmetal medical devices (surfaces sterilization only) and</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

<b>Feature</b>	<b>Modified Vis-U-All Low Temperature Sterilization Pouch (proposed) K222440</b>	<b>Vis-U-All Low Temperature Sterilization Pouch (K183401)</b>	<b>Comparison</b>
	restricted spaces, such as the hinged portion of forceps and scissors.  <u>STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle</u> Medical devices including: <ul style="list-style-type: none"> <li>• most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with <math>\geq 1</math> mm ID and <math>\leq 875</math> mm in length</li> <li>• accessory devices that are normally connected to a flexible endoscope during use</li> <li>• flexible endoscopes without lumens</li> </ul>	instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.  <u>STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle</u> Medical devices including: <ul style="list-style-type: none"> <li>• most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with <math>\geq 1</math> mm ID and <math>\leq 875</math> mm in length</li> <li>• accessory devices that are normally connected to a flexible endoscope during use</li> <li>• flexible endoscopes without lumens</li> </ul>	
Device Features	Chevron end of pouches for ease of opening Chemical process indicator for EO	Chevron end of pouches for ease of opening Chemical process indicator for EO	Identical
Maintenance of Sterility	1 year	1 year	Identical
Materials of Construction	Tyvek and plastic	Tyvek and plastic	Identical
Types	Self Seal, Heat Seal, Tubing	Self Seal, Heat Seal, Tubing	Identical

**6. Summary of Performance Testing**

**Table 3** summarizes the testing of the Vis-U-All Low Temperature Sterilization Pouches/Tubing to demonstrate that the proposed pouch is qualified for use in V-PRO Low temperature Sterilization Systems.

**Table 3.** Performance Test Summary

<b>Test</b>	<b>Result</b>	<b>Conclusion</b>
<b>½ Cycle Verification of Mated Surfaces</b>	Sterile efficacy was demonstrated for mated surfaces under worst case conditions in the V-PRO Sterilizer cycles.	PASS
<b>½ Cycle Efficacy</b>	The standard injection weight resulted in all sterile results within the validation load used to qualify each sterilizer cycle.	PASS
<b>Simulated Use Test</b>	Simulated use testing verified the ability of the sterilizer cycles to sterilize medical devices under worst-case processing conditions.	PASS
<b>In Use Test</b>	The in use investigation demonstrated the ability of the V-PRO Sterilizer cycles to sterilize patient-soiled, clinically-cleaned, medical instruments.	PASS

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION**  
**K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing**  
**PRO-LITE Sterilization Trays**

---

**7. Conclusion**

Based on the intended use, technological characteristics and non-clinical performance data, the subject device is as safe, as effective and performs as well as or better than the legally marketed predicate device (K183401), Class II (21 CFR 880.6850), product code FRG.

**510(k) Summary  
For  
PRO-LITE™ Sterilization Tray**

**Sponsor Facility**

STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060  
Phone: (440) 354-2600  
Fax No: (440) 357-9198

Contact: Anthony Piotrkowski  
Director, Regulatory Affairs

Telephone: (440) 392-7437  
Fax No: (440) 357-9198  
e-mail: [tony\\_piotrkowski@steris.com](mailto:tony_piotrkowski@steris.com)

Submission Date: September 7, 2022

Premarket Notification Number: K222440

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION**  
**K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing**  
**PRO-LITE Sterilization Tray**

**7. Device Name**

Trade Name: PRO-LITE Sterilization Tray  
 Common/usual Name: Sterilization Trays, cassettes and other accessories  
 Classification Name: Sterilization Wrap  
 Classification: 21 CFR 880.6850  
 Class: II  
 Product Code: KCT, FRG

**8. Predicate Device**

PRO-LITE Sterilization Tray, K183402

**9. Description of Device**

The PRO-LITE Sterilization Trays contain medical devices for sterilization in the V-PRO Low Temperature Sterilization Systems and STERRAD Sterilizers. The trays are available in various sizes, outlined in **Table 4**, to accommodate the loads to be processed. The proposed trays are identical in design to the predicate Sterilization Tray (K183402) and are composed of a base and a lid. The lid includes clamping mechanisms designed to secure the lid onto the base. There are numerous holes in the base and lid for sterilant penetration. The tray is categorized as a cassette and requires complete enclosure in a legally-marketed sterilization wrap or pouch to maintain sterility of the enclosed devices. Both the base and the lid for the proposed tray are made of a mineral-filled polypropylene material.

**Table 4.** External Dimensions of Tray Line

Model	Description (in)	Model	Description (in)	Model	Description (in)
VP0045	13 x 4.5 x 2.25	VP0048	13 x 7.75 x 2.25	VP0051	12 x 11.75 x 4
VP0046	19 x 4.5 x 2.25	VP0049	19 x 7.75 x 2.25	VP0052	25 x 11.75 x 4
VP0047	25 x 4.5 x 2.25	VP0050	27 x 7.75 x 4	VP0053	25 x 14 x 5

Optional instrument organizers are provided as accessories to the trays and allow stabilization of various cylindrical medical devices during processing. **Table 5** lists the organizer sizes. The organizers are either “blank” and are used to partition the tray or have a device holding portion into which the devices are inserted. At the organizer base is a flapped groove that is used to position the organizer within a PRO-LITE Sterilization Tray.

**Table 5.** Instrument Organizer Model Numbers

Model	Description	Model	Description
VP0054	Blank, Tall	VP0055	Blank, Short
VP0063	3 mm - 7mm, Tall	VP0067	3 mm - 7 mm, Short
VP0064	7 mm - 11mm, Tall	VP0068	7 mm -11 mm, Short
VP0065	11 mm - 15 mm, Tall	VP0069	11 mm - 15 mm, Short
VP0066	15 mm - 19mm, Tall	VP0070	15 mm – 19 mm, Short

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

---

Optional sterilization mats are provided as accessories to the trays. The mats, which are used to cushion and stabilize devices placed into the trays, are available in sizes as shown in **Table 6** to fit the nine PRO-LITE Sterilization Trays. The mats are a diamond grid design with “fingers” that extend from each corner of the diamond and at the midpoint of each diamond side. The fingers cushion and stabilize instruments, helping to prevent the instruments from freely moving in the tray during packaging, sterilization and storage. The cushioning and stabilization qualities help protect delicate instruments placed into the trays.

**Table 6.** Silicone Mat Model Numbers

<b>Model</b>	<b>Description (in)</b>	<b>Model</b>	<b>Description (in)</b>	<b>Model</b>	<b>Description (in)</b>
VP0071	13 x 4.5	VP0074	13 x 7.75	VP0077	12 x 11.75
VP0072	19 x 4.5	VP0075	19 x 7.75	VP0078	25 x 11.75
VP0073	25 x 4.5	VP0076	27 x 7.75	VP0079	25 x 14

The purpose of this submission is to expand claims for the use of these tray models in the following sterilizer cycles:

- V-PRO 1, 1 Plus, maX and maX 2 Lumen Cycle
- V-PRO 1 Plus, maX and maX 2 Non Lumen Cycle
- V-PRO maX and maX 2 Flexible Cycle
- V-PRO maX 2 Fast Non Lumen Cycle
- V-PRO 60 and s2 Lumen, Non Lumen and Flexible Cycles
- V-PRO s2 Fast Cycle

**10. Intended Use/ Indications for Use**

The tray indications for use are included in Table 7 which compares the proposed device to the predicate.

Instrument organizers are optional accessories intended to stabilize cylindrical medical instruments within the PRO-LITE Sterilization Trays.

<b>Model</b>	<b>Description</b>	<b>Model</b>	<b>Description</b>
VP0054	Blank, Tall	VP0055	Blank, Short
VP0063	3 mm - 7 mm, Tall	VP0067	3 mm - 7 mm, Short
VP0064	7 mm - 11 mm, Tall	VP0068	7 mm - 11 mm, Short
VP0065	11 mm - 15 mm, Tall	VP0069	11 mm - 15 mm, Short
VP0066	15 mm - 19 mm, Tall	VP0070	15 mm - 19 mm, Short

Sterilization mats are optional accessories intended to cushion and stabilize instruments within the PRO-LITE Sterilization Trays.

<b>Model</b>	<b>Description (in)</b>	<b>Model</b>	<b>Description (in)</b>	<b>Model</b>	<b>Description (in)</b>
VP0071	13 x 4.5	VP0074	13 x 7.75	VP0077	12 x 11.75
VP0072	19 x 4.5	VP0075	19 x 7.75	VP0078	25 x 11.75
VP0073	25 x 4.5	VP0076	27 x 7.75	VP0079	25 x 14

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

**11. Summary of Technical Characteristics**

The proposed PRO-LITE sterilization trays, sterilization mats and instrument organizers are identical in composition to the claimed predicate device. The technical characteristics are summarized below in **Table 7**.

**Table 7. Summary of Tray Physical Description and Technological Properties**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
Intended Use / Indications for Use	<p>The PRO-LITE Sterilization Trays are used to contain medical devices for sterilization in the following Cycles:</p> <ul style="list-style-type: none"> <li>• Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO Low Temperature Sterilization Systems</li> <li>• Default Cycle of the STERRAD®* 100S Sterilizer</li> <li>• Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers</li> <li>• Standard, Flex Scope, Express and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers</li> </ul> <p>*STERRAD and ALLClear are trademarks of Advanced Sterilization Products</p> <p>Prior to placing in the Sterilizer, the trays must either be:</p> <ul style="list-style-type: none"> <li>• wrapped with a legally marketed sterilization wrap for use in the Sterilizers listed above</li> </ul> <p>or</p> <ul style="list-style-type: none"> <li>• placed inside a legally marketed pouch for enclosing trays in the Sterilizers listed above.</li> </ul> <p>The PRO-LITE Sterilization Trays are not intended to maintain sterility; they are intended to be used in conjunction with a validated, FDA-cleared sterilization wrap or pouch in order to maintain sterility of the enclosed medical instruments.</p>	<p>The PRO-LITE Sterilization Trays are used to contain medical devices for sterilization in the following Cycles:</p> <ul style="list-style-type: none"> <li>• Lumen, Non Lumen, Flexible, Fast Non Lumen and Fast Cycles of the V-PRO Low Temperature Sterilization Systems</li> <li>• Default Cycle of the STERRAD®* 100S Sterilizer</li> <li>• Standard and Advanced Cycles of the STERRAD NX and NX with ALLClear Technology Sterilizers</li> <li>• Standard, Flex Scope, Express and DUO Cycles of the STERRAD 100NX and 100NX with ALLClear Technology Sterilizers</li> </ul> <p>*STERRAD and ALLClear are trademarks of Advanced Sterilization Products</p> <p>Prior to placing in the Sterilizer, the trays must either be:</p> <ul style="list-style-type: none"> <li>• wrapped with a legally marketed sterilization wrap for use in the Sterilizers listed above</li> </ul> <p>or</p> <ul style="list-style-type: none"> <li>• placed inside a legally marketed pouch for enclosing trays in the sterilizers listed above</li> </ul> <p>The PRO-LITE Sterilization Trays are not intended to maintain sterility; they are intended to be used in conjunction with a validated, FDA-cleared sterilization wrap or pouch in order to maintain sterility of the enclosed medical instruments.</p>	<p>Identical except for:</p> <ul style="list-style-type: none"> <li>• Clarification of lumen claims for V-PRO lumen cycles</li> <li>• Clarification of endoscope claims for V-PRO flexible cycles</li> </ul> <p>Removal of “stainless steel or titanium” from V-PRO Non Lumen cycles.</p>



**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052</p> <p><u>V-PRO 60 and s2 Lumen Cycle:</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes with the following configurations:               <ul style="list-style-type: none"> <li>○ Single or dual channeled devices with stainless steel lumens that are:                   <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID x <math>\leq 542</math> mm in length</li> </ul> </li> <li>○ Triple channeled devices with stainless steel lumens that are:                   <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p>or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 60 and s2 Non Lumen Cycle:</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 and s2 Flexible Cycle:</u> <u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p>	<p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052</p> <p><u>V-PRO 60 and s2 Lumen Cycle:</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes with the following configurations:               <ul style="list-style-type: none"> <li>○ Single or dual lumen devices                   <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ Triple lumen devices                   <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p>or</p> <ul style="list-style-type: none"> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> <p><u>V-PRO 60 and s2 Non Lumen Cycle:</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO 60 and s2 Flexible Cycle:</u> <u>Load 1:</u> One flexible surgical endoscope or bronchoscope with a light cord (if not integral to endoscope) and mat without any additional load. The flexible endoscope may be a:</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<ul style="list-style-type: none"> <li>• Single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes with the following configurations:</p> <ul style="list-style-type: none"> <li>• <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>• <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> <li>• <math>\geq 1.8</math> mm ID and <math>\leq 542</math> mm in length</li> </ul> <p>Intended Sterilization Cycle and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049</p> <p><u>V-PRO s2 Fast Cycle:</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes, and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps and scissors.</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations: <ul style="list-style-type: none"> <li>○ Single or dual channeled devices with stainless steel lumens <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 542</math> mm in length</li> </ul> </li> <li>○ Triple channeled devices with stainless steel lumens that are either <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li>• Single or dual lumen device with lumens that are <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> </ul> <p><u>Load 2:</u> Non-lumened devices including non-lumened rigid, semi-rigid, and flexible endoscopes and non-lumened devices with diffusion-restricted areas such as the hinged portion of forceps or scissors. Medical devices, including rigid and semi-rigid endoscopes with the following configurations:</p> <ul style="list-style-type: none"> <li>• <math>\geq 2.0</math> mm ID and <math>\leq 400</math> mm in length</li> <li>• <math>\geq 0.76</math> mm ID and <math>\leq 233</math> mm in length</li> <li>• <math>\geq 1.0</math> mm ID and <math>\leq 254</math> mm in length</li> </ul> <p><u>V-PRO s2 Fast Cycle:</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Medical devices (including single, dual and triple channeled rigid and semi-rigid endoscopes) with the following configurations: <ul style="list-style-type: none"> <li>▪ Single or dual channeled devices with stainless steel lumens that are <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ Triple channeled devices with stainless steel lumens that are either <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052, VP0053</p> <p><u>V-PRO 1, 1 Plus, maX, and maX 2 Lumen Cycle:</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Medical devices, including single, dual or triple channeled stainless steel lumens that are: <ul style="list-style-type: none"> <li>○ <math>\geq 0.77</math> mm ID and <math>\leq 527</math> mm in length</li> <li>○ <math>\geq 0.8</math> mm ID and <math>\leq 542</math> mm in length</li> <li>○ <math>\geq 0.48</math> mm ID and <math>\leq 100</math> mm in length</li> </ul> </li> <li>• Medical devices with dead end stainless steel lumens that are <math>\geq 1.3</math> mm ID and <math>\leq 73</math> mm in length</li> <li>• Devices with rigid non-metallic lumens (such as those used in endoscope sheaths, take-apart forceps and trocars) that are: <ul style="list-style-type: none"> <li>○ <math>\geq 3</math> mm ID and <math>\leq 298</math> mm in length</li> <li>○ <math>\geq 4</math> mm ID and <math>\leq 424</math> mm in length</li> </ul> </li> </ul> <p><u>V-PRO 1 Plus, maX, and maX 2 Non Lumen Cycle:</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</p> <p><u>V-PRO maX, and maX 2 Flexible Cycle:</u></p> <p><u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.</p>	<p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0050, VP0051, VP0052, VP0053</p> <p><u>V-PRO 1, 1 Plus, maX, and maX 2 Lumen Cycle:</u></p> <ul style="list-style-type: none"> <li>• Instruments with diffusion-restricted spaces such as the hinged portions of forceps and scissors</li> <li>• Non-lumened devices including non-lumened rigid and semi-rigid endoscopes</li> <li>• Medical devices, including single, dual and triple channeled rigid and semi-rigid endoscopes, with the following configurations: <ul style="list-style-type: none"> <li>○ single or dual lumen devices <ul style="list-style-type: none"> <li>▪ <math>\geq 0.77</math> mm ID and <math>\leq 410</math> mm in length</li> </ul> </li> <li>○ triple lumen devices <ul style="list-style-type: none"> <li>▪ <math>\geq 1.2</math> mm ID and <math>\leq 275</math> mm in length</li> <li>▪ <math>\geq 1.8</math> mm ID and <math>\leq 310</math> mm in length</li> <li>▪ <math>\geq 2.8</math> mm ID and <math>\leq 317</math> mm in length</li> </ul> </li> </ul> </li> </ul> <p><u>V-PRO 1 Plus, maX, and maX 2 Non Lumen Cycle:</u></p> <p>Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p><u>V-PRO maX, and maX 2 Flexible Cycle:</u></p> <p><u>Load 1:</u> Single or dual lumen surgical flexible endoscopes (such as those used in ENT, Urology and Surgical Care) and bronchoscopes with a light cord (if not integral to endoscope) and mat with no additional load.</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<p>The flexible endoscopes may contain either a single or dual lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</p> <p><u>Load 2:</u></p> <ul style="list-style-type: none"> <li>• Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors</li> <li>• Single, dual or triple channel stainless steel lumens that are <math>\geq 0.48</math> mm ID and <math>\leq 100</math> mm in length</li> </ul> <p><u>V-PRO maX 2 Fast Non Lumen Cycle:</u> Non-lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with diffusion-restricted spaces such as the hinged portion of forceps and scissors.</p> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052</p> <p><u>STERRAD 100S Default Cycle:</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Metal and nonmetal lumened instruments with</p> <ul style="list-style-type: none"> <li>• <math>\geq 6</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 125</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 250</math> mm in length</li> <li>• <math>\geq 3</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049</p> <p><u>STERRAD NX and NX with ALLClear Technology Standard Cycle:</u></p>	<p>The flexible endoscopes may contain either:</p> <ul style="list-style-type: none"> <li>• a single lumen that is <math>\geq 1</math> mm ID and <math>\leq 1050</math> mm in length</li> <li>• or two lumens with: <ul style="list-style-type: none"> <li>○ one lumen that is <math>\geq 1</math> mm ID and <math>\leq 990</math> mm in length</li> <li>○ and the other lumen that is <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>Load 2:</u> Non-lumened instruments including instruments with diffusion-restricted areas such as the hinged portion of forceps or scissors.</p> <p><u>V-PRO maX 2 Fast Non Lumen Cycle:</u> Non-Lumened devices including non-lumened rigid, semi-rigid and flexible endoscopes and non-lumened devices with stainless steel or titanium diffusion-restricted spaces such as the hinged portion of forceps and scissors</p> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0047, VP0048, VP0049, VP0050, VP0051, VP0052</p> <p><u>STERRAD 100S Default Cycle:</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Metal and nonmetal lumened instruments with</p> <ul style="list-style-type: none"> <li>• <math>\geq 6</math> mm ID and <math>\leq 310</math> mm in length</li> </ul> <p>Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 125</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 250</math> mm in length</li> <li>• <math>\geq 3</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<p>Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 150</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAND NX and NX with ALLClear Technology Advanced Cycle:</u> Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>○ a single stainless steel lumen with: <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> </li> <li>○ single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049, VP0051, VP0052, VP0053</p> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle:</u> Metal and nonmetal medical devices including instruments with have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 0.7</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle:</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices, including most flexible endoscopes, with:</p>	<p><u>STERRAD NX and NX with ALLClear Technology Standard Cycle:</u> Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <ul style="list-style-type: none"> <li>• <math>\geq 1</math> mm ID and <math>\leq 150</math> mm in length</li> <li>• <math>\geq 2</math> mm ID and <math>\leq 400</math> mm in length</li> </ul> <p><u>STERRAND NX and NX with ALLClear Technology Advanced Cycle:</u> Metal and non-metal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>○ a single stainless steel lumen with: <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 500</math> mm in length</li> </ul> </li> <li>○ single channel polyethylene and Teflon (polytetrafluoroethylene) <ul style="list-style-type: none"> <li>○ <math>\geq 1</math> mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p>Intended Sterilization Cycles and Intended Tray Load for Tray Models: VP0045, VP0046, VP0048, VP0049, VP0051, VP0052, VP0053</p> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Standard Cycle:</u> Metal and nonmetal medical devices including instruments with have diffusion-restricted spaces, such as the hinged portion of forceps and scissors. Medical devices with a single stainless steel lumen with:</p> <p><math>\geq 0.7</math> mm ID and <math>\leq 500</math> mm in length</p> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Flex Scope Cycle:</u> Metal and nonmetal medical devices including instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p>	

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

Feature	PRO-LITE Sterilization Tray (proposed) K222440	PRO-LITE Sterilization Tray (K183402)	Comparison
	<ul style="list-style-type: none"> <li>Single channel polyethylene and Teflon (polytetrafluoroethylene)                             <ul style="list-style-type: none"> <li>o <math>\geq 1</math>mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Express Cycle:</u> Metal and nonmetal devices surfaces and instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle:</u> Medical devices including:</p> <ul style="list-style-type: none"> <li>most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with <math>\geq 1</math> mm ID and <math>\leq 875</math> mm in length</li> <li>accessory devices that are normally connected to a flexible endoscope during use</li> <li>flexible endoscopes without lumens</li> </ul>	<p>Medical devices, including most flexible endoscopes, with:</p> <ul style="list-style-type: none"> <li>Single channel polyethylene and Teflon (polytetrafluoroethylene)                             <ul style="list-style-type: none"> <li>o <math>\geq 1</math>mm ID and <math>\leq 850</math> mm in length</li> </ul> </li> </ul> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Express Cycle:</u> Metal and nonmetal devices surfaces and instruments which have diffusion-restricted spaces, such as the hinged portion of forceps and scissors.</p> <p><u>STERRAD 100NX and 100NX with ALLClear Technology Duo Cycle:</u> Medical devices including:</p> <ul style="list-style-type: none"> <li>most flexible endoscopes with a single channel of polyethylene and Teflon (polytetrafluoroethylene) with <math>\geq 1</math> mm ID and <math>\leq 875</math> mm in length</li> <li>accessory devices that are normally connected to a flexible endoscope during use</li> <li>flexible endoscopes without lumens</li> </ul>	
Vent to Volume Ratio	All trays are the same: 0.135 in <sup>-1</sup>	All trays are the same: 0.135 in <sup>-1</sup>	Identical
Tray Composition	Mineral-filled polypropylene, stainless steel	Mineral-filled polypropylene, stainless steel	Identical
Instrument Organizer Composition	Medical Grade Silicone, UPS grade VI	Medical Grade Silicone, UPS grade VI	Identical
Mat Composition	Medical Grade Silicone, UPS grade VI	Medical Grade Silicone, UPS grade VI	Identical

**12. Summary of Non-clinical Tests**

Performance testing to demonstrate substantial equivalence to the predicate has been completed and is summarized in **Table 8** below.

**Table 8.** Performance Test Summary

Test	Result	Conclusion
<b>½ Cycle Verification of Mated Surfaces</b>	Sterile efficacy was demonstrated for mated surfaces under worst case conditions in the V-PRO Sterilizer cycles.	PASS
<b>½ Cycle Efficacy</b>	The standard injection weight resulted in all sterile results within the validation load used to qualify each sterilizer cycle.	PASS

**STERIS SPECIAL 510(k) PREMARKET NOTIFICATION  
K222440 Vis-U-All Low Temperature Sterilization Pouches/Tubing  
PRO-LITE Sterilization Trays**

<b>Test</b>	<b>Result</b>	<b>Conclusion</b>
<b>Simulated Use Test</b>	Simulated use testing verified the ability of the sterilizer cycles to sterilize medical devices under worst-case processing conditions.	PASS
<b>In Use Test</b>	The in use investigation demonstrated the ability of the V-PRO Sterilizer cycles to sterilize patient-soiled, clinically-cleaned, medical instruments.	PASS

**7. Conclusion**

Based on the intended use, technological characteristics and non-clinical performance data, the subject device is as safe, as effective and performs as well as or better than the legally marketed predicate device (K183402), Class II (21 CFR 880.6850), product code KCT.