



December 8, 2020

Medentika GmbH
% Jennifer Jackson
Director, Regulatory Affairs
Straumann USA, LLC
60 Minuteman Road
Andover, Massachusetts 01810

Re: K201711

Trade/Device Name: CARES[®] Screw - Retained Bars and Bridges (SRBB)

Regulation Number: 21 CFR 872.3630

Regulation Name: Endosseous Dental Implant Abutment

Regulatory Class: Class II

Product Code: NHA

Dated: November 5, 2020

Received: November 9, 2020

Dear Jennifer Jackson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Andrew Steen
Assistant Director
DHT1B: Division of Dental Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K201711

Device Name:

CARES® Screw - Retained Bars and Bridges (SRBB)

Indications for Use (Describe)

CARES® Screw-retained Bridges and Bars are indicated for use as bars and bridges that attach to implants to provide support for prosthetic reconstructions such as bridges and overdentures. The final processed products have the purpose of restoring chewing function. CARES® Screw-retained Bridges and Bars are indicated for screw-retained restorations.

CARES® Screw - Retained Bars and Bridges are designed to interface with the following dental implant systems (Implant System Compatibility / Series / Implant diameter (mm) / Platform diameter (mm)):

Nobel Biocare Replace Select / E-Series / Diameters 3.5, 4.3, 5.0, 6.0 / Platform 3.5, 4.3, 5.0, 6.0
Dentsply Implants - ASTRA TECH OsseoSpeed EV / EV-Series / 3.0, 3.6, 4.2, 4.8, 5.4 / Platform 3.0, 3.6, 4.2, 4.8, 5.4
Nobel Biocare NobelActive / F-Series / Diameter 3.0, 3.5, 4.3, 5.0 / Platform 3.0, 3.5, 3.9 (4.3), 3.9 (5.0)
Neodent - Grand Morse / GM Series / Diameters 3.5, 3.75, 4.0, 4.3, 5.0, 6.0 / Platform 3.0
Biomet 3i - Certain / H-Series / Diameter 3.25, 4.0, 5.0 / Platform 3.4, 4.1, 5.0
Biomet 3i - External Hex / I-Series / Diameter 3.25, 3.75, 4.0, 5.0 / Platform 3.4, 4.1, 5.0
Nobel Biocare - Brånemark System / K-Series / Diameter 3.3, 3.75, 4.0, 5.0 / Platform 3.5, 4.1, 4.1, 5.1
Zimmer Dental Tapered Screw-vent / R-Series / Diameter 3.3, 3.7, 4.1, 4.7, 6.0 / Platform 3.5, 4.5, 5.7
Dentsply Implants - ASTRA TECH OsseoSpeed TX / S-Series / Diameter 3.5, 4.0, 4.5, 5.0 / Platform 3.5, 4.0, 4.5, 5.0
Dentsply Implants - XiVE S / T-Series / Diameter 3.4, 3.8, 4.5, 5.5 / Platform 3.4, 3.8, 4.5, 5.5

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug
Administration
Office of Chief Information Officer Paperwork Reduction
Act (PRA) Staff PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K201711 – Traditional 510(k)
CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

5 510(k) Summary

5.1 Submitter's Contact Information

Submitter: Straumann USA, LLC
60 Minuteman Road
Andover, MA 01810
Registration No.: 1222315 Owner/Operator No.: 9005052

On the behalf of:

Medentika GmbH
Hammweg 8-10
76549 Hügelsheim, Germany
Registration No.: 3008770646 Owner/Operator No.: 10034561

Contact Person: Jennifer M. Jackson, MS
Director of Regulatory Affairs and Quality
Phone Number: +1-978-747-2509
Email: jennifer.jackson@straumann.com

Prepared By & Alternate Contact: Olivier Russo
Regulatory Affairs and Compliance Manager
Institut Straumann AG
Phone number: +41 61 965 1260

Date Prepared: December 8, 2020

5.2 Name of the Device

Trade Names: CARES® Screw - Retained Bars and Bridges (SRBB)

Common Name: CARES® Screw - Retained Bars and Bridges

Classification Name: Abutment, Implant, Dental, Endosseous

Regulation Number: 21 CFR 872.3630

Device Classification: II

Product Code(s): NHA

Classification Panel: Dental

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

5.3 Predicate Device(s)

Primary Predicate:

- K190097 – Straumann CARES® Screw-Retained Bridges and Bars

Reference Devices:

- K170838 – Medentika CAD/CAM TiBases
- K180536 – Neodent GM (Grand Morse) implant Line
- K120414 – OsseoSpeed™ Plus
- K180564 – MRI Safety Information Labeling Change
- K150203 – Medentika CAD/CAM Abutments

5.4 Device Description

The CARES® Screw Retain Bars and Bridges, referenced in Table 1, are used for the restoration of different dental implants systems with different endosteal diameters, lengths and platforms. The bars and bridges presented in the premarket notification submission (identified as “SRBB” for Screw Retained Bridges and Bars) are designed to interface with different implant connections. They allow for individual customization regarding function and esthetics. They attach directly to dental implants. The devices are intended to be finished into a bridge or overdenture using standard dental laboratory techniques and materials. Screw Retained Bars and Bridges (SRBB) devices facilitate customization to meet the functional and esthetic requirements of the individual patient. They are patient-specific medical devices, i.e. they are designed by a dental professional (clinician or dental technician) and fabricated by Medentika specifically for an individual patient. SRBB devices are designed via Computer Aided Design (CAD). After importing a scan of the patient model, Commercial Off-The-Shelf (COTS) Software includes the ability to generate digital restoration models incorporating the subject devices. The digital restoration model is transferred to the milling center where the restoration is produced using Computer Aided Manufacturing (CAM)-techniques. The devices are made entirely of titanium Grade 5 according to ASTM F136, *Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications*.

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

Reference	Product Description	BRIDGE / BAR SYSTEM
CR 12000	CARES® s-ret bridge, TAV, 2 elements	BRIDGE
CR 12001	CARES® s-ret bridge, TAV, 3 elements	
CR 12002	CARES® s-ret bridge, TAV, 4 elements	
CR 12003	CARES® s-ret bridge, TAV, 5 elements	
CR 12004	CARES® s-ret bridge, TAV, 6 elements	
CR 12005	CARES® s-ret bridge, TAV, 7 elements	
CR 12006	CARES® s-ret bridge, TAV, 8 elements	
CR 12007	CARES® s-ret bridge, TAV, 9 elements	
CR 12008	CARES® s-ret bridge, TAV, 10 elements	
CR 12009	CARES® s-ret bridge, TAV, 11 elements	
CR 12010	CARES® s-ret bridge, TAV, 12 elements	
CR 12011	CARES® s-ret bridge, TAV, 13 elements	
CR 12012	CARES® s-ret bridge, TAV, 14 elements	
CR 12013	CARES® s-ret bridge, TAV, 15 elements	
CR 12014	CARES® s-ret bridge, TAV, 16 elements	
CR 12015	CARES® Basic Bar, TAV, 2 Impl.	BASIC BAR
CR 12016	CARES® Basic Bar, TAV, 3 Impl.	
CR 12017	CARES® Basic Bar, TAV, 4 Impl.	
CR 12018	CARES® Basic Bar, TAV, 5 Impl.	
CR 12019	CARES® Basic Bar, TAV, 6 Impl.	
CR 12020	CARES® Basic Bar, TAV, 7 Impl.	
CR 12021	CARES® Basic Bar, TAV, 8 Impl.	
CR 12022	CARES® Basic Bar, TAV, 9 Impl.	
CR 12023	CARES® Basic Bar, TAV, 10 Impl.	
CR 12024	CARES® Advanced Bar, TAV, 2 Impl.	ADVANCED BAR
CR 12025	CARES® Advanced Bar, TAV, 3 Impl.	
CR 12026	CARES® Advanced Bar, TAV, 4 Impl.	
CR 12027	CARES® Advanced Bar, TAV, 5 Impl.	
CR 12028	CARES® Advanced Bar, TAV, 6 Impl.	
CR 12029	CARES® Advanced Bar, TAV, 7 Impl.	
CR 12030	CARES® Advanced Bar, TAV, 8 Impl.	
CR 12031	CARES® Advanced Bar, TAV, 9 Impl.	
CR 12032	CARES® Advanced Bar, TAV, 10 Impl.	
CR 12033	CARES® Complex Bar, TAV, 2 Impl.	COMPLEX BAR
CR 12034	CARES® Complex Bar, TAV, 3 Impl.	
CR 12035	CARES® Complex Bar, TAV, 4 Impl.	
CR 12036	CARES® Complex Bar, TAV, 5 Impl.	
CR 12037	CARES® Complex Bar, TAV, 6 Impl.	
CR 12038	CARES® Complex Bar, TAV, 7 Impl.	
CR 12039	CARES® Complex Bar, TAV, 8 Impl.	
CR 12040	CARES® Complex Bar, TAV, 9 Impl.	
CR 12041	CARES® Complex Bar, TAV, 10 Impl.	

Table 1 – Listing of devices

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

5.5 Intended Use

CARES® Screw-retained Bars and Bridges are prosthetic components directly connected to the endosseous dental implant as an aid in prosthetic rehabilitations.

5.6 Indications for Use

CARES® Screw-retained Bridges and Bars are indicated for use as bars and bridges that attach to implants to provide support for prosthetic reconstructions such as bridges and overdentures. The final processed products have the purpose of restoring chewing function. CARES® Screw-retained Bridges and Bars are indicated for screw-retained restorations.

CARES® Screw - Retained Bars and Bridges are designed to interface with the following dental implant systems (Implant System Compatibility / Series / Implant diameter (mm) / Platform diameter (mm)):

Nobel Biocare Replace Select / E-Series / Diameters 3.5, 4.3, 5.0, 6.0 / Platform 3.5, 4.3, 5.0, 6.0

Dentsply Implants - ASTRA TECH OsseoSpeed EV / EV-Series / 3.0, 3.6, 4.2, 4.8, 5.4 / Platform 3.0, 3.6, 4.2, 4.8, 5.4

Nobel Biocare NobelActive / F-Series / Diameter 3.0, 3.5, 4.3, 5.0 / Platform 3.0, 3.5, 3.9 (4.3), 3.9 (5.0)

Neodent - Grand Morse / GM Series / Diameters 3.5, 3.75, 4.0, 4.3, 5.0, 6.0 / Platform 3.0

Biomet 3i - Certain / H-Series / Diameter 3.25, 4.0, 5.0 / Platform 3.4, 4.1, 5.0

Biomet 3i – External Hex / I-Series / Diameter 3.25, 3.75, 4.0, 5.0 / Platform 3.4, 4.1, 5.0

Nobel Biocare - Brånemark System / K-Series / Diameter 3.3, 3.75, 4.0, 5.0 / Platform 3.5, 4.1, 4.1, 5.1

Zimmer Dental Tapered Screw-vent / R-Series / Diameter 3.3, 3.7, 4.1, 4.7, 6.0 / Platform 3.5, 4.5, 5.7

Dentsply Implants - ASTRA TECH OsseoSpeed TX / S-Series / Diameter 3.5, 4.0, 4.5, 5.0 / Platform 3.5, 4.0, 4.5, 5.0

Dentsply Implants – XiVE S / T- Series / Diameter 3.4, 3.8, 4.5, 5.5 / Platform 3.4, 3.8, 4.5, 5.5

5.7 Technological Characteristics

The technological characteristics of the subject devices are compared to the primary predicate and reference devices in Table 2. The Indications for Use Statement for the subject devices has been modified from the primary predicate (K190097) to include the implant-to-abutment compatibilities.

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

The implant-to-abutment compatibilities (E-, EV-, F-, H-, I-, K-, R-, S-, and T-Series) are identical to the reference device, excluding the 3.0 mm implant diameter/3.0 mm platform diameter DENTSPLY Implants EV Series (K120414). Reverse Engineering Analysis for the addition of this new compatibility was conducted using OEM implant bodies, OEM abutments, and OEM abutment screws. Reverse engineering was not needed for the addition of the implant-to-abutment compatibilities for the GM-Series, because the implant-to-abutment connection dimensions were shared as part of a business partnership with the OEM company. The device materials of the subject device are identical to the reference device K170838.

Please note that the MRI compatibility reference device K180564 – MRI Safety Information Labeling Change is included for the MRI only. K150203 is included for reference to dynamic fatigue testing (this is the primary predicate to K170838).

The design step defined in the subject device allows to control the same design limits implemented in the primary predicate (K190097). The design of the CARES® Screw - Retained Bars and Bridges Framework is performed by a trained technician using a Commercial Off-The-Shelf (COTS) Software. The design parameters are identical between the subject devices and the primary predicate (K190097) and are respected during the design of the restoration. Once the design is finalized, the design file is digitally transferred to a Straumann milling center holding the appropriate Establishment Registration. The prosthetic restoration is then manufactured according 21 CFR 820 Quality System Regulations. After production, all CARES® Screw - Retained Bars and Bridges undergoes a rigid Quality Control. If the result is correct, production sends the device to the customer. The manufacturing process for the subject devices is identical to the primary predicate (K190097).

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

Feature	Subject Device	Predicate Device	Reference Device
	<i>CARES® Screw - Retained Bars and Bridges</i>	<i>K190097 - Straumann CARES® Screw-Retained Bridges and Bars</i>	<i>K170838 – Medentika CAD/CAM TiBases</i>
Indications for use	<p>CARES® Screw-retained Bridges and Bars are indicated for use as bars and bridges that attach to implants to provide support for prosthetic reconstructions such as bridges and overdentures. The final processed products have the purpose of restoring chewing function. CARES® Screw-retained Bridges and Bars are indicated for screw-retained restorations.</p> <p>CARES® Screw - Retained Bars and Bridges are designed to interface with the following dental implant systems (Implant System Compatibility / Series / Implant diameter (mm) / Platform diameter (mm)):</p> <p>Nobel Biocare Replace Select / E-Series / Diameters 3.5, 4.3, 5.0, 6.0 / Platform 3.5, 4.3, 5.0, 6.0</p> <p>Dentsply Implants - ASTRA TECH OsseoSpeed EV / EV-Series / 3.0, 3.6, 4.2, 4.8, 5.4 / Platform 3.0, 3.6, 4.2, 4.8, 5.4</p> <p>Nobel Biocare NobelActive / F-Series / Diameter 3.0, 3.5, 4.3, 5.0 / Platform 3.0, 3.5, 3.9 (4.3), 3.9 (5.0)</p> <p>Neodent - Grand Morse / GM Series / Diameters 3.5, 3.75, 4.0, 4.3, 5.0, 6.0 / Platform 3.0</p> <p>Biomet 3i - Certain / H-Series / Diameter 3.25, 4.0, 5.0 / Platform 3.4, 4.1, 5.0</p> <p>Biomet 3i – External Hex / I-Series / Diameter 3.25, 3.75, 4.0, 5.0 / Platform 3.4, 4.1, 5.0</p> <p>Nobel Biocare - Brånemark System / K-Series / Diameter 3.3, 3.75, 4.0, 5.0 / Platform 3.5, 4.1, 4.1, 5.1</p> <p>Zimmer Dental Tapered Screw-vent / R-Series / Diameter 3.3, 3.7, 4.1, 4.7, 6.0 / Platform 3.5, 4.5, 5.7</p> <p>Dentsply Implants - ASTRA TECH OsseoSpeed TX / S-Series / Diameter 3.5, 4.0, 4.5, 5.0 / Platform 3.5, 4.0, 4.5, 5.0</p> <p>Dentsply Implants – XiVE S / T- Series / Diameter 3.4, 3.8, 4.5, 5.5 / Platform 3.4, 3.8, 4.5, 5.5</p>	<p>Straumann® CARES® Screw-retained Bridges and Bars are indicated for use as bars and bridges that attach to implants to provide support for prosthetic reconstructions such as bridges and overdentures. The final processed products have the purpose of restoring chewing function. Straumann® CARES® Screw-retained Bridges and Bars are indicated for Screw-retained restorations. Straumann® CARES® Screw-retained Bridges and Bars are designed to interface with the Bone Level (BL), Tissue Level (TL), and BLX implants of the Straumann Dental Implant System (SDIS).</p>	<p>Medentika TiBase CAD/CAM Abutments are intended for use with dental implants as a support for single or multiple tooth prostheses in the maxilla or mandible of a partially or fully edentulous patient.</p> <p>(Implant System Compatibility/ Series / Implant diameter (mm) / Platform Diameter (mm)):</p> <p>Nobel Biocare Replace™ Select / E-Series / 3.5, 4.3, 5.0, 6.0 / 3.5, 4.3, 5.0, 6.0</p> <p>Dentsply Implants/ASTRA TECH OsseoSpeed EV/ EV-Series / 3.6, 4.2, 4.8, 5.4 / 3.6, 4.2, 4.8, 5.4</p> <p>Nobel Biocare NobelActive™ / F-Series / 3.5, 4.3, 5.0 / 3.5, 3.9 (4.3), 3.9 (5.0)</p> <p>Biomet 3i Osseotite® Certain® / H-Series / 3.25, 4.0, 5.0 / 3.4, 4.1, 5.0</p> <p>Biomet 3i Osseotite® / I-Series / 3.25, 3.75, 4.0, 5.0 / 3.4, 4.1, 5.0</p> <p>Nobel Biocare Brånemark / K-Series / 3.3, 3.75, 4.0, 5.0 / 3.5, 4.1, 4.1, 5.1</p> <p>Straumann Bone Level / L-Series / 3.3, 4.1, 4.8 / 3.3, 4.1, 4.8</p> <p>Straumann Standard / N-Series / 3.3, 4.1, 4.8 / 3.5 (NNC), 4.8, 6.5</p> <p>Zimmer Tapered Screw-Vent® / R-Series / 3.3, 3.7, 4.1, 4.7, 6.0 / 3.5, 4.5, 5.7</p> <p>Astra Tech OsseoSpeed™ / S-Series / 3.5, 4.0, 4.5, 5.0 / 3.5, 4.0, 4.5, 5.0</p> <p>Dentsply Friadent® Frialit/XiVE® / T-Series / 3.4, 3.8, 4.5, 5.5 / 3.4, 3.8, 4.5, 5.5</p> <p>Dentsply Friadent® Ankylos® / Y-Series / 3.5, 4.5, 5.5, 7.0 / 3.5, 4.5, 5.5, 7.0</p> <p>Medentika TiBase is intended for use with the Straumann® CARES® System. All digitally designed copings and/or crowns are intended to be sent to Straumann for manufacture at a validated milling center.</p>

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

Feature	Subject Device	Predicate Device	Reference Device
	<i>CARES® Screw - Retained Bars and Bridges</i>	<i>K190097 - Straumann CARES® Screw-Retained Bridges and Bars</i>	<i>K170838 – Medentika CAD/CAM TiBases</i>
Material of the restoration	Ti6Al4V, Titanium grade 5, conforming ASTM F 136	<u>Restorations:</u> Cobalt Chromium Alloy (CoCr) Titanium Grade 4	Ti6Al4V, Titanium grade 5, conforming ASTM F 136
Material of screws	Ti6Al4V, Titanium grade 5, conforming ASTM F 136	<u>Screws:</u> Titanium-Aluminum-Niobium alloy (Ti-6Al-7Nb)	Ti6Al4V, Titanium grade 5, conforming ASTM F 136
Interface	Non engaging	Non engaging	Engaging
Restoration Types Supported	Bridges from 2 units to 16 units (full-arch) Bars from 2 units to 10 units	Bridges from 2 units to 16 units (full-arch) Bars from 2 units to 10 units	Single restorations
Design Workflow	CAD- Commercial Off-The-Shelf (COTS) Software	CAD - CARES Visual Software	Straumann CARES System
Manufacturing Workflow	Straumann Milling Center	Straumann Milling Center	Straumann Milling Center
Design Software	Manual design by a trained dental technician under the Straumann Quality System using Commercial Off-The-Shelf (COTS) Software and based on defined design limit by work instruction and training	Straumann CARES Visual Software	Straumann CARES Visual Software
Design Limits for Patient Specific Component	Critical geometry parameters are controlled during CAD design step based on defined design limit by work instruction and training. After manufacturing the devices undergo rigid quality controls.	Critical geometry parameters are enforced by CARES Visual Software limits	Critical geometry parameters are enforced by CARES Visual Software limits
Sterility	Delivered non-sterile	Delivered non-sterile	Delivered non-sterile
Sterilization by end user	Yes	Yes	Yes
Method of sterilization	Moist heat (steam autoclave)	Moist heat (steam autoclave)	Moist heat (steam autoclave)

K201711 – Traditional 510(k)
CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

Feature	Subject Device	Predicate Device	Reference Device
	<i>CARES® Screw - Retained Bars and Bridges</i>	<i>K190097 - Straumann CARES® Screw-Retained Bridges and Bars</i>	<i>K170838 – Medentika CAD/CAM TiBases</i>
Mode of Action	Screw-retained	Screw-retained	Screw-retained
Reusable	No	No	No
Packaging	Box package for the milled restoration and the screws	Box package for the restoration and the screws	The abutment in reclosable poly bag and the corresponding screws are individually packed in a plastic bag, co-packaged in a craft board box

Table 2 – Technological Characteristics

K201711 – Traditional 510(k)

CARES® Screw - Retained Bars and Bridges (SRBB)

510(k) Summary

5.8 Performance Testing

The implant to abutment connections for the E-, EV-, F-, H-, I-, K-, R-, S-, and T-Series are identical to the implant to abutment compatibilities of the reference devices (K170838), except for the addition of the 3.0 mm implant diameter/3.0 mm platform diameter DENTSPLY Implants EV Series (K120414), for which Reverse Engineering Analysis was conducted using OEM implant bodies, OEM abutments, and OEM abutment screws. For the addition of the new line of abutments (GM-Series), the implant to abutment connection dimensions were shared as part of a business partnership.

The change in design workflow from the CARES Visual Software to COTS software does not require additional testing to support substantial equivalence as the design parameters are identical between the subject devices and the primary predicate (K190097) and are respected during the design of the restoration by a trained dental technician under the Straumann Quality System to ensure the pre-defined design limits have been followed.

Dynamic fatigue testing was conducted according to the FDA guidance document “*Guidance for Industry and FDA Staff – Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Abutments*” and demonstrated that the CARES® Screw - Retained Bars and Bridges are equivalent to the predicate and reference devices.

The subject device materials and manufacturing process are identical to the primary predicate (K190097) and reference (K170838) devices, therefore, no new issues regarding biocompatibility were raised.

The sterilization process for the CARES® Screw - Retained Bars and Bridges, as recommended in the labeling, was validated according to applicable recommendations in the FDA guidance document “*Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling, issued on March 17, 2015*”, and conforming to the ISO 17665-1 and ISO/TS 17665-2.

5.9 Conclusion

The documentation submitted in this premarket notification demonstrates the CARES® Screw - Retained Bars and Bridges are substantially equivalent to the primary predicate and reference devices.