



March 24, 2022

Rizhao Sanqi Medical & Health Articles Co., Ltd.  
% Ivy Wang  
Technical Manager  
Shanghai Sungo Management Consulting Co., Ltd.  
Room 1309, Dongfang Building, 1500# Century Ave.,  
Shanghai, Shanghai 200122  
China

Re: K202903  
Trade/Device Name: Surgical Face Mask  
Regulation Number: 21 CFR 878.4040  
Regulation Name: Surgical Apparel  
Regulatory Class: Class II  
Product Code: FXX  
Dated: December 24, 2021  
Received: December 29, 2021

Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K202903

Device Name  
Surgical Face Mask

### Indications for Use (Describe)

The Surgical Face Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. When worn properly, these face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.

Level 1 Surgical Face Mask models: SQ-1001

Level 2 Surgical Face Mask models: SQ-2001, SQ-2001H, SQ-2004, SQ-2004H

Level 3 Surgical Face Mask models: SQ-3001, SQ-3001H

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(K) Summary

**K202903**

Date of Summary prepared: 2021-12-24

### A. Applicant:

Rizhao Sanqi Medical & Health Articles Co., Ltd.

Address: Heshan Industry park, Donggang District, Rizhao, Shandong Province, China

Contact person: Julia Yu

Title: Sales Manager

Tel: +86-633-7663737

Fax: +86-633-7653718

Email: 3qmedical@3qcn.cc

Submission Correspondent:

Primary contact: Ms. Ivy Wang

Shanghai SUNGO Management Consulting Co., Ltd.

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: [haiyu.wang@sungoglobal.com](mailto:haiyu.wang@sungoglobal.com)

Secondary contact: Mr. Raymond Luo

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: [fda.sungo@gmail.com](mailto:fda.sungo@gmail.com)

### B. Device:

Trade Name: SURGICAL FACE MASK

Common Name: SURGICAL MASK

Model:

Table 1 – Surgical face mask model numbers

Mask Style	Ear loops	Tie-on	Color
Level 1	SQ-1001	-	Blue, white, pink, green, yellow
Level 2	SQ-2001	SQ-2001H	
Level 2 with visor	SQ-2004	SQ-2004H	
Level 3	SQ-3001	SQ-3001H	

### Regulatory Information

Classification Name: Surgical Face Mask

Classification: Class II

Product code: FXX

Regulation Number: 878.4040

Review Panel: Surgical Apparel

**C. Predicate device:**

K160269	Surgical Face Masks (Ear loops and Tie-on)	SAN-M PACKAGE CO., LTD.
---------	--	-------------------------

**D. Indications for use of the device:**

The Surgical Face Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. When worn properly, these face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.

**E. Device Description:**

The Surgical Face Masks are 3-layer or 4-layer, flat-folded masks constructed of nonwoven polypropylene materials. The mask is provided with ear loops (Spandex) or ties (polypropylene). A malleable nose clip is placed in the layers of facemask for comfort and individualized fit. **The surgical face mask will be provided in blue, white, pink, green, yellow and with the option for a visor.** The surgical face masks are single-use, disposable devices, provided non-sterile.

**F. Technological Comparison with predicate device**

Table 2 General Comparison

Device	Subject Device			Predicate Device			Result
<b>510K #</b>	K202903			K160269			
<b>Manufacturer</b>	Rizhao Sanqi Medical & Health Articles Co., Ltd.			SAN-M PACKAGE CO., LTD.			
<b>Model Name</b>	SURGICAL FACE MASK			Surgical Face Masks (Ear loops and Tie-on)			Similar
<b>Classification</b>	Class II Device, FXX (21 CFR878.4040)			Class II Device, FXX (21 CFR878.4040)			Same
<b>Intend use</b>	The Surgical Face Mask is intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. When worn properly, these face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile. This is a single use, disposable device(s), provided non-sterile.			The surgical face masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids, and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single-use, disposable device, provided non-sterile.			Same
<b>Materials</b>	Level 1	Level 2	Level 3	Level 1	Level 2	Level 3	-
<b>Outer layer</b>	Spunbond Polypropylene			Polypropylene			Same
<b>Inner layer</b>	Spunbond Polypropylene			Polypropylene			Same
<b>Filter layer</b>	Melt-blown	1.Melt-blown		1. Polypropylene spunbond			Similar

	Polypropylene	Polypropylene 2.Microporous Film	2. Polypropylene meltblown	
<b>Nose wire</b>	Steel coated by polypropylene		Polyethylene coated steel wire	Different
<b>Ear loops</b>	Spandex		Polyester, polyurethane	Different
<b>Tie-on</b>	Spunbond Polypropylene		Polypropylene spunbond or polyester spunbond	Similar
<b>Design Features</b>	Ear Loops, Tie-on		Ear Loops, Tie-on	Same
<b>Offer with visor</b>	Yes		Yes	Same
<b>Mask style</b>	3 Flat Pleated	4 Flat Pleated	4 Flat Pleated	Similar
<b>Color</b>	Blue, white, pink, green, yellow		white or blue	Different
<b>Dimension (Length)</b>	175±5mm		175±5mm	Same
<b>Dimension (Width)</b>	95±5mm		95±5mm	Same
<b>OTC use</b>	Yes		Yes	Same
<b>Sterility</b>	Non-Sterile		Non-Sterile	Same
<b>Use</b>	Single Use, Disposable		Single Use, Disposable	Same
<b>Biocompatibility</b>	Non-cytotoxic, Non-sensitizing, non-irritating		Non-cytotoxic, Non-sensitizing, non-irritating	Same
<b>Performance Testing (ASTM 2100)</b>	Level 1	Level 2	Level 3	-
<b>Fluid Resistance</b>	Meet ASTM F1862-17		Meet ASTM F1862-13	similar
<b>Particulate Filtration Efficiency</b>	Meet ASTM F2299-17		Meet ASTM F2299-03	Similar
<b>Bacterial Filtration Efficiency</b>	Meet ASTM F2101-19		Meet ASTM F2101-14	Similar
<b>Differential Pressure</b>	Meet EN 14683: 2019, Annex C		Meet MIL-M36945C	Different
<b>Flammability</b>	Meet 16 CFR 1610		Meet 16 CFR 1610	Similar
<b>Biocompatibility</b>	Non-cytotoxic, non-sensitizing, non-irritating		Non-cytotoxic, non-sensitizing, non-irritating	Same

### G. Non-Clinical Test Conclusion

Non-clinical tests were conducted to verify that the proposed device met all design specifications as was same to the predicate device. The test results demonstrated that the proposed device complies with the following standards and the requirements stated in the Guidance for Industry and FDA Staff: *Surgical Masks – Premarket Notification [510(k)] Submission* issued on March 5, 2004:

- ISO 10993-5: 2009 Biological Evaluation of Medical Devices -- Part 5: Tests For In Vitro Cytotoxicity
- ISO 10993-10: 2010 Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin

Sensitization

- ASTM F2100, Standard Specification for Performance of Materials Used In Medical Face Masks
- ASTM F1862, Standard Test Method for Resistance of Medical Face Masks To Penetration by Synthetic Blood (Horizontal Projection of Fixed Volume At A Known Velocity);
- EN 14683, Medical Face Masks—Requirements and Test Methods;
- ASTM F2101, Standard Test Method for Evaluating the Bacterial Filtration Efficiency (BFE) Of Medical Face Mask Materials, Using A Biological Aerosol of Staphylococcus Aureus;
- ASTM F2299, Standard test method for determining the initial efficiency of materials used in medical face masks to penetration by particulates using latex spheres;
- 16 CFR 1610, Standard for the Flammability of clothing textiles;

Mask Style	Ear loops	Tie-on	Color
Level 1	SQ-1001	-	Blue, white, pink, green, yellow
Level 2	SQ-2001	SQ-2001H	
Level 2 with visor	SQ-2004	SQ-2004H	
Level 3	SQ-3001	SQ-3001H	

For level 1: SQ-1001

Item	Purpose	Acceptance Criteria	Result
		Level 1	3 nonconsecutive lots tested
<b>Fluid Resistance Performance ASTM F1862</b>	The purpose of the performance testing is to demonstrate the functionality of the subject device.	29 out of 32 pass at 80 mmHg	PASS
<b>Particulate Filtration Efficiency ASTM F2299</b>		≥ 95%	PASS
<b>Bacterial Filtration Efficiency ASTM F2101</b>		≥ 95%	PASS
<b>Differential Pressure (Delta P) EN 14683 Annex C</b>		<5.0 mmH <sub>2</sub> O/cm <sup>2</sup>	PASS
<b>Flammability 16 CFR 1610</b>		Class 1	PASS

For: Blue, white, pink, green, yellow

Item	Purpose	Acceptance Criteria	Result
<b>Cytotoxicity</b>	The purpose of the testing is to demonstrate the safety of the subject device.	Non-Cytotoxic	Under the conditions of the study, the device is non-cytotoxic.
<b>Irritation</b>		Non-Irritating	Under the conditions of the study, the device is non-irritating.
<b>Sensitization</b>		Non-Sensitizing	Under the conditions of the study, the device is non-sensitizing

For level 2: SQ-2001 and SQ-2001H

Item	Purpose	Acceptance Criteria	Result
		Level 2	3 nonconsecutive lots tested
<b>Fluid Resistance Performance ASTM F1862</b>	The purpose of the performance testing is to demonstrate the functionality of the subject device.	29 out of 32 pass at 120 mmHg	PASS
<b>Particulate Filtration Efficiency ASTM F2299</b>		$\geq 98\%$	PASS
<b>Bacterial Filtration Efficiency ASTM F2101</b>		$\geq 98\%$	PASS
<b>Differential Pressure (Delta P) EN 14683 Annex C</b>		$<6.0 \text{ mmH}_2\text{O}/\text{cm}^2$	PASS
<b>Flammability 16 CFR 1610</b>		Class 1	PASS

For: Blue, white, pink, green, yellow

Item	Purpose	Acceptance Criteria	Result
<b>Cytotoxicity</b>	The purpose of the testing is to demonstrate the safety of the subject device.	Non-Cytotoxic	Under the conditions of the study, the device is non-cytotoxic.
<b>Irritation</b>		Non-Irritating	Under the conditions of the study, the device is non-irritating.
<b>Sensitization</b>		Non-Sensitizing	Under the conditions of the study, the device is non-sensitizing

For Level 2 with visor: SQ-2004 and SQ-2004H

Item	Purpose	Acceptance Criteria	Result
		Level 2	3 nonconsecutive lots tested
<b>Fluid Resistance Performance ASTM F1862</b>	The purpose of the performance testing is to demonstrate the functionality of the subject device.	29 out of 32 pass at 120 mmHg	PASS
<b>Particulate Filtration Efficiency ASTM F2299</b>		$\geq 98\%$	PASS
<b>Bacterial Filtration Efficiency ASTM F2101</b>		$\geq 98\%$	PASS



Rizhao Sanqi Medical & Health Articles Co., Ltd.  
Heshan Industry park, Donggang District, Rizhao, Shandong Province, China

<b>Differential Pressure (Delta P) EN 14683 Annex C</b>		<6.0 mmH <sub>2</sub> O/cm <sup>2</sup>	PASS
<b>Flammability 16 CFR 1610</b>		Class 1	PASS

For: Blue, white, pink, green, yellow

Item	Purpose	Acceptance Criteria	Result
<b>Cytotoxicity</b>	The purpose of the testing is to demonstrate the safety of the subject device.	Non-Cytotoxic	Under the conditions of the study, the device is non-cytotoxic.
<b>Irritation</b>		Non-Irritating	Under the conditions of the study, the device is non-irritating.
<b>Sensitization</b>		Non-Sensitizing	Under the conditions of the study, the device is non-sensitizing

For level 3: SQ-3001 and SQ-3001H

Item	Purpose	Acceptance Criteria	Result
		Level 3	3 nonconsecutive lots tested
<b>Fluid Resistance Performance ASTM F1862</b>	The purpose of the performance testing is to demonstrate the functionality of the subject device.	29 out of 32 pass at 160 mmHg	PASS
<b>Particulate Filtration Efficiency ASTM F2299</b>		≥ 98%	PASS
<b>Bacterial Filtration Efficiency ASTM F2101</b>		≥ 98%	PASS
<b>Differential Pressure (Delta P) EN 14683 Annex C</b>		<6.0 mmH <sub>2</sub> O/cm <sup>2</sup>	PASS
<b>Flammability 16 CFR 1610</b>		Class 1	PASS

For: Blue, white, pink, green, yellow

Item	Purpose	Acceptance Criteria	Result
<b>Cytotoxicity</b>	The purpose of the testing is to demonstrate the safety of the subject device.	Non-Cytotoxic	Under the conditions of the study, the device is non-cytotoxic.
<b>Irritation</b>		Non-Irritating	Under the conditions of the study, the device is non-irritating.
<b>Sensitization</b>		Non-Sensitizing	Under the conditions of the study, the device is non-sensitizing

#### H. Clinical Test Conclusion

No clinical study is included in this submission.

## **I. Conclusion**

The conclusions drawn from the nonclinical tests demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device K160269.