



April 22, 2022

Fisher & Paykel Healthcare Ltd  
Reena Daken  
Regulatory Affairs Market Manager  
15 Maurice Paykel Place, East Tamaki  
Auckland, 2013  
New Zealand

Re: K212031

Trade/Device Name: F&P 850 AirSpiral Adult NIV and NHF Circuit Kit  
Regulation Number: 21 CFR 868.5450  
Regulation Name: Respiratory Gas Humidifier  
Regulatory Class: Class II  
Product Code: BTT  
Dated: March 22, 2022  
Received: March 22, 2022

Dear Reena Daken:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Brandon Blakely, Ph.D.  
Respiratory Devices Team  
DHT1C: Division of Sleep Disordered  
Breathing, Respiratory and  
Anesthesia Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K212031

Device Name  
F&P 850 AirSpiral Adult NIV and NHF Circuit Kit

Indications for Use (Describe)

For the delivery of heated, humidified breathing gases to spontaneously breathing adult patients. This breathing set is suitable for use with Fisher & Paykel Healthcare MR850 Humidifiers in hospital and long-term care environments.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# 510(k) Summary

As Required by 21 CFR 807.92

## I. SUBMITTER

<b>Company Name and Address</b>	Fisher & Paykel Healthcare Limited 15 Maurice Paykel Place East Tamaki Auckland 2013, New Zealand Telephone: +64 9 574 0100
<b>Prepared and Submitted by</b>	Nicholas Yap Regulatory Affairs Specialist
<b>Contact Person</b>	Reena Daken Regulatory Affairs Manager Telephone: +64 9 574 0100 Email: <a href="mailto:reena.daken@fphcare.co.nz">reena.daken@fphcare.co.nz</a>
<b>Date Prepared</b>	22 April 2022

## II. DEVICE

<b>Device Name</b>	F&P 850 AirSpiral Adult NIV and NHF Circuit Kit
<b>Common/Usual Name</b>	Heated Breathing Tube
<b>Classification Name</b>	Respiratory gas humidifier
<b>Regulatory Class</b>	Class II (21 CFR §868.5450)
<b>Product Code</b>	BTT

## III. PREDICATE DEVICE

<b>Predicate Device</b>	F&P AirSpiral Heated Breathing Tube	K162553
<b>Secondary Predicate Device</b>	Bipap Vision Ventilatory Support	K982454
<b>Reference Device(s)</b>	F&P RT380 Adult Evaqua 2 Dual Heated Breathing Circuit	K122432
	<ul style="list-style-type: none"><li>Used to support claims of substantial equivalence with respect to performance and system compatibility</li></ul>	

## **DEVICE DESCRIPTION**

The F&P 850 AirSpiral Adult NIV and NHF Circuit Kit (850A61) is designed to provide a gas conduit between respiratory support equipment and a patient interface such as a mask, nasal prongs or tracheostomy interface for a spontaneously breathing Adult patient.

The Inspiratory Limb is an extruded bubble tube design 1.6 metres in length, and is intended to connect the gas path from humidification chamber to the patient interface. Additionally, temperature sensor clips and a gown clip are included to aid in cable management and to affix the circuit to a gown or sheet. The subject device is single use, prescription only and is provided in a non-sterile state.

The Disposable Exhalation Port is an optional accessory for NIV applications. It features a perforated ventilation hole and a 22mm/15mm coaxial taper connector for connecting to patient interfaces and a 22mm tapered male connector for connecting to the patient end of the Inspiratory Limb. It also features a capped port for connecting to the Pressure Line.

## **IV. INDICATIONS FOR USE**

For the delivery of heated, humidified breathing gases to spontaneously breathing adult patients. This breathing set is suitable for use with Fisher & Paykel Healthcare MR850 Humidifiers in hospital and long-term care environments.

## V. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

A comparison of intended use, features, and performance specifications demonstrate that the subject device is substantially equivalent to the predicate device.

Parameter	Subject Device F&P 850A61	Predicate Device F&P AirSpiral Heated Breathing Tube (K162553)	Comments
<b>Classification</b>			
<b>Manufacturer</b>	Fisher & Paykel Healthcare Ltd.	Fisher & Paykel Healthcare Ltd.	Identical
<b>Device Regulation</b>	Class II, Regulation: 21 CFR 868.5270	Class II, Regulation: 21 CFR 868.5270	
<b>Product Code</b>	BTT	BTT	
<b>Classification Panel</b>	Anesthesiology	Anesthesiology	
<b>Intended Use and Indications for Use</b>			
<b>Intended Use</b>	Heated breathing tube for delivery of humidified respiratory gases	Heated breathing tube for delivery of humidified respiratory gases	Identical
<b>Indications for Use</b>	For the delivery of heated, humidified breathing gases to spontaneously breathing adult patients. This breathing set is suitable for use with Fisher & Paykel Healthcare MR850 Humidifiers in hospital and long-term care environments.	Heated breathing tube for delivery of humidified respiratory gases. For use with AIRVO and AIRVO2 Series humidifiers in hospitals and long-term care facilities. For use at flows from 2 to 60 L/min depending on the patient interface.	Equivalent
<b>Availability</b>	Prescription use (Part 21 CFR 801 Subpart D)	Prescription use (Part 21 CFR 801 Subpart D)	Identical
<b>Patient Population</b>	Spontaneously breathing adult patients requiring flows between 10 – 120 L/min depending on the required therapy and patient interface.	Spontaneously breathing patient requiring flows between 2 – 60 L/min.	Equivalent

Parameter	Subject Device F&P 850A61	Predicate Device F&P AirSpiral Heated Breathing Tube (K162553)	Comments
<b>Operating Environment</b>	Hospital and long-term care environments	Hospital and long-term care facilities	Identical
<b>Reusability</b>	Single use	Single use	Identical
<b>Sterility</b>	Device not provided sterile	Device not provided sterile	Identical
<b>Life Supporting or Life Sustaining</b>	No	No	Identical
<b>Operation and Safety Features</b>			
<b>Principle of Operation</b>	Connects to a humidifier. Powered by the humidifier to maintain heat as a conduit for humidified heated respiratory gases.  Connects the humidifier to the patient interface.	Connects to a humidifier. Powered by the humidifier to maintain heat as a conduit for humidified heated respiratory gases.  Connects the humidifier to the patient interface.	Identical
<b>Tube Dimensions</b>	<b>Length:</b> 1.6 m <b>Internal Diameter:</b> 17 mm	<b>Length:</b> 1.8 m <b>Internal Diameter:</b> 13.1 mm	Equivalent
<b>Design of the Tube Construction</b>	Heated tube to transport humidified gas between the humidification chamber and the Patient End  Two-spiral wall tubing design: <ul style="list-style-type: none"> <li>Insulating spiral made of Spiral-wound bubble tubing</li> </ul> Heating spiral encapsulating heater wire	Heated tube to transport humidified gas between the Unit End and the Patient End  Two-spiral wall tubing design: <ul style="list-style-type: none"> <li>Insulating spiral made of Spiral-wound bubble tubing</li> </ul> Heating spiral encapsulating heater wire	Identical
<b>Temperature Sensor</b>	A removable temperature and flow sensor are used at the humidification chamber outlet port which is an accessory to the F&P MR850 humidifier.	An integrated temperature sensor at the patient end of the tube eliminates the need for external probes, cables or adaptors.	Similar

Parameter	Subject Device F&P 850A61	Predicate Device F&P AirSpiral Heated Breathing Tube (K162553)	Comments
	A removable temperature sensor is used at the patient end port which is an accessory to the F&P MR850 humidifier.		
<b>Heater Wire Lumen Design</b>	Consists of filament conductor. Position in tube: Double helix wound moulded into bead of tube wall. Filament insulation: Filament is moulded into the bead. The bead provides thermal insulation.	Consists of filament conductor. Position in tube: Double helix wound moulded into bead of tube wall. Filament insulation: Filament is moulded into the bead. The bead provides thermal insulation.	Identical
<b>Clips</b>	Allows the breathing tube to be attached to clothing or bedding and to position the tube to suit the clinician and patient	Allows the breathing tube to be attached to clothing or bedding and to position the tube to suit the clinician and patient	Identical
<b>Humidifier Compatibility</b>	F&P MR850 Respiratory Humidifier	F&P AIRVO / AIRVO2 Series Humidifier F&P myAIRVO / myAIRVO2 Series Humidifier	Similar
<b>Water Chamber Compatibility</b>	MR290V autofill chamber <ul style="list-style-type: none"> <li>• Single patient use only</li> </ul>	900PT290E auto fill chamber, <ul style="list-style-type: none"> <li>• Single patient use only</li> </ul> HC360 manual fill chamber, <ul style="list-style-type: none"> <li>• Single patient reuse</li> </ul>	Identical
<b>Patient Interface Compatibility</b>	<b>HF Indication:</b> <ul style="list-style-type: none"> <li>• Interface with 22 mm male connector (ISO 5356) e.g. OPT944 Nasal Cannula (K162553)</li> </ul> <b>NIV Indication:</b>	<b>HF Indication:</b> <ul style="list-style-type: none"> <li>• OPT942, OPT944, OPT946 Nasal Cannula</li> <li>• OPT970 Tracheostomy Direct Connection</li> <li>• OPT980 Mask Interface Adapter</li> </ul>	Equivalent



Parameter	Subject Device F&P 850A61	Predicate Device F&P AirSpiral Heated Breathing Tube (K162553)	Comments
	<ul style="list-style-type: none"> <li>Non-vented mask with 22 mm male connector (ISO 5356) and exhalation port e.g. RT045 Mask (K170367)</li> <li>Vented mask with 22 mm female connector (ISO 5366) e.g. RT047 Mask (K191624)</li> </ul> All are single patient use only	<ul style="list-style-type: none"> <li>OPT316, OPT318 Junior Nasal Cannula</li> </ul> All are single patient use only	
<b>Performance</b>			
<b>Flow Range</b>	10 – 120 L/min  <b>HF (MR850 Invasive Mode):</b> <ul style="list-style-type: none"> <li>10 – 60 L/min</li> </ul> <b>NIV (MR850 Mask Mode):</b> <ul style="list-style-type: none"> <li>10 – 120 L/min</li> </ul>	2 – 60 L/min  <b>Default Mode:</b> <ul style="list-style-type: none"> <li>10 – 60 L/min</li> </ul> <b>Junior Mode:</b> <ul style="list-style-type: none"> <li>2 – 5 L/min</li> </ul>	The expanded flow rate range of the subject device is supported by the reference device (K122432), which is also indicated for use up to 120 L/min when used with MR850.
<b>Shelf Life</b>	2 Years	5 Years	Both devices were preconditioned using equivalent methods to simulate shelf life and performance was verified thereafter.
<b>Useful Life</b>	Single patient use – 14 days in Hospital / Long-term care facilities	Single patient use – 14 days in Hospital / Long-term care facilities	Identical
<b>Humidity Delivery</b>	Noninvasive Mode: > 12 mg/L Invasive Mode: > 33 mg/L	Noninvasive indications: >12 mg/L <ul style="list-style-type: none"> <li>Default and Junior modes</li> </ul> Invasive indications: > 33 mg/L <ul style="list-style-type: none"> <li>Default and Junior modes</li> </ul>	Identical

Technical Characteristic	Secondary Predicate Device Philips Disposable Exhalation Port (K982454)	Subject Device F&P Disposable Exhalation Port	Comments
<b>Indications for Use</b>	This item is supplied as a component in a complete breathing system and supports the respective indications for use.	This item is supplied as a component in the circuit kit and supports the respective indications for use.	Identical
<b>Availability</b>	Prescription only	Prescription only	
<b>Patient Population</b>	Adult	Adult	
<b>Patient State</b>	Spontaneously breathing	Spontaneously breathing	
<b>Patient Monitoring</b>	Appropriate patient monitoring	Appropriate patient monitoring	
<b>Operating Environment</b>	Hospital and long-term care	Hospital and long-term care	
<b>User Group</b>	Respiratory therapists, sleep lab technicians or nurses under the supervision or direction of a physician	Under the supervision of trained medical personnel	
<b>Reusability</b>	Single Patient Use	Single Patient Use	
<b>Useful Life</b>	-	14 days	
<b>Technical Specifications and Features</b>			
<b>Connection to interface</b>	ISO 5356-1 Conical connectors (22mm male and 15mm female)	ISO 5356-1 Conical connectors (22mm male and 15mm female)	Identical
<b>Connection to circuit</b>	ISO 5356-1 Conical connectors (22mm male)	ISO 5356-1 Conical connectors (22mm male)	
<b>Connection to pressure port</b>	Compatible with 1/8" pressure line and port cap	Compatible with 1/8" pressure line and port cap	
<b>Inhalation Path Diameter (min)</b>	15.2 mm	15.2 mm	
<b>Exhalation Port functional performance</b>	>15 L/min @ 4 cmH2O	>15 L/min @ 4 cmH2O	Equivalent
<b>Pressure port functional performance</b>	3 mL air volume	3.5 mL air volume	Equivalent
<b>Inhalation Path functional performance</b>	17% Tidal volume discrepancy	13% Tidal volume discrepancy	Equivalent
<b>Materials</b>			
<b>Materials</b>	Polycarbonate	Polycarbonate	Identical

## VI. PERFORMANCE DATA

### Summary of non-clinical tests

The F&P 850A61 has been tested to applicable requirements of the following standards:

Standards and Designation Number	Standards Title
ISO 80601-2-74:2017	Medical electrical equipment – Part 2-74: Particular requirements for basic safety and essential performance of respiratory humidifying equipment
IEC 60601-1:2005+AMD1:2012	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2:2014	Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests
AIM 7351731 Rev 2.0:2017	Medical Electrical Equipment and System Electromagnetic Immunity Test for Exposure to Radio Frequency Identification Readers
ISO 5356-1:2015	Anaesthetic and respiratory equipment – Conical connectors: Part 1: Cones and sockets
ISO 5367:2014	Anaesthetic and respiratory equipment – Breathing sets and connectors
IEC 60601-1-6 ed 3.1:2013	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability
IEC 62366-1:2015	Medical devices — Part 1: Application of usability engineering to medical devices
ISO 10993-1:2018	Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process
ISO 10993-3:2014	Biological evaluation of medical devices — Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
ISO 10993-18:2020	Biological evaluation of medical devices — Part 18: Chemical characterization of medical device materials within a risk management process
ISO 18562-1:2017	Biocompatibility evaluation of breathing gas pathways in healthcare applications – Part 1: Evaluation and testing within a risk management process
ISTA 3A	Packaged Products for Parcel Delivery System Shipment 70 kg (150 lb) or Less

In addition to the standards listed above, the following testing has been completed to demonstrate substantial equivalence:

**Exhalation Port Equivalence:**

A comparison of the flow rates through the exhalation port for the subject and predicate DEP was tested. Testing showed there is minimal difference in flow rate between the predicate and subject DEP at any of the pressures tested.

**Pressure Port Equivalence:**

A comparison of the inspired volume required for the inspiration of a breath was performed using the subject and predicate DEP.

The comparative testing showed the tidal volumes of the subject DEP within 1 mL of the predicate DEP and can be considered to have an equivalent sensitivity to a breath being taken by the patient.

**Inhalation Path Equivalence:**

A comparison of the tidal volumes delivered to the patient in relation to the volumes set on the ventilator showed the configurations using the subject DEP have a smaller tidal volume difference than those using the predicate DEP, and hence allow equivalent accuracy of delivery of therapy to the patient.

**Human Factors Evaluation**

A Human Factors and Usability Engineering validation study was conducted, and documentation was provided as per FDA's Guidance for Industry and FDA Staff, "Applying Human Factors and Usability Engineering to Medical Devices".

Validation study participants included 15 experienced, US Licensed Respiratory Therapists, and the report demonstrates that the critical tasks of the 850A61 breathing set can be completed by experienced respiratory therapists.

## **VII. CONCLUSIONS**

The F&P 850A61 is substantially equivalent to the predicate device based on intended use, patient population, comparison of the technological characteristics and performance specifications. In addition, the conclusions drawn from the non-clinical tests demonstrate that the device is substantially equivalent to the legally marketed predicate device.