



February 24, 2021

Sirius Medical Systems B.V.  
Bram Schermers  
Clinical Application Specialist  
High Tech Campus 41  
Eindhoven, North Brabant 5656AE  
Netherlands

Re: K200734

Trade/Device Name: Sirius Pintuition Seed, Sirius Pintuition Probe, Sirius Pintuition Base Unit  
Regulation Number: 21 CFR 878.4300  
Regulation Name: Implantable Clip  
Regulatory Class: Class II  
Product Code: PBY  
Dated: January 18, 2021  
Received: January 25, 2021

Dear Bram Schermers:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Cindy Chowdhury, Ph.D., M.B.A.  
Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K200734

Device Name

Sirius Pintuition Localization System (Sirius Pintuition Seed and Sirius Pintuition Detector)

Indications for Use (Describe)

The Sirius Pintuition Seed is intended to be placed percutaneously in the breast to mark temporarily (< 30 days) a lumpectomy site intended for surgical removal. Using imaging guidance (such as ultrasound or radiography) or aided by non-imaging guidance (Sirius Pintuition Detector) the Sirius Pintuition Seed is located and surgically removed with the target tissue.

The Sirius Pintuition Detector is intended for the non-imaging detection and localization of the Sirius Pintuition Seed that has been implanted in a lumpectomy site intended for surgical removal.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.


This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

	<b>PREMARKET NOTIFICATION 510(k)</b>		
	<i>Meta</i>	<i>Identification</i>	<i>Version</i>
		001015	2.0
	<i>Classification</i>		<i>Page</i>
	PUBLIC		15 of 121
<i>Title</i>			
Sirius Pintuition Localization System (Sirius Pintuition Seed and Sirius Pintuition Detector)			

## 5 510(k) Summary

### 5.1 Submitter Information

**Submitter's name:** Sirius Medical Systems B.V.  
**Address:** High Tech Campus 41  
5656 AE  
Eindhoven  
The Netherlands  
**Contact Person:** Bram Schermers  
Product Owner  
**Telephone:** 0031 6 2011 6299  
**E-mail:** [bram.schermers@sirius-medical.com](mailto:bram.schermers@sirius-medical.com)  
**Date summary prepared:** March 9<sup>th</sup>, 2020

### 5.2 Device Information

**Trade name:** Sirius Pintuition Seed and Sirius Pintuition Detector  
**Common name / device:** Temporary Tissue Marker  
**Regulation description:** Implantable Clip  
**Regulation number:** 21 CFR 878.4300  
**Regulatory Class:** Class II  
**Review Panel:** General & Plastic Surgery  
**Product Code:** PBY

### 5.3 Predicate Device

Endomag Magseed and Sentimag (K163541)


### 5.4 Device Description

The Sirius Pintuition Seed and Sirius Pintuition Detector are part of the Sirius Pintuition Localization System.

The Sirius Pintuition Seed is a small (1.65 x 5mm) Titanium tissue marker that is intended to be placed percutaneously in the breast for temporary (<30 days) marking of a lumpectomy site intended for surgical removal. The device is supplied single-use, sterile and pre-loaded within its delivery needle (12cm or 20cm length).

The Pintuition Detector is designed to detect the presence and proximity of the implanted Pintuition Seed. It consists of a mains-powered, table-top Pintuition Base Unit, and a cable-connected, reusable Pintuition Probe. Using the Pintuition Probe, a user may use the Pintuition Detector prior to and during breast surgery to plan the surgical approach and guide surgery. The location of the seed is fed back to the user using audible and visual cues (distance in mm).

The principle of operation is magnetism, the Pintuition Seed is associated with a magnetic field which the Pintuition Detector utilizes to determine the location of the Pintuition Seed.

	<b>PREMARKET NOTIFICATION 510(k)</b>		
	<i>Meta</i>	<i>Identification</i>	<i>Version</i>
		001015	2.0
	<i>Classification</i>		<i>Page</i>
	PUBLIC		16 of 121
<i>Title</i>			
Sirius Pintuition Localization System (Sirius Pintuition Seed and Sirius Pintuition Detector)			


## 5.5 Intended Use

The Sirius Pintuition Seed is intended to be placed percutaneously in the breast to mark temporarily (< 30 days) a lumpectomy site intended for surgical removal. Using imaging guidance (such as ultrasound or radiography) or aided by non-imaging guidance (Sirius Pintuition Detector) the Sirius Pintuition Seed is located and surgically removed with the target tissue.


The Sirius Pintuition Detector is intended for the non-imaging detection and localization of the Sirius Pintuition Seed that has been implanted in a lumpectomy site intended for surgical removal.

## 5.6 Summary of Technological Characteristics

Elements of Comparison	Sirius Pintuition System	Predicate: Magseed Magnetic Marker System	Comparison
<b>SYSTEM</b>			
<b>510(k) ID</b>	[TBD]	K163541	N/A
<b>Regulation Number</b>	§878.4300	§878.4300	Same
<b>Regulation Description</b>	Implantable Clip	Implantable Clip	Same
<b>Regulatory Class</b>	Class II	Class II	Same
<b>Product Code</b>	PBY	PBY	Same
<b>Intended use</b>	Temporary marking of a breast lumpectomy site for surgical removal	Temporary marking of a breast lumpectomy site for surgical removal	Same
<b>Indications for use</b>	<p>The Sirius Pintuition Seed is intended to be placed percutaneously in the breast to mark temporarily (&lt; 30 days) a lumpectomy site intended for surgical removal. Using imaging guidance (such as ultrasound or radiography) or aided by non-imaging guidance (Sirius Pintuition Detector) the Sirius Pintuition Seed is located and surgically removed with the target tissue.</p> <p>The Sirius Pintuition Detector is intended for the non-imaging detection and localization of the Sirius Pintuition Seed that has been implanted in a lumpectomy site intended for surgical removal.</p>	<p>The Endomag Magseed Magnetic Marker is intended to be placed percutaneously in the breast to mark temporarily (&lt; 30 days) a lumpectomy site intended for surgical removal. Using imaging guidance (such as ultrasound or radiography) or aided by non-imaging guidance (Endomag Sentimag® System) the Endomag Magseed Magnetic Marker is located and surgically removed with the target tissue.</p> <p>The Endomag Sentimag® System is intended for the non-imaging detection and localization of the "Endomag Magseed Magnetic Marker" that has been implanted in a lumpectomy site intended for surgical removal</p>	Same
<b>Type of Use</b>	Prescription Use	Prescription Use	Same
<b>Anatomical Locations</b>	Breast	Breast	Same

	<b>PREMARKET NOTIFICATION 510(k)</b>		
	<i>Meta</i>	<i>Identification</i>	<i>Version</i>
		001015	2.0
	<i>Classification</i>		<i>Page</i>
	PUBLIC		17 of 121
<i>Title</i>			
Sirius Pintuition Localization System (Sirius Pintuition Seed and Sirius Pintuition Detector)			

Elements of Comparison	Sirius Pintuition System	Predicate: Magseed Magnetic Marker System	Comparison
<b>Technological Characteristics</b>	<p>The Sirius Pintuition System utilizes non-imaging technologies that are comprised of a console that incorporates electronics and a simple user interface, plus a probe handpiece.</p> <p>A location marker (Pintuition Seed) is placed percutaneously in situ at the clinical target site by a delivery system and then the detector handpiece is used for the intraoperative detection and localization of the implanted marker.</p> <p>The handpiece is connected by a flexible cable to a console unit that provides the user with a visual indication of the presence and proximity of the marker.</p>	<p>The Magseed Magnetic Marker System utilizes non-imaging technologies that are comprised of a console that incorporates electronics and a simple user interface, plus a probe handpiece.</p> <p>A location marker (Magseed) is placed percutaneously in situ at the clinical target site by a delivery system and then the detector handpiece is used for the intraoperative detection and localization of the implanted marker.</p> <p>The handpiece is connected by a flexible cable to a console unit that provides the user with a visual indication of the presence and proximity of the marker.</p>	Same
<b>DETECTOR</b>			
<b>Probe type</b>	Handheld, flexible, cord-connected, reusable	Handheld, flexible, cord-connected, reusable	Same
<b>Probe tissue contacting material</b>	Poly Ether Ether Ketone (PEEK)	Poly Ether Ether Ketone (PEEK)	Same
<b>User Feedback</b>	Real-time visual and audible	Real-time visual and audible	Same
<b>Sensing Depth</b>	0-50 mm	0-30 mm	Similar. Sirius Pintuition System has improved sensing depth.
<b>SEED/MARKER</b>			
<b>Seed/Marker Materials</b>	Commercially Pure Titanium Grade II (Tissue-contacting) Neodymium magnet (Internal)	Surgical Grade Stainless Steel	Similar. Both biocompatible exterior. Internal material different due to different use of magnetism.
<b>Seed/Marker diameter</b>	1.65mm	1.0mm	Similar. Pintuition Seed is slightly larger to enable greater sensing depth.
<b>Seed/Marker length</b>	5.20mm	5.0mm	Similar
<b>Sterility</b>	Ethylene Oxide	Ethylene Oxide	Same
<b>Visibility</b>	X-ray, Ultrasound	X-ray, Ultrasound	Same
<b>DELIVERY DEVICE</b>			
<b>Type</b>	Preloaded, single-use, needle implanter	Preloaded, single-use, needle implanter	Same
<b>Material</b>	304 Stainless Steel	304 Stainless Steel	Same

	<b>PREMARKET NOTIFICATION 510(k)</b>		
	<i>Meta</i>	<i>Identification</i>	<i>Version</i>
		001015	2.0
	<i>Classification</i>		<i>Page</i>
	PUBLIC		18 of 121
<i>Title</i>			
Sirius Pintuition Localization System (Sirius Pintuition Seed and Sirius Pintuition Detector)			

Elements of Comparison	Sirius Pintuition System	Predicate: Magseed Magnetic Marker System	Comparison
<b>Delivery device gauge</b>	14G	18G	Similar. Pintuition gauge is slightly larger to enable implantation of larger Seed.

## 5.7 Summary of Non-Clinical Performance Data

Testing was conducted to evaluate and characterize the safety and performance of the Sirius Pintuition Localization System. Pre-clinical testing included:

- Design verification
- System accuracy and range verification
- Biocompatibility evaluation
- MR safety testing
- Sterilization validation
- Packaging validation
- Shelf life validation
- Electrical safety testing

## 5.8 Summary of Clinical Performance Data

An analysis of available data was conducted to evaluate and characterize the clinical safety and performance of the Sirius Pintuition Localization System. The clinical data support the safety and effectiveness of the device:

- Clinical Evaluation, including clinical safety and performance data with the actual device, a previous version of the device and an extensive evaluation of available literature data pertaining to the previous version device, the predicate device and additional benchmark devices.

## 5.9 Conclusion

The Sirius Pintuition Localization System has the same Intended Use as the predicate device. The different technological characteristics have not led to additional questions of safety or effectiveness. The data presented in the remainder of this submission show substantial equivalence of the Sirius Pintuition Localization System with the predicate device.