



February 3, 2021

Guilin HBM Sanitary Protections, Inc
% Shelley Li
Director
Shanghai Landlink Medical Information Technology Co., Ltd.
Room 703, 705, Baohua International Plaza,
West Guangzhong Road 555, Jingan
Shanghai, 200071
China

Re: K202627

Trade/Device Name: Medical Face Masks
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical Apparel
Regulatory Class: Class II
Product Code: FXX
Dated: November 9, 2020
Received: November 13, 2020

Dear Shelley Li:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For CAPT Elizabeth Claverie, M.S.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K202627

Device Name
Medical Face Masks

Indications for Use (Describe)

The Medical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary: K202627

I. Submitter

Guilin HBM Sanitary Protections, Inc
No.1-2, Shuijing East Road, Economic and Technological Development Area, Guilin
541805, China
Contact person: Liu Xin
Position: General Manager Assistant
Tel.: +86-13761713489
E-mail: 724095871@qq.com

Preparation date: Jan.30, 2021

Submission Correspondent

Ms. Shelley Li
Tel.: +86-021-80317637
Shanghai Landlink Medical Information Technology Co., Ltd.
E-mail: shelley.li@landlink-healthcare.com

II. Proposed Device

| | |
|-----------------------|--------------------|
| Trade Name of Device: | Medical Face Masks |
| Common name: | Surgical Mask |
| Regulation Number: | 21 CFR 878.4040 |
| Regulatory Class: | Class II |
| Product code: | FXX |
| Review Panel | General Hospital |

III. Predicate Devices

| | |
|-----------------|--|
| 510(k) Number: | K153496 |
| Trade name: | Disposable Surgical Face Mask |
| Common name: | Surgical Mask |
| Classification: | Class II |
| Product Code: | FXX |
| Manufacturer | Xiantao Rayxin Medical Product Co., Ltd. |

IV. Device description

The Medical Face Masks are Flat Pleated type mask, utilizing Ear Loops way for wearing, and they all have Nose Piece design for fitting the face mask around the nose. The Medical Face Masks are manufactured with three layers. The outer layer is made of spun-bonded polypropylene (PP) non-woven fabric. The middle layer with filtration function is made of melt blown polypropylene (PP) non-woven fabric. The inner layer contact with face is made of spun-bonded polypropylene (PP) non-woven fabric.

The Medical Face Masks are single use, disposable device, provided non-sterile.

V. Indication for use

The Medical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile.

VI. Comparison of technological characteristics with the predicate devices

Table 1 Comparison of characteristics

| Item | Proposed device (K202627) | Predicate device (K153496) | Comparison |
|--------------------|--|--|-------------------|
| Product name | Medial Face Masks | Disposable Surgical Face Mask | Similar |
| Product Code | FXX | FXX | Same |
| Regulation No. | 21 CFR 878.4040 | 21 CFR 878.4040 | Same |
| Class | Class II | Class II | Same |
| Mask style | Flat-pleated, ear loop, 3 layers | Flat-pleated, ear loop, tie-on, 3 Layers | Similar |
| Indication for use | The Medical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. | The Disposable Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and | Same |

| | | | | |
|---|--------------|---|---|------------------------|
| | | These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile. | particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluids. This is a single use, disposable device(s), provided non-sterile. | |
| Materials | Inner layer | Spun-bond polypropylene | Spun-bond polypropylene | Same |
| | Middle layer | Melt blown polypropylene filter | Melt blown polypropylene filter | Same |
| | Outer layer | Spun-bond polypropylene | Spun-bond polypropylene | Same |
| | Ear loop | 85% Spandex + 15% Polyester | Polyester | Different ¹ |
| | Nose piece | Iron wire covered polypropylene | Malleable aluminum wire | Different ¹ |
| Color | | Blue | Blue | Same |
| Length | | 17.5cm | 17.5±1cm | Same |
| Width | | 9.5cm | 9.5±1cm | Same |
| OTC use | | Yes | Yes | Same |
| Sterile | | Non-sterile | Non-sterile | Same |
| Single for use | | Yes | Yes | Same |
| Biocompatibility | | Confirm to the requirements of ISO 10993 series standards | Confirm to the requirements of ISO 10993 series standards | Same |
| Fluid Resistance Performance ASTM F1862 | | Pass at 120mmHg | pass at 120mmHg | Similar |

| | | | |
|---|--|--|------|
| Particulate Filtration Efficiency ASTM F2299 | Average 99.82% at 0.1µm | Average 99.74% at 0.1µm | |
| Bacterial Filtration Efficiency ASTM F2101 | Average 99.78% | Average 99.4% | |
| Differential Pressure | Average 3.3 mmH ₂ O/cm ² | Average 2.7 mmH ₂ O/cm ² | |
| Flammability 16 CFR 1610 | Class I Non Flammable | Class 1 Non Flammable | |
| In Vitro Cytotoxicity ISO 10993-5 | Under the conditions of this study the device is non-cytotoxic | Under the conditions of this study the device is non-cytotoxic | Same |
| Skin Irritation ISO 10993-10 | Under the conditions of this study the device is non-irritating | Under the conditions of this study the device is non-irritating | Same |
| Skin Sensitization ISO 10993-10 | Under the conditions of this study the device is non-sensitizing | Under the conditions of this study the device is non-sensitizing | Same |

¹ The difference in the materials does not raise additional questions for safety and effectiveness of the device. The biocompatibility evaluation test of the subject devices have been performed on the final finished device which includes all construction materials and color additives. The test results shows pass the requirements.

VII. Summary of Non-Clinical Testing

Non-clinical performance tests were conducted to verify that the proposed device met all design specifications. The below table shows the test results of test article, which demonstrated that the proposed device complies with the related standards:

| Methodology/ Standard | Purpose | Acceptance Criteria | Results |
|--------------------------|---------------------------------|---------------------------------|-------------------------|
| ASTM F1862M-17 | Fluid Resistance Performance | 29 out of 32 pass at 120mmHg | 29 out of 32 pass at |

| | | | |
|-----------------------|-----------------------------------|--------------------------------------|--|
| | | | 120mmHg |
| ASTM F2299 | Particulate Filtration Efficiency | ≥98% | Average 99.82% at 0.1µm |
| ASTM F2101-19 | Bacterial Filtration Efficiency | ≥98% | Average 99.78% |
| EN 14683:2019 Annex C | Differential Pressure | <6mmH ₂ O/cm ² | Average 3.3 mmH ₂ O/cm ² |
| 16 CFR 1610 | Flammability | Class I Non Flammable | Meet Class I |

VIII. Clinical Testing

No clinical study is included in this submission.

IX. Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the subject device is as safe as effective, and performs as well as or better than the legally marketed predicate device K153496.