



July 29, 2021

Largent Health, LLC  
% Elizabeth O’Keeffe  
Regulatory Affairs Director  
Secure BioMed Evaluations  
7828 Hickory Flat Highway Suite 120  
Woodstock, Georgia 30188

Re: K210115

Trade/Device Name: FiteBac CC OrthoSeal  
Regulation Number: 21 CFR 872.3750  
Regulation Name: Bracket adhesive resin and tooth conditioner  
Regulatory Class: Class II  
Product Code: DYH  
Dated: May 28, 2021  
Received: June 1, 2021

Dear Elizabeth O’Keeffe:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Michael E. Adjodha, M.ChE.  
Assistant Director  
DHT1B: Division of Dental and ENT Devices  
OHT1: Office of Ophthalmic, Anesthesia,  
Respiratory, ENT and Dental Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES

Food and Drug Administration

**Indications for Use**

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023

See PRA Statement on last page.

510(k) Number (if known)

K210115

Device Name

FiteBac® CC OrthoSeal

Indications for Use (Describe)

FiteBac® CC OrthoSeal with K18 QAMS is a light cure primer intended to prepare the tooth surface prior to bonding orthodontic appliances to etched enamel.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services

Food and Drug Administration

Office of Chief Information Officer

Paperwork Reduction Act (PRA) Staff

*PRAStaff@fda.hhs.gov*

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) SUMMARY: FiteBac® CC OrthoSeal (K210115)

**Company Name:** Largent Health, LLC  
3698 Largent Way NW, Suite 101  
Marietta, GA 30064  
Phone: 770-218-6221

**Date Prepared:** July 28, 2021

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR § 807.92.

### 1. **General Information:**

**Applicant Name:**

Largent Health, LLC  
3698 Largent Way NW, Suite 101  
Marietta, GA 30064  
Phone: 770-218-6221

**Establishment Registration Number:** 3015709927

### 2. **Contact Person:**

Elizabeth O’Keeffe, Ph.D.  
Secure BioMed Evaluations  
7828 Hickory Flat Hwy, Suite 120  
Woodstock, GA 30188  
770-837-2681

[Regulatory@SecureBME.com](mailto:Regulatory@SecureBME.com)

**Secondary Contact:**

Patricia D. Jones, VP of Regulatory  
Secure BioMed Evaluations  
7828 Hickory Flat Highway, Suite 120  
Woodstock, GA 30188  
770-837-2681

[Regulatory@SecureBME.com](mailto:Regulatory@SecureBME.com)

### 3. **Device Name and Classification:**

<b>Trade Name:</b>	FiteBac® CC OrthoSeal
<b>Classification Name:</b>	Bracket adhesive resin and tooth conditioner
<b>Common Name:</b>	Adhesive, Bracket and Tooth Conditioner
<b>Classification Panel:</b>	Dental
<b>Regulation Number:</b>	21 CFR §872.3750
<b>Device Class:</b>	II
<b>Product Codes:</b>	DYH

**4. Primary Predicate Device:**

**Trade Name:** Ultradent Opal® Seal™  
**510(k) Clearance:** K090355  
**Clearance Date:** April 29, 2009  
**Classification Name:** Bracket adhesive resin and tooth conditioner  
**Common Name:** Adhesive, Bracket and Tooth Conditioner  
**Classification Panel:** Dental  
**Regulation Number:** 21 CFR §872.3750  
**Device Class:** II  
**Product Codes:** DYH

**Predicate Device:**

**Trade Name:** Transbond Light Cured Orthodontic Adhesive (Transbond XT Primer)  
**510(k) Clearance:** K880393  
**Clearance Date:** March 31, 1988  
**Classification Name:** Bracket adhesive resin and tooth conditioner  
**Common Name:** Adhesive, Bracket and Tooth Conditioner  
**Classification Panel:** Dental  
**Regulation Number:** 21 CFR §872.3750  
**Device Class:** II  
**Product Codes:** DYH

**Reference Device:**

**Trade Name:** Orthodontic Acrylic 2  
**510(k) Clearance:** K163482  
**Clearance Date:** May 31, 2017  
**Classification Name:** Denture relining, repairing, or rebasing resin.  
**Classification Panel:** Dental  
**Regulation Number:** 21 CFR §872.3760  
**Device Class:** II  
**Product Codes:** EBI

**5. Device Description:**

The FiteBac® CC OrthoSeal is a 5% K18 Quaternary Ammonium Methacrylate (QAMS) functionalized, light cured orthodontic sealant used to prepare etched enamel for orthodontic bonding.

**Technological Characteristics:**

The subject device has substantially equivalent technological characteristics as the predicate devices and/or the reference device in terms of principles of operation, intended use, material performance, and biocompatibility.

**6. Indications for Use:**

FiteBac® CC OrthoSeal with K18 QAMS is a light cure primer intended to prepare the tooth surface prior to bonding orthodontic appliances to etched enamel.

**7. Substantial Equivalence:**

FiteBac® CC OrthoSeal has the same intended use, principles of operation and substantially equivalent technological characteristics as Ultradent’s Opal® Seal™ (K090355). Although FiteBac® CC OrthoSeal differs from the predicate devices in the presence of K18 QAMS, FiteBac® CC OrthoSeal and the reference device Lang Orthodontic Acrylic 2 both contain this same K18 QAMS compound. Both the subject device and the predicate devices have been shown to prime teeth surfaces without negatively effecting the bond strength of orthodontic adhesives. FiteBac® CC OrthoSeal is as safe and effective as the predicate and reference comparator devices currently cleared for marketing in the United States. Additional biocompatibility and non-clinical testing demonstrate FiteBac® CC OrthoSeal is as safe and effective as these devices and does not raise additional questions of safety and effectiveness.

**8. Comparison of Technological Characteristics with the Predicate devices:**

A comparison of the subject device to the predicate devices and reference device is shown in the following table.

Characteristics	Subject Device FiteBac® CC OrthoSeal	Primary Predicate Ultradent Opal® Seal™ K090355	Predicate Device Transbond Light Cured Orthodontic Adhesive (Transbond XT Primer) K0880393	*Reference Device Lang Orthodontic Acrylic 2 K163482	Comparison Results
<b>Indication for Use Statement</b>	FiteBac® CC OrthoSeal with K18 QAMS is a light cure primer intended to prepare the tooth surface prior to bonding orthodontic appliances to etched enamel.	Opal® Seal is a light cure primer that is used when bonding orthodontic appliances to etched enamel.	N/A	N/A	Comparable
<b>Principle of Operation</b>	Methacrylate Sealant	Methacrylate Sealant	Methacrylate Sealant	Methacrylate-based material	Equivalent to predicates
<b>Available Configurations</b>	Liquid	Liquid in a single barrel syringe	Liquid	Liquid and Powder	Equivalent to predicate devices
<b>Size</b>	6 mL	1.2- or 3-mL Syringe	6 mL	Multiple Sizes	Equivalent to secondary predicate device
<b>Package</b>	Polyethylene Bottles	NA	Polyethylene Bottles	Polyethylene Bottles	Equivalent to secondary predicate device

Characteristics	Subject Device FiteBac® CC OrthoSeal	Primary Predicate Ultradent Opal® Seal™ K090355	Predicate Device Transbond Light Cured Orthodontic Adhesive (Transbond XT Primer) K0880393	*Reference Device Lang Orthodontic Acrylic 2 K163482	Comparison Results
Biocompatibility	Yes - see table 12.2.1	Yes	Yes	Yes	Equivalent to predicate devices
CC OrthoSeal Light Curing Specification Comparison with Ultradent Opal® Seal™ (K090355)					
Characteristics	Subject Device FiteBac® CC OrthoSeal K210115		Primary Predicate Ultradent Opal® Seal™ K090355		
Time to Cure	10 seconds	20 seconds	20 seconds	10 seconds	
Wavelength	< 600 mW/cm <sup>2</sup>	≥ 600 mW/cm <sup>2</sup>	< 600 mW/cm <sup>2</sup>	≥ 600 mW/cm <sup>2</sup>	

**\*Reference Device Comparison Rationale:** The reference device Lang Orthodontic Acrylic 2 (K163482) was chosen as a technological reference due to both the subject device and the reference device containing K18 QAMS. See relevant data provided in MAF2174.

### 9. Non-clinical Testing – Bench Study Comparison

The subject device has mechanical properties substantially equivalent to commercially available devices with the same intended uses. The following characteristics were evaluated:

#### Performance Testing:

- Bond Strength Testing

#### Non-clinical Testing – Biocompatibility:

- Sensitization
- Cytotoxicity
- Irritation

#### Testing Summary:

Results of all conducted testing was found acceptable and do not raise any new issues of safety or effectiveness.

**10. Performance Standards:**

FiteBac® CC OrthoSeal complies with the following performance standards:

Recognition #	Area	Title of Standard	Reference Number and Date	Standards Development Organization
2-258	Biocompatibility	Biological Evaluation of Medical Devices – Part 1: Evaluation in Testing Within a Risk Management Process	10993-1:2018	ISO
2-245	Biocompatibility	Biological Evaluation of Medical Devices – Part 5: Tests for In Vitro Cytotoxicity	10993-5:2009	ISO
2-174	Biocompatibility	Biological Evaluation of Medical Devices – Part 10: Tests for Skin Irritation and Delayed Hypersensitivity	0993-10:2010	ISO

**11. Conclusion as to Substantial Equivalence:**

The FiteBac® CC OrthoSeal is comparable to listed Predicate Devices and any differences do not raise new questions of safety and effectiveness. When compared to marketed devices, performance data demonstrates that the FiteBac® CC OrthoSeal is at least as safe and effective as the Predicate Devices listed in this submission.