



February 7, 2020

Abbott Point of Care Inc.
Susan Tibedo
Director, Regulatory Affairs
400 College Road East
Princeton, NJ 08540

Re: K183688

Trade/Device Name: i-STAT CHEM8+ cartridge with the i-STAT 1 System
Regulation Number: 21 CFR 862.1665
Regulation Name: Sodium Test System
Regulatory Class: Class II
Product Code: JGS, CDS, CEM, CGZ
Dated: January 9, 2020
Received: January 10, 2020

Dear Susan Tibedo:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Kellie Kelm, Ph.D.
Acting Director
Division of Chemistry
and Toxicology Devices
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K183688

Device Name

i-STAT CHEM8+ cartridge with the i-STAT 1 System

Indications for Use (Describe)

The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the in vitro quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings.

Sodium measurements are used for monitoring electrolyte imbalances.

Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.

Chloride measurements are primarily used in the diagnosis, monitoring, and treatment of electrolyte and metabolic disorders including, but not limited to, cystic fibrosis, diabetic acidosis, and hydration disorders.

Blood urea nitrogen measurements are used for the diagnosis, monitoring, and treatment of certain renal and metabolic diseases.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

The information in this 510(k) summary is being submitted in accordance with the requirements of 21 CFR 807.92.

1. Submitter Information

Owner Abbott Point of Care Inc.
400 College Road East
Princeton, NJ 08540

Contact Primary: Susan Tibedo
Director Regulatory Affairs
susan.tibedo@abbott.com
Phone: 609-454-9360

Secondary: Maria Figueroa
Manager Regulatory Affairs
maria.l.figueroa@abbott.com
Phone: 609-454-9271

Date Prepared February 5, 2020

510(k) Number k183688

2. Device Information

Proprietary Name i-STAT CHEM8+ cartridge with i-STAT 1 System

Common Name Chemistry test, analyzer, handheld

Product code	Device Classification name	Regulation Number	Class	Panel
JGS	Electrode, Ion Specific, Sodium	862.1665	II	Clinical Chemistry
CEM	Electrode, Ion Specific, Potassium	862.1600	II	Clinical Chemistry
CGZ	Electrode, Ion Specific Chloride	862.1170	II	Clinical Chemistry
CDS	Electrode, Ion Specific, Urea Nitrogen	862.1770	II	Clinical Chemistry

3. Predicate Device

Proprietary Name SYNCHRON Systems Sodium Reagent on UniCel DxC
600/800 SYNCHRON Clinical System

 SYNCHRON Systems Potassium Reagent on UniCel DxC
600/800 SYNCHRON Clinical System

 SYNCHRON Systems Chloride Reagent on UniCel DxC
600/800 SYNCHRON Clinical System

 SYNCHRON Systems BUN Reagent on UniCel DxC 600/800
SYNCHRON Clinical System

510(k) Number K042291

Product code	Device Classification name	Regulation Number	Class	Panel
JGS	Electrode, Ion Specific, Sodium	862.1665	II	Clinical Chemistry
CEM	Electrode, Ion Specific, Potassium	862.1600	II	Clinical Chemistry
CGZ	Electrode, Ion Specific Chloride	862.1170	II	Clinical Chemistry
LFP	Conductivity Rate, Urea Nitrogen	862.1770	II	Clinical Chemistry

4. Device Description

The i-STAT CHEM8+ test cartridge contains test reagents to analyze whole blood at the point of care or in the clinical laboratory for sodium, potassium, chloride and blood urea nitrogen. The test is contained in a single-use, disposable cartridge. Cartridges require two to three drops of whole blood which are typically applied to the cartridge using a transfer device.

The i-STAT 1 Analyzer is a handheld, *in vitro* diagnostic analytical device designed to run only i-STAT test cartridges. The instrument interacts with the cartridge to move fluid across the sensors and generate a quantitative result (within approximately 2 minutes).

The i-STAT 1 System is comprised of the i-STAT 1 analyzer, the i-STAT test cartridges and accessories (i-STAT 1 Downloader/Recharger, electronic simulator and portable printer). The system is designed for use by trained medical professionals at the patient point of care or in the clinical laboratory and is for prescription use only.

5. Intended Use Statement

The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the *in vitro* quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings.

Sodium measurements are used for monitoring electrolyte imbalances.

Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.

Chloride measurements are primarily used in the diagnosis, monitoring, and treatment of electrolyte and metabolic disorders including, but not limited to, cystic fibrosis, diabetic acidosis, and hydration disorders.

Blood urea nitrogen measurements are used for the diagnosis, monitoring, and treatment of certain renal and metabolic diseases.

6. Summary Comparison of Technological Characteristics

Similarities and Differences: System (Test and Instrument) for Sodium		
Feature or Characteristic	Predicate SYNCHRON Systems Sodium Reagent on UniCel DxC 600/800 SYNCHRON Clinical System (K042291)	Candidate Sodium Test with i-STAT 1 System
Intended Use	<p>The sodium test systems are intended for the quantitative determination of sodium concentration in human serum, plasma or urine.</p> <p>Sodium measurements are used in the diagnosis and treatment of aldosteronism (excessive secretion of hormone aldosterone), diabetes insipidus, adrenal hypertension, Addison's disease, dehydration, inappropriate antidiuretic hormone secretion, or other diseases involving electrolyte imbalance.</p>	<p>The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the <i>in vitro</i> quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings.</p> <p>Sodium measurements are used for monitoring electrolyte imbalances.</p>
Reportable Range	100-200 mmol/L (Serum or Plasma) 10-300 mmol/L (Urine)	Sodium: 100-180 mmol/L (mEq/L)
Sample Type	Serum, plasma, urine	Arterial or venous whole blood
Sample Volume	0.5 mL (500 µL)	95 µL
Preparation	Sample tubes prepared and then processed within analyzer	Ready to Use
Traceability	NIST SRM 919	NIST SRM 956
Calibration	Must be conducted every 24 hours and with each new reagent	1-point on-board contained within the cartridge
Time to Test (Sample Stability)	<p>Serum or plasma: Within 8 hours at room temperature, or up to 48 hours if stored at +2°C to +8°C</p> <p>Urine: Within 2 hours of collection</p>	Heparinized samples: within 30 minutes of collection
Principle of Measurement	Ion selective electrode	Ion selective electrode
Reagent Format	Reagent handling system, stored within analyzer	Cartridge
Reagent Storage and Stability	Room temperature	2°C to 8°C (35-46°F)
Analyzer Type	Floor Model	Handheld

Similarities and Differences: System (Test and Instrument) for Potassium		
Feature or Characteristic	Predicate SYNCHRON Systems Potassium Reagent on UniCel DxC 600/800 SYNCHRON Clinical System (K042291)	Candidate Potassium Test with i-STAT 1 System
Intended Use	<p>The potassium test systems are intended the quantitative determination of potassium concentration in human serum, plasma or urine.</p> <p>Potassium measurements are used in the diagnosis and treatment of hypokalemia (metabolic alkalosis, metabolic acidosis or the absence of acid-base disturbances), hyperkalemia (over administration of potassium, acidosis, or crush injuries), renal failure, Addison's disease or other diseases involving electrolyte imbalance.</p>	<p>The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the <i>in vitro</i> quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings</p> <p>Potassium measurements are used in the diagnosis and monitoring of diseases and clinical conditions that manifest high and low potassium levels.</p>
Reportable Range	1-15.0 mmol/L (Serum or Plasma)	2.0 – 9.0 mmol/L (mEq/L)
Sample Type	Serum, plasma, urine	Arterial or venous whole blood
Sample Volume	0.5 mL (500 µL)	95 µL
Sample Preparation	Sample tubes prepared and then processed within analyzer	Ready to Use
Traceability	NIST SRM 918	NIST SRM 956
Calibration	Must be conducted every 24 hours and with each new reagent	1-point on-board contained within the cartridge
Time to Test (Sample Stability)	<p>Serum or plasma: Within 8 hours at room temperature, or up to 48 hours if stored at +2°C to +8°C</p> <p>Urine: Within 2 hours of collection</p>	<p>Non-anticoagulated samples: within 3 minutes of collection</p> <p>Heparinized samples: within 30 minutes of collection</p>
Principle of Measurement	Ion selective electrode	Ion selective electrode
Reagent Format	Reagent handling system, stored within analyzer	Cartridge
Reagent Storage and Stability	Room temperature	2°C to 8°C (35-46°F)
Analyzer Type	Floor Model	Handheld

Similarities and Differences: System (Test and Instrument) for Chloride		
Feature or Characteristic	Predicate SYNCHRON Systems Chloride Reagent on UniCel DxC 600/800 SYNCHRON Clinical System (K042291)	Candidate Chloride Test with i-STAT 1 System
Intended Use	<p>The chloride test systems are intended for the quantitative determination of chloride concentration in human serum, plasma, urine or cerebrospinal fluid (CSF).</p> <p>Chloride measurements are used in the diagnosis and treatment of electrolyte and metabolic disorders such as cystic fibrosis and diabetic acidosis.</p>	<p>The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the <i>in vitro</i> quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings.</p> <p>Chloride measurements are primarily used in the diagnosis, monitoring, and treatment of electrolyte and metabolic disorders including, but not limited to, cystic fibrosis, diabetic acidosis, and hydration disorders.</p>
Reportable Range	50-200 mmol/L (Serum or Plasma)	65-140 mmol/L (mEq/L)
Sample Type	Serum, plasma, urine	Arterial or, venous whole blood
Sample Volume	0.5 mL (500 µL)	95 µL
Sample Preparation	Sample tubes prepared and then processed within analyzer	Ready to Use
Traceability	NIST SRM 918/919	NIST SRM 956
Calibration	Must be conducted every 24 hours and with each new reagent	1-point on-board contained within the cartridge
Time to Test (Sample Stability)	<p>Serum or plasma: Within 8 hours at room temperature, or up to 48 hours if stored at +2°C to +8°C</p> <p>Urine: Within 2 hours of collection</p>	<p>Non-anticoagulated samples: within 3 minutes of collection</p> <p>Heparinized samples: within 30 minutes of collection</p>
Principle of Measurement	Ion selective electrode	Ion selective electrode
Reagent Format	Reagent handling system, stored within analyzer	Cartridge
Reagent Storage and Stability	Room temperature	2°C to 8°C (35-46°F)
Analyzer Type	Floor Model	Handheld

Similarities and Differences: System (Test and Instrument) for Blood Urea Nitrogen (BUN)		
Feature or Characteristic	Predicate SYNCHRON Systems BUN Reagent on UniCel DxC 600/800 SYNCHRON Clinical System (K042291)	Candidate BUN Test with i-STAT 1 System
Intended Use	The blood urea nitrogen (BUN) test systems are intended for the quantitative determination of urea nitrogen or urea concentration in human serum, plasma or urine. Urea nitrogen or urea measurements are used in the diagnosis and treatment of certain renal and metabolic diseases.	The i-STAT CHEM8+ cartridge with the i-STAT 1 System is intended for use in the <i>in vitro</i> quantification of sodium, potassium, chloride and blood urea nitrogen in arterial or venous whole blood in point of care or clinical laboratory settings Blood urea nitrogen measurements are used for the diagnosis, monitoring, and treatment of certain renal and metabolic diseases.
Reportable Range	1-150 mg/dL (Serum or Plasma) 10-1500 mg/dL (Urine)	3-140 mg/dL
Sample Type	Serum, plasma, urine	Arterial or venous whole blood
Sample Volume	0.5 mL (500 µL)	95 µL
Sample Preparation	Sample tubes prepared and then processed within analyzer	Ready to Use
Traceability	NIST SRM 912	NIST SRM 909
Calibration	Must be conducted every 24 hours and with each new reagent	1-point on-board contained within the cartridge
Time to Test (Sample Stability)	Serum or plasma: Within 8 hours at room temperature, or up to 48 hours at +2°C to +8°C Urine: Within 2 hours of collection	Non-anticoagulated samples: within 3 minutes of collection Heparinized samples: within 30 minutes of collection
Principle of Measurement	Conductivity	Ion selective electrode
Reagent Format	Reagent handling system, stored within analyzer	Cartridge
Reagent Storage and Stability	2°C to 8°C (35-46°F)	2°C to 8°C (35-46°F)
Analyzer Type	Floor Model	Handheld

7. Performance Characteristics

Analytical Performance

a. Precision

Precision 20 days (Aqueous Materials)

The precision of the i-STAT Potassium, Chloride, and BUN Tests on the i-STAT 1 Analyzer were evaluated using 5 levels of aqueous materials. This 20-day multi-day precision testing was based on CLSI document EP05-A3: *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. The study

was conducted using multiple instruments and one test cartridge lot over 20 days at one site. Total precision ('within-laboratory', S_T), within-run, (S_r), between-run, (S_{rr}) and between-day, (S_{dd}) were estimated for each level. The results of the 20-day precision study are shown in Table 1.

Table 1: 20-day Precision Study Results (i-STAT 1 Analyzer)											
i-STAT Test	Fluid Level	N	Mean	Total		Within-run		Between-run		Between-day	
				S_T	CV_T (%)	S_r	CV_r (%)	S_{rr}	CV_{rr} (%)	S_{dd}	CV_{dd} (%)
Sodium (mmol/L)	CV L1	81	100.0	0.27	0.3	0.25	0.3	0.07	0.1	0.07	0.1
	CV L2	81	121.6	0.35	0.3	0.33	0.3	0.09	0.1	0.10	0.1
	CV L3	81	134.8	0.31	0.2	0.28	0.2	0.09	0.1	0.10	0.1
	CV L4	80	160.4	0.41	0.3	0.39	0.2	0.10	0.1	0.10	0.1
	CV L5	80	178	0.42	0.2	0.40	0.2	0.10	0.1	0.11	0.1
Potassium (mmol/L)	CV L1	81	2.07	0.006	0.3	0.006	0.3	0.001	0.05	0.002	0.1
	CV L2	81	2.83	0.011	0.4	0.011	0.4	0.003	0.1	0.002	0.1
	CV L3	81	3.69	0.010	0.3	0.008	0.2	0.004	0.1	0.003	0.1
	CV L4	80	6.17	0.018	0.3	0.017	0.3	0.002	0.03	0.007	0.1
	CV L5	80	7.75	0.033	0.4	0.027	0.3	0.017	0.2	0.009	0.1
Chloride (mmol/L)	CV L1	81	70.9	0.47	0.7	0.43	0.6	0.12	0.2	0.12	0.2
	CV L2	81	76.2	0.53	0.7	0.50	0.7	0.13	0.2	0.13	0.2
	CV L3	81	89.2	0.33	0.4	0.29	0.3	0.08	0.1	0.12	0.1
	CV L4	80	107.9	0.43	0.4	0.40	0.4	0.11	0.1	0.12	0.1
	CV L5	80	122.3	0.48	0.4	0.44	0.4	0.15	0.1	0.12	0.1
BUN (mg/dL)	CV L1	81	107.3	0.91	0.8	0.81	0.8	0.41	0.4	0.04	0.04
	CV L2	81	59.7	0.92	1.5	0.86	1.4	0.23	0.4	0.20	0.3
	CV L3	81	10.5	0.12	1.1	0.11	1.0	0.05	0.5	0.03	0.3
	CV L4	80	8.1	0.18	2.2	0.17	2.1	0.03	0.4	0.05	0.6
	CV L5	80	4.1	0.15	3.7	0.14	3.4	0.04	1.0	0.03	0.7

Precision (Whole Blood)

The whole blood precision of the i-STAT Sodium, Potassium, Chloride, and BUN Tests on i-STAT 1 Analyzer were evaluated using venous whole blood (native or altered) samples targeted to three levels within the test reportable range.

One test cartridge lot was used across 3 point of care sites. At each site, each sample was tested 3 times on each of 7 i-STAT 1 Analyzers (total of 21 test results per sample per instrument). The results of the whole blood precision using the i-STAT 1 Analyzer are shown in **Table 2**.

Table 2: Whole Blood Precision Results – i-STAT 1 Wireless Analyzer										
i-STAT Test	Level	Site	N	Mean	Within-Analyzer		Total			
					SD	%CV	SD	SD 95% CI	%CV	%CV 95% CI
Sodium (mmol/L)	≤ 134	01	21	110.3	0.49	0.4	0.49	0.39 to 0.78	0.4	0.4 to 0.7
		02	21	123.9	0.38	0.3	0.38	0.31 to 0.57	0.3	0.3 to 0.5
		03	21	108.3	0.44	0.4	0.49	0.39 to 0.84	0.5	0.4 to 0.8
	135-145	01	21	139.2	0.58	0.4	0.58	0.48 to 0.83	0.4	0.3 to 0.6
		01	21	138.0	0.22	0.2	0.22	0.17 to 0.35	0.2	0.1 to 0.3
		01	21	136.5	0.65	0.5	0.68	0.54 to 1.13	0.5	0.4 to 0.8
		02	20	138.3	0.45	0.3	0.45	0.36 to 0.73	0.3	0.3 to 0.5
		02	14	141.8	0.46	0.3	0.46	0.38 to 0.81	0.3	0.3 to 0.6
		02	21	141.9	0.31	0.2	0.31	0.25 to 0.48	0.2	0.2 to 0.3
		02	21	140.0	0.38	0.3	0.38	0.31 to 0.63	0.3	0.2 to 0.5
		02	21	138.9	0.44	0.3	0.48	0.38 to 0.83	0.3	0.3 to 0.6
		02	21	138.3	0.53	0.4	0.53	0.44 to 0.76	0.4	0.3 to 0.5
		03	21	140.9	0.31	0.2	0.36	0.29 to 0.64	0.3	0.2 to 0.5
		03	21	141.0	0.22	0.2	0.22	0.17 to 0.35	0.2	0.1 to 0.2
		03	21	142.1	0.31	0.2	0.31	0.25 to 0.48	0.2	0.2 to 0.3
	03	21	143.0	0.22	0.2	0.22	0.17 to 0.35	0.2	0.1 to 0.2	
	03	21	139.8	0.44	0.3	0.44	0.35 to 0.71	0.3	0.3 to 0.5	
	≥ 146	01	21	150.1	0.62	0.4	0.62	0.51 to 0.88	0.4	0.3 to 0.6
02		21	163.2	0.49	0.3	0.49	0.40 to 0.70	0.3	0.2 to 0.4	
03		21	150.1	0.38	0.3	0.38	0.31 to 0.57	0.3	0.2 to 0.4	
Potassium (mmol/L)	2.75-3.25	1	21	2.80	0.022	0.8	0.022	0.017 to 0.035	0.8	0.6 to 1.3
		2	21	2.80	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0
		3	21	3.05	0.058	1.9	0.058	0.048 to 0.083	1.9	1.6 to 2.7
	>3.25 - <5.55	1	21	4.30	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0
		1	21	4.00	0.022	0.6	0.022	0.017 to 0.035	0.6	0.4 to 0.9
		1	21	3.98	0.049	1.2	0.049	0.040 to 0.070	1.2	1.0 to 1.8
		2	20	4.95	0.059	1.2	0.059	0.049 to 0.085	1.2	1.0 to 1.7
		2	14	4.20	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0
		2	21	4.13	0.053	1.3	0.053	0.044 to 0.076	1.3	1.1 to 1.8
		2	21	4.24	0.058	1.4	0.058	0.048 to 0.082	1.4	1.1 to 1.9
		2	21	4.29	0.038	0.9	0.038	0.031 to 0.057	0.9	0.7 to 1.3
		2	21	4.10	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0
		3	21	4.80	0.031	0.6	0.032	0.025 to 0.052	0.7	0.5 to 1.1
		3	21	3.97	0.053	1.3	0.053	0.044 to 0.078	1.3	1.1 to 2.0
		3	21	4.00	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0
	3	21	3.40	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0	
	5.55 - 6.05	1	21	5.80	0.022	0.4	0.022	0.017 to 0.035	0.4	0.3 to 0.6
		2	21	5.85	0.058	1.0	0.058	0.048 to 0.083	1.0	0.8 to 1.4
3		21	5.71	0.038	0.7	0.038	0.031 to 0.057	0.7	0.5 to 1.0	
7.25 - 7.75	1	21	7.60	0.000	0.0	0.000	0.000 to 0.000	0.0	0.0 to 0.0	
	2	21	7.64	0.053	0.7	0.053	0.044 to 0.081	0.7	0.6 to 1.1	
	3	21	7.73	0.053	0.7	0.053	0.044 to 0.078	0.7	0.6 to 1.0	

Table 2: Whole Blood Precision Results – i-STAT 1 Wireless Analyzer										
i-STAT Test	Level	Site	N	Mean	Within-Analyzer		Total			
					SD	%CV	SD	SD 95% CI	%CV	%CV 95% CI
Chloride (mmol/L)	<80	01	21	77.0	0.53	0.7	0.53	0.44 to 0.80	0.7	0.6 to 1.0
		02	21	77.4	0.95	1.2	0.95	0.78 to 1.40	1.2	1.0 to 1.8
		03	21	76.9	0.62	0.8	0.66	0.53 to 1.11	0.9	0.7 to 1.4
	90-112	01	21	102.0	0.62	0.6	0.62	0.51 to 0.89	0.6	0.5 to 0.9
		01	21	97.0	0.31	0.3	0.32	0.25 to 0.52	0.3	0.3 to 0.5
		01	21	102.3	0.53	0.5	0.53	0.44 to 0.78	0.5	0.4 to 0.8
		02	21	100.2	0.32	0.3	0.37	0.29 to 0.66	0.4	0.3 to 0.7
		02	21	101.3	0.53	0.5	0.53	0.44 to 0.90	0.5	0.4 to 0.9
		02	21	103.5	0.38	0.4	0.52	0.40 to 0.99	0.5	0.4 to 1.0
		02	21	102.1	0.31	0.3	0.31	0.25 to 0.48	0.3	0.2 to 0.5
		02	21	101.4	0.44	0.4	0.51	0.41 to 0.91	0.5	0.4 to 0.9
		02	21	101.0	0.22	0.2	0.22	0.17 to 0.35	0.2	0.2 to 0.3
		03	21	104.0	0.49	0.5	0.50	0.40 to 0.82	0.5	0.4 to 0.8
		03	21	104.0	0.62	0.6	0.62	0.50 to 0.95	0.6	0.5 to 0.9
		03	21	103.1	0.49	0.5	0.54	0.43 to 0.94	0.5	0.4 to 0.9
	03	21	101.0	0.31	0.3	0.32	0.25 to 0.52	0.3	0.2 to 0.5	
	03	21	99.8	0.44	0.4	0.44	0.36 to 0.64	0.4	0.4 to 0.6	
	>120	01	21	126.1	0.58	0.5	0.58	0.47 to 0.86	0.5	0.4 to 0.7
02		21	123.8	0.69	0.6	0.69	0.57 to 0.98	0.6	0.5 to 0.8	
03		21	123.2	0.44	0.4	0.55	0.43 to 1.00	0.4	0.3 to 0.8	
BUN (mg/dL)	< 10	01	21	5.0	0.00	0.0	0.00	0.00 to 0.00	0.0	0.0 to 0.0
		01	21	5.5	0.49	8.9	0.51	0.41 to 0.86	9.4	7.5 to 15.6
		02	21	7.0	0.00	0.0	0.00	0.00 to 0.00	0.0	0.0 to 0.0
		03	21	6.9	0.31	4.5	0.36	0.29 to 0.64	5.3	4.2 to 9.3
	10-25	01	21	14.0	0.00	0.0	0.00	0.00 to 0.00	0.0	0.0 to 0.0
		02	20	18.0	0.23	1.3	0.23	0.18 to 0.37	1.3	1.0 to 2.1
		02	14	23.5	0.60	2.5	0.60	0.50 to 1.00	2.5	2.1 to 4.3
		02	21	20.9	0.38	1.8	0.38	0.31 to 0.57	1.8	1.5 to 2.7
		02	21	10.0	0.00	0.0	0.00	0.00 to 0.00	0.0	0.0 to 0.0
		03	21	11.0	0.22	2.0	0.22	0.17 to 0.35	2.0	1.5 to 3.2
		03	21	14.0	0.00	0.0	0.00	0.00 to 0.00	0.0	0.0 to 0.0
	03	21	13.9	0.31	2.2	0.31	0.25 to 0.48	2.2	1.8 to 3.5	
	25-50	01	21	38.0	0.44	1.1	0.50	0.40 to 0.88	1.3	1.1 to 2.3
		02	21	46.0	1.11	2.4	1.12	0.89 to 1.81	2.4	1.9 to 3.9
		03	21	27.8	0.62	2.2	0.71	0.56 to 1.24	2.6	2.0 to 4.5
	> 110	01	21	111.8	2.82	2.5	2.82	2.31 to 4.20	2.5	2.1 to 3.8
		02	21	125.0	1.72	1.4	1.97	1.56 to 3.44	1.6	1.2 to 2.8
		03	21	118.6	1.83	1.5	1.83	1.50 to 2.71	1.5	1.3 to 2.3

b. Linearity

The study was designed based on CLSI EP06-A: *Evaluation of the linearity of quantitative measurement procedures.*

The linearity of the i-STAT Sodium, Potassium, Chloride, and BUN tests on the i-STAT 1 Analyzer was evaluated by preparing whole blood samples of varying analyte levels that spanned the reportable range of the tests. The i-STAT Sodium, Potassium, Chloride, and BUN tests demonstrated linearity over the reportable range as shown in **Table 3.** Regression summary of the Sodium, Potassium, Chloride, and BUN response

versus the concentration of the whole blood samples of varying analyte levels is also provided in **Table 3**.

i-STAT Test	Reportable Range	Range Tested	Slope	Intercept	R ²
Sodium (mmol/L)	100 - 180	89 – 205	1.077	-5.244	0.9994
Potassium (mmol/L)	2.0 – 9.0	1.7 – 10.4	1.109	-0.136	0.9999
Chloride (mmol/L)	65 – 140	76 – 158 WB 54 – 157 Plasma	0.867 0.969	7.353 2.130	0.9983 0.9987
BUN (mg/dL)	3 – 140	2 - 143	0.946	0.209	0.9950

c. Limit of Quantitation (LoQ)

The study was based on the CLSI EP17-A2: *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline–Second Edition*.

The LoQs of the i-STAT Sodium, Potassium and BUN tests were evaluated on the i-STAT 1 Analyzer using whole blood that was altered to low sodium (<100 mmol/L), potassium (< 2.0 mmol/L), or BUN (<3 mg/dL) concentrations and two test cartridge lots. For Chloride, plasma samples that were altered to <65 mmol/L were evaluated as whole blood samples altered to achieve a chloride concentration <65 mmol/L have been found to hemolyze. The LoQs for the i-STAT Sodium Potassium, Chloride, and BUN tests on the i-STAT 1 Analyzer were determined to be 91 mmol/L, 1.5 mmol/L, 56 mmol/L, and 1 mg/dL, respectively.

d. Interference

The interference performance of the i-STAT Sodium, Potassium, Chloride, and BUN tests on the i-STAT 1 System was evaluated using whole blood samples based on CLSI EP07 ED3: *Interference Testing in Clinical Chemistry, Third Edition*. The effect of each substance was evaluated by comparing the performance of a control sample, spiked with blank solvent solution, with the test results from a sample spiked with the potentially interfering substance at the toxic/pathological concentration based on CLSI EP37 ED1: *Supplemental Tables for Interference Testing in Clinical Chemistry, First Edition*, as applicable. A substance was identified as an interferent if the difference between the control and test samples was outside of the allowed error (Ea) for the i-STAT test.

Tables 4, 5, 6 and 7 contain the list of potentially interfering substances tested for the i-STAT Sodium, Potassium, Chloride, and BUN tests and the interference results.

i-STAT Sodium Test

Table 4: Summary of Substance Tested and Interference Results for the i-STAT Sodium test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
Acetaminophen	1.03	15.6	No	
N-Acetyl-L-Cysteine	0.92	15.0	No	
Acetylsalicylic Acid	0.167	3.0	No	
Ammonium Chloride	2.0 ^a	10.7	No	
Ascorbic Acid	0.298	5.25	No	
Bilirubin	0.684	40	No	
Lithium Bromide	37.5 ^a	325.69	No	
β-Hydroxybutyric Acid	6.0 ^a	62.47	No	
Calcium Chloride	5	20	No	
Hemoglobin	10 g/L	1000	No	
Sodium Heparin	3.30 U/mL	330 U/dL	No	
Ibuprofen	1.06	21.9	No	
Lithium Lactate	10	90	No	
Lithium Chloride	3.2 ^a	13.6	No	
Lithium Salicylate	0.207	2.86	No	
Magnesium Chloride	4.1	10	No	
Sodium Thiosulfate	16.7 ^a	264.04	Yes	Increased results ≥ 3.1 mmol/L
Triglyceride	16.94	1500	No	
Uric Acid	1.4	23.5	No	
Cholesterol	10.3	400	No	

i-STAT Potassium Test

Table 5: Summary of Substance Tested and Interference Results for the i-STAT Potassium test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
Acetaminophen	1.03	15.6	No	
N-Acetyl-L-Cysteine	0.92	15.0	No	
Ammonium Chloride	2.0 ^a	10.7	No	
Ascorbic Acid	0.298	5.25	No	
Benzalkonium Chloride	0.03	1.13	No	
β-Hydroxybutyric Acid	6.0 ^a	62.47	No	
Bilirubin	0.684	40	No	
Lithium Bromide	37.5 ^a	325.69	No	
Calcium Chloride	5.0	20	No	
Hemoglobin	10 g/L	1000	No	
Lithium Lactate	10	90	No	
Lithium Chloride	3.2 ^a	13.6	No	
Magnesium Chloride	4.1	10	No	
Lithium Salicylate	0.207	2.86	No	

^a No CLSI EP37 test concentration available. The molecular weight of the substance tested was used to convert the test concentration from mmol/L to mg/dL. The molecular weight of each substance could vary depending on the form chosen.

Table 5: Summary of Substance Tested and Interference Results for the i-STAT Potassium test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
Sodium Thiosulfate	16.7 ^{aa}	264.04	No	
Triglyceride	16.94	1500	No	
Cholesterol	10.3	400	No	

i-STAT Chloride Test

Table 6: Summary of Substance Tested and Interference Results for the i-STAT Chloride test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
Acetaminophen	1.03	15.6	No	
N-Acetyl-L-Cysteine	0.92	15.0	No	
L-Ascorbic Acid	0.298	5.25	No	
Bicarbonate	35.0 ^a	294	No	
Bilirubin	0.684	40	No	
Lithium Bromide	37.5	325.69	Yes	Increased results ≥ 2.4 mmol/L
β-Hydroxybutyric Acid	6.0 ^a	62.47	No	
Hemoglobin	10 g/L	1000	No	
Sodium Iodide	2.99 ^a	44.82	No	
Lithium Lactate	10	90	No	
Sodium Oxalate	0.09	1.206	No	
Lithium Salicylate	0.207	2.86	No	
Lithium Thiocyanate	0.898	5.22	No	
Sodium Thiosulfate	16.7 ^a	264.04	Yes	Increased results ≥ 4.19 mmol/L
Triglyceride	16.94	1500	No	
Cholesterol	10.3	400	No	

i-STAT BUN Test

Table 7: Summary of Substance Tested and Interference Results for i-STAT BUN test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
Acetaminophen	1.03	15.6	No	
N-Acetyl-L-Cysteine	0.92	15.0	No	
L-Ascorbic Acid	0.298	5.25	No	
β-Hydroxybutyric Acid	6.0 ^a	62.47	No	
Bilirubin	0.684	40	No	
Lithium Bromide	37.5 ^a	325.69	No	
Hemoglobin	10 g/L	1000	No	
Hydroxyurea	0.405	3.08	No	
Lithium Lactate	10	90	No	

^a No CLSI EP37 test concentration available. The molecular weight of the substance tested was used to convert the test concentration from mmol/L to mg/dL. The molecular weight of each substance could vary depending on the form chosen.

Table 7: Summary of Substance Tested and Interference Results for i-STAT BUN test				
Substance	Test Concentration		Interference (Yes/No)	Interference Results
	mmol/L	mg/dL		
pH	8.0 pH units	N/A	No	
Lithium Salicylate	0.207	2.86	No	
Lithium Thiocyanate	0.898	5.22	No	
Sodium Thiosulfate	16.7 ^a	264.04	No	
Triglyceride	16.94	1500	Yes	Increased results ≥ 10.2 mmol/L
Cholesterol	10.3	400	No	

Comparison Study

e. Method Comparison with Predicate Device

Method comparison was demonstrated in a study comparing the i-STAT Sodium Potassium, Chloride, and BUN tests performance on the i-STAT 1 System to the Beckman DxC. The study was based on CLSI guideline EP09c-ED3. Venous and arterial blood specimens were evaluated and analyzed on the i-STAT 1 Analyzer against plasma specimens on the Beckman DxC. A Passing-Bablok linear regression analysis was performed using the first replicate result from the i-STAT 1 Analyzer versus the singlicate result of the comparative method.

Table 8: Method Comparison Results				
i-STAT Test	N	Slope	Intercept	r
Sodium	187	1.00	2.00	0.96
Potassium	189	1.00	0.00	0.99
Chloride	176	1.000	0.000	0.96
Bun	184	0.940	1.675	0.99

8. Conclusion

The results of these studies demonstrate that performance of the i-STAT CHEM8+ Sodium, Potassium, Chloride and BUN test with the i-STAT 1 System are substantially equivalent to the comparative method.