



May 28, 2021

Medos International SARL  
% Ariell Joiner, Ph.D.  
Senior Regulatory Affairs Program Lead  
CERENOVUS  
6303 Blue Lagoon Drive, Suite 315  
Miami, Florida 33126

Re: K211344

Trade/Device Name: MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3  
XSFT, and GALAXY G3 Mini Microcoil Delivery Systems

Regulation Number: 21 CFR 882.5950

Regulation Name: Neurovascular Embolization Device

Regulatory Class: Class II

Product Code: HCG, KRD

Dated: April 30, 2021

Received: May 3, 2021

Dear Dr. Ariell Joiner:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Naira Muradyan, Ph.D.  
Assistant Director  
DHT5A: Division of Neurosurgical,  
Neurointerventional  
and Neurodiagnostic Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K211344

Device Name  
MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT, and GALAXY G3 Mini Microcoil Delivery Systems

### Indications for Use (Describe)

MICRUSFRAME, DELTAFILL, and DELTAXSFT Microcoil Delivery Systems are intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and are also intended for arterial and venous embolizations in the peripheral vasculature.

The GALAXY G3 Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.

The GALAXY G3 XSFT Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms.

The GALAXY G3 Mini Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 510(k) Summary – K211344

**I. Submitter** Medos International SARL  
Chemin-Blanc 38  
2400 Le Locle, Switzerland

Tel: +1 (908) 249-0182

Contact Person: Ariell Joiner  
Email: [ajoiner@its.jnj.com](mailto:ajoiner@its.jnj.com)

Date Prepared: April 30, 2021

### II. Devices

Table 1. Device	
<b>Device Proprietary Name</b>	MICRUSFRAME , DELTAFILL, DELTAXSFT , GALAXY G3, GALAXY G3 XSFT, and GALAXY G3 Mini Microcoil Delivery Systems
<b>Common or Usual Name</b>	Device, Neurovascular Embolization & Vascular, for Promoting Embolization
<b>Classification Name</b>	Class II 21 CFR 882.5950 – Device, Neurovascular Embolization Class II 21 CFR 870.3300 – Device, Vascular, For Promoting Embolization
<b>Regulatory Classification</b>	II
<b>Product Codes</b>	HCG, KR D

### III. Predicate Devices

The predicate devices are listed below in **Table 2**.

Table 2: Predicate Devices			
510(k) Number	Date Cleared	Name	Manufacturer
K171747	Jul 14, 2017	MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3 FILL, and GALAXY G3 XSFT Microcoil Delivery Systems	Medos International SARL*
K171862	Sep 25, 2017	GALAXY G3 Mini Microcoil Delivery System	
K150319	Jun 12, 2015	MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, and GALAXY G3 XSFT Microcoil Delivery Systems	

\*510(k)s were previously held by Codman & Shurtleff, Inc. The 510(k)s were transferred to Medos International SARL, which is the current holder of the 510(k)s. These devices are marketed on behalf of Medos International SARL under the CERENOVUS brand.


*Continued on next page*

## 510(k) Summary, Continued

### IV. Device Description

The MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT, and GALAXY G3 Mini Microcoil Delivery Systems consist of three components: a Microcoil System, a connecting cable, and a Detachment Control Box (DCB). Each component is sold separately.

As shown in **Figure 1**, the Microcoil System consists of a microcoil attached to a Device Positioning Unit (DPU).



The Microcoil System is packaged in an introducer sheath designed to protect the coil in the packaging dispenser and to provide support for introducing the coil into the microcatheter. The microcoil is the implantable segment of the device and is detached from the Device Positioning Unit (DPU) using the Detachment Control System (Detachment Control Box and connecting cable).

*For the MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, and GALAXY G3 XSFT microcoils:*

- The microcoil is fabricated from a platinum alloy wire. The wire is wound into a primary coil which may contain either a polypropylene suture (SR) or an absorbable polymer suture and then formed into a secondary shape. The secondary shape may be spherical, complex, or helical.

*For the GALAXY G3 Mini microcoils:*

- The microcoil is fabricated from a platinum alloy wire. The wire is wound into a primary coil which contains a polypropylene suture (SR) and then formed into a secondary shape. The secondary shape is complex.

*For all Microcoil Delivery Systems:*

- The DPU is a variable stiffness wire and has a radiopaque marker band located three (3) cm from its distal end. The Device Positioning Unit includes five (5) fluoro saver markers on the proximal section of the shaft. The markers are intended to indicate when the tip of the microcoil is approaching the tip of the microcatheter. When the distal-most marker reaches the proximal end of the Rotating Hemostatic Valve (RHV) on the microcatheter, the tip of the coil is approaching the tip of the microcatheter and fluoroscopy should be used to guide further coil insertion.
- The introducer sheath has three main components: an introducer tip, a translucent introducer body, and a re-sheathing tool.

*Continued on next page*

---

## 510(k) Summary, Continued

---

### IV. Device Description, continued

The ENPOWER Detachment Control Box (DCB) provides the energy necessary to allow for a thermo-mechanical detachment of the microcoil from the DPU. The connecting cable delivers the energy necessary to detach the embolic coil from the Microcoil System's detachment zone. The connecting cable is connected between the Microcoil System's hub connector on the DPU and the output connector on the DCB.

- The connecting cables may be one of two types: one with a remote detach button (the ENPOWER Control Cable) catalog no. ECB000182-00, or one without a detach button (standard connecting cable) catalog no. CCB00157-00.
- The ENPOWER Detachment Control Box, catalog no. DCB2000500, works with the ENPOWER Control Cable and with the standard connecting cable.

---

### V. Indications for Use

MICRUSFRAME, DELTAFILL, and DELTAXSFT Microcoil Delivery Systems are intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and are also intended for arterial and venous embolizations in the peripheral vasculature.

The GALAXY G3 Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.

The GALAXY G3 XSFT Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms.

The GALAXY G3 Mini Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.

---

### VI. Comparison of Technological Characteristics with Predicate Device

The devices in this submission include minor design changes (reduced distal outer diameter and taper transition) only to the Core Wire of the Device Positioning Unit. There are no modifications to the material of the Core Wire, to other elements of the Device Positioning Unit, to the components or materials of the microcoil, to the introducer, or to the ENPOWER Detachment Control System.

Endovascular coil embolization is the technological principle for both the subject and predicate devices. This technology is based on placing embolic coils in the neurovascular or peripheral vasculature to reduce or block blood flow. The subject devices and predicate devices are based on the same overall technological characteristics as shown in **Table 3-1** and **Table 3-2**.

---

*Continued on next page*

## 510(k) Summary, Continued

VI.  
 Comparison of  
 Technological  
 Characteristics  
 with Predicate  
 Device,  
 continued

<b>Table 3-1. Technological Characteristics of the Predicate and Subject Device</b>		
<b>Description</b>	<b>Predicate Devices: MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT Microcoil Delivery Systems (K150319, K171747)</b>	<b>This Submission: MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT Microcoil Delivery Systems</b>
<b>Indications for Use</b>	<p>MICRUSFRAME, DELTAFILL, and DELTAXSFT Microcoil Delivery Systems are intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and are also intended for arterial and venous embolizations in the peripheral vasculature.</p> <p>The GALAXY G3 Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.</p> <p>The GALAXY G3 XSFT Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms.</p>	Same as Predicate
<b>Microcoil</b>		
<b>Microcoil Material</b>	Platinum/Tungsten	Same as Predicate
<b>Microcoil Primary Wind</b>	Triangular or Cylindrical	Same as Predicate
<b>Microcoil Secondary Shape</b>	Complex, Helical, or Spherical	Same as Predicate
<b>Microcoil Stretch-Resistant</b>	PGA= Polyglycolic Acid Suture PP= Polypropylene Suture	Same as Predicate
<b>Primary Coil Wind Outer Diameter (OD)</b>	0.009” – 0.015”	Same as Predicate
<b>Secondary Shape OD Ranges</b>	1.5mm – 24mm	Same as Predicate
<b>Microcoil Length Ranges</b>	1cm - 60cm	Same as Predicate
<b>Delivery System</b>		
<b>Delivery System Type</b>	Wire Shaft with radiopaque marker	Same as Predicate
<b>Delivery System Introducer Sheath</b>	HDPE Introducer	Same as Predicate
<b>Delivery System Resheathing Tool</b>	Nylon 12	Same as Predicate
<b>Introducer Tip Wall Thickness</b>	0.0103”	Same as Predicate
<b>Introducer Sheath Length</b>	120cm and 81cm	Same as Predicate

*Continued on next page*

## 510(k) Summary, Continued

VI.  
 Comparison of  
 Technological  
 Characteristics  
 with Predicate  
 Device,  
 continued

<b>Table 3-1. Technological Characteristics of the Predicate and Subject Device, continued</b>		
<b>Description</b>	<b>Predicate Devices: MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT Microcoil Delivery Systems (K150319, K171747)</b>	<b>This Submission: MICRUSFRAME, DELTAFILL, DELTAXSFT, GALAXY G3, GALAXY G3 XSFT Microcoil Delivery Systems</b>
<b>Delivery System, continued</b>		
<b>Device Positioning Unit (DPU) Delivery System Length</b>	190cm ± 5cm	Same as Predicate
<b>Device Positioning Unit Diameter</b>	0.0159"	Same as Predicate
<b>Fluoroscopy Saver Markers</b>	Five Markers Located on the Proximal Section of the Shaft	Same as Predicate
<b>Fluoro Saver Marker Microcatheter Compatibility</b>	150cm Length	Same as Predicate
<b>Mechanism of Detachment</b>	Connection to Microcoil System: Uses Connecting Cable or ENPOWER Control Cable	Same as Predicate
	Detachment: Thermo-Mechanical System uses the ENPOWER Detachment Control Box (DCB) with ENPOWER Control Cable or Connecting Cable	
<b>Sterilization, Shelf Life, and Packaging</b>		
<b>Sterilization Method</b>	Electron Beam Radiation	Same as Predicate
<b>Shelf Life</b>	3 years	Same as Predicate
<b>Packaging</b>	Packaged in a plastic hoop and enclosed in a pouch	Same as Predicate

<b>Table 3-2. Technological Characteristics of the Predicate and Subject Device</b>		
<b>Description</b>	<b>Predicate Device: GALAXY G3 Mini Microcoil Delivery System (K171862)</b>	<b>This Submission: GALAXY G3 Mini Microcoil Delivery System</b>
<b>Indications for Use</b>	The GALAXY G3 Mini Microcoil Delivery System is intended for endovascular embolization of intracranial aneurysms, other neurovascular abnormalities such as arteriovenous malformations and arteriovenous fistulae, and is also intended for arterial and venous embolizations in the peripheral vasculature.	Same as Predicate

Continued on next page



## 510(k) Summary, Continued

VI.  
 Comparison of  
 Technological  
 Characteristics  
 with Predicate  
 Device,  
 continued

<b>Table 3-2. Technological Characteristics of the Predicate and Subject Device, continued</b>		
<b>Description</b>	<b>Predicate Device: GALAXY G3 Mini Microcoil Delivery System (K171862)</b>	<b>This Submission: GALAXY G3 Mini Microcoil Delivery System</b>
<b>Microcoil</b>		
<b>Microcoil Material</b>	Platinum/Tungsten	Same as Predicate
<b>Microcoil Primary Wind</b>	Cylindrical	Same as Predicate
<b>Microcoil Secondary Shape</b>	Complex	Same as Predicate
<b>Microcoil Stretch-Resistant</b>	PP= Polypropylene Suture	Same as Predicate
<b>Primary Coil Wind Outer Diameter (OD)</b>	0.009"	Same as Predicate
<b>Secondary Shape OD Ranges</b>	1.0mm – 3.0mm	Same as Predicate
<b>Microcoil Length Ranges</b>	1cm – 8cm	Same as Predicate
<b>Delivery System</b>		
<b>Delivery System Type</b>	Wire Shaft with radiopaque marker	Same as Predicate
<b>Delivery System Introducer Sheath</b>	HDPE Introducer	Same as Predicate
<b>Delivery System Resheathing Tool</b>	Nylon 12	Same as Predicate
<b>Introducer Tip Wall Thickness</b>	0.0103"	Same as Predicate
<b>Introducer Sheath Length</b>	81cm and 120cm	Same as Predicate
<b>Device Positioning Unit (DPU) Delivery System Length</b>	190cm ± 5cm	Same as Predicate
<b>Device Positioning Unit Diameter</b>	0.0159"	Same as Predicate
<b>Fluoroscopy Saver Markers</b>	Five Markers Located on the Proximal Section of the Shaft	Same as Predicate
<b>Fluoro Saver Marker Microcatheter Compatibility</b>	150cm Length	Same as Predicate
<b>Mechanism of Detachment</b>	Connection to Microcoil System: Uses Connecting Cable or ENPOWER Control Cable	Same as Predicate
	Detachment: Thermo-Mechanical System uses the ENPOWER Detachment Control Box (DCB) with ENPOWER Control Cable or Connecting Cable	
<b>Sterilization, Shelf Life, and Packaging</b>		
<b>Sterilization Method</b>	Ethylene Oxide (EtO)	Same as Predicate
<b>Shelf Life</b>	3 years	Same as Predicate
<b>Packaging</b>	Packaged in a plastic hoop and enclosed in a pouch	Same as Predicate

Continued on next page

## 510(k) Summary, Continued

### VII. Non-Clinical Performance Data

Table 4: Verification Testing		
Test	Test Method Summary	Results
Tracking Force (delivery)	The purpose of the Track Force test was to evaluate the force it takes to deliver the subject device through a microcatheter; utilizing the system Catheter Performance Simulation System (CPSS).	Pass Samples passed the established acceptance criteria

Table 5: Validation Testing		
Test	Test Method Summary	Results
Simulated Use: Tracking Force (delivery)	The purpose of the Tracking Force test was to evaluate the ease of deliverability of the coils through the microcatheter in a clinically relevant bench top model.	Pass Samples passed the established acceptance criteria
Simulated Use: Microcatheter Stability	The purpose of the Microcatheter Stability test was to evaluate the amount of movement of the microcatheter tip from its starting position in a clinically relevant bench top model.	Pass Samples passed the established acceptance criteria

#### Animal Testing

N/A – No animal studies were required as appropriate verification and validation of the modifications was achieved based on the bench testing.

#### Shelf Life Testing

N/A – Changes did not impact the shelf-life of the product.

#### Biocompatibility Testing

N/A – Changes did not impact biocompatibility.

#### Sterilization

N/A – Changes did not impact sterilization.

### VIII. Clinical Performance Data

Clinical studies were not required as appropriate verification and validation of the minor design modifications was achieved based on the bench testing.

### IX. Conclusion

The minor design modifications made to the Core Wire of the Device Positioning Unit do not alter the intended use or indications for use of the subject devices. The overall technological characteristics of the subject and predicate devices remain the same. The risk assessment and successful verification and validation testing raises no new questions regarding the safety and effectiveness of the devices. Therefore, the modified devices are substantially equivalent to their respective predicate devices.

*End of 510(k) Summary*