

This patient guide is for those who have significant mitral regurgitation (MR) and are considered at prohibitive risk for mitral valve surgery by your Heart Team.

The information in this booklet will help you understand more about a less invasive procedure called transcatheter mitral valve repair.

Be sure to ask your Heart Team to explain all of your treatment options and the possible risks and benefits of each.





Table of Contents

The Heart3-4
Mitral Valve Repair Options
Deciding on the Appropriate Treatment Option For You
The Edwards PASCAL Precision System – Implant7–8
The Transcatheter Valve Repair Procedure9–10
After Your Transcatheter Heart Valve
Repair Procedure11–12
Clinical Data13–14
Risks of the Transcatheter Valve Repair Procedure15–10
Warnings and Precautions 17–18


Edwards Lifesciences is the global leader in patient focused medical innovations for structural heart disease and has been helping critically ill patients for over 60 years. Driven by a passion to help patients, Edwards works to improve outcomes and enhance the lives of patients.

The Edwards PASCAL Precision system is designed to repair the mitral valve and help reduce mitral regurgitation.

The Heart



Heart with Degenerative Mitral Regurgitation



Mitral Valve Repair Options

Understanding Your Treatment Options

If you have mitral regurgitation and your doctor has evaluated you as prohibitive risk for open heart surgery, transcatheter valve repair may be an option for you. Only a specialized Heart Team can determine which treatment option is appropriate for you.

There are a few different ways to treat degenerative mitral regurgitation. Two common methods include treating with medical therapy, or repairing the mitral valve in order to prevent backward leakage and help blood flow properly.

Transcatheter Mitral Valve Repair


Deciding on the Appropriate Treatment Option for You

What Are the Benefits of Treating Mitral Regurgitation?


reating mitral regurgitation is key to maintaining your heart health, and could make a difference for your quality of life.

What Is the Best Treatment Option for You?

A specialized doctor on a Heart
Team will evaluate you for treatment
options. They will consider factors
about your health to decide the
most appropriate treatment
option for you.



The Edwards PASCAL Precision System – Implant



The Transcatheter Valve Repair Procedure

This section describes what happens during transcatheter mitral valve repair. It is intended as a general overview, and your experience may be different. Please talk to your doctor for more information about what you should expect.




Your doctor will make a small cut in your groin area to

Your doctor will then insert a tube-like device (delivery catheter) through the small cut and direct it to your mitral valve under imaging guidance. The PASCAL implant is attached to the tip of the delivery catheter.

Once the implant is properly positioned in your mitral valve, your doctor will grasp the two leaflets of your mitral

elease the implant from the delivery catheter and remove the delivery catheter from





After Your Transcatheter Heart Valve Repair Procedure

What Happens After the Transcatheter Mitral Valve Repair Procedure?

Transcatheter Mitral Valve Repair Follow-Up Visits

Regular check-ups with your doctor are very important. Y your doctor to have your heart valve checked at 30 days and up to 5 years after your procedure.

However, call or see your doctor whenever you have questions or concerns about your health.



linic u an

Edwards mplant.

s of your ve repair. ical or

Edwards

ur doctor



SAMPLE PATIENT

Implanting Physician
SAMPLE PHYSICIAN

Hospital
SAMPLE HOSPITAL
CITY STATE, COUNTRY ZIP CODE
Serial Model

For more information on your implant card, please go to Edwards.com




Clinical Data

Edwards Transcatheter Valve Repair Clinical Data

If you were to undergo mitral valve repair with the PASCAL Precision system, the risks you could expect are shown in the following table. These risks were based on a clinical study of the PASCAL Precision system.

Transcatheter Mitral Repair with the PASCAL Precision System			
Major Complications			
Death from any cause	2 out of 100	5 out of 100	
Death from a heart related cause	1 out of 100	1 out of 100	
Stroke	1 out of 100	1 out of 100	
Heart attack	1 out of 100	1 out of 100	
Kidney failure requiring dialysis	1 out of 100	1 out of 100	
Severe bleeding	3 out of 100	4 out of 100	
Unplanned repeat procedure or surgery	1 out of 100	2 out of 100	



Risks of the Transcatheter Valve Repair Procedure

What Are the Risks?

- Stroke
- Serious bleeding
- Unplanned repeat procedure or surgery



.


- X-ray contrast media
- Have an active infection of the mitral valve endocarditis)
- Have mitral regurgitation caused by rheumatic disease
- Have evidence of blood clots in the heart or veins leading to the heart



- Abnormal heart rhythms or cardiac arrest, which may require a pacemaker
- Abnormal low or high blood pressure
- Allergic reaction to anesthetic, contrast, heparin, Nitinol (Nickel and Titanium) and/or other medications
- Aneurysm or pseudoaneurysm
- Bleeding, stomach bleeding, hemolysis, or decreased blood count, which may require transfusion
- Blood clots in the legs (Deep Vein Thrombosis)
- Blood clots, particles, catheter fragments or air in the blood vessels, lungs, body or brain
- · Cardiogenic shock
- · Chest pain
- Damage or puncture of the heart or blood vessels that may require surgery
- Damage, injury to, narrowing, or tearing of the mitral valve or other valve structures
- Damage to the swallowing passage (esophagus), with possible puncture or narrowing
- · Dislodgement of a previous implant
- Failure to retrieve any PASCAL Precision system components

- · Fever or infection, including of the heart valve
- Fluid or blood around the heart or lungs
- Heart attack
- Implant deterioration (wear, tear, fracture or other), malposition, clotting, movement or embolization
- Kidney failure
- Lab values that are not normal
- Nerve injury, paralysis or neurological symptoms, including problems with movement or walking
- Organ failure, including heart failure
- Pain
- Respiratory compromise that may require prolonged need for a respirator
- Shortness of breath, fainting or dizziness, nausea, swelling, weakness, diminished exercise ability
- Skin burn, injury or tissue changes due to exposure to X-rays
- Single leaflet device attachment (SLDA)
- Vascular injury or trauma, including decreased blood flow, dissection or occlusion
- · Worsening of valvular insufficiency
- · Wound healing infection or slow healing

Warnings and Precautions



Precautions Prior to Use

Precautions After Use

- .
- .

For More Information about the Edwards Transcatheter Valve Repair Procedure:

Toll-free phone in the USA: Online: Mail:

888-713-1564

Email address:

 $patient_support@edwards.com$

