

December 1, 2020

Hi-Care Thai Gloves Co. Ltd. % Manoj Zacharias Consultant Liberty Management Group Ltd. 75 Executive Dr. STE 114 Aurora, Illinois 60504

Re: K202377

Trade/Device Name: Latex Examination Powder Free Gloves

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LYY Dated: October 31, 2020 Received: November 3, 2020

Dear Manoj Zacharias:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For: CAPT Elizabeth Claverie, M.S.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2020 See PRA Statement below.

K202377				
Device Name Latex Examination Powder Free Gloves				
Indications for Use (Describe)				
Latex Examination Powder Free Gloves are disposable device intended for medical purpose that is worn on the examiner's hand to prevent contamination between patient and examiner				
Type of Use (Select one or both, as applicable)				
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)			
CONTINUE ON A SEPARA	ATE PAGE IF NEEDED.			

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510K SUMMARY

K202377

As required by 21CFR§807.92(c)

A. APPLICANT INFORMATION

Applicant Hi-Care Thai Gloves Co. Ltd.

Address 199Moo.11, T.Banpru, A.Hatyai, Songkhla 90250, Thailand.

Phone +66-92 225 8472 Fax +66-74-291044

E-mail daniel@hicarethai.com

Contact Person Mr. Daniel John

Designation International Operations Manager,

Contact Number +66-92 225 8472 Contact Email daniel@hicarethai.com Date Submitted August 5th, 2020

B. DEVICE IDENTIFICATION

Name of the device Latex Examination Powder Free Gloves

Proprietary or trade name Palm Care Latex Examination Powder Free Gloves

Common or usual name

Latex Examination Powder Free Gloves

Classification name Patient Examination Glove

Device Classification Class I Product Code LYY

Regulation Number 21 CFR 880.6250 Review Panel General Hospital

C. PREDICATE DEVICE

Predicate device JR Medic Powder Free Latex Examination Gloves

Predicate 510(K) K192329 Regulatory Class Class I Product code LYY

D. DESCRIPTION OF THE DEVICE

Latex Examination Powder Free Gloves are manufactured to meet all the current specifications listed under the ASTM Specification D 3578-05 (Reapproved 2015), Standard Specification for Rubber Examination Gloves. They are made from Natural Rubber Latex. These gloves are natural in color (no color is added) and are powder free.

E. INDICATIONS FOR USE/INTENDED USE OF THE DEVICE:

Latex Examination Powder Free Gloves are disposable devices intended for medical purpose that are worn on the examiner's hand to prevent contamination between patient and examiner.

F. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS BETWEEN PREDICATE AND SUBJECT DEVICES

CHARACTERSTICS	STANDARDS	DEVIC	COMPARISON		
		PREDICATE	CUR		
510(k) Number		K192329	K202377		
Name of device		JR MEDIC Blue Latex Examination Powder Free Gloves	Latex Examination Powder Free Gloves		
Dimensions- Length	ASTMD3578-05	Length $> 230 \text{ mm}$	Length > 230 mm		
	(Reapproved		Size	Average	Similar
	2015)		X-Small	303	
			Small	304	
			Medium	304	
			Large	305	
Dimensions- Width	ASTMD3578-05 Width Min 95+/- (Reapproved 10 mm (for medium size)	Width Min 95+/-10 mm (for medium size)		Similar	
		`	Size	Average	
	•	,	X-Small	75	
			Small	84	
			Medium	94	
			Large	105	
Physical Properties-	ASTMD3578-05	Before Ageing		Ageing	Similar
Tensile Strength	(Reapproved	Tensile Strength	Tensile Streng	th > 18 Mpa	
	2015)	> 18 Mpa_	Size	Actual value	
		After Ageing	X-Small	22.13	
		Tensile Strength > 14 Mpa	Small	22.20	
		> 14 Mpa	Medium	22.25	
			Large	22.28	
			After Ageing Tensile Strength > 14 Mpa		Similar
			Size	Actual value	
			X-Small	18.57	
			Small	18.64	
			Medium	18.70	
			Large	18.74	

CHARACTERSTICS	STANDARDS	DEVICE 1	COMPARISON		
		PREDICATE	CU	JRRENT	
510(K) Number		K192329	K202377		1
Physical Properties- Ultimate Elongation	ASTMD3578-05 (Reapproved 2015)	Before Ageing Ultimate Elongation > 650%	Before Ageing Ultimate Elongation > 650% Size Actual value		Similar
		After Ageing	X-Small	867	-
		Ultimate Elongation	Small	871	-
		>500%	Medium	874	_
			Large	877	_
				er Ageing	1
				ongation > 500%	
			Size	Actual value	1
			X-Small	845	-
			Small	848	
			Medium	854	1
			Large	860	1
Thickness	ASTMD3578-05	Palm > 0.08 mm		> 0.08 mm	Similar
111101111000	(Reapproved 2015)	Finger > 0.08 mm		> 0.08 mm	
	, 11	S	Size	Palm Finger	
				(Actual (Actual	
				value) value)	
			X-Small	0.16 0.21	
			Small	0.16 0.21	
			Medium	0.16 0.21	
			Large	0.16 0.21	
Powder Free Residue	ASTMDD 3578-10	≤2 mg/glove		ng/glove	Similar
	(Reapproved2015)		Size	Residual	
				powder content	
			X-Small	(mg/glove)	
				0.21	-
			Small	0.21	-
			Medium	0.22	1
	D : 01:	TT 1 .1 11.1	Large	0.22	σ.
	Primary Skin	Under the condition	Under the condition of study, not		Same
	Irritation-ISO 10993-10:2010(E)	of study, not an irritant	an irritant		
	Dermal Sensitization-	Under the conditions	Under the cor	nditions of the	Same
Biocompatibility	ISO 10993-	of the study, not a	study, not a se		
	10:2010(E)	sensitizer	,		
	In vitro	Under the conditions	Under the conditions of the study, non-cytotoxic		Same
	cytotoxicity	of the study,			
	ISO10993-5	non-cytotoxic			
	:2009(E) Material mediated	Under the	Tindan the same difference Cale		G.
	Pyrogenicity	conditions of the	Under the conditions of the study, non-pyrogenic		Same
	ISO 10993-	study non pyrogenic	study, non-py	10501110	
	11:2017(E) / USP				
	41<151>				

CHARACTERSTICS	STANDARDS	DEVICE P	COMPARISON	
		PREDICATE	PREDICATE CURRENT	
510(K) Number		K192329	K202377	
	Bacterial Endotoxin	No data available	<20EU/pair of gloves	
Biocompatibility	test USP 42<85>			
	Acute Systemic Toxicity Test ISO 10993-11:2017(E)	Under the condition of study, not systemic toxic	Under the conditions of study, the device extracts do not pose a systemic toxicity concern	same
Water Tight (1000 ml)	ASTM D5151-06 (Reapproved 2015)	Passes AQL-1.5	Passes AQL-1.5	Same
Intended use		JR MEDIC Blue Latex Examination Powder Free Gloves are disposable devices intended for medical purpose that are won on the examiner's hand to prevent contamination between patient and examiner.	Latex Examination Powder Free Gloves are disposable devices intended for medical purpose that are won on the examiner's hand to prevent contamination between patient and examiner.	Same
Material	-	Natural Latex	Natural Latex	Identical
Color	-	Blue	Natural (No color is added)	different
Texture	-	Finger Texture	Finger texture	Identical
Size	ASTMD 3578-5 (Reapproved 2015)	Small, Medium, Large & X Large	X Small, Small, Medium, Large,	Similar
Single Use	Medical Glove Guidance Manual - Labeling	Single Use	Single Use	Same
Sterile/non sterile	-	Non sterile	Non sterile	Same
Powder/Powder free	-	Powder free	Powder free	Same
Label and Labeling	FDA Label requirements	Meets FDA's label and labeling requirements	Meets FDA's label and labeling requirements	Same
Manufacturer(s)	-	JR Engineering & Medical Technologies (M) SDN.BHD. Malaysia.	Hi-Care Thai Gloves Co. Ltd.	

G. SUMMARY OF NON-CLINICAL PERFORMANCEDATA

TEST METHOD	PURPOSE	ACCEPTANCE	RESULT
A CTM D2579 05 (D	To determine	CRITERIA Min 230 mm for all sizes	X-Small:-303 mm
ASTM D3578-05 (Reapproved 2015) Standard Specification for Rubber	the length of the	Min 230 mm for all sizes	Small:-304mm
Examination Gloves	gloves		Medium:-304mm
Examination Gloves	gioves		Large:-305mm
ASTM D3578-05 (Reapproved 2015)	To determine the	X-Small:-70+/-10 mm	X-Small:-75 mm
Standard Specification for Rubber	width of the gloves	Small:-80+/-10mm	Small:-84 mm
Examination Gloves		Medium:-95+/-10 mm	Medium:-94mm
		Large:-111+/-10 mm	Large:-105mm

TEST METHOD	PURPOSE	ACCEPTANCE CRITERIA	RESULT			
ASTM D3578-05 (Reapproved	To determine the	Palm 0.08 mm min	Size	Palm	Finger	
2015) Standard Specification	thickness of the	Finger 0.08 mm min	X-Small	0.16mm	0.21mm	
For Rubber Examination	gloves	for all sizes	Small	0.16mm	0.21mm 0.21mm	
Gloves			Medium	0.16mm	0.21mm 0.21mm	
			Large	0.16mm	0.21mm 0.21mm	
		Before Ageing	Size	Before	After	
		Tensile Strength	Size	ageing	ageing	
	To Determine the	18Mpa Min for all		agemg	ugenig	
	physical properties-	sizes	X-Small	22.13Mpa	18.57Mpa	
	Tensile strength	After Ageing Tensile	Small	22.20Mpa	18.64Mpa	
	8	Strength 14Mpa Min	Medium	22.25Mpa	18.70Mpa	
		for all sizes	Large	22.28Mpa	18.74Mpa	
		Before Ageing	Size	Before	After	
ASTM D3578-05 (Reapproved		Ultimate elongation		ageing	ageing	
2015) Standard Specification	To Determine the	650% Min for all				
for Rubber Examination	physical properties-	sizes After Ageing	X-Small	867%	845%	
Gloves	Ultimate Elongation	Ultimate Elongation	Small	871%	848%	
		500% Min for all	Medium	874%	854%	
		sizes	Large	877%	860%	
			Size	Before	NA	
	To Determine the		*** 0 11	ageing		
	physical properties-	Before Ageing 5.5	X-Small	5.1 Mpa		
	stress at 500%	Mpa Max for all sizes	Small	5.1Mpa		
	Elongation		Medium Large	5.2 Mpa		
ASTM D5151-06 (Reapproved	To determine the	AQL 2.5		5.2Mpa asses AQL 1.	5	
2015) Standard Test Method	holes in the gloves	AQL 2.3	Gioves Pa	isses AQL 1.	3	
for Detection of Holes in	noies in the gloves					
Medical Gloves						
ASTM D6124-06 (Reapproved	To determine the	2 Mg/Glove Max	Size	Residual	Powder	
2017) Standard Test Method	residual powder in the	<i>G</i>	Content			
for Residual Powder on	gloves		X-Small			
Medical Gloves			Small	0.21 mg/g	glove	
			Medium	0.22 mg/g	glove	
			Large	0.22 mg/g		
ASTM D 5712-95 (Re	To determine the	$200 \ \mu g / \ dm^2 \ Max$	Size	Extractal	Extractable Protein	
approved 2010) Standard Test	extractable protein in	for all sizes		content		
Method for the Analysis of	the gloves.		X-Small	50 μg/ dm		
Protein in Natural Rubber			Small	50 μg/ dm		
			Medium	50 μg/dm		
			Large	50 μg/ dm	l²	

H. SUMMARY OF CLINICAL PERFORMANCEDATA

Not applicable.

I. CONCLUSION

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission Latex Examination Powder Free Gloves (K202377) is as safe, as effective, and performs as well as or better than the legally marketed predicate device K192329.