

December 2, 2020

Hi-Care Thai Gloves Co. Ltd. % Manoj Zacharias Consultant Liberty Management Group Ltd. 75 Executive Dr. STE 114 Aurora, Illinois 60504

Re: K202384

Trade/Device Name: Blue Nitrile Examination Gloves Powder Free

Regulation Number: 21 CFR 880.6250

Regulation Name: Non-Powdered Patient Examination Glove

Regulatory Class: Class I, reserved

Product Code: LZA
Dated: October 31, 2020
Received: November 3, 2020

Dear Manoj Zacharias:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to https://www.fda.gov/medical-device-problems.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance) and CDRH Learn (https://www.fda.gov/training-and-continuing-education/cdrh-learn). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For: CAPT Elizabeth Claverie, M.S.
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

Indications for Use

510(k) Number (if known)

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

K202384			
Device Name Blue Nitrile Examination Gloves Powder Free			
Indications for Use (Describe) Blue Nitrile Examination Gloves Powder Free is disposable devices intended for medical purpose that are worn on the examiner's hand to prevent contamination between patient and examiner.			
Type of Use (Select one or both, as applicable)			
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)		

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY K202384

Prepared following 21CFR§807.92(C)

A. APPLICANT INFORMATION

510(K) Owner's Name Hi-Care Thai Gloves Co. Ltd.

Address 199Moo.11, T.Banpru, A.Hatyai, Songkhla 90250,

Thailand.

Phone +66-92 225 8472 Fax +66-74-291044

E-mail daniel@hicarethai.com

Contact Person Mr. Daniel John

Designation International Operations Manager,

Contact Number +66-92 225 8472 Contact Email daniel@hicarethai.com Date 510(K) summary prepared August 14th, 2020

B. DEVICE IDENTIFICATION

Name of the device Blue Nitrile Examination Gloves Powder Free Product proprietary or trade name Palm Care Blue Nitrile Examination Gloves Powder

Free

Common or usual name Patient Examination Gloves
Classification name Patient Examination Gloves

Device Classification Class-1 Product Code LZA

Regulation Number 21 CFR 880.6250 Review Panel General Hospital

C. PREDICATE DEVICE

Predicate Device JR Engineering & Medical Technologies (M)

SDN.BHD.

510(k) Number K192333

Regulatory Class 1 Product code LZA

D. DESCRIPTION OF THE DEVICE:

Blue Nitrile Examination Gloves Powder Free are Class I patient examination gloves bearing the product code Nitrile - LZA (21CFR880.6250). The gloves are made from acrylonitrile-butadiene copolymer dispersion. These gloves are blue in color and are powder free.

E. INTENDED USE OF THE DEVICE:

Blue Nitrile Examination Gloves Powder Free is disposable devices intended for medical purpose that are worn on the examiner's hand to prevent contamination between patient and examiner.

F. SUMMARY OF THE TECHNOLOGICAL CHARACTERISTICS OF THE DEVICE COMPARED TO THE PREDICATE DEVICE

CHARACTERSTICS	STANDARDS	DEVICE PERFO	REMARKS	
		PREDICATE	SUBJECT	-
510(K) Number		K192333	K202384	
Name of device		JR MEDIC Blue Nitrile Examination Gloves Powder free	Blue Nitrile Examination Gloves Powder free	
Dimensions	ASTMD 6319-10 (Reapproved 2015)	Length Min 230 m Width Min 95+/-10mm (for medium size)	Length Min 230 mm Width Min 95+/mm (for medium size)	Same
Physical Properties	ASTMD 6319-10 (Reapproved 2015)	Before Aging Tensile Strength min 14 Mpa Ultimate Elongation Min 500% After Aging Tensile Strength min 14 Mpa Ultimate Elongation Min 400%	Before Aging Tensile Strength min 14 Mpa Ultimate Elongation Min 500% After Aging Tensile Strength min 14 Mpa Ultimate Elongation Ultimate Elongation Min 400%	Same
Thickness	ASTMD 6319-10 (Reapproved 2015)	Palm min 0.05 mm Finger min 0.05 mm	Palm min 0.05 mm Finger min 0.05 mm	Same
Powder Residue	ASTMD 6319-10	≤2 mg/glove	≤2 mg/glove	Similar
	Primary Skin Irritation- ISO 10993- 10:2010(E)	Under the condition of study not an irritant	Under the condition of study not an irritant	Same
Biocompatibility	Dermal Sensitization- ISO 10993-10:2010(E)	Under the conditions of the study not a sensitizer	Under the conditions of the study not a sensitizer	Same
	In vitro cytotoxicity ISO10993-5 :2009(E)	Under the conditions of the study, cytotoxic	Under the conditions of the study, non- cytotoxic	different
	Acute Systemic Toxicity Test ISO 10993-11:2017(E)	Under the condition of study not systemic toxic	Under the conditions of study, the device extracts do not pose a systemic toxicity concern	Same
	Material Mediated Pyrogenicity ISO 10993-11:2017(E) / USP 41<151>	Under the conditions of the study non pyrogenic	Under the conditions of the study non pyrogenic	Same
	Bacterial Endotoxin test USP 42<85>	No data available	<20EU/pair of gloves	

CHARACTERSTICS	STANDARDS	DEVICE PER	FORMANCE	REMARKS
		PREDICATE	CURRENT	
Water Tight (1000 ml)	ASTM D 5151-06	Passes AQL-1.5	Passes AQL-1.5	Similar
Intended use		JR MEDIC Blue Nitrile Examination Gloves Powder-free is disposable devices intended for medical purpose that are won on the examiner's hand to prevent contamination between patient and examiner.	Blue Nitrile Examination Gloves Powder free is disposable devices intended for medical purpose that are won on the examiner's hand to prevent contamination between patient and examiner.	Similar
Material	ASTMD 6910-10 (Reapproved 2015)	Nitrile (NBR)	Nitrile (NBR)	Same
Color	-	Blue	Blue	Same
Texture	-	Finger Texture	Finger texture	Same
Size	ASTMD 6319-10 (Reapproved 2015)	Extra Small, Small, Medium, Large, Extra Large	Extra Small, Small, Medium, Large, Extra Large	Same
Single Use	Medical Glove Guidance Manual - Labeling	Single Use	Single Use	Same
Manufacturer(s)	-	JR Engineering & Medical Technologies (M) SDN.BHD. Malaysia.	Hi-Care Thai Gloves Co. Ltd. Thailand.	

There are no significant differences between the two products and are identical in terms of intended use, materials, design, manufacturing methods. Both devices meet the ASTM standard D6913-10(Reapproved 2015).

G. NON-CLINICAL TESTING SUMMARY

PERFORMANCE DATA

Test Method	Purpose	Acceptance Criteria	Result
ASTM D6319-10 (Reapproved	To determine the length	Min 230 mm for all sizes	X-Small:-305 mm
2015) Standard Specification for	of the gloves		Small:-305 mm
Nitrile Examination Gloves for			Medium:-305mm
Medical Application			Large:-306mm
			X-Large:-306mm
ASTM D6319-10 (Reapproved	To determine the width	X-small: -70+/-10 mm	X-small-74 mm
2015) Standard Specification for	of the gloves	Small: -80+/-10 mm	Small: -84 mm
Nitrile Examination Gloves for	_	Medium: -95+/-10mm	Medium: -94 mm
Medical Application		Large: -110+/-10 mm	Large: -105 mm
		X-Large: -120+/-10 mm	X-Large: -115 mm

TEST METHOD	PURPOSE	ACCEPTANCE	RESULT		
		CRITERIA			
ASTM D6319-10 (Reapproved	To determine	Palm 0.05 mm min	Size	Palm	Finger
2015) Standard Specification for	the thickness	Finger 0.05 mm min for	X-Small	0.10mm	0.13mm
Nitrile Examination Gloves for	of the gloves	all sizes	Small	0.10mm	0.13mm
Medical Application			Medium	0.10mm	0.14mm
			Large	0.10mm	0.13mm
			X-Large	0.10mm	0.14mm
ASTM D5151-06 (Reapproved	To determine	AQL 1.5	Gloves Passes AQL 1.5		
2015) Standard Test Method for	the holes in the				
Detection of Holes in Medical	gloves				
Gloves					
ASTM D6124-06 (Reapproved	To determine	2 Mg/Glove Max	Size	Residual Powder Content	
2017) Standard Test Method for	the residual		X-small	0.21mg/glove	
Residual Powder on Medical	powder in the		Small	0.21mg/glove	
Gloves	gloves		Medium	0.22 mg/glo	
			Large	0.22 mg/glove	
			X-Large	Large 0.22 mg/glove	
ASTM D6319-10	To Determine	Before Ageing	Size	Before	After
(Reapproved 2015) Standard	the physical	Tensile Strength 14Mpa		ageing	ageing
Specification for Nitrile	properties-	Minimal for all sizes	X-Small	18.43Mpa	17.49Mpa
Examination Gloves for Medical	Tensile strength	After Ageing Tensile	Small	18.49Mpa	17.51Mpa
Application		Strength 14Mpa	Medium	18.54Mpa	17.65Mpa
		Minimal for all sizes	Large	18.57Mpa	17.67Mpa
			X-Large	18.64Mpa	17.71Mpa
	To Determine	Before Ageing	Size	Before	After
	the physical	Ultimate Elongation		ageing	ageing
	properties-	500% Min for all sizes	X-Small	685%	654%
	Ultimate	After Ageing Ultimate	Small	688%	658%
	Elongation	Elongation 400% Min	Medium	692%	661%
		for all sizes	Large	695%	665%
			X-Large	698%	669%

H. CLINICAL TESTING SUMMARY

Not applicable - Clinical data is not needed for gloves or for most devices cleared by the 510(k) process.

I. CONCLUSION

The conclusions drawn from the non-clinical test demonstrate that the subject device K202384 in 510(K) submission, Blue Nitrile Examination Gloves Powder Free is as safe, as effective, and performs as well as or better than the legally marketed predicate device K192333.