



December 29, 2020

Stryker Neurovascular
Kathy Nguyen
Staff Regulatory Affairs Specialist
47900 Bayside Parkway
Fremont, California 94538

Re: K203219

Trade/Device Name: Trevo XP ProVue Retriever, Trevo NXT ProVue Retriever

Regulation Number: 21 CFR 882.5600

Regulation Name: Neurovascular Mechanical Thrombectomy Device for Acute Ischemic Stroke Treatment

Regulatory Class: Class II

Product Code: POL, NRY

Dated: October 30, 2020

Received: November 2, 2020

Dear Kathy Nguyen:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Naira Muradyan, Ph.D.
Assistant Director
DHT5A: Division of Neurosurgical,
Neurointerventional
and Neurodiagnostic Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K203219

Device Name
Trevor XP ProVue Retriever and Trevor NXT ProVue Retriever

Indications for Use (Describe)

1. The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should start within 6 hours of symptom onset.
2. The Trevor Retriever is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.
3. The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid artery (ICA) or middle cerebral artery (MCA)-MI segments with smaller core infarcts (0-50 cc for age < 80 years, 0-20 cc for age ≥ 80 years). Endovascular therapy with the device should start within 6-24 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Introduction:

According to the requirements of 21 CFR 807.92, the following information provides sufficient details to understand the basis for a determination of substantial equivalence.

Submitter Name, Address and Content:

Submitter: Stryker Neurovascular
47900 Bayside Parkway
Fremont, CA 94538-6515
(FDA Registration Number: 3008853977)

Contact: **Kathy Nguyen**
Staff Regulatory Affairs Specialist
Phone: 510-413-2366
Fax: 510-413-2588
Email: Kathy.Nguyen@stryker.com

Date Prepared: December 22, 2020

Device Name and Classification:

Trade/Proprietary Name: Trevo XP ProVue Retriever and Trevo NXT™ ProVue Retriever

Common Name: Trevo Retriever

Classification Name: Neurovascular Mechanical Thrombectomy Device for Acute Ischemic Stroke Treatment, 21 CFR 882.5600, Class II
Percutaneous Catheter, 21 CFR 870.1250 – Class II

Product Code: POL, NRY

Legally Marketed Predicate Device

Name of Predicate Device	Name of Manufacturer	510(k) Number
Trevo XP ProVue Retriever	Stryker Neurovascular	K190779
Trevo NXT™ ProVue Retriever	Stryker Neurovascular	K200117

Device Description

The Trevo Retriever consists of a flexible, tapered core wire with a shaped section at the distal end. Platinum markers at the distal end allow fluoroscopic visualization. In addition, the shaped section is also radiopaque. Retriever dimensions are indicated on product label. The Retriever delivery wire has a hydrophilic coating on the distal 101cm length to reduce friction during use. The Retriever has a shaft marker to indicate proximity of Retriever tip relative to Microcatheter tip. A Torque Device is provided with the Retriever to facilitate manipulation and retrieval. The Retriever comes preloaded in an insertion tool to introduce the Retriever into a Microcatheter.

Indications for Use

1. The Trevo Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should start within 6 hours of symptom onset.
2. The Trevo Retriever is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.
3. The Trevo Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid

artery (ICA) or middle cerebral artery (MCA)-M1 segments with smaller core infarcts (0-50 cc for age < 80 years, 0-20 cc for age ≥ 80 years). Endovascular therapy with the device should start within 6-24 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy.

Technological Characteristics and Product Feature Comparison

Stryker Neurovascular has demonstrated the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever is substantially equivalent to the Predicate device, Trevo XP ProVue Retriever (**K190779**) and Trevo NXT ProVue Retriever (**K20117**) based on the same indications for use, device design, materials, manufacturing, packaging and sterilization methods.

Comparison of Technological Characteristics with the Primary Predicate Device:

The intended use remains as a neurovascular mechanical thrombectomy device for acute ischemic stroke treatment used in the treatment of acute ischemic stroke to improve clinical outcomes. The technological characteristics and principles of operation remain unchanged.

The Subject device is substantially equivalent to the Predicate device cleared under **K190779 & K200117** based on the following:

- Same indications for use
- Same materials and manufacturing processes
- Same device design and technology
- Same biocompatibility information
- Same materials and processes for packaging
- Same sterilization method and process for devices

A comparison of the Subject device with the Predicate device is summarized in **Table 1** (Trevo XP ProVue Retriever) and **Table 2** (Trevo NXT ProVue Retriever) below.

Table 1. Product Feature Comparison of Subject Device to Predicate Device

Feature	<u>Predicate Device</u> Trevor XP ProVue Retriever (K190779)	<u>Subject Device</u> Trevor XP ProVue Retriever
Intended Use	Neurovascular mechanical thrombectomy device for acute ischemic stroke treatment used in the treatment of acute ischemic stroke to improve clinical outcomes.	Same
Indications for Use	<ol style="list-style-type: none"> 1. The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should start within 6 hours of symptom onset. 2. The Trevor Retriever is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. 3. The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid artery (ICA) or middle cerebral artery (MCA)-M1 segments with smaller core infarcts (0-50 cc for age < 80 years, 0-20 cc for age ≥ 80 years). Endovascular therapy with the device should start within 6-24 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy. 	Same
Regulation Number/ Name/ Class/ Product Code	<p>21 CFR 882.5600, Neurovascular Mechanical Thrombectomy Device for Acute Ischemic Stroke Treatment, Class II, POL</p> <p>21 CFR 870.1250, Percutaneous Catheter, Class II, NRY</p>	Same
Target Population	Patients experiencing acute ischemic stroke	Same
Anatomical Sites	Neurovasculature	Same

Feature	<p align="center">Predicate Device Trevo XP ProVue Retriever (K190779)</p>	<p align="center">Subject Device Trevo XP ProVue Retriever</p>																										
Device Description	<p>The Retriever consists of a flexible, tapered core wire with a shaped section at the distal end. Platinum markers at the distal end allow fluoroscopic visualization. In addition, the shaped section is also radiopaque. Retriever dimensions are indicated on product label. The Retriever has a hydrophilic coating to reduce friction during use. The Retriever has a shaft marker to indicate proximity of Retriever tip relative to Microcatheter tip. A torque device is provided with the Retriever to facilitate manipulation. The torque device is used to lock the core wire to the microcatheter during the procedure. Locking of the torque device to the wire allows the microcatheter and Retriever to be retracted as a system during clot retrieval. An insertion tool is provided to introduce the Retriever into a Microcatheter. The Insertion Tool is a sheath in which the Retriever comes preloaded. Once half the retriever's length is inserted into the microcatheter, the insertion tool is removed. Retrievers have a modified proximal end that permits attachment of the Abbott Vascular DOC Guide Wire Extension (REF 22260). Joining Guide Wire Extension to Retriever facilitates removal or exchange of a catheter while maintaining Retriever position in anatomy. After exchange has been completed, the extension can be detached.</p>	Same																										
Principle of Operation	<p>The Trevo Retriever is delivered to the thrombus using a microcatheter. The Microcatheter is then retracted to deploy the shaped section of the Retriever. If using an Aspiration Catheter, remove Microcatheter. Advance Aspiration Catheter over proximal section of Retriever while applying aspiration using a 60 mL syringe or an aspiration pump. The Retriever with Microcatheter or Aspiration Catheter are pulled back to capture the thrombus. The Retriever, thrombus, and Microcatheter or Aspiration Catheter are removed as a unit from the body.</p>	Same																										
Procedural Steps Aspiration Source	Syringe, Aspiration pump	Same																										
Sizes	3x20mm 4x20mm 4x30mm 6x25mm	Same																										
Accessory Devices	Insertion tool and torque device provided within product package	Same																										
Compatibility	<table border="1" data-bbox="505 1356 980 1829"> <thead> <tr> <th data-bbox="505 1356 623 1549">Retriever Size</th> <th data-bbox="623 1356 695 1549">Trevo Pro14 Microcatheter</th> <th data-bbox="695 1356 797 1549">Trevo Pro18 Microcatheter</th> <th data-bbox="797 1356 889 1549">Excelsior[®] XT-27[™] Microcatheters (150cm x 6cm straight REF XT275081)</th> <th data-bbox="889 1356 980 1549">Recommended Minimum Vessel ID (mm)</th> </tr> </thead> <tbody> <tr> <td data-bbox="505 1549 623 1591">Trevo XP ProVue</td> <td align="center" data-bbox="623 1549 695 1591">√</td> <td align="center" data-bbox="695 1549 797 1591">√</td> <td data-bbox="797 1549 889 1591"></td> <td align="center" data-bbox="889 1549 980 1829" rowspan="5">2.5</td> </tr> <tr> <td data-bbox="505 1591 623 1633">Trevo XP ProVue</td> <td data-bbox="623 1591 695 1633"></td> <td align="center" data-bbox="695 1591 797 1633">√</td> <td data-bbox="797 1591 889 1633"></td> </tr> <tr> <td data-bbox="505 1633 623 1686">Trevo ProVue 4X20mm</td> <td data-bbox="623 1633 695 1686"></td> <td align="center" data-bbox="695 1633 797 1686">√</td> <td data-bbox="797 1633 889 1686"></td> </tr> <tr> <td data-bbox="505 1686 623 1759">Trevo XP ProVue 4X30mm</td> <td data-bbox="623 1686 695 1759"></td> <td align="center" data-bbox="695 1686 797 1759">√</td> <td align="center" data-bbox="797 1686 889 1759">√</td> </tr> <tr> <td data-bbox="505 1759 623 1829">Trevo XP ProVue 6X25mm</td> <td data-bbox="623 1759 695 1829"></td> <td data-bbox="695 1759 797 1829"></td> <td align="center" data-bbox="797 1759 889 1829">√</td> </tr> </tbody> </table>	Retriever Size	Trevo Pro14 Microcatheter	Trevo Pro18 Microcatheter	Excelsior [®] XT-27 [™] Microcatheters (150cm x 6cm straight REF XT275081)	Recommended Minimum Vessel ID (mm)	Trevo XP ProVue	√	√		2.5	Trevo XP ProVue		√		Trevo ProVue 4X20mm		√		Trevo XP ProVue 4X30mm		√	√	Trevo XP ProVue 6X25mm			√	Same
Retriever Size	Trevo Pro14 Microcatheter	Trevo Pro18 Microcatheter	Excelsior [®] XT-27 [™] Microcatheters (150cm x 6cm straight REF XT275081)	Recommended Minimum Vessel ID (mm)																								
Trevo XP ProVue	√	√		2.5																								
Trevo XP ProVue		√																										
Trevo ProVue 4X20mm		√																										
Trevo XP ProVue 4X30mm		√	√																									
Trevo XP ProVue 6X25mm			√																									
Core Wire	Nitinol (nickel titanium alloy)	Same																										


Feature	Predicate Device Trevo XP ProVue Retriever (K190779)	Subject Device Trevo XP ProVue Retriever										
Shaped Section	Nitinol	Same										
Distal Coil	Platinum/Tungsten	Same										
Shaped Section Radiopaque Wire	Platinum/Tungsten	Same										
Outer Jacket	Not applicable.	Not applicable.										
Mid Coil	304 Stainless Steel	Same										
Proximal Coil / Proximal Support / Length	304 Stainless Steel	Same										
Solder	Gold/Tin	Same										
Hydrophilic Coating	Sodium hyaluronate mixture, L578 formulation	Same										
Dimensional Drawing	<p>3x20, 4x20mm Retriever sizes:</p> <p>4x30, 6x25mm Retriever sizes:</p>	Same										
Core Wire Dimension	Flatten Thickness 0.0030" +/- 0.0003"	Flatten Thickness 0.0040" +/- 0.0003"										
Overall Length (A)	180, 190cm	Same										
Total Shaped Section Length (nominal) (B)	32, 36, 40, 44mm	Same										
Full Diameter Length (C) 4x30 & 6x25	20, 25, 30mm	Same										
Active Shaped Section Length (C) 3x20 & 4x20 (D) 4x30 & 6x25	<table border="1" data-bbox="505 1541 932 1730"> <thead> <tr> <th>Retriever Size (mm)</th> <th>Dimension C (mm)</th> </tr> </thead> <tbody> <tr> <td>3x20</td> <td>NA</td> </tr> <tr> <td>4x20</td> <td>NA</td> </tr> <tr> <td>4x30</td> <td>35</td> </tr> <tr> <td>6x25</td> <td>30</td> </tr> </tbody> </table>	Retriever Size (mm)	Dimension C (mm)	3x20	NA	4x20	NA	4x30	35	6x25	30	Same
Retriever Size (mm)	Dimension C (mm)											
3x20	NA											
4x20	NA											
4x30	35											
6x25	30											
Shaped Section Diameter (nominal) (D) 3x20 & 4x20 (E) 4x30 & 6x25	3, 4, 6mm	Same										

Feature	Predicate Device Trevor XP ProVue Retriever (K190779)	Subject Device Trevor XP ProVue Retriever												
Proximal Core Wire Diameter (E) 3x20 & 4x20 (F) 4x30 & 6x25	<table border="1"> <thead> <tr> <th>Retriever Size (mm)</th> <th>Dimension E (inches)</th> </tr> </thead> <tbody> <tr> <td>3x20</td> <td>0.015</td> </tr> <tr> <td>4x20</td> <td>0.018</td> </tr> </tbody> </table> <table border="1"> <thead> <tr> <th>Retriever Size (mm)</th> <th>Dimension F (inches)</th> </tr> </thead> <tbody> <tr> <td>4x30</td> <td>0.018</td> </tr> <tr> <td>6x25</td> <td>0.018</td> </tr> </tbody> </table>	Retriever Size (mm)	Dimension E (inches)	3x20	0.015	4x20	0.018	Retriever Size (mm)	Dimension F (inches)	4x30	0.018	6x25	0.018	Same
Retriever Size (mm)	Dimension E (inches)													
3x20	0.015													
4x20	0.018													
Retriever Size (mm)	Dimension F (inches)													
4x30	0.018													
6x25	0.018													
Materials and Configuration	Polyethylene Hoop, polycarbonate mounting card, Tyvek/Film Pouch, HDPE Tubing Clips, Chipboard carton	Same												
Sterilization Method	100% EtO	Same												
How Supplied	Sterile/Single Use	Same												

Table 2. Product Feature Comparison of Subject Device to Predicate Device

Feature	Predicate Device Trevor NXT ProVue Retriever (K200117)	Subject Device Trevor NXT ProVue Retriever
Regulation Number	<ul style="list-style-type: none"> 21 CFR 882.5600 21 CFR 870.1250 	Same
Regulation Name	<ul style="list-style-type: none"> Neurovascular Mechanical Thrombectomy Device for Acute Ischemic Stroke Treatment Percutaneous Catheter 	Same
Classification	Class II	Same
Product Code	<ul style="list-style-type: none"> POL NRY 	Same
Intended Use	Neurovascular mechanical thrombectomy device for acute ischemic stroke treatment used in the treatment of acute ischemic stroke to improve clinical outcomes.	Same
Indications for Use	<ol style="list-style-type: none"> The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion, and smaller core infarcts who have first received intravenous tissue plasminogen activator (IV t-PA). Endovascular therapy with the device should start within 6 hours of symptom onset. The Trevor Retriever is intended to restore blood flow in the neurovasculature by removing thrombus in patients experiencing ischemic stroke within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment. The Trevor Retriever is indicated for use to restore blood flow in the neurovasculature by removing thrombus 	Same

Feature	<u>Predicate Device</u> Trevor NXT ProVue Retriever (K200117)	<u>Subject Device</u> Trevor NXT ProVue Retriever
	for the treatment of acute ischemic stroke to reduce disability in patients with a persistent, proximal anterior circulation, large vessel occlusion of the internal carotid artery (ICA) or middle cerebral artery (MCA)-M1 segments with smaller core infarcts (0-50 cc for age < 80 years, 0-20 cc for age ≥ 80 years). Endovascular therapy with the device should start within 6-24 hours of time last seen well in patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy.	
Target Population	Patients experiencing acute ischemic stroke	Same
Anatomical Sites	Neurovasculature	Same
TECHNOLOGICAL CHARACTERISTICS		
Device Description	The Retriever consists of a flexible, tapered core wire with a shaped section at the distal end. Platinum markers at the distal end allow fluoroscopic visualization. In addition, the shaped section is also radiopaque. Retriever dimensions are indicated on product label. The Retriever delivery wire has a hydrophilic coating on the distal 101cm to reduce friction during use. The Retriever has a shaft marker to indicate proximity of Retriever tip relative to Microcatheter tip. A torque device is provided with the Retriever to facilitate manipulation and retrieval. The Retriever comes preloaded in an insertion tool to introduce the Retriever into a Microcatheter.	Same
Principle of Operation	The Trevor Retriever is delivered to the thrombus using a microcatheter. The Microcatheter is then retracted to deploy the shaped section of the Retriever. If using an Aspiration Catheter, remove Microcatheter. Advance Aspiration Catheter over proximal section of Retriever while applying aspiration using a 60 mL syringe or an aspiration pump. The Retriever with Microcatheter or Aspiration Catheter are pulled back to capture the thrombus. The Retriever, thrombus, and Microcatheter or Aspiration Catheter are removed as a unit from the body.	Same
Procedural Steps Aspiration Source	Syringe, Aspiration pump	Same
Accessory Devices	Insertion tool and torque device provided within product package	Same
Sizes	3x32mm 4x28mm 4x41mm 6x37mm	Same

Feature	Predicate Device Trevo NXT ProVue Retriever (K200117)	Subject Device Trevo NXT ProVue Retriever																											
Compatibility	<table border="1"> <thead> <tr> <th data-bbox="505 300 594 499">Retriever Size</th> <th data-bbox="594 300 662 499">Trevo Pro14 Microcatheter</th> <th data-bbox="662 300 730 499">Trevo Trak 21 Microcatheter*</th> <th data-bbox="730 300 799 499">Trevo Pro18 Microcatheter</th> <th data-bbox="799 300 914 499">Excelsior[®] XT-27 Microcatheter (REF XT275081)</th> <th data-bbox="914 300 995 499">Recommended Minimum Vessel ID (mm)</th> </tr> </thead> <tbody> <tr> <td data-bbox="505 499 594 548">Trevo NXT 3x32</td> <td data-bbox="594 499 662 548">√</td> <td data-bbox="662 499 730 548">√</td> <td data-bbox="730 499 799 548">√</td> <td data-bbox="799 499 914 548">√</td> <td data-bbox="914 499 995 709" rowspan="4">2.5</td> </tr> <tr> <td data-bbox="505 548 594 596">Trevo NXT 4x28</td> <td data-bbox="594 548 662 596"></td> <td data-bbox="662 548 730 596">√</td> <td data-bbox="730 548 799 596">√</td> <td data-bbox="799 548 914 596">√</td> </tr> <tr> <td data-bbox="505 596 594 644">Trevo NXT 4x41</td> <td data-bbox="594 596 662 644"></td> <td data-bbox="662 596 730 644">√</td> <td data-bbox="730 596 799 644">√</td> <td data-bbox="799 596 914 644">√</td> </tr> <tr> <td data-bbox="505 644 594 709">Trevo NXT 6x37</td> <td data-bbox="594 644 662 709"></td> <td data-bbox="662 644 730 709">√</td> <td data-bbox="730 644 799 709">√</td> <td data-bbox="799 644 914 709">√</td> </tr> </tbody> </table>	Retriever Size	Trevo Pro14 Microcatheter	Trevo Trak 21 Microcatheter*	Trevo Pro18 Microcatheter	Excelsior [®] XT-27 Microcatheter (REF XT275081)	Recommended Minimum Vessel ID (mm)	Trevo NXT 3x32	√	√	√	√	2.5	Trevo NXT 4x28		√	√	√	Trevo NXT 4x41		√	√	√	Trevo NXT 6x37		√	√	√	Same
Retriever Size	Trevo Pro14 Microcatheter	Trevo Trak 21 Microcatheter*	Trevo Pro18 Microcatheter	Excelsior [®] XT-27 Microcatheter (REF XT275081)	Recommended Minimum Vessel ID (mm)																								
Trevo NXT 3x32	√	√	√	√	2.5																								
Trevo NXT 4x28		√	√	√																									
Trevo NXT 4x41		√	√	√																									
Trevo NXT 6x37		√	√	√																									
	Balloon Guide Catheters and Aspiration Catheters (commercially available aspiration catheters with minimum inner diameter 0.058 inches (1.47mm)) are recommended for use during thrombus removal procedures.	Same																											
MATERIALS																													
Core Wire	Nitinol (nickel titanium alloy)	Same																											
Shaped Section	Nitinol	Same																											
Distal Coil	Platinum/Tungsten	Same																											
Shaped Section Radiopaque Wire	Platinum/Tungsten	Same																											
Mid Coil	304 Stainless Steel	Same																											
Proximal Coil	Pebax	Same																											
Solder	Gold/Tin	Same																											
Hydrophilic Coating	Sodium hyaluronate mixture, T070 formulation	Same																											
DIMENSIONS																													
Dimensional Drawing	All Retriever sizes: 	Same																											
Core Wire Dimensions	Flatten Thickness 0.0030" +/- 0.0003"	Flatten Thickness 0.0040" +/- 0.0003"																											
Overall Length (A)	200cm	Same																											
Total Shaped Section Length (B)	32, 36, 40, 44mm	Same																											

Feature	<u>Predicate Device</u> Trevo NXT ProVue Retriever (K200117)	<u>Subject Device</u> Trevo NXT ProVue Retriever
Full Diameter Length (C)	21, 25, 30, 35mm	Same
Shaped Section Diameter (D)	3, 4, 6mm	Same
Delivery Core Wire Outer Diameter (E)	0.015, 0.019”	Same
Cell Coverage Length (F)	28, 32, 37, 41mm	Same
PACKAGING		
Materials and Configuration	Polyethylene Hoop, polycarbonate mounting card, Tyvek/Film Pouch, Chipboard carton	Same
Sterilization Method	100% EtO	Same
How Supplied	Sterile/Single Use	Same

Risk Assessment

Risk assessment of the modifications has been conducted in accordance with EN ISO 14971. Stryker Neurovascular has determined the modifications to the Predicate device raise no new questions of safety or effectiveness.

Results of verification testing are appropriate for use in determining that the Trevo Retriever devices are substantially equivalent to the Predicate device. Furthermore, the modifications associated with the Subject device did not result in the identification of any new failure modes nor were there any changes to existing failure modes, including no change to severity or occurrence, and, therefore, no change to overall residual risk.

Testing Summary

There are no changes to the device intended use or indications for use statement. Other than the modification of the core wire dimension, there are no changes in the device materials, manufacturing, packaging and sterilization methods; therefore, biocompatibility data, bench performance data, sterilization and stability data from the Predicate device (**K10779 & K200117**) are directly applicable. Relevant testing data supporting the Subject device are summarized as follows.

Performance Data – Bench Testing

Testing on the Subject device was conducted for core wire dimension change. The results of design verification testing conducted on the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever demonstrates that it is substantially equivalent to the legally marketed Predicate device. The design verification bench testing are summarized in **Table 3** below.

Table 3- Performance Data – Bench Testing

Test	Test Method Summary	Conclusions
Insertion Tool Use, Dimensions	<p><u>Purpose:</u> Evaluate device integrity after use of insertion tool</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.
Retriever / Microcatheter Deliverability (Track Test- First Push)	<p><u>Purpose:</u> To determine the trackability of retriever</p> <p><u>Method:</u> Measure the force required to track a retriever through a microcatheter placed in a clinically relevant vasculature model.</p>	Deliverability force met acceptance criteria.
Retriever In-Vivo Resheathability	<p><u>Purpose:</u> Evaluate device integrity after re-sheathing into microcatheter</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.
Retractability of Retriever into Balloon Guide Catheter	<p><u>Purpose:</u> Evaluate device integrity after retracting retriever into Balloon Guide Catheter</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.
Retriever Shaft Kink Resistance	<p><u>Purpose:</u> Evaluate device integrity after tracking through a tortuous anatomy</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.

Test	Test Method Summary	Conclusions
Reloadability into Insertion Tool	<p><u>Purpose:</u> Evaluate device integrity after reloading retriever into insertion tool</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.
Retriever Shaped Section Torque/Tensile Durability	<p><u>Purpose:</u> To determine durability of shaped section after multiple torque and tensile loading cycles.</p> <p><u>Method:</u> Apply two clockwise (CW) or counter clockwise (CCW) revolutions and load proximal 50% of exposed shaped section to 1.7 lbf five times, then pull to failure. Device shall withstand six cycles at 1.7 lbf before fracture.</p>	Torque/tensile durability met acceptance criteria.
Retriever Proximal Curl Resistance	<p><u>Purpose:</u> To determine if curl formation is present after the required revolutions.</p> <p><u>Method:</u> Apply 5 clockwise (CW) or 5 counter clockwise (CCW) revolutions, the core wire shall not form a curl or fracture.</p>	Curl resistance met acceptance criteria.
Retriever Mid Joint Tensile Strength	<p><u>Purpose:</u> To determine the mid joint tensile strength of a retriever</p> <p><u>Method:</u> Identify mid solder joint(s) and subject to tensile testing. The device shall withstand 1.7 lbf minimum before wire or solder joint failure.</p>	Mid Joint Tensile met acceptance criteria.
Retriever platinum wire and joint durability	<p><u>Purpose:</u> To determine the durability of platinum wire weaves and solder joints of a retriever</p> <p><u>Method:</u> Identify Mid-Joint and position tubing with Retriever inside. Wrap tubing containing shaped section around Dowel Pin in one layer until the entire shaped section is wrapped around the pin, unwrap and repeat a total of 6 times. Unsheath device and examine the platinum wire and solder joint areas of Retriever.</p>	Platinum wire and joint met acceptance criteria.
Retriever platinum wire joint tensile strength	<p><u>Purpose:</u> To determine the platinum wire solder joint tensile strength of a retriever</p> <p><u>Method:</u> Free Platinum wires by cutting platinum wire and then unweaving from Nitinol struts. The joint is subjected to tensile testing. The joint shall withstand 0.1 lbf minimum before wire or solder joint failure</p>	Platinum wire solder joint met acceptance criteria.

Test	Test Method Summary	Conclusions
ISO Fracture	<p><u>Purpose:</u> For testing for fracture of endovascular retrievers (called retriever herein) based on EN ISO 11070:2014, Sterile single-use intravascular introducers, dilators and guidewires, Annex F, Test method for fracture of guidewires.</p> <p><u>Method:</u> The retriever wire is wound around a cylindrical former for at least eight complete turns, then unwound and examined for fractures. The diameter of the cylindrical former is less than or equal to ten times the maximum outside diameter (OD) of the wire.</p>	Fracture resistance met acceptance criteria.
ISO Flexure	<p><u>Purpose:</u> For testing for resistance of guide wires, delivery wires, or introducer wires to damage by flexing based on EN ISO 11070:2014, Annex G: Sterile, single-use intravascular catheter introducers, Test for resistance of guide wires to damage by flexing.</p> <p><u>Method:</u> The retriever delivery wire is subjected to 20 cycles of repeated reverse bending and straightening, then examined for damage and flaking of lamination.</p>	Flexure resistance met acceptance criteria.
Retrievability of Retriever with Intermediate Catheter	<p><u>Purpose:</u> Evaluate device integrity after retrieval of Retriever with Intermediate Catheter</p> <p><u>Method:</u> Device is prepared per IFU and inspected prior to use. The device is subjected to 3 runs of simulated use in a clinically relevant vasculature model and inspected after each run.</p>	Device integrity met acceptance criteria.

Performance Data – Animal, Clinical

No animal or clinical studies were conducted as there is no change to the indications for use or the fundamental scientific technology. Substantial equivalence of the Subject device has been established to the Predicate devices through the results of bench testing.

Shelf Life Testing

Shelf life testing previously conducted for the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever were leveraged to support the changes to the device and can be found in **K190779** and **K200117**. Shelf life testing was not performed on the Subject device since there was no impact to device material, or safety and efficacy as a result of the core wire design modification. As with the Predicate device, the Subject device is labeled with a 2-year shelf life.

Sterilization

Sterilization evaluation previously conducted for the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever were leveraged to support the labeling changes to the device and can be found in **K190779** and **K200117**. The Trevo XP ProVue Retrievers and Trevo NXT ProVue Retrievers are sterilized with 100% Ethylene Oxide and provided sterile. A sterility assurance level (SAL) of 10^{-6} has been demonstrated. The Trevo XP ProVue Retrievers and Trevo NXT ProVue Retrievers meet EO residuals per EN ISO 10993-7 for limited contacting, externally communicated devices. The Trevo NXT ProVue Retrievers are for single use only.

Biocompatibility

Biocompatibility testing previously conducted for the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever was leveraged to support the core wire design modification to the device and can be found in **K190779** and **K200117**. The results of biocompatibility testing, and biological safety evaluation of the Trevo XP ProVue Retriever and Trevo NXT ProVue Retrievers demonstrate that the devices meet biological safety requirements per ISO 10993-1 for externally communicating medical devices with circulating blood contact for less than 24 hours. The Trevo XP ProVue Retriever and Trevo NXT ProVue Retrievers are considered to have no residual risk of biological hazards. Also, the devices and its packaging do not contain detectable latex. Therefore, the Trevo XP ProVue Retriever and Trevo NXT ProVue Retriever devices, accessories, and primary packaging are considered biocompatible for their intended use.

Conclusion

Based on the conclusions drawn from risk assessment and the bench testing results summarized above, the Subject device demonstrates substantial equivalence to the legally marketed Predicate device (**K190779** and **K200117**).