

November 10, 2021

Shandong Shangwei Medical Products Co.,Ltd
% Chu Xiaoan
Official Correspondent
Beijing Easy-Link Company
Rm. F302 Bldg., 41, Jing Cheng Ya Ju,
Courtyard 6 of Southern Dou Ge Zhuang, Chaoyung District
Beijing, 100121
China

Re: K210779

Trade/Device Name: Powder Free Nitrile Patient Examination Gloves, Blue Color Regulation Number: 21 CFR 880.6250 Regulation Name: Non-powdered patient examination glove Regulatory Class: Class I, reserved Product Code: LZA Dated: September 17, 2021 Received: October 19, 2021

Dear Chu Xiaoan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <u>https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems</u>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</u>) and CDRH Learn (<u>https://www.fda.gov/training-and-continuing-education/cdrh-learn</u>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<u>https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</u>) for more information or contact DICE by email (<u>DICE@fda.hhs.gov</u>) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD Assistant Director DHT4B: Division of Infection Control and Plastic Surgery Devices OHT4: Office of Surgical and Infection Control Devices Office of Product Evaluation and Quality Center for Devices and Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

# **Indications for Use**

510(k) Number (if known)

#### K210779

Device Name

Powder Free Nitrile Patient Examination Gloves, Blue Color

Indications for Use (Describe)

Powder Free Nitrile Patient Examination Gloves, Blue Color is a non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff *PRAStaff@fda.hhs.gov* 

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Form Approved: OMB No. 0910-0120 Expiration Date: 06/30/2020 See PRA Statement below.

#### 510(K) Summary

## K210779

This summary of 510(k) is being submitted in accordance with 21 CFR 807.

Premarket Notification [510(k)] Summary

## 1.0 Submitter:

Submitter's name: Shandong Shangwei Medical Products Co.,Ltd Submitter's address: North Road, Fumin Avenue,Qinghe Street, Caoxian County, Heze City, Shandong Province, 274400,P.R. China Name of contact person: Ms. Li Hua Phone number: 0086-530-2069711 Date of preparation: 2021-11-08

## 2.0 Name of the Device

Proprietary/Trade name: Powder Free Nitrile Patient Examination Gloves, Blue Color Common Name: Patient Examination gloves Classification Name: Non-powdered Patient examination glove Device Classification: I Regulation: 21 CFR 880.6250 Panel: General Hospital Product Code: LZA

## 3.0 Predicate device

Device Name: Powder Free Nitrile Patient Examination Glove, Blue Color Company name: Tangshan Zhonghong Pulin Plastic Co., Ltd. 510(K) Number: K120970

#### **4.0 Device Description:**

The proposed device is Powder Free Nitrile Examination Gloves. The proposed device is blue. The proposed device is non-sterile.

## 5.0 Indications for Use Statement:

Powder Free Nitrile Patient Examination Gloves, Blue Color is a non-sterile disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

reenhological Characteristic Comparison.					
Features &	Predicate Device	Subject Device	Result of		
Description	(K120970)	(K210779)	Comparison		
Product name	Powder Free Nitrile Patient	Powder Free Nitrile Patient	Same		
	Examination Glove, Blue Color	Examination Gloves, Blue Color			
Regulation Number	21CFR880.6250	21CFR880.6250	Same		
Product Code	LZA	LZA	Same		
Color	Blue	Blue	Same		
Size	Small/ Medium/	Small/ Medium/	Same		
	Large/X large	Large/X large			
Indications for Use	Powder Free Nitrile Patient Examination Glove, Blue Color is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.	Powder Free Nitrile Patient Examination Gloves, Blue Color is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.	Same		

#### 6.0 Technological Characteristic Comparison:

Device Description	Meets ASTM D6319-10		Meets ASTM D6319-10		Same
and Specifications	10000 10110 20017-10		(Reapproved 2015)		2
DimensionsLength ILS-2 AQL4.0 (ASTMD 6319-10)	≥230mm min		232 mm min for all sizes		Similar
Dimensions	Small	70-90 mm	Small	76-90 mm	Similar
Width	Medium	85-105mm	Medium	87-102 mm	
IL S-2 AQL4.0	Large	100-120mm	Large	108-119mm	
(ASTM D6319-10)	X large	110-130 mm	X large	115-128 mm	
Dimensions	Finger 0.05mm min		Thickness (mm) min.		Similar
Thickness	Palm 0.05mm min.	-	Finger 0.08		
IL S-2 AQL4.0			Palm 0.08		
(ASTM D6319-10)					
Physical Properties	Before aging/after aging Tensile		Before Aging Elongation (%):		Similar
IL S-2 AQL4.0	Strength≥ 14MPa	00	550-610 After Aging Elongation (%):		
(ASTM D D6319-	8 -				
10)	Before aging Elong	gation $\geq 500\%$	450-570		
,	After aging Elongat				
			Before Aging Tensil	e Strength	
			(MPa): 18-25		
			After Aging Tensile	Strength	
			(MPa): 17-22		
Freedom from	Meets		1) Inspection Level I AQL2.5, and		Similar
Pinholes	• 21 CFR 800.20		Accept/Reject c		
	• ASTM D6319-1	10	2) Water leakage to		
Inspection			noncompliance	is allowed.	
Level I					
AQL2.5			1) (1 1 1	5	Similar
Residual Powder	h - 1	1	1) Checked or	1	Similar
(ASTM D 6124-	below 2mg of residual		sub-samples (N=5).		
06(Reaffirmation 2011))	powder		2) Result as following: Mean: 0.1mg/pcs		
2011))			wiean. 0.1mg/pes		
Materials used to	Nitrile		Nitrile		Same
fabricate the devices	1 (10110		1 110000		2
Single Patient Use	Single Patient Use		Single Patient Use		Same
U	C		C		
Biocompatibility	Under the condition	s of this study,	Under the conditions	s of this study,	Similar
	the test article was a non- irritant or non- sensitizer (ISO 10993-		the test article was a non- irritant or non- sensitizer (ISO 10993-10:		
	10:2002/Amd.1:		Third Edition 2010-0	08-01)	
			~ • • •		Diff
			Cytotoxicity study me		Different
			5 Third edition 2009-06-01		D'CC
	N/A		Under the conditions		Different
			device extracts do no		
			systemic toxicity con	ncern (ISO	
Labolina	-Powder Free		10993-11:2017) -Powder Free		Same
Labeling	-Powder Free -Patient Examinatio	n	-Powder Free -Patient Examination	Glove	Same
	Glove	11	-Single Use Only		
	-Single Use Only		- Manufactured For:		
	- Manufactured For		- Manufactured For: - Lot		
	- Infantifictured For	•	- Lot -Blue color		
	-Blue color		- Non sterile		
	- Non sterile				
					L

## 7.0 Discussion of Non-clinical and Performance Testing

Non-clinical tests were conducted to verify that the proposed device met all design specifications. The test results demonstrated that the proposed device complies with the following standards:

ISO 10993-5:2009 Biological Evaluation of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity

ISO 10993-10:2010 Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.

ISO 10993-11:2017 Biological evaluation of medical devices - Part 11: Tests for systemic toxicity

ASTM D6124-06 (Reapproved 2017), Standard Test Method for Residual Powder on Medical Gloves

ASTM D5151-19, Standard Test Method for Detection of Holes in Medical Gloves.

ASTM D6319-10(Reapproved 2015), Standard Specification for Nitrile Examination Gloves for Medical Application.

Test Methodology	Purpose	Acceptance Criteria		Results
ASTM D 6319-		Length	≥230mm	
06(Reapproved		Width	Small	76-90 mm
2015).			Medium	87-102 mm
	Dimension		Large	108-119mm
			X large	115-128 mm
		Thickness	Fingertip	≥0.08mm
			Palm	≥0.08mm
ASTM D 6319- 06(Reapproved 2015).	Physical	Tensile strength (Before & After aging)	≥14MPa	17-25
	Properties	Before aging Elongation	≥500%	550-610
		After aging Elongation	≥400%	450-570
<ul> <li>21 CFR 800.20</li> <li>ASTM D 6319- 06(Reapproved 2015).</li> <li>ASTM D5151-19</li> </ul>	Freedom from pinholes	Water leakage test: Inspection Level I, AQL2.5, and Accept/Reject criteria of 10/11.		5 noncompliance is allowed. Pass
<ul> <li>ASTM D6319- 10(Reapproved 2015)</li> <li>ASTM D6124-06 (Reapproved 2017),</li> </ul>	Powder Residual	Meets <2mg/glove		Mean: 0.1mg/pcs Pass
Primary Skin Irritation in rabbits ISO 10993-10: 2010- 08-01	Biocompatibil ity	Under the conditions of the study, the subject device is not a primary skin irritant.		Passes
Dermal sensitization in the guinea pig ISO 10993-10: 2010-08- 01		Under the conditions of the study, the subject device is not a skin sensitizer.		Passes

The test article was added to L929 cells measured by MTT assay ISO 10993-5: 2009	Under the conditions of this study, the test article was non- cytotoxicity to L-929 cells.	Under the conditions of the study, cytotoxic. Additional testing was performed to determine if this was a systemic toxicity concern.
Acute Systemic Toxicity Systemic injection in mice ISO 10993-11:2017	Under the conditions of study the device extracts do not pose a systemic toxicity concern	Pass

## 8.0 Discussion of Clinical and Performance Testing

Clinical testing is not needed for this device.

### 9.0 Conclusion

The conclusions drawn from the nonclinical tests demonstrate that the proposed device is as safe, as effective, and performs as well as or better than the legally marketed predicated device.