



September 30, 2022

Suzhou Borage Medical Technology Co., Ltd.
% Grace Liu
Consultant
Shenzhen Joyantech Consulting Co. Ltd
1713A, 17th Floor, Block A, Zhongguan Times Square,
Nanshan District
Shenzhen, Guangdong 518000
China

Re: K222046
Trade/Device Name: Unscented Menstrual Tampon
Regulation Number: 21 CFR§ 884.5470
Regulation Name: Unscented Menstrual Tampon
Regulatory Class: II
Product Code: HEB
Dated: June 13, 2022
Received: July 11, 2022

Dear Grace Liu:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jason R. Roberts, Ph.D.
Assistant Director
DHT3B: Division of Reproductive,
Gynecology and Urology Devices
OHT3: Office of GastroRenal, ObGyn,
General Hospital and Urology Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K222046

Device Name
Unscented Menstrual Tampon

Indications for Use (Describe)

The unscented menstrual tampon is intended for insertion into the vagina for the absorption of menstrual or other vaginal discharge.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary - K222046

1. Contact Details

1.1 Applicant information

| | |
|-----------------------|--|
| Applicant Name | Suzhou Borage Medical Technology Co., Ltd. |
| Address | No.9 Shiheshan Road, Dongshan Town, Wuzhong District, Suzhou City, Jiangsu Province, 215107, China |
| Contact person | Zhou Nana |
| Phone No. | +86-512-66280008 |
| E-mail | znn@borage.com.cn |
| Date Prepared | 2022-09-27 |

1.2 Submission Correspondent

| | |
|--|---|
|  <p>Phone No.</p> <p>Contact person</p> <p>Contact person's e-mail</p> <p>Website</p> | Shenzhen Joyantech Consulting Co., Ltd |
| | 1713A, 17th Floor, Block A, Zhongguan Times Square, Nanshan District, Shenzhen, Guangdong Province, China |
| | +86-755-86069197 |
| | Grace Liu; Field Fu; |
| | grace@cefda.com ; field@cefda.com |
| | http://www.cefda.com |

2. Device Information

| | |
|----------------------------|------------------------------|
| Trade name | Unscented Menstrual Tampon |
| Common name | Unscented Menstrual Tampon |
| Classification | II |
| Classification name | Tampon, Menstrual, Unscented |
| Product code | HEB |
| Regulation No. | 21 CFR 884.5470 |

3. Legally Marketed Predicate Device

| | |
|----------------------|---|
| Trade Name | Tosama 100% Organic Cotton Menstrual Tampon |
| 510(k) Number | K15117 |
| Product Code | HEB |
| Manufacturer | TOSAMA, d.o.o. |

The predicate device has not been subject to a design-related recall.

4. Device Description

The proposed device is a traditional unscented menstrual tampon, and it has two types (i.e. digital tampon and applicator tampon). Each device consists of a tampon, including an absorbent pledget ("absorbent core") completely surrounded by overwrap ("security veil") and a removal string ("withdrawal cord"), and an applicator (only for the applicator tampon). The tampon is of the traditional cylindrical, bullet-like shape. The applicator

has a smooth, rounded tip to ease insertion.

The proposed device is provided in 4 absorbencies: light ($\leq 6g$), regular (6~9g), super (9~12g) and super plus (12~15g). Each device is individually wrapped and packaged in sealed multi-unit containers for retail sale. And it is provided non-sterile and for single use only.

5. Intended Use/Indication for Use

The unscented menstrual tampon is intended for insertion into the vagina for the absorption of menstrual or other vaginal discharge.

6. Technological Characteristics Comparison

Table 1 Technological Characteristics Comparison Table

| Comparison item | Proposed Device (K222046) | Predicate Device (K15117) | Comment |
|-----------------------------------|--|---|-----------|
| Manufacturer | Suzhou Borage Medical Technology Co., Ltd. | TOSAMA, d.o.o. | None |
| Product Name | Unscented Menstrual Tampons | Tosama 100% Organic Cotton Menstrual Tampon | None |
| Product Code | HEB | HEB | Same |
| Regulation Number | 21 CFR § 884.5470 | 21 CFR § 884.5470 | Same |
| Classification | Class II | Class II | Same |
| Intended Use/ Indications for Use | The unscented menstrual tampon is intended for insertion into the vagina for the absorption of menstrual or other vaginal discharge. | The Tosama 100% Organic Cotton Menstrual Tampon is intended for insertion into the vagina for the absorption of menstrual or other vaginal discharge. | Same |
| Single Use | Yes | Yes | Same |
| Sterility | Non-sterile | Non-sterile | Same |
| Design | Tampon with cylindrical shape and bullet-like tip. Applicator with smooth, rounded tip. | Tampon with cylindrical shape and bullet-like tip. Applicator with smooth, rounded tip. | Same |
| Absorbency | | | |
| Light | $\leq 6g$ | N/A | Different |
| Regular | 6~9g | 6~9g | |
| Super | 9~12g | 9~12g | |
| Super plus | 12~15g | 12~15g | |
| Product dimensions | | | |
| Digital tampon | | | Different |
| Light | Pledget length | (40~50) mm | |
| | | | N/A |

| | | | | |
|-------------------|-----------------------|----------------|---------------------|-----------|
| | Pledget diameter | (9.5~12.5) mm | | |
| | Removal string length | (125~165) mm | | |
| Regular | Pledget length | (40~50) mm | (43.0~46.0) mm | |
| | Pledget diameter | (10.5~13.5) mm | (11.0~12.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| Super | Pledget length | (45~55) mm | (48.0~51.0) mm | |
| | Pledget diameter | (11.5~14.5) mm | (12.0~13.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| Super plus | Pledget length | (45~55) mm | (48.0~51.0) mm | |
| | Pledget diameter | (13.3~16.3) mm | (14.0~15.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| Applicator tampon | | | | |
| Light | Pledget length | (40~50) mm | N/A | |
| | Pledget diameter | (9.5~12.5) mm | | |
| | Removal string length | (125~165) mm | | |
| | Applicator length | 125 mm | | |
| | Applicator diameter | 15 mm | | |
| Regular | Pledget length | (40~50) mm | (43.0~46.0) mm | |
| | Pledget diameter | (10.5~13.5) mm | (11.0~12.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| | Applicator length | 125 mm | 120 mm | |
| | Applicator diameter | 15 mm | 13 mm | |
| Super | Pledget length | (45~55) mm | (43.0~46.0) mm | |
| | Pledget diameter | (11.5~14.5) mm | (12.0~13.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| | Applicator length | 126 mm | 120 mm | |
| | Applicator diameter | 16.5 mm | 16 mm | |
| Super plus | Pledget length | (45~55) mm | (43.0~46.0) mm | |
| | Pledget diameter | (13.3~16.3) mm | (14.0~15.0) mm | |
| | Removal string length | (125~165) mm | (130~160) mm | |
| | Applicator length | 126 mm | 120 mm | |
| | Applicator diameter | 16.5 mm | 18 mm | |
| Component | Pledget | 100% Viscose | 100% Organic Cotton | Different |

| | | | | |
|---|--------------------------------------|--|-------------------------------|-----------|
| Materials | Overwrap | Polyethylene and Polyethyleneterephthalate | 100% Organic Cotton | |
| | Removal string | Polyester and Viscose | 100% Organic Cotton | |
| | Applicator | Polyethylene and Polypropylene | TPO | |
| Additives and Finishing Agents | Anti-wicking agent of removal string | Paraffin-based resin | Paraffin-based resin | Same |
| | Finishing agent of pledget | Polyoxyethylene stearate | Not available | Different |
| Complies with ISO 10993-1 | | Yes | Yes | Same |
| Compliance with microbiology requirements of FDA Guidance for Tampons | | Yes | Yes | Same |
| Labeling | | Complied with 21 CFR part 801 | Complied with 21 CFR part 801 | Same |

The proposed device has the similar indication for use as the predicate device as well as comparable technical characteristics, and the differences don't raise any additional questions for safety and effectiveness. Therefore, the subject device has the same intended use as the predicate.

7. Summary of Non-clinical Testing

Non-clinical testing was conducted to verify that the proposed device can meet all design specifications and perform similarly to the predicate device. The following tests were conducted.

Performance Testing

The following performance characteristics were assessed in accordance with the FDA guidance document "*Guidance for Industry and FDA Staff - Menstrual Tampons and Pads: Information for Premarket Notification Submissions (510(k)s)*" issued on July 27, 2005.

- Dimensions
- Absorbency
- Removal string strength
- Fiber shedding
- Tampon integrity
- Chemical residues

All samples met the predefined acceptance criteria.

Biocompatibility Testing

Biocompatibility studies were performed in accordance with the FDA guidance document “*Guidance for Industry and Food and Drug Administration Staff - Use of International Standard ISO 10993-1, ‘Biological Evaluation of Medical Devices - Part 1: Evaluation and testing within a risk management process’*” issued on September 4, 2020 as follows:

- In vitro cytotoxicity test per ISO 10993-5:2009
- Skin sensitization test per ISO 10993-10:2010
- Vaginal irritation test per ISO 10993-10:2010
- Acute systemic toxicity test per ISO 10993-11:2017 (tampon only)

All tests were performed on the tampon and applicator separately. The results demonstrated that the proposed device is non-cytotoxic, non-irritating, non-sensitizing, and non-systemically toxic.

Microbiology Testing

Per the FDA guidance document “*Guidance for Industry and FDA Staff - Menstrual Tampons and Pads: Information for Premarket Notification Submissions (510(k)s)*” issued on July 27, 2005, the following microbiology testing was conducted on the final, finished form of the proposed device, and the test results showed that the proposed device does not:

- enhance the growth of *Staphylococcus aureus*;
- increase the production of Toxic Shock Syndrome Toxin-1 (TSST-1);
- alter the growth of normal vaginal microflora.

8. Clinical Testing

No clinical study is included in this submission.

9. Conclusions

The nonclinical tests demonstrate that the proposed device is as safe and effective, as the legally marketed device (K15117). Therefore, the subject device is substantially equivalent to the predicate.