



March 27, 2020

Diadent Group International
Kab Sun, Lee
Quality Assurance Manager
16, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu
Cheongu-si, 28161 KOREA

Re: K190091
Trade/Device Name: Diapaste
Regulation Number: 21 CFR 872.3820
Regulation Name: Root Canal Filling Resin
Regulatory Class: Class II
Product Code: KIF
Dated: February 19, 2020
Received: February 20, 2020

Dear Kab Sun, Lee:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for Srinivas "Nandu" Nandkumar, Ph.D.
Director
DHT1B: Division of Dental Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K190091

Device Name
DIAPASTE

Indications for Use (Describe)

Aqueous ointment material that temporarily fills the root canal for the following indications:

- Apexification
- Temporary root filling
- Root canal filling for primary teeth
- Vital pulpotomy
- Temporary pulp capping

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Diadent Group International

510(k) Summary

1 Application Information

Date Prepared:	March 26, 2020
Company Name and Address:	DiaDent Group International 16, Osongsaengmyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, Republic of Korea
Contact Person:	Kab Sun, Lee Quality Assurance Manager Phone: +82-43-266-2315 FAX: +82-43-235-2315 Email: diadent32@diadent.co.kr

2 Device Information

Device Type:	Root Canal Filling Resin
Regulation Description:	Root Canal Filling Resin
Review Panel:	Dental
Regulation Number:	21 CFR 872.3820
Product Code:	KIF
Device Class:	II
Device Name:	DIAPASTE

3 Predicate Devices

The legally marketed devices to which substantial equivalence is being claimed are:

	Primary Predicate Device	Reference Predicate Device
510(k) Number:	K032605	K060365
Applicant:	Meta Biomed Co., Ltd.	Ivoclar Vivadent
Device Name:	Metapaste	ApexCal
Regulation Number:	21 CFR 872.3820	21 CFR 872.3250
Product Code:	KIF	EJK
Device Class:	II	II

4 Device Configuration

Each model configurations of subject device are described as following:

Model Name	Contents
DiaPaste Type A	1 syringe(2g), 20 disposable tips, 1silicone cap
DiaPaste Type B	1 syringe(2g), 10 disposable tips, 1silicone cap
DiaPaste Type C	1 syringe(2g), 4 disposable tips, 1silicone cap
DiaPaste Refill Kit	1 syringe(2g)
DiaPaste Intro Kit	1 syringe(0.5g), 4 disposable tips 1silicone cap

5 Device Description

Diapaste is a calcium hydroxide paste with barium sulfate, used as a temporary root canal filling material.

Diadent Group International

6 Intended Use/Indications for Use

Aqueous ointment material that temporarily fills the root canal for the following indications:

- Apexification
- Temporary root filling
- Root canal filling for primary teeth
- Vital pulpotomy
- Temporary pulp capping

7 Comparison Table and Discussion

This device compares to the legally marketed devices as follows:

	Subject Device	Primary Predicate Device	Reference Predicate Device	
	DiaPaste	Metapaste	ApexCal	Discussion
Regulation Number:	21 CFR 872.3820	21 CFR 872.3820	21 CFR 872.3250	
510(k) Number	K190091	K032605	K060365	-
Indication for Use	Aqueous ointment material that temporarily fills the root canal for the following indications: <ul style="list-style-type: none"> - Apexification - Temporary root filling - Root canal filling for primary teeth - Vital pulpotomy - Temporary pulp capping 	Metapaste is a biocompatible root canal sealer used for the temporary filling of root canals after endodontic surgery. Metapaste can be used on its own and for vital pulpectomies in deciduous teeth. Metapaste is intended for use by qualified healthcare personnel trained in its use.	ApexCAL is a Calcium Hydroxide paste that has a creamy consistency and is suitable for several indications including: <ul style="list-style-type: none"> * Temporary disinfectant dressings in the obturation of root canals; * Indirect pulp capping or management of deep caries lesions; or * Direct pulp capping. 	
Composition	<ul style="list-style-type: none"> - Calcium hydroxide - Calcium carbonate - Zinc oxide - Barium sulfate - Water - Polyethylene glycol 400 - Polysorbate (Tween80) 	Calcium Hydroxide Barium Sulfate Poly propylene glycol	<ul style="list-style-type: none"> - Calcium hydroxide - Bismuth carbonate - Zinc oxide - Polyethylene glycol - glycerin - Water 	
Period of Use	Temporary (remains in the body for 29 days or less)	Temporary	Temporary	equivalent
Physical properties	Conformed to ISO 6876 <ul style="list-style-type: none"> - Flowability - Film Thickness - Radio-opacity 	Conformed to ISO 6876		equivalent
Biocompatibility	Biocompatible <ul style="list-style-type: none"> - Cytotoxicity - Sensitization - Acute Systemic toxicity - Oral Mucosa irritation - Genotoxicity- Bacterial Reverse Mutation - Genotoxicity- Micronucleus test 	Biocompatible Freedom from toxicity per ISO/TR 7405 Agar diffusion test	Biocompatible	equivalent
Package Contents	<ul style="list-style-type: none"> •Syringe •Disposable Tip •Silicone Cap 	<ul style="list-style-type: none"> •2.2g paste in a syringe •Disposable tips •One ring rotator for the direction control of the tip 	<ul style="list-style-type: none"> • 2 syringes 2.5g each • 15 application tips 	-

Diadent Group International

Differences

- Indication For Use

DiaPaste (K190091)	Primary Predicate Device (Metapaste, K032605)	Discussion
<p>Aqueous ointment material that temporarily fills the root canal for the following indications:</p> <ul style="list-style-type: none"> - Apexification - Temporary root filling - Root canal filling for primary teeth • Vital pulpotomy - Temporary pulp capping 	<p>Metapaste is a biocompatible root canal sealer used for the temporary filling of root canals after endodontic surgery.</p> <p>Metapaste can be used on its own and for vital pulpectomies in deciduous teeth.</p> <p>Metapaste is intended for use by qualified healthcare personnel trained in its use.</p>	<p>Diapaste and Primary predicate device (Metapaste) can be used for apexification & pulpectomy technique for root canal treatment as temporary root canal filling materials. The Indication for Use of Subject and predicate devices has same contents. In the DiaPaste's Indication for Use, there are contents of Pulpotomy and Temporary pulp capping that is not found in Indication of use of Metapaste. It is a technique commonly used of temporary root canal fillings and as such does not affect the equivalence of Indication of subject and predicate devices.</p>

-Composition

Subject Device	Primary Predicate Device (K032605)	Reference Predicate Device(K060365)	Discussion
<ul style="list-style-type: none"> - Calcium hydroxide - Calcium carbonate - Zinc oxide - Barium sulfate (Radiopaque) - Water - Polyethylene glycol 400 - Polysorbate (Tween80) 	<p>Calcium Hydroxide Barium Sulfate (Radiopaque) Poly propylene glycol</p>	<ul style="list-style-type: none"> - Calcium hydroxide - Bismuth carbonate - Zinc oxide - Polyethylene glycol - glycerin - Water 	<p>The main ingredient of subject and predicate devices is Calcium Hydroxide. and the functions of other additives are similar (base, Radiopaque, pH controller, Stabilizer). Also, the biocompatibility of subject device was confirmed by biological safety study.</p>

Diadent Group International**8. Non-Clinical performance data**

This device has demonstrated conformance with non-clinical performance requirements through evaluation and testing in accordance with the following harmonized standards:

-ISO 6876	Root Canal Sealing materials
-ISO 7405	Dentistry - Evaluation Of Biocompatibility Of Medical Devices Used In Dentistry [Including: Amendment 1 (2013)]
-ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
-ISO 10993-3	Biological Evaluation of Medical Devices – Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity
-ISO 10993-5	Biological evaluation of medical devices - Part 5. Tests for in vitro Cytotoxicity
-ISO 10993-10	Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization
-ISO 10993-11	10993-11 – Biological Evaluation of Medical Devices – Part 11: Tests for Systemic Toxicity

The following Performance tests were conducted.

- Flowability
- Film thickness
- Radiopacity
- Biocompatibility test

Performance test (Physical properties):

Test	Requirement	Result	Pass/Fail
Flowability	Each disc to have a diameter of not less than 17mm	21mm (mean value)	Pass
Film thickness	No more than 50µm	45µm (mean value)	Pass
Radiopacity	Equivalent to not less than 3mm of Aluminum	Not less than 3mm of Aluminum (approx. 5.3mm)	Pass

9. Conclusion

Based on the above information and all data provided in this submission, including comparison of intended uses, technological characteristics discussion of differences of subject and predicate devices, show that the subject device and the predicate device have similar technical characteristic and chemical composition.

It is demonstrated that the subject device and the legally marketed devices identified in this submission are substantially equivalent.