



December 13, 2021

Svelte Medical Systems, Inc.
Mark Kielek
Vice President Quality Assurance and Regulatory Affairs
675 Central Avenue, Suite 2
New Providence, New Jersey 07974

Re: P210014

Trade/Device Name: SLENDER Sirolimus-Eluting Coronary Stent Integrated Delivery System and
DIRECT Sirolimus-Eluting Coronary Stent Rapid Exchange Delivery System

Product Code: NIQ

Filed: April 1, 2021

Amended: September 14, 2021

Dear Mark Kielek:

The Center for Devices and Radiological Health (CDRH) of the Food and Drug Administration (FDA) has completed its review of your premarket approval application (PMA) for the SLENDER Sirolimus-Eluting Coronary Stent Integrated Delivery System and DIRECT Sirolimus-Eluting Coronary Stent Rapid Exchange Delivery System. The SLENDER Sirolimus-Eluting Coronary Stent Integrated Delivery System is indicated for improving coronary artery luminal diameter in patients with symptomatic heart disease due to atherosclerotic lesions ≤ 24 mm in length in native coronary arteries with ≥ 2.25 mm to ≤ 4.00 mm reference vessel diameters, using direct stenting or pre-dilatation interventional techniques. The DIRECT Sirolimus-Eluting Coronary Stent Rapid Exchange Delivery System is indicated for improving coronary artery luminal diameter in patients with symptomatic heart disease due to atherosclerotic lesions ≤ 34 mm in length in native coronary arteries with ≥ 2.25 mm to ≤ 4.00 mm reference vessel diameters, using direct stenting or pre-dilatation interventional techniques. We are pleased to inform you that the PMA is approved. You may begin commercial distribution of the device in accordance with the conditions of approval described below. Although this letter refers to your product as a device, please be aware that some approved products may instead be combination products. The Premarket Approval Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfPMA/pma.cfm> identifies combination product submissions.

The sale and distribution of this device are restricted to prescription use in accordance with 21 CFR 801.109 and under section 515(d)(1)(B)(ii) of the Federal Food, Drug, and Cosmetic Act (the act). The device is further restricted under section 515(d)(1)(B)(ii) of the act insofar as the labeling must specify the specific training or experience practitioners need in order to use the device. FDA has determined that these restrictions on sale and distribution are necessary to provide reasonable assurance of the safety and effectiveness of the device. Your device is therefore a restricted device subject to the requirements in sections 502(q) and (r) of the act, in addition to the many other FDA requirements governing the manufacture, distribution, and marketing of devices.

Expiration dating for this device has been established and approved at 24 months at controlled room temperature.

Continued approval of the PMA is contingent upon the submission of periodic reports, required under 21 CFR 814.84, at intervals of one year (unless otherwise specified) from the date of approval of the original PMA. This report, identified as "Annual Report" and bearing the applicable PMA reference number, should be submitted to the address below. The Annual Report should indicate the beginning and ending date of the period covered by the report and should include the information required by 21 CFR 814.84.

In addition to the above, and in order to provide continued reasonable assurance of the safety and effectiveness of the PMA device, the Annual Report must include, separately for each model number (if applicable), the number of devices sold and distributed during the reporting period, including those distributed to distributors. The distribution data will serve as a denominator and provide necessary context for FDA to ascertain the frequency and prevalence of adverse events, as FDA evaluates the continued safety and effectiveness of the device.

You have agreed to provide the following non-clinical information in a separate report, annually, which may be followed by a PMA supplement where applicable.

1. Long-term drug stability studies will be completed on five finished product batches representing the commercial process each year, with product codes rotating in a 4 year cycle, as agreed upon in your response to FDA's Major Deficiency Letter issued June 30, 2021. All batches for these studies will be stored at Long Term Conditions of $25^{\circ}\text{C} \pm 2^{\circ}\text{C}/60\% \text{RH} \pm 5\%$, per ICH Q1A(R2). Testing for all studies will occur at 0, 3, 6, 9, 12, 18, 24, and 36 months.

Be advised that failure to comply with any post-approval requirement, including test protocol, sampling size, sampling plan, and acceptance criteria, constitutes grounds for FDA withdrawal of approval of the PMA in accordance with 21 CFR 814.82(c) and 814.46(a)(2).

In addition to the Annual Report requirements, you must provide the following data in post-approval study (PAS) reports for each PAS listed below.

You must obtain approval of your PAS protocol(s) within 60 days from the date of this order. Within 30 days of your receipt of this letter, you must submit a PMA supplement that includes a complete protocol of your post-approval study described below. Your PMA supplement should be clearly labeled as a "PMA Post-Approval Study Protocol" as noted below and submitted to the address below. Please reference the PMA number above to facilitate processing. If there are multiple protocols being finalized after PMA approval, please submit each protocol as a separate PMA supplement.

2. PMA Post-Approval Study – Continued Follow-Up of OPTIMIZE Clinical Study. The OPTIMIZE Clinical Study (G160227/S004) is a single-blind, randomized, active-control, multi-center clinical study which enrolled 1,645 subjects. The OPTIMIZE Clinical Study was designed to compare the safety and efficacy of the Svelte Sirolimus-Eluting Coronary Stent Integrated Delivery System (Svelte DES-IDS) and Svelte Sirolimus-Eluting Coronary Stent Rapid Exchange Delivery System (Svelte DES-RX) to the commercially available Abbott Vascular XIENCE or Boston Scientific Promus Drug-Eluting Coronary Stents (control DES) through 5 years post-index procedure. The

primary endpoint is Target Lesion Failure (TLF) at 12 months post-procedure, defined as cardiac death, Target Vessel Myocardial Infarction (TVMI, including Q wave and non-Q wave) or clinically-driven Target Lesion Revascularization (TLR) by percutaneous or surgical methods. You must collect and report clinical outcomes to FDA through 5 years post-procedure on patients enrolled in the OPTIMIZE Clinical Study.

3. PMA Post-Approval Study – Svelte Post-Approval Study. The Svelte PAS is a prospective, multicenter, non-randomized study intended to monitor and evaluate the safety and efficacy outcomes of the Svelte DES post-PMA approval of the SLENDER IDS and DIRECT RX systems in a real world setting. The study will enroll approximately 500 subjects with coronary artery disease (CAD) at up to 50 clinical sites, with at least one-half of clinical sites and study subjects from the United States. The primary endpoint for all study subjects enrolled in the Svelte PAS is percentage of subjects with Target Lesion Failure (TLF) at 12 months post-procedure, defined as cardiac death, non-fatal Target Vessel Myocardial Infarction (TVMI, including Q wave and non-Q wave) or clinically-driven Target Lesion Revascularization (TLR) by percutaneous or surgical methods. All subjects will require pre-procedure enrollment to ensure pre- and post-procedural cardiac biomarkers are collected. A central core lab will be used in the assessment of cardiac biomarkers. Follow-up contacts post-procedure will be made for clinical assessment at 30 days and 1, 2 and 3 years. There are additionally several secondary endpoints. You must collect and report clinical outcomes to FDA through at least 3 years post-procedure on patients enrolled in the Svelte PAS.

From the time of Svelte PAS study protocol approval, you must meet the following timelines for

- First subject enrolled within 6 months
- 20% of subjects enrolled within 12 months
- 50% of subjects enrolled within 18 months
- 100% of subjects enrolled within 24 months
- Submission of Final study report: 3 months from study completion (i.e. last subject, last follow-up date)

In addition, you must submit separate periodic reports on the progress of the Svelte PAS as follows:

- PAS Progress Reports every six (6) months until subject enrollment has been completed, and annually thereafter.
- If any enrollment milestones are not met, you must begin submitting quarterly enrollment status reports (i.e., every 3 months), in addition to your periodic (6-months) PAS Progress Reports, until FDA notifies you otherwise.

For all other condition of approval studies, you must submit separate PAS Progress Reports for each study, every six (6) months for the first two (years) and annually thereafter, unless otherwise specified by FDA.

Each PAS report should be submitted to the address below identified as a "PMA Post-Approval Study Report" in accordance with how the study is identified above and bearing the applicable PMA reference number.

Be advised that failure to comply with any post-approval requirement, including initiation, enrollment, periodic reporting, pre- and post-procedure cardiac biomarker draws, and use of a central core lab, constitutes grounds for FDA withdrawal of approval of the PMA in accordance with 21 CFR 814.82(c) and 814.46(a)(2).

Be advised that the failure to conduct any such study in compliance with the good clinical laboratory practices in 21 CFR part 58 (if a non-clinical study subject to part 58) or the institutional review board regulations in 21 CFR part 56 and the informed consent regulations in 21 CFR part 50 (if a clinical study involving human subjects) may be grounds for FDA withdrawal of approval of the PMA in accordance with 21 CFR 814.46(a)(3)-(4).

Be advised that protocol information, interim and final results will be published on the Post Approval Study Webpage https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfPMA/pma_pas.cfm.

In addition, the results from any post approval study should be included in the labeling as these data become available. Any updated labeling must be submitted to FDA in the form of a PMA Supplement. For more information on post-approval studies, see the FDA guidance document entitled, "Procedures for Handling Post-Approval Studies Imposed by PMA Order" (<https://www.fda.gov/media/71327/download>).

This is a reminder that as of September 24, 2014, class III devices are subject to certain provisions of the final Unique Device Identification (UDI) rule. These provisions include the requirement to provide a UDI on the device label and packages (21 CFR 801.20), format dates on the device label in accordance with 21 CFR 801.18, and submit data to the Global Unique Device Identification Database (GUDID) (21 CFR 830 Subpart E). Additionally, 21 CFR 814.84 (b)(4) requires PMA annual reports submitted after September 24, 2014, to identify each device identifier currently in use for the subject device, and the device identifiers for devices that have been discontinued since the previous periodic report. It is not necessary to identify any device identifier discontinued prior to December 23, 2013. Combination Products may also be subject to UDI requirements (see 21 CFR 801.30). For more information on these requirements, please see the UDI website, <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-udi-system>.

Before making any change affecting the safety or effectiveness of the PMA device, you must submit a PMA supplement or an alternate submission (30-day notice) in accordance with 21 CFR 814.39. All PMA supplements and alternate submissions (30-day notice) must comply with the applicable requirements in 21 CFR 814.39. For more information, please refer to the FDA guidance document entitled, "Modifications to Devices Subject to Premarket Approval (PMA) - The PMA Supplement Decision-Making Process" <https://www.fda.gov/media/81431/download>.

You are reminded that many FDA requirements govern the manufacture, distribution, and marketing of devices. For example, in accordance with the Medical Device Reporting (MDR) regulation, 21 CFR 803.50 and 21 CFR 803.52 for devices or post-marketing safety reporting (21 CFR 4, Subpart B) for combination products, you are required to report adverse events for this device. Manufacturers of medical devices, including in vitro diagnostic devices, are required to report to FDA no later than 30 calendar days after the day they receive or otherwise becomes aware of information, from any source, that reasonably suggests that one of their marketed devices:

1. May have caused or contributed to a death or serious injury; or
2. Has malfunctioned and such device or similar device marketed by the manufacturer would be likely to cause or contribute to a death or serious injury if the malfunction were to recur.

Additional information on MDR, including how, when, and where to report, is available at <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems> and on combination product post-marketing safety reporting is available at (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>).

In accordance with the recall requirements specified in 21 CFR 806.10 for devices or the post-marketing safety reporting requirements (21 CFR 4, Subpart B) for combination products, you are required to submit a written report to FDA of any correction or removal of this device initiated by you to: (1) reduce a risk to health posed by the device; or (2) remedy a violation of the act caused by the device which may present a risk to health, with certain exceptions specified in 21 CFR 806.10(a)(2). Additional information on recalls is available at

<https://www.fda.gov/safety/recalls-market-withdrawals-safety-alerts/industry-guidance-recalls>.

CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading. CDRH will notify the public of its decision to approve your PMA by making available, among other information, a summary of the safety and effectiveness data upon which the approval is based. The information can be found on the FDA CDRH Internet Home Page located at

<https://www.fda.gov/medical-devices/device-approvals-denials-and-clearances/pma-approvals>. Written requests for this information can also be made to the Food and Drug Administration, Dockets Management Branch, (HFA-305), 5630 Fishers Lane, Rm. 1061, Rockville, MD 20852. The written request should include the PMA number or docket number. Within 30 days from the date that this information is placed on the Internet, any interested person may seek review of this decision by submitting a petition for review under section 515(g) of the act and requesting either a hearing or review by an independent advisory committee. FDA may, for good cause, extend this 30-day filing period.

Failure to comply with any post-approval requirement constitutes a ground for withdrawal of approval of a PMA. The introduction or delivery for introduction into interstate commerce of a device that is not in compliance with its conditions of approval is a violation of law.

You are reminded that, as soon as possible and before commercial distribution of your device, you must submit an amendment to this PMA submission with a copy of all final labeling. Final labeling that is identical to the labeling approved in draft form will not routinely be reviewed by FDA staff when accompanied by a cover letter stating that the final labeling is identical to the labeling approved in draft form. If the final labeling is not identical, any changes from the final draft labeling should be highlighted and explained in the amendment.

All required documents should be submitted, unless otherwise specified, to the address below and should reference the above PMA number to facilitate processing.

U.S. Food and Drug Administration
Center for Devices and Radiological Health
Document Control Center - WO66-G609
10903 New Hampshire Avenue
Silver Spring, MD 20993-0002

If you have any questions concerning this approval order, please contact Chelsea Virgile, PhD at 301-796-3125 or Chelsea.Virgile@fda.hhs.gov.

Sincerely,

Brian D. Pullin -S

Brian Pullin

Director

DHT2C: Division of Coronary

and Peripheral Intervention Devices

OHT2: Office of Cardiovascular Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health