



February 20, 2020

de Götzen S.r.l
% Dario Bandiera
Quality Manager
via Roma, 45
21057 Olgiate Olona, Varese
ITALY

Re: K191719
Trade/Device Name: X-MIND Trium
Regulation Number: 21 CFR 892.1750
Regulation Name: Computed tomography x-ray system
Regulatory Class: Class II
Product Code: OAS, MUH
Dated: January 10, 2020
Received: January 22, 2020

Dear Dario Bandiera:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For

Thalia T. Mills, Ph.D.
Director
Division of Radiological Health
OHT7: Office of In Vitro Diagnostics
and Radiological Health
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K191719

Device Name
X-MIND trium

Indications for Use (Describe)

Indications for Use (Describe)

X-MIND trium is a digital panoramic, cephalometric and tomographic extra-oral X-ray system, indicated for use in:

- producing panoramic X-ray images for diagnostic examination of dentition (teeth), jaws and oral structures;
- producing radiographs of maxillofacial region and parts of the skull for cephalometric examination, if equipped with CEPH arm;
- producing radiographs of hands and wrists for carpus examination, if equipped with CEPH arm;
- producing tomographic images of the oral and maxillofacial region, for diagnostic examination of dentition (teeth), jaws, oral structures and some cranial bones, if equipped with CBCT option.

From a clinical point of view, X-MIND trium can be applied for the following medical indications:

- Generic dentistry
- Dental implantology
- Dental surgery
- Maxillo-facial surgery
- Cephalometric analysis
- Carpus radiology

The target patient population includes adults and pediatric patients from 5 years old [~21 kg (46 lb); 113 cm (44.5 in) standing height]; anyway the sustainability to X-ray exposure must be evaluated by surgeons, dentists and qualified and authorized physicians.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.


This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov


"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | K191719 |

Index

| | |
|--|----|
| 510 (k) Summary | 2 |
| I. SUBMITTER..... | 2 |
| II. DEVICE..... | 2 |
| III. PREDICATE DEVICE | 3 |
| IV. DEVICE OVERVIEW..... | 3 |
| V. MAIN CHANGES TO THE PREVIOUS CLEARED DEVICE ARE | 4 |
| VI. OVERVIEW OF THE SYSTEM..... | 8 |
| VII. OVERVIEW OF THE SCANNER..... | 10 |
| VIII. EQUIPMENT CONFIGURATIONS | 13 |
| IX. EXAMS | 15 |
| X. Setup configuration | 15 |
| Wall mounted | 15 |
| Self-standing..... | 17 |
| XI. ACCESSORIES | 18 |
| XIII. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE...21 | |
| XIV. PERFORMANCE, SAFETY, EMC TESTS | 55 |
| XV. CONCLUSIONS..... | 57 |

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 1 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | K191719 |

510 (k) Summary

The summary of this 510(k) is being submitted in accordance with the requirements of 21 CFR Part 807.92.

I. SUBMITTER

| | |
|----------------|--|
| Owner's name: | de Götzen S.r.l. – ACTEON Group |
| Address: | via Roma, 45 – 21057 Olgiate Olona (VA), Italy |
| Tel. | +39 0331 376760 |
| Fax | +39 0331 376763 |
| Dario Bandiera | dario.bandiera@acteongroup.com |
| Alvise Reither | alvise.reither@acteongroup.com |
| Date: | 10/01/2020 |


Table 1: Submitter information

II. DEVICE

| | |
|-----------------------|---|
| Name of the device: | X-MIND trium |
| Common or Usual name: | Dental panoramic, cephalometric and CBCT X-ray system |
| Classification name: | Computed tomography X-ray system (21 CFR 892.1750) |
| Regulatory class: | II |
| Product Code: | OAS, MUH |

Table 2: device

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 2 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

III. PREDICATE DEVICE

Legally marketed device to which equivalence is claimed is:

| PREDICATE DEVICE | |
|-------------------------|----------------------------------|
| Device name | X-MIND trium |
| Manufacturer | de Götzen S.r.l. – ACTEON Group |
| Device product code | OAS |
| Subsequent product code | MUH |
| Regulation number | 892.1750 |
| Regulation name | Computed tomography x-ray system |
| Clearance date | November 15, 2016 |
| 510(k) number | K160166 |

Table 3: Predicate device information

This predicate has not been subjected to a design-related recall.

No reference devices were used in this submission.

IV. DEVICE OVERVIEW

The modified X-MIND trium has the same indication for use of X-MIND trium (as identified in the labelling) and it is based on the same fundamental scientific technology.


It is a digital panoramic, cephalometric and tomographic extra-oral X-ray system, indicated for use in:

- producing panoramic X-ray images for diagnostic examination of dentition (teeth), jaws and oral structures;
- producing radiographs of maxillofacial region and parts of the skull for cephalometric examination, if equipped with CEPH arm;
- producing radiographs of hands and wrists for carpus examination, if equipped with CEPH arm;
- producing tomographic images of the oral and maxillofacial region, for diagnostic examination of dentition (teeth), jaws, oral structures and some cranial bones, if equipped with CBCT option.

From a clinical point of view, X-MIND trium can be applied for the following medical indications:

- Generic dentistry
- Dental implantology

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 3 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

- Dental surgery
- Maxillo-facial surgery
- Cephalometric analysis
- Carpus radiology

The target patient population includes adults and pediatric patients from 5 years old [~21 kg (46 lb); 113 cm (44.5 in) standing height]; anyway, the sustainability to X-ray exposure must be evaluated by surgeons, dentists and qualified and authorized physicians.

V. MAIN CHANGES TO THE PREVIOUS CLEARED DEVICE ARE

Main changes to the previous cleared device are addition of a new accessory: new sliding cassette called X-MIND trium TLD cassette including a new model of CBCT sensor (C12903D-40 - Hamamatsu) that can be used for both panoramic and CBCT examination; introduction, as alternative components, of CEPH and PAN sensors with wider admitted range of power supply value: (from 4.9 – 5.1 VDC to 8 – 12 VDC).

Since CBCT and PAN sensor usually operates on a different SID (Source to Image Distance), a new movement of the cassette has been added, to give to the device the possibility to modify the SID.

The possibility to configure the SID allows also the possibility to perform special exam for child: reduced SID distance that means possible reduction of x-ray dose to the patient.





Figure 1 addition of TLD Sliding cassette

The lateral movement of the cassette is obtained with the same motor used in the cleared version of X-MIND trium.

For the movement we use a new linear guide with ball bearing runners.

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 4 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

The SID movements uses the same motor as the one installed in the handle rest area of the patient (satellite group), and two linear guides with ball bearing runners as the cleared version of TRIUM

The cassette is made of a steel plate with the new CBCT sensor and relevant shielding lead on the back.




Figure 2: new TLD cassette scheme

The new CBCT sensor has nominal dimensions of 15x12 cm, while the sensor used in cleared version of X-MIND trium is 12x12.

The 15 cm height is required to allow panoramic examinations.

The new sensor C12903D-40 has 8V power supply input while the previous was 5V. Similarly, the new SID movements require a related driver. For that reason, a new AUXCBCT board has been designed with code: PE5541_AUXCBCTPLUS.

The overall insulation diagram of the whole device has not changed in respect of the cleared X-MIND trium, apart from the absence of the contacts for PAN and CEPH sensors that aren't anymore detachable (improvement)

The relevant insulation of the board PE5541_AUXCBCTPLUS_A01 is the same as the PE4016_AUXCBCT_A04, already installed (and certified) on the V1 version of the TRIUM.

Due to the possibility to reduce the SID, special exams for child has been introduced and at the same time we have developed countermeasures to avoid that an adult is positioned during a special child exam.

The new workflow on the software introduces a new step: the doctor uses a tool provided with the machine to check if the user fits the proper dimension for the examination with the reduced SID.

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 5 of 58 |


| | | |
|---|--------------------------------|---------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 |
| | | 002_k Summary |



Figure 3: use of head measurement tool


X-MIND trium and the X-MIND trium TLD are provided with a new functionality, the Low Dose functionality.

The Low Dose functionality allows the user to perform the CBCT exams by using the Low Dose protocols for all kind of patients (MAN, WOMAN and CHILD) as an alternative of the Standard Dose protocols.

The selection of a Low Dose protocol is performed in the AIS interface during the specification of the exam settings (see Figure 4).

A Low Dose protocol applies at least the half of the anodic current than the corresponding Standard Dose protocol, so ensuring a radiation dose saving for the patient of about 50% (refer to the declared DAP in the AIS interface, Figure 5); the image quality by the Low Dose protocols ensures the same diagnostic value as the corresponding Standard Dose protocols (Basing on the bench test results on image quality and radiation dose obtained).

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 6 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |





Figure 4 - AIS interface for the definition of the Low Dose protocol

| Acquisition Parameters | Standard | Low Dose |
|----------------------------|----------|----------|
| Current (mA) | 8 | 4 |
| Voltage (KV) | 90 | 45 |
| Current Time Product (mAs) | 55.9 | 27.9 |
| DAP (mGycm ²) | 902.7 | 451.35 |

Figure 5 - Acquisition parameters and dose declaration between Standard Dose and Low Dose protocols

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 7 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

VI. OVERVIEW OF THE SYSTEM








Figure 6: overview of the system

| | | |
|--------|---|---|
| Part 1 |  | - 1 Gb, 8-wire Ethernet cable cat. 6>250mhz, 27 AWG, 10 m long - Operator's Workstation |
| Part 2 |  | - Control cable, 10 wires, 5 twisted pairs 10 m long aluminium shielded - X-MIND trium Remote Control - X-MIND trium Light |
| Part 3 |  | - Power cord, 3G1,5 AWG 16, 300 / 500 V, 10 m long -Connection to the mains by means of clamps and screws (Equipment installed permanently in compliance with IEC 60601-1) |

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 8 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| | | |
|--------|---|---|
| Part 4 |  | <ul style="list-style-type: none"> - X-MIND trium main connector consisting of 3 modules Communication, control and power supply - Moving part of the X-MIND trium connector - Fixed part of the X-MIND trium connector |
|--------|---|---|

Table 4: details of the system

The system consists of:

- X-MIND trium SCANNER
- OPERATOR'S WORKSTATION

The operator can control all aspects of the medical device by means of the workstation that can be supplied by de Götzen.

The workstation must have installed with the following modules:

- **AIS** equipment management software + **AIS 2DApp** diagnostic analysis
- **AIS 3DApp** dataset display software CBCT

Communication with the medical device occurs by means of a sturdy and reliable Ethernet protocol.

The workstation allows the operator to perform the following procedures:

- Calibrations of the medical device
- Acquisition parameters setting
- Getting an acquisition
- Image visualization and post processing
- Database management
- Periodic quality Tests


- X-MIND trium REMOTE CONTROL

The X-MIND trium remote control must be in a safe place protected against radiations, in compliance with the local standards in force concerning ionising radiation protection.

The X-MIND trium remote control allows the operator to activate or deactivate X-ray emission from the control room. This consists of two switches, one for exposure and one for emergency, which control device operation.

- X-MIND trium Light

| | | |
|------------|----------------------|--------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 9 of 58 |

| | | |
|---|--|-----------------------------------|
|  | <p style="text-align: center;">X-MIND trium 510(k) Summary</p> | <p>VOL. 005 002_k Summary</p> |
| <p>Via Roma, 45 21057 Olgiate Olona (VA)</p> | | |

The X-MIND trium Light is an indicator light that warns that X-ray emission is in progress.

VII. OVERVIEW OF THE SCANNER





Figure 7: overview of the scanner


| | | |
|-------------------|-----------------------------|----------------------|
| <p>Rev. n. 2</p> | <p>X-MIND trium 510 (k)</p> | <p>VOL. 005</p> |
| <p>10/01/2020</p> | | <p>Page 10 of 58</p> |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

The X-MIND trium scanner consists of the following parts:

| | | |
|---|------------------------|--|
| 1 | Control panel | The control panel provides an intuitive overview of the system to move the mobile column, turn the positioning lasers on and off and activate the X-MIND trium medical device. |
| 2 | X-ray generator | The X-ray assembly is the source of the X-ray beam during the rotation of the U-Arm. The pulsed beam is modelled by an automatic collimator, whereas the electronic control ensures stability and accuracy of the selected loading factors (exposure time, kVp and anodic current). An aluminium filter is used to harden the beam and remove low-energy ionising radiations, thereby obtaining suitable radiation quality whilst reducing its dose absorbed by the patient. |
| 3 | Sliding body | The sliding body is the mobile part of the column that supports U-arm. |
| 4 | U-arm | The U-Arm supports the Detector sliding group and the X-ray generator. This is the rotating part of the medical device, which moves around the patient during the image acquisition phase. |
| 5 | Detector sliding group | It contains detectors that allow acquire images. This cassette can be provided in the following versions: <ul style="list-style-type: none"> • PAN Only cassette; • Standard PAN/ CBCT Cassette; • TRUE LOW DOSE cassette that permit to perform both PAN and CBCT exams. (* only with X-MIND trium TLD version) |
| 6 | CBCT detector | This flat panel detector is indicated for use in generating radiographic images of the maxillo-facial region, more specifically it is dedicated for CBCT acquisitions. in X-MIND trium TLD version there is only one image detector, for PAN and CBCT exam; the equipment automatically exposes in front of X-Ray source the right portion of detector, depending on the selected exam; |
| 7 | PAN detector | This flat panel detector is indicated for use in generating radiographic images of the maxillo-facial region, more specifically it is dedicated for PAN acquisitions. |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 11 of 58 |


| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| | | |
|----|-----------------|--|
| | | * in X-MIND trium TLD version there is only one image detector, for PAN and CBCT exam; the equipment automatically exposes in front of X-Ray source the right portion of detector, depending on the selected exam; |
| 7 | TLD CASSETTE* | * only for X-MIND trium TLD version Sliding cassette that lets to activate true low dose exams. Some models of X-MIND trium can be upgraded in the field to TLD cassette configuration. |
| 8 | Patient support | The patient support allows stabilising and immobilising the patient. It is equipped with instruments to move the patient's head in all directions and find the best superimposition between the patient's anatomy and the Field of View (FOV). |
| 9 | F group | It is the whole assembled mobile group of the device. It is the moving part of the medical device, which adapts the acquisition geometry to the patient's anatomy and stance (sitting or standing). It supports the U-Arm and head support. |
| 10 | Column | The fixed column supports the entire structure of the medical device. This contains the motor that raises and lowers the mobile column, the main control panel of the system and the main switch. In the wall installation version, the fixed column is equipped with an anchor plate; in the free-standing version, on the other hand, the medical device is equipped with a footrest. <i>CEPH</i> arm extension (not represented in Overview of the medical device) |

Table 5: scanner parts

The ceph arm extension can be positioned both on the right or left side of the vertical column, with the patient positioning aiming device (ear and nasion rest) for CEPH exams;

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 12 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

the CEPH secondary collimator, translating during X Rays (to stay aligned with the X-Ray beam emerging from the tubehead on overhead carriage);
the CEPH image detector, translating during X Rays (to stay aligned with the X-Ray beam emerging from the secondary collimator on the CEPH arm itself)

VIII. EQUIPMENT CONFIGURATIONS

X-MIND trium equipment can be sold in these different configurations:

PAN only

the equipment can carry out uniquely exams of the PAN group (PAN, TMJ, Sinus);
the image detector is not movable from PAN bay;
all the exams are carried out with the same extension in height (146mm height on image receptor), primary collimator being of fixed dimension;
the equipment can be upgraded in the field to more advanced configurations


PAN / CBCT

the equipment can carry out uniquely exams of the PAN group (PAN, TMJ, Sinus) and of the CBCT group;
there are two image detectors, one for PAN and one for CBCT exam; the equipment automatically exposes in front of X-Ray source the right detector, depending on the selected exam;
the PAN exams can be carried out with different extension in height (up to 146mm height on image receptor), primary collimator being automatically adjusted in height, depending on exam selected by the user;
the CBCT exams can be carried out at different FOV dimensions;
the equipment can be upgraded in the field to PAN/CEPH/CBCT configuration by simple procedure, adding CEPH arm with CEPH image detector replacing the existing PAN image detector;
the equipment has a touch screen panel as main user interface;

PAN / CEPH

the equipment can carry out uniquely exams of the PAN group (PAN, TMJ, Sinus) and of the CEPH group (AP/PA, LL, carpus);
the image detector is unique, both for PAN and CEPH exams, and must be manually moved from PAN bay to CEPH bay and viceversa, depending on the selected exam;
the exams can be carried out with different extension in height (up to 146mm height on image receptor for PAN, up to 220mm height on image receptor for CEPH), primary collimator being automatically adjusted in height depending on exam selected by the user;

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 13 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

the equipment can be upgraded in the field to PAN/CBCT/CEPH configuration by simple procedure, adding CBCT image detector and bite block holder with vertical movement;

the equipment mounts PAN detector sliding unity.

PAN / CBCT / CEPH (PAN Full)

the equipment can carry out all feasible exams of the PAN group (PAN, TMJ, Sinus), of the CEPH group (AP/PA, LL, carpus) and of the CBCT group;

there are two image detectors, one for PAN/CEPH and one for CBCT exam; to carry out PAN or CEPH exam the PAN/CEPH image detector must be manually moved from PAN bay to CEPH bay and viceversa;

to carry out PAN or CBCT exam the equipment automatically exposes in front of X-Ray source the right detector, depending on the selected exam;

the PAN and CEPH exams can be carried out with different extension in height (up to 146mm height on image receptor for PAN, up to 220mm height on image receptor for CEPH), primary collimator being automatically adjusted in height depending on exam selected by the user;

the CBCT exams can be carried at different FOV dimensions;

the equipment has a touch screen panel as main user interface;

PAN / CBCT TLD version

The equipment can carry out both exams of the PAN group (PAN, TMJ, Sinus) and of the CBCT group; There is only one image detector, for PAN and CBCT exam; the equipment automatically exposes in front of X-Ray source the right portion of detector, depending on the selected exam;

PAN exams can be carried out with different extension in height (up to 146 mm height on image receptor), since the collimator height can be adjusted depending on the exam selected by the user;

CBCT exams can be carried out at different FOV dimensions;

The equipment can be upgraded in the field to CEPH configuration by simple procedure, adding CEPH arm with CEPH image detector.

PAN / CBCT / CEPH TLD version


The equipment can carry out all feasible exams of the PAN group (PAN, TMJ, Sinus), of the CEPH group (AP/PA, LL, carpus) and of the CBCT group;

To carry out PAN or CEPH or CBCT exam the equipment automatically exposes in front of X-Ray source the right detector, depending on the selected exam;

PAN and CEPH exams can be carried out with different extension in height (up to 146 mm height on image receptor for PAN, up to 220 mm height on image receptor for CEPH), since the primary collimator can be automatically adjusted in height depending on exam selected by the user;

the image detector is not movable from CEPH bay;

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 14 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

CBCT exams can be carried out at different FOV dimensions;

IX. EXAMS

In its different configurations X-MIND trium equipment carries out the following exams (also said projections):

PAN group

- Standard panorama
- Child special SID panorama (only with TLD cassette)
- TMJ
- Bitewing
- Sectorial
- Maxillary sinuses

CBCT type

- Full view
- Extended view (TBC)
- Stitched (TBC)
- Child special SID CBCT (only with TLD cassette)

CEPHALOGRAPHIC type


- Frontal
- Lateral
 - o Child lateral (reduced height and width) (non c'è nel manuale)
 - o Complete lateral (non c'è nel manuale)
- Hand acquisition (special support needed)

X. Setup configuration

Wall mounted

For this type of installation, the column is provided with two dedicated plates to fasten the equipment to the wall.

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 15 of 58 |

| | | |
|---|--|-----------------------------------|
|  | <p style="text-align: center;">X-MIND trium 510(k) Summary</p> | <p>VOL. 005 002_k Summary</p> |
| <p>Via Roma, 45 21057 Olgiate Olona (VA)</p> | | |





Figure 8: wall mounted configuration

| | |
|---|---------------------|
| 1 | 2 expansion anchors |
| 2 | 2 expansion anchors |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 16 of 58 |

| | | |
|---|--|-----------------------------------|
|  | <p style="text-align: center;">X-MIND trium 510(k) Summary</p> | <p>VOL. 005 002_k Summary</p> |
| <p>Via Roma, 45 21057 Olgiate Olona (VA)</p> | | |

Self-standing

The self-standing accessory allows to install X-MIND trium without structurally fix it to the wall; the self-standing alone is designed to give full structural support to X-MIND trium for installation where the walls do not guarantee such support.

Although self-standing is sufficient to provide alone the needed structural support to X-MIND trium, it is anyway highly recommended to fix the X-MIND trium also to the wall, to reduce as much as possible the vibrations during the exams; this wall fixation has not safety implications, so it is responsibility of the customer to use it or not.





Figure 9: wall mounted base

| | | |
|-------------------|-----------------------------|----------------------|
| <p>Rev. n. 2</p> | <p>X-MIND trium 510 (k)</p> | <p>VOL. 005</p> |
| <p>10/01/2020</p> | | <p>Page 17 of 58</p> |

| | | |
|---|--------------------------------|---------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 |
| | | 002_k Summary |

XI. ACCESSORIES









| | |
|------------------------------|--|
| XMT SELF STANDING BASE (FWT) |  |
| XMT REMOTE BUTTON |  |
| EXTERNAL LIGHT 100 - 240 V |  |
| XMT CARPUS SUPPORT |  |
| XMt UPG CEPH |  |
| XMt UPG 3D |  |
| XMt TRUE LOW DOSE CASSETTE |  |

Table 6: accessories


| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 18 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

XII. CONFIGURATIONS



| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 19 of 58 |

| | | |
|---|--|-----------------------------------|
|  | <p>X-MIND trium 510(k) Summary</p> | <p>VOL. 005 002_k Summary</p> |
| <p>Via Roma, 45 21057 Olgiate Olona (VA)</p> | | |

TRUE
LOW DOSE




X-MIND trium TLD Pan 3D



X-MIND trium TLD Pan Ceph 3D

| | | |
|------|---|---|
| Pan | ● | ● |
| 3D | ● | ● |
| Ceph | ○ | ● |


| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 20 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

XIII. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE


| Characteristic | Predicate | Proposed | Change Id |
|--------------------|--|--|-----------|
| | X-MIND trium | X-MIND trium | |
| Manufacturer | de Götzen S.r.l. – Acteon Group Via Roma, 45 21057 Olgiate Olona VA - Italia | de Götzen S.r.l. – Acteon Group Via Roma, 45 21057 Olgiate Olona VA - Italia | same |
| 510 (k) number | K160166 | | |
| Indication for use | <p>X-MIND trium is a digital panoramic, cephalometric and tomographic extra-oral X-ray system, indicated for use in:</p> <ul style="list-style-type: none"> - producing panoramic X-ray images for diagnostic examination of dentition (teeth), jaws and oral structures; - producing radiographs of maxillofacial region and parts of the skull for cephalometric examination, if equipped with CEPH arm; - producing radiographs of hands and wrists for carpus examination, if equipped with CEPH arm; | <p>X-MIND trium is a digital panoramic, cephalometric and tomographic extra-oral X-ray system, indicated for use in:</p> <ul style="list-style-type: none"> - producing panoramic X-ray images for diagnostic examination of dentition (teeth), jaws and oral structures; - producing radiographs of maxillofacial region and parts of the skull for cephalometric examination, if equipped with CEPH arm; - producing radiographs of hands and wrists for carpus examination, if equipped with CEPH arm; | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 21 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| Characteristic | Predicate | Proposed | Change Id |
|----------------------------|--|--|-----------|
| | X-MIND trium | X-MIND trium | |
| | - producing tomographic images of the oral and maxillofacial region, for diagnostic examination of dentition (teeth), jaws, oral structures and some cranial bones, if equipped with CBCT option. | - producing tomographic images of the oral and maxillofacial region, for diagnostic examination of dentition (teeth), jaws, oral structures and some cranial bones, if equipped with CBCT option. | |
| Performance specifications | 1) Panoramic 2) Cephalometric (optional) 3) CBCT (optional) | 1) Panoramic 2) Cephalometric (optional) 3) CBCT (optional) | same |
| Exam mode | 1) PAN: - standard panoramic - child panoramic - TMJ - bitewing - sectorial panoramic with improved orthogonality - maxillary sinuses 2) CEPH: - frontal (AP/PA) - lateral (LL) - hand acquisition – carpus (special support needed) | 1) PAN: - standard panoramic - child panoramic - TMJ - bitewing - sectorial panoramic with improved orthogonality - maxillary sinuses 2) CEPH: - frontal (AP/PA) - lateral (LL) - hand acquisition – carpus (special support needed) | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 22 of 58 |

| | | |
|---|--|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| Characteristic | Predicate | Proposed | Change Id |
|---|--|--|-----------|
| | X-MIND trium | X-MIND trium | |
| | 3) CBCT: - full view - extended view | 3) CBCT: - full view - extended view | |
| Type of power supply | Single-phase alternating current power supply | Single-phase alternating current power supply | same |
| Nominal voltage | 100 V - 240V | 100 V - 240V | same |
| Maximum voltage variation | ±10% | ±10% | same |
| Nominal current | 7 A (@ 240 V) 15 A (@ 100 V) | 7 A (@ 240 V) 15 A (@ 100 V) | same |
| Frequency | 50 - 60 Hz | 50 - 60 Hz | same |
| Operation mode | Continuous operation with intermittent X-ray loading | Continuous operation with intermittent X-ray loading | same |
| Protection against electrical shock | Class I | Class I | same |
| Degree of protection against electrical shock | Type B | Type B | same |
| Emission | Group 1, Class A | Group 1, Class A | same |
| Manufacturer of the X-ray tube | Skand X | Skand X | same |
| X-ray tube model | CEI OPX/105 | CEI OPX/105 | same |
| Nominal High Voltage | 60 - 85 kV for PAN/CEPH | 60 - 85 kV for PAN/CEPH | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 23 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| Characteristic | Predicate | Proposed | Change Id |
|-----------------------------|--|--|-----------|
| | X-MIND trium | X-MIND trium | |
| | 80 - 90 kV for CBCT | 80 - 90 kV for CBCT | |
| Anode heat storage capacity | 30 kJ | 30 kJ | same |
| Nominal focal spot diameter | 0.5 mm | 0.5 mm | same |
| Anode material | Tungsten | Tungsten | same |
| Target Angle | 5° | 5° | same |
| Anodic current | 4 - 10 mA for PAN/CEPH 4 - 10 mA for CBCT | 4 - 10 mA for PAN/CEPH 4 - 10 mA for CBCT | same |
| Total filtration | 2.8 mm Al for PAN/CEPH at 85 kV 7.0 mm Al for CBCT at 90 kV | 2.8 mm Al for PAN/CEPH at 85 kV 7.0 mm Al for CBCT at 90 kV | same |
| HVL | 3.4 mm Al for PAN/CEPH at 85 kV 5.2 mm Al for CBCT at 90 kV | 3.4 mm Al for PAN/CEPH at 85 kV 5.2 mm Al for CBCT at 90 kV | same |
| High Voltage | kV ± 10% | kV ± 10% | same |
| Current in the tube | mA ± 20% | mA ± 20% | same |
| Exposure time seconds | Seconds ± (5% + 50ms) | Seconds ± (5% + 50ms) | same |
| Detector technology | CMOS | CMOS | same |
| Detector manufacturer | Hamamatsu | Hamamatsu | same |


| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 24 of 58 |

| | | |
|---|--|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| Characteristic | Predicate | Proposed | Change Id |
|---------------------------------|--|---|-----------|
| | X-MIND trium | X-MIND trium | |
| Detector model | PAN: C10500D-42 | PAN: C10500D-70 | 1 |
| | CBCT: C12280D-40 | PAN/CBCT: C12903D-40G (TLD) | 3 |
| | CEPH: C10502D-42 | CEPH: C10502D-70 | 2 |
| Pixel size | 100 µm | 100 µm | same |
| Typical image size | PAN: 1480 px x 2658 px (C10500D-42) CEPH: 2232 px x 2673 px (C10502D-42) | PAN: 1480 px x 2658 px (C10500D-70) | 1 |
| | | PAN: 1400 px x 2658 px (C12903D-40G) | 3 |
| | | CEPH: 2232 px x 2673 px (C10502D-70) | 2 |
| Active area | PAN: 148 mm x 6 mm (C10500D-42) CBCT: 121.6 mm x 123.2 mm (C12280D-40) CEPH: 220 mm x 6 mm (C10502D-42) | PAN: 148 mm x 6 mm (C10500D-70) | 1 |
| | | PAN: 140 mm x 6 mm (C12903D-40G) | 3 |
| | | CBCT: 139.2 mm x 121.2 mm (C12903D-40G) | 3 |
| | | CEPH: 223,2 mm x 6 mm (C10502D-70) | 2 |
| Theoretical detector | 4.5 lp/mm | 4.5 lp/mm | 1, 2 |
| Actual in imaging ¹ | 4.5 lp/mm | 4.5 lp/mm | 1, 2 |
| Image resolution ² : | PAN: high contrast resolution: 3.1 lp/mm | PAN: high contrast resolution: 3.1 lp/mm | 1 |
| | PAN: low contrast resolution: Four holes are visible | PAN: low contrast resolution: Four holes are visible | 1 |


² Resolution data are taken from results of bench test performed on both devices (See section 15 Bench tests)

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 25 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| Characteristic | Predicate | Proposed | Change Id |
|------------------------------------|--|--|-----------|
| | X-MIND trium | X-MIND trium | |
| | CEPH: high contrast resolution: 3.1 lp/mm CEPH: low contrast resolution Four holes are visible | CEPH: high contrast resolution: 3.1 lp/mm CEPH: low contrast resolution Four holes are visible | 2 2 |
| X-ray exposure time | PAN: 16.8 s | PAN: 16.8 s | 1 |
| | CEPH:18 s | CEPH:18 s | 2 |
| Detector technology | Flat panel CMOS | Flat panel CMOS | same |
| Detector manufacturer | Hamamatsu | Hamamatsu | same |
| Detector model | C12280D-40 | C12903D-40G | 3 |
| Image voxel size | 75 µm (minimum slice thickness) | 75 µm (minimum slice thickness) | 3 |
| FOV (diameter x height) | 40 mm x 40 mm | 40 mm x 40 mm | 3, 4 |
| | 60 mm x 60 mm | 60 mm x 60 mm | |
| | 80 mm x 80 mm | 80 mm x 80 mm | |
| | 80 mm x 90 mm | 80 mm x 90 mm | |
| | 110 mm x 80 mm | 110 mm x 80 mm | |
| Detector resolution ² : | High contrast resolution: V10%: 2.07 lp/mm: | High contrast resolution: V10%: 2.81 lp/mm: | 3 |
| | High contrast resolution: V50%: 0.72 lp/mm | High contrast resolution: V50%: 0.95 lp/mm | |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 26 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| Characteristic | Predicate | Proposed | Change Id |
|----------------------------------|--|--|--------------|
| | | X-MIND trium | X-MIND trium |
| Rotation | 360° | 360° | same |
| Scanning time | 12 s to 30 s maximum | 12 s to 30 s maximum | same |
| X-ray real exposure time | 6s (Normal quality) 7.2s (Medium quality) 9s (High quality) | 6s (Normal quality) 7.2s (Medium quality) 9s (High quality) | same |
| Source Image Distance (SID) | PAN: 550mm CBCT: 633mm | PAN: 550mm (adult) PAN: 490mm (special small child) CBCT: 610 (adult) CBCT: 588 (special small child) | 4 |
| Protocol types (CBCT exams only) | Standard Dose protocols | Standard Dose protocols Low Dose (LD) protocols | 7 |
| Dimensions | Footprint: PAN/CBCT: max 1100 mm x 1505 mm With CEPH: max 1715 mm x 1500 mm Height: 2350 mm | Footprint: PAN/CBCT: max 1100 mm x 1505 mm With CEPH: max 1715 mm x 1500 mm Height: 2350 mm | same |
| Weight | PAN: 230 kg PAN/CBCT: 240 kg PAN/CBCT With CEPH: 280 kg | Wall mount: PAN: 230 kg PAN/CBCT: 240 kg PAN/CBCT With CEPH: 280 kg Self standing: PAN: kg 280 | 6 |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 27 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| Characteristic | Predicate | Proposed | Change Id |
|--------------------------|------------------------------------|--|--------------|
| | | X-MIND trium | X-MIND trium |
| | | PAN/CBCT: kg 290 PAN/CBCT With CEPH: kg 330 | |
| Mechanical configuration | Wall mount | Wall mount Self standing | 6 |
| Patient position | Standing, seated or on wheel-chair | Standing, seated or on wheel-chair | same |

Table 7 Comparison between proposed and predicate devices.

Both X-MIND trium and X-MIND trium TLD are digital panoramic, cephalometric and tomographic extra-oral X-ray systems. For analysis of patients' images, X-MIND trium can be provided with a dedicated workstation which has the following modules installed:

- ACTEON IMAGING SUITE (AIS) + 2D diagnostic analysis
- AIS 3D app for CBCT dataset only (clearance K173041)

The 2D diagnostic analysis and the AIS 3D app (CBCT) provide a complete set of tools for visualization and post-processing of 2D radiographs (such as Panoramic, Cephalometric and Carpus radiographs) and 3D CBCT datasets.


AIS 3D app also includes special features for computer aided dental implant planning.

X-MIND trium TLD introduce a new SID movement (source to image distance) and the possibility to perform both CBCT and PAN exam with the same flat panel. In order to support his new movement and cassette a new version of AUXCBCT board, AUXCBCT Plus, has been developed with the related firmware,

X-MIND trium software is described in detail in Section 13 of this submission.

The following summary table shows a comparison of the functionalities for image processing, enhancing and visualization available in X-MIND trium cleared version and new version that includes LD functionality

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 28 of 58 |


| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

This comparison is to provide evidence of substantial equivalence between the tools and means made available by both software module.

Please note that 3D functionalities are now covered by clearance K173041


| Software characteristic | X-MIND trium | X-MIND trium | Change Id |
|--|---|---|-----------|
| | 2D analysis + AIS 3D app | 2D analysis + AIS 3D app | |
| Supported Image format | DICOM, JPEG, BMP, TIFF, GIFF, PNG, DGI (proprietary format) | DICOM, JPEG, BMP, TIFF, GIFF, PNG, DGI (proprietary format) | same |
| Functionalities and features for 2D images: single 3D tomographic sections, panoramic and cephalometric radiographs | | | |
| Basic image visualization features | + Zoom and pan + Rotation + Flipping / mirroring, vertical / horizontal | + Zoom and pan + Rotation + Flipping / mirroring, vertical / horizontal | same |
| Measuring function | + Measuring distances + Measuring angles + Adding annotation + Viewing Pixel profile | + Measuring distances + Measuring angles + Adding annotation + Viewing Pixel profile | same |
| Visualization filters | + Gamma filter + Sharp filter + Equalize filter + Emboss filter + Reverse filter | + Gamma filter + Sharp filter + Equalize filter + Emboss filter + Reverse filter | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 29 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| | | | |
|--|---|---|------|
| Pseudocolour analysis LUTs | + Gray + Inverse Spectrum | + Gray + Inverse Spectrum | same |
| Windowing | + Contrast-brightness set + Lower-Upper level + SHS optimal windowing | + Contrast-brightness set + Lower-Upper level + SHS optimal windowing | same |
| Reporting | Included | Included | same |
| Functionalities and features for CBCT 3D dataset | | | |
| Basic image visualization features | + Zoom and pan + Rotation + Flipping / mirroring, vertical / horizontal | + Zoom and pan + Rotation + Flipping / mirroring, vertical / horizontal | same |
| Measuring function | + Measuring distances + Measuring angles + Adding annotation + Viewing Voxel profile + Measuring polygon + Measuring volumes + Measuring distance on 3D | + Measuring distances + Measuring angles + Adding annotation + Viewing Voxel profile + Measuring polygon + Measuring volumes + Measuring distance on 3D | same |
| Processing filters | + Sharp filter | + Sharp filter | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 30 of 58 |

| | | |
|---|--|---------------|
|  | X-MIND trium 510(k) Summary | VOL. 005 |
| | | 002_k Summary |
| Via Roma, 45 21057 Olgiate Olona (VA) | | |


| | | | |
|--|--|---|------|
| | + Median filter + SHARP filter: denoiser and edge enhancer + STAR filter: streak and metal artifacts reduction | + Median filter + SHARP filter: denoiser and edge enhancer + STAR filter: streak and metal artifacts reduction | |
| Multi Planar Reconstruction (MPR) | + Axial plane + Coronal plane + Sagittal plane | + Axial plane + Coronal plane + Sagittal plane | same |
| Oblique MPR | Included | Included | same |
| Single and multiple cross sections visualization | Included + Setting number of sections + Setting distance of sections | Included + Setting number of sections + Setting distance of sections | same |
| Windowing | + Contrast-brightness set + Lower-Upper level + Editing preset | + Contrast-brightness set + Lower-Upper level + Editing preset | same |
| Display mode | + MPR (voxel intensity) + Ray Sum + Maximum intensity projection (MIP) + Minimum intensity projection (minMIP) + SCOUT | + MPR (voxel intensity) + Ray Sum + Maximum intensity projection (MIP) + Minimum intensity projection (minMIP) | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 31 of 58 |

| | | |
|---|--|---------------|
|  | X-MIND trium 510(k) Summary | VOL. 005 |
| | | 002_k Summary |
| Via Roma, 45 21057 Olgiate Olona (VA) | | |

| | | | |
|--|---|---|------|
| | + Volume thickness selection (for MIP and Ray Sum) | + SCOUT + Volume thickness selection (for MIP and Ray Sum) | |
| 3D volume rendering | + Preset template + Editing of user-defined template | + Preset template + Editing of user-defined template | same |
| Surface Rendering | + Surface generation + STL export | + Surface generation + STL export | same |
| Curve Planar Reformation (CPR) for dental analysis | + Dentascan wizard | + Dentascan wizard | same |
| Advanced 3D visualization tools | + 3D clipping: 3D cut by orthogonal planes + 3D thin slab rendering + 3D sculpting: 3D cut by user-defined planes + 3D flying mode: 3D virtual endoscopy | + 3D clipping: 3D cut by orthogonal planes + 3D thin slab rendering + 3D sculpting: 3D cut by user-defined planes + 3D flying mode: 3D virtual endoscopy | same |
| Advanced 3D functions | + Implant planning wizard + Nerve tracing + Importing 3D STL object in patient volume + STL implant library + Measuring bone density at implant sites | + Implant planning wizard + Nerve tracing + Importing 3D STL object in patient volume + STL implant library | same |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 32 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| | | | |
|---------------------------------|---|--|------|
| | + Bone/tissue removal | + Measuring bone density at implant sites + Bone/tissue removal | |
| General purpose functionalities | | | |
| Screenshot capture function | Image capture (AIS 3D app) 2D analysis module (Export) | Image capture (AIS 3D app) 2D analysis module (Export) | same |

Table 8: comparison between Acteon Imaging Suite and predicate device software.

The above tables show that the two devices have the same intended use.

We have conducted the risk analysis and the necessary verification and validation activities to demonstrate that the design outputs of the modified device meet the design input requirements according to our design and development procedure in compliance with 21 CFR § 820.30 and design control guidance for medical device manufacturers.

The impact of the modification on the device and its components on the risk has been evaluated Re-applying Fault Tree Analysis (FTA) and identifying additional events that can generate undesired consequences (top events) originally identified with Preliminary Hazard Analysis (PHA).


Change in detectors have been verified following “Guidance for the Submission of 510(k)’s for Solid State X-ray Imaging Devices”

Bench testing was used to evaluate the trade-off between the image quality and radiation dose of the images obtained with dose reduction due to addition of two main aspects:

- 1) TLD cassette (when using special exam for child)
- 2) Low Dose functionality, CBCT mode only

Based on the bench test results on image quality and radiation dose obtained with X-MIND trium cleared device and X-MIND trium proposed device (with TLD cassette and separately with Low Dose functionality), we can assume that the clinical performance obtained with the above changes is equal or better than the ones of the predicate device; these results, confirm the equivalence of X-MIND trium proposed device to X-MIND trium predicate device.


| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 33 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| Id | Device Change | Risks | Verification/Validation Method(s) | Acceptance Criteria | Summary of results ³ |
|----|-------------------|--|--|--|---------------------------------|
| 1 | New PAN detector | 21A. Ionising radiation hazard Harm: biological effects due to X-ray ionising radiation due to event 4) Unwanted x-ray exposure due to poor image quality (inadequate for diagnostic purpose) | Guidance for the Submission of 510(k)'s for Solid State X-ray Imaging Devices" | Same technology DQE values for new detector are in the range ± 0.1 on the whole frequency range (1 - 4 LP/mm) | Same PASS |
| 2 | New CEPH detector | 21A. Ionising radiation hazard Harm: biological effects due to X-ray ionising radiation due to event 4) | Guidance for the Submission of 510(k)'s for Solid State X-ray Imaging Devices" | Same technology DQE values for new detector are in the range ± 0.1 on the whole frequency range (1 - 4 LP/mm) | Same PASS |


³ Complete results are available in substantial equivalence discussion section 11 and in specific reports of V&V and bench tests available in sections 15, 16

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 34 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 | |
|--------|-----------------------|---|--|---|---|---|
| | | Unwanted x-ray exposure due to poor image quality (inadequate for diagnostic purpose) | | | | |
| 3 | New PAN-CBCT detector | New detector | 21A. Ionising radiation hazard Harm: biological effects due to X-ray ionising radiation due to event 4) Unwanted x-ray exposure due to poor image quality (inadequate for diagnostic purpose) | Bench test according to same protocol used in X-MIND trium Submission K160166 | Same or superior trade-off between dose and image quality (i.e. System Indicator (SI) in CBCT) $S_{DAP} \geq S_{DAP \text{ predicate}}$ $S_{DAP} = \frac{CNR}{DAP \left(\frac{1}{2 \cdot V_{50\%}} \right)^2}$ CNR is the contrast to noise ratio (overaged value from 5 measurements) DAP is the dose area product declared in IFU V50% is the resolution indicator at 50% modulation | <u>Woman Medium Standard quality</u> FOV 80X80-80X90 $S_{DAP \text{ TLD}} = > S_{DAP \text{ pred}}$ PASS <hr/> <u>Woman Medium Standard quality</u> FOV 110X80 $S_{DAP \text{ TLD}} >$ PASS <hr/> <u>Child medium</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 35 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|--|---|
| | | | | | <u>Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Adult</u> $S_{DAP\ TLD} >$ $S_{DAP\ pred}$ PASS |
| | | | | | <u>Child medium Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Child</u> $S_{DAP\ TLD} >$ $S_{DAP\ pred}$ PASS |
| | | | | Voxel Noise [a.u.] with TLD cassette including new PAN-CBCT sensor ($VN_{TLD} \leq$ Voxel Noise [a.u.] predicate (VN) | <u>Woman Medium Standard quality</u> FOV <u>80X90-</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 36 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|---------------------|--|
| | | | | | VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Woman Medium Standard quality</u> FOV <u>110X80</u> VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Child medium Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Adult</u> VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Child medium Standard d quality</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 37 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|---|--|
| | | | | | <u>FOV</u> <u>80X80-</u> <u>80X90</u> <u>SID Child</u> $VN_{TLD} < VN$ <u>PASS</u> |
| | | | | NPS _{0D} [a.u. ²] with TLD cassette including new PAN-CBCT sensor $\leq NPS_{0D}$ [a.u. ²] pedicate | <u>Woman</u> <u>Medium</u> <u>Standard</u> <u>quality</u> <u>FOV</u> <u>80X90-</u> $NPS_{0D\ TLD} < NPS_{0D}$ <u>PASS</u> <u>Woman</u> <u>Medium</u> <u>Standard</u> <u>quality</u> <u>FOV</u> <u>110X80</u> $NPS_{0D\ TLD} < NPS_{0D}$ <u>PASS</u> |
| | | | | | <u>Child</u> <u>medium</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 38 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---|--|---|
| | | | | | <u>Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Adult</u> NPS _{0D TLD} < NPS _{0D} <u>PASS</u> |
| | | | | | <u>Child medium Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID child</u> NPS _{0D TLD} < NPS _{0D} <u>PASS</u> |
| | | | IEC 60601-2-63 203.8.5.3 Correspondence | c) along each of the two axes of the IMAGE RECEPTION AREA, the | <u>SID Adult edge ≤ 20mm PASS</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 39 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---|---|--|
| | | | between X-RAY FIELD and EFFECTIVE IMAGE RECEPTION AREA (rectangular area) | edges of the X-RAY FIELD shall not exceed the corresponding edges of the EFFECTIVE IMAGE RECEPTION AREA by more than 2 cm or 3 % of the indicated FOCAL SPOT TO IMAGE RECEPTOR DISTANCE when the IMAGE RECEPTION PLANE is normal to the X-RAY BEAM AXIS, whichever the larger; | <u>Sid Child</u> edge ≤ 20mm PASS |
| | | | | d) the sum of the discrepancies on both axes shall not exceed 3 cm or 4 % of the indicated FOCAL SPOT TO IMAGE RECEPTOR DISTANCE, whichever the larger | < 30mm PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 40 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|--|--|
| | | | | e) along the axis of the IMAGE RECEPTION AREA that is parallel to the direction of the scanning, the X-RAY FIELD shall not exceed the EFFECTIVE IMAGE RECEPTION more than 1 mm on each side. | ≤ 1 PASS |
| | | | | f) along the axis of the IMAGE RECEPTION AREA that is perpendicular to the direction of the scanning the X-RAY FIELD shall not exceed the EFFECTIVE IMAGE RECEPTION AREA | =0 PASS |
| | | | 21 CFR 1020.31 f) 4) | field does not exceed each dimension of the image receptor by more than 2 percent of the SID, | <u>PAN</u> ≤ 2% PASS |
| | | | | | <u>PAN reduced SID</u> ≤ 2% PASS |
| | | | | | <u>CBCT</u> ≤ 2% PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 41 of 58 |

| | | |
|---|--|---------------------------|
|  | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | Via Roma, 45 21057 Olgiate Olona (VA) | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|--------------------------|--|---------------------------------------|---|--|
| | | | | | <u>CBCT</u> <u>reduced</u> <u>SID</u> ≤2% PASS |
| | | | | align the center of the x-ray field with the center of the image receptor to within 2 percent of the SID, | <u>Only for</u> <u>PAN and</u> <u>CEPH</u> PASS |
| | Addition of FOV | No variation in intended use | N.A. | N..A. | N.A. |
| 4 | Addition of TLD cassette | New TLD cassette 3D. Cutting or severing, trapping, crushing, impact, instability. 2) Fall of suspended and moving masses | ES60601-1: 2005/(R)2012 and A1:2012, | Safety factor ≥2.5 | k >>2.5 respect to local collapse PASS |
| | New Child special exam | 29A Unintended use of the device due to | IEC 60601-1-6 | The criteria to pass usability verification are that all the comments are positive. | PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 42 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---|---|---|---|--|
| | | human factors Event 6) Misplace adult patients during child exam Inadequate arrangement of controls and indicators | | | |
| | New tool Child special exam | 4C. Contact allergy due to Inappropriate materials, not biocompatible | Same material, same plastic injection process of chin rest, bite block ISO 10993-5: 2009 ISO 10993-10: 2010 | Cytotoxicity: No cytotoxic Hypersensitivity: Not sensitizing – | PASS PASS |
| | True Low dose (TLD) (reduction of SID) | 21A. Ionising radiation hazard Harm: biological effects due to X-ray | Bench test PAN mode according to same protocol used in X-MIND trium Submission | Absence of artifacts Equal Hight contrast resolution | <u>Both verified</u> <u>PASS</u> <u>Both 3.1 lp/mm</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 43 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|---|--|---|--|
| | | ionising radiation due to event 4) Unwanted x-ray exposure due to poor image quality (inadequate for diagnostic purpose) | K160166 | | <u>PASS</u> |
| | | | | Equal Low contrast resolution | <u>Both 4 holes visible</u> <u>PASS</u> |
| | | | | Declared DAP PAN CHILD NEDIUM SID ≤predicate | <u>38.77<52.7</u> <u>PASS</u> |
| | | | Bench test CBCT mode according to same protocol used in X-MIND trium Submission K160166 $S_{DAP} \geq S_{DAP \text{ predicate}}$ $S_{DAP} = \frac{CNR}{DAP \left(\frac{1}{2 \cdot V_{50\%}} \right)^2} \cdot CNR$ is the contrast to noise | Same or superior trade-off between dose and image quality (i.e. System Indicator (SI) in CBCT) <u>Woman Medium Standard quality</u> <u>FOV 80X80-80X90</u> $S_{DAP \text{ TLD}} > S_{DAP \text{ pred}}$ <u>PASS</u> | |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 44 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|--|---|
| | | | | ratio (overaged value from 5 measurements) DAP is the dose area product declared in IFU V50% is the resolution indicator at 50% modulation | <u>Woman</u> <u>Medium</u> <u>Standard</u> <u>quality</u> <u>FOV</u> <u>110X80</u> $S_{DAP} TLD >$ $S_{DAP} pred$ PASS |
| | | | | | <u>Child</u> <u>medium</u> <u>Standard</u> <u>d quality</u> <u>FOV</u> <u>80X80-</u> <u>80X90</u> <u>SID Adult</u> $S_{DAP} TLD =$ $> S_{DAP}$ pred PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 45 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|--|--|
| | | | | | <u>Child medium Standard d quality</u> <u>FOV 80X80- 80X90</u> <u>SID Child</u> $S_{DAP\ TLD} >$ $S_{DAP\ pred}$ PASS |
| | | | | Voxel Noise [a.u.] with TLD cassette including new PAN-CBCT sensor ($VN_{TLD} \leq$ Voxel Noise [a.u.] predicate (VN) | <u>Woman Medium Standard quality</u> <u>FOV 80X90-</u> $VN_{TLD} <$ VN PASS |
| | | | | | <u>Woman Medium Standard quality</u> <u>FOV 110X80 110X90</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 46 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|---------------------|--|
| | | | | | VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Child medium Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Adult</u> VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Child medium Standard d quality</u> FOV <u>80X80- 80X90</u> <u>SID Child</u> VN _{TLD} < VN <u>PASS</u> |
| | | | | | <u>Woman Medium</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 47 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|---|---|
| | | | | NPS _{0D} [a.u. ²] with TLD cassette including new PAN-CBCT sensor ≤ NPS _{0D} [a.u. ²] predicate | <u>Standard quality</u> <u>FOV 80X90 90X90- NPS_{0D} TLD < NPS_{0D} PASS</u> <u>Woman Medium Standard quality</u> <u>FOV 110X80 110X90 NPS_{0D} TLD < NPS_{0D} PASS</u> <u>Child medium Standard d quality</u> <u>FOV 80X80- 80X90 SID Adult</u> |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 48 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|------------------------|--|---|---|
| | | | | | NPS _{0D TLD} < NPS _{0D} <u>PASS</u> <u>Child medium Standard d quality</u> <u>FOV 80X80- 80X90</u> <u>SID Child</u> NPS _{0D TLD} < NPS _{0D} <u>PASS</u> |
| | | Detectability study | As requested by FDA RAI 23th August 2019 | Same pathological or clinically relevant conditions identifiable both in the STD protocol scan and in the TLD protocol one | PASS All pathologi cal or clinically relevant condition s identified with STD protocol scan have been identifies also with |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 49 of 58 |

| | | |
|---|--|--|
|  | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | Via Roma, 45 21057 Olgiate Olona (VA) |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|--|--|---|------------------------|---|
| | | | | | the LD protocol one Note Pos. 45 "worn out dental neck" not identified nor in standard mode nor in TLD configura tion |
| | Software (new moveme nt managed by firmware) | PEMS 1.6 Unintended movement of the U- arm, patient support, SID regulation axis or mobile column /patient | Same protocol used protocol used in X- MIND trium Submission K160166 | All test shall be PASS | All test PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 50 of 58 |

| | | |
|---|--|-------------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| I d | Device Change | | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|--|--|--|--|---|---|
| 5 | Change of beam limitin g device | Change of aperture dimensio ns | 5A Biological effects due to X-ray ionising radiation. 2) Unwanted X-ray due to unsuitable or fault shielding (causing leakage radiation), filtration or collimation | IEC 60601-2-63 203.8.5.3 Correspondence between X-RAY FIELD and EFFECTIVE IMAGE RECEPTION AREA (rectangular area) | c) along each of the two axes of the IMAGE RECEPTION AREA, the edges of the X-RAY FIELD shall not exceed the corresponding edges of the EFFECTIVE IMAGE RECEPTION AREA by more than 2 cm or 3 % of the indicated FOCAL SPOT TO IMAGE RECEPTOR DISTANCE when the IMAGE RECEPTION PLANE is normal to the X-RAY BEAM AXIS, whichever the larger; | <u>SID Adult</u> edge < 20mm PASS |
| | | | | 21 CFR 1020.31 f) 4) | field does not exceed each dimension of the image receptor by more than 2 percent of the SID, | <u>PAN</u> < 2% PASS |
| | | | | | | <u>PAN reduced SID</u> < 2% PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 51 of 58 |

| | | |
|---|--|---------------------------|
|  | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | Via Roma, 45 21057 Olgiate Olona (VA) | |


| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|--|---|--|---|--|
| | | | | | <u>CBCT</u> < 2% PASS |
| | | | | | <u>CBCT</u> <u>reduced</u> <u>SID</u> < 2% PASS |
| | | | | align the center of the x-ray field with the center of the image receptor to within 2 percent of the SID, | PASS Only for PAN and CEPH |
| | | | 21 CFR 1020.30 (k) | Same geometry as the previous version | PASS |
| 6 | Addition of self standing version self standing base | Hazardous situations: 3D . Cutting or severing, trapping, crushing, impact, instability. Event (8) Unbalancing of the self-standing device due to base mechanical | IEC 60601-1 9.4.2 instability overbalance | At 10° in all directions no overbalance | PASS |
| | | | 9.8.2 tensile safety factor | Safety factor ≥ 4 | >4 with respect to yielding (local plastic strains) >4 with respect to local collapse |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 52 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |


| Id | Device Change | Risks | Verification/Validation Method(s) | Acceptance Criteria | Summary of results |
|----|----------------------------|--|---|--|--|
| | | failure or base-column connection mechanical failure. | | | PASS |
| 7 | The Low Dose functionality | X-ray: AIS 5 14 Biological effects due to X-ray ionising radiation Event 2: Overexposure of X-ray due to the need of exposure to be repeated due to fault in AIS software functions - New acquisition | Bench test according to same protocol used in X-MIND trium Submission K160166 | Same or superior trade-off between dose and image quality (i.e. System Indicator (SI) in CBCT) $S_{DAP} \text{ low dose} \geq S_{DAP} \text{ standard protocol}$ $S_{DAP} = \frac{CNR}{DAP \left(\frac{1}{2 \cdot V_{50\%}} \right)^2}$ CNR (average on 5 measurements) DAP (declared in IFU) | <u>Man</u> $S_{DAP} \text{ low dose} > S_{DAP} \text{ standard protocol}$ PASS <u>Woman</u> $S_{DAP} \text{ low dose} > S_{DAP} \text{ standard protocol}$ PASS <u>Child</u> $S_{DAP} \text{ low dose} > S_{DAP} \text{ standard protocol}$ PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 53 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---------------------------------------|--|--|
| | | | | $\Delta\%_{VN} \leq 10\%$ $\Delta\%_{VN} = \frac{VNld - VNs}{100 * VNsd}$ <i>VNld</i> = Voxel Noise [a.u.] low dose protocol <i>VNsd</i> = Voxel Noise [a.u.] standard protocol | <u>Man</u> $\Delta_{VN}\% < 10\%$ PASS <u>Woman</u> $\Delta_{VN}\% < 10\%$ PASS <u>Child</u> $\Delta_{VN}\% < 10\%$ PASS |
| | | | | $\Delta\%_{NPS} \leq 30\%$ $\Delta\%_{NPS} = \frac{NPSld - NPSsd}{100 * NPSsd}$ <i>NPSld</i> = NPS^{OD} [a.u. ²] low dose protocol <i>NPSsd</i> = NPS^{OD} [a.u. ²] standard dose protocoll | <u>Man</u> $\Delta_{NPS} < 30\%$ PASS <u>Woman</u> $\Delta_{NPS} < 30\%$ PASS <u>Child</u> $\Delta_{NPS} < 30\%$ PASS |

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 54 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

| I d | Device Change | Risks | Verification/Validat ion Method(s) | Acceptance Criteria | Summary of results 3 |
|--------|---------------|-------|---|--|--|
| | | | Detectability study as requested by FDA RAI 23th August 2019 | Same pathological or clinically relevant conditions identifiable both in the STD protocol scan and in the LD protocol one | PASS All pathologi cal or clinically relevant condition s have been identified both in the STD protocol scan and in the LD protocol one |

Table 9: change traceability matrix

Details on Risk Management file, verification and validation are available in section 16

XIV. PERFORMANCE, SAFETY, EMC TESTS

The following performance data were provided in support of the substantial equivalence determination.


Electrical safety and electromagnetic compatibility (EMC)

Electrical safety and EMC testing were conducted on X-MIND trium.

The performance tests on X-MIND trium V2 have been conducted by a Nationally Recognized Testing Laboratory (NRTL) in order to verify:

- compliance with general requirements for basic safety and essential performance of medical electrical equipment, dental intra-oral X-ray equipment, diagnostic X-ray equipment, programmable electrical medical systems, high-voltage generators of diagnostic X-ray generators, X-ray source assemblies and X-ray tube assemblies for medical diagnosis;
- compliance with usability requirements;

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 55 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

- compliance with electromagnetic compatibility requirements.

The following table illustrates the standards to which X-MIND trium V2 complies, compared to those related to X-MIND trium

| Proposed device | Predicate device |
|---------------------------------|-------------------------------------|
| X-MIND trium | X-MIND trium |
| IEC 60601-1 3 rd ed. | IEC 60601-1 3.1 rd ed. |
| ANSI/AAMI ES60601-1: 2005 | ES60601-1: 2005/(R)2012 and A1:2012 |
| CAN/CSA C22.2 No. 60601-1:08 | CAN/CSA C22.2 No. 60601-1:14 |
| IEC 60601-1-2: 2007 | IEC 60601-1-2: 2014 |
| IEC 60601-1-3: 2008 | IEC 60601-1-3: 2013 |
| IEC 60601-2-63: 2012 | IEC 60601-2-63: 2017 |
| IEC 60601-1-6: 2010 | IEC 60601-1-6: 2013 |
| IEC 62366: 2007 | IEC 62366: 2016 |
| ISO 14971 ⁴ | ISO 14971 |

Table 10: comparison between the standard applied

All the test reports for X-MIND trium are available in Section 16 of this submission.


X-MIND trium has been tested according to approved verification protocols to assure its conformity to the following parts of USA Code of Federal Regulations relating to PERFORMANCE STANDARDS FOR IONIZING RADIATION EMITTING PRODUCTS USA

21 CFR §1020.30 Diagnostic x-ray systems and their major components.

21 CFR §1020.31 Radiographic equipment.

21 CFR §1020.33 Computed tomography (CT) equipment.

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 56 of 58 |

| | | |
|---|--------------------------------|---------------------------|
|  Via Roma, 45 21057 Olgiate Olona (VA) | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| | | |

Software Verification and Validation Testing

Software verification and validation testing were conducted and documentation was provided as recommended by FDA's Guidance for Industry and FDA Staff, "Guidance for the content of premarket submissions for software contained in medical devices". The software for this device was considered a "moderate" level of concern, since a failure or latent flow could indirectly result in minor injury to the patient or operator through incorrect or delayed information or through the action of a care provider.

Performance testing - Bench

Bench testing was used to evaluate the trade-off between the image quality and radiation dose of the images obtained with dose reduction due to addition of two main aspects:

- TLD cassette (when using special exam for child)
- Low Dose functionality, CBCT mode only

The technical bench tests issued to evaluate the performance and safety of the TLD (True Low Dose) cassette, reduced SIDs and the LD (Low Dose) functionality have been integrated by a clinical detectability analysis that provided evidence that the overall diagnostic quality of the CBCT datasets and dental panoramic images obtained by using the TLD or the LD functionality is basically the same as the images obtained in standard scan protocols.

The detectability analysis consisted in submitting the radiographs of cadavers to a set of intended users (radiologist, dentist, oral surgeon), who evaluated the visibility of anatomical structures and features of interests in the different scanning conditions.


Based on the bench test results on image quality and radiation dose and on the results of the detectability study obtained with X-MIND trium cleared device and X-MIND trium proposed device (with TLD cassette and separately with Low Dose functionality), we can assume that the clinical performance obtained with the above changes is equal or better than the ones of the predicate device; these results, confirmed the equivalence of X-MIND trium proposed device to X-MIND trium predicate device.

XV. [FDA guidance documents](#)

Here below is shown the list of FDA guidance documents that X-MIND trium conforms to or that we referred during the development

- FDA Guidance for Industry and FDA Staff Format for Traditional and Abbreviated 510(k)s August 2005

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 57 of 58 |

| | | |
|---|--|---------------------------|
|  | X-MIND trium 510(k) Summary | VOL. 005 002_k Summary |
| Via Roma, 45 21057 Olgiate Olona (VA) | | |

- FDA Guidance for the Submission of 510(k)'s for Solid State X-ray Imaging Devices September 2016
- FDA Guideline Use of International Standard ISO 10993-1,"Biological evaluation of medical devices Part 1: Evaluation and testing June 2016
- FDA Guideline Recommended Content and format of complete Test reports for NON-ClinicalBench performance Testing in 510K submission April 2019.
- FDA Guidance for "Pediatric Information for X-ray Imaging Device Premarket Notifications" dated November 2017
- FDA Guidance for the content of premarket submissions for software contained in medical devices May 2005
- FDA Guidance Content of premarket submission for management of cybersecurity October, 2014

I. CONCLUSIONS

Based on a comparison of intended use, indications, constructions, construction materials, principal of operation, features, technical data and performances including diagnostic quality of the acquired images,

X-MIND trium including changes listed above is safe and effective to perform its intended use as well as substantially equivalent predicate device X-MIND trium (already cleared version)

| | | |
|------------|----------------------|---------------|
| Rev. n. 2 | X-MIND trium 510 (k) | VOL. 005 |
| 10/01/2020 | | Page 58 of 58 |