



May 27, 2020

Rhythmink International, LLC  
Gabriel Orsinger  
Vice President of Engineering and R&D  
1140 First Street South  
Columbia, South Carolina 29209

Re: K200984

Trade/Device Name: Guardian Needle Electrode  
Regulation Number: 21 CFR 882.1350  
Regulation Name: Needle Electrode  
Regulatory Class: Class II  
Product Code: GXZ  
Dated: April 9, 2020  
Received: April 14, 2020

Dear Gabriel Orsinger:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

for  
Jay Gupta  
Assistant Director  
DHT5A: Division of Neurosurgical,  
Neurointerventional  
and Neurodiagnostic Devices  
OHT5: Office of Neurological  
and Physical Medicine Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

Device Name

Guardian Needle™ Electrode

Indications for Use (Describe)

RhythmLink International Subdermal Needle Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals. Examples include: Electromyography (EMG), Electroencephalography (EEG), and Nerve potential signals. The electrodes are sterile and for single patient use.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

## 510(k) Summary – K200984

<b>807.92(a)(1) Submitter Information</b>	Rhythmink International, LLC 1140 First Street South Columbia, SC 29209 Phone: 803-252-1222 FDA Registration #: 1067162
<b>Official Correspondent</b>	Gabriel Orsinger, PhD Vice President of Engineering and R&D Email: gorsinger@rhythmink.com Phone: 803-365-9664
<b>Summary Date</b>	May 27, 2020
<b>807.92(a)(2) Device Identification</b>	<b>Device Trade Name:</b> Guardian Needle™ Electrode <b>Classification Name:</b> Needle Electrode <b>Product Code:</b> GXZ <b>Classification:</b> 21 CFR 882.1350 Class II <b>Classification Panel:</b> Neurology
<b>807.92(a)(3) Predicate Device</b>	<b>Device Trade Name:</b> Subdermal Needle Electrodes <b>Classification Name:</b> Needle Electrode <b>510(k) Number:</b> K022914 <b>Product Code:</b> GXZ <b>Classification:</b> 21 CFR 882.1350 Class II <b>Classification Panel:</b> Neurology <b>Manufacturer:</b> Rhythmink International, LLC 1140 First Street South Columbia, SC 29209
<b>807.92(a)(4) Device Description</b>	Guardian Needle™ Electrodes are a sterile, single-use device. The electrodes are applied during the study of biopotentials such as electromyography (EMG), electroencephalography (EEG), nerve conduction and stimulation/response. The electrodes are invasive as they are placed subcutaneously or in contact with nerve or muscle tissue. The needle is housed inside a sheath until application and secured to the patient during use with adhesive tape.
<b>807.92(a)(5) Intended Use</b>	Rhythmink International Guardian Needle™ Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals. Examples include: Electromyography (EMG), Electroencephalography (EEG) and Nerve potential signals. The electrodes are sterile and for single patient use only
<b>807.92(a)(6) Technological Characteristics</b>	The technological characteristics of Guardian Needle™ Electrodes are identical to the predicate device (K022914), with the addition of a needle sheath and adhesive tape, both of which have been assessed to not change the fundamental scientific technology, intended use, performance, or safety and effectiveness of the predicate device (reference Substantial Equivalence of

Technological Characteristics table, below). The test methods were identical to those used to assess the predicate device.

**807.92(b)(1)  
Summary of Non-  
Clinical Tests**

Guardian Needle™ Electrodes are substantially equivalent in technology, safety, and effectiveness as the predicate device (K022914), as demonstrated by the test results.

The Guardian Needle™ Electrode was assessed by applying the TIR 28:2009 Guidance, from which it was determined that the subject device does not challenge the Submitter's Master Challenge Device and can thus fall under the previously validated EtO sterilization cycle. Reevaluation of the EtO residuals was performed to confirm low levels of EtO and Ethylene Chlorohydrin (ECH) after a 24-hour aeration time below the maximum limits described in ANSI/AAMI/ISO 10993-7:2008/(R) 2012.

The biocompatibility evaluation was conducted within the risk management framework and in compliance with ISO 10993 standards. This biocompatibility evaluation establishes the biological safety for the Guardian Needle™ Electrode with intact skin contact for a limited ( $\leq 24$  hours) duration

Functional performance equivalency was determined by electrical and adhesive benchtop testing, as follows:

- Electrical Continuity
- Adhesion Testing

Benchtop performance testing passed predetermined acceptance criteria, demonstrating Guardian Needle™ Electrodes are equivalent to the predicate device in functionality, safety, and effectiveness.

**807.92(b)(2)  
Clinical Tests**

No Clinical Tests were conducted as referenced in 21 CFR 807.92(b)(2).

**807.92(b)(3)  
Clinical Summary**

No Clinical Tests were conducted as referenced in 21 CFR 807.92(b)(3).

**Substantial Equivalence Table**

<b>Characteristic</b>	<b>Subject Device: Guardian Needle™ Electrodes</b>	<b>Predicate Device: Subdermal Needle Electrodes</b>	<b>Substantial Equivalence (SE)<sup>1</sup></b>
<b>510(k) Number</b>	K200984	K022914 Subdermal Needle Electrodes	
<b>Manufacturer</b>	Rhythmink International, LLC	Rhythmink International, LLC	SE
<b>Device Class</b>	Class II	Class II	SE
<b>Product Code</b>	GXZ	GXZ	SE
<b>Regulatory Name</b>	Subdermal Needle	Subdermal Needle	SE
<b>Device Type</b>	Subdermal Needle Electrodes	Subdermal Needle Electrodes	SE
<b>Regulation #</b>	21 CFR 882.1350	21 CFR 882.1350	SE
<b>Intended Use</b>	Rhythmink International Subdermal Needle Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals. Examples include: Electromyography (EMG), Electroencephalography (EEG), and Nerve potential signals. The electrodes are sterile and for single patient use.	Rhythmink International Subdermal Needle Electrodes are intended for use with recording, monitoring and stimulation equipment for the purpose of recording of biopotential signals. Examples include: Electromyography (EMG), Electroencephalography (EEG), and Nerve potential signals. The electrodes are sterile and for single patient use.	SE
<b>Anatomical Site(s)</b>	Subdermal muscle or nerve tissue	Subdermal muscle or nerve tissue	SE
<b>Environment usage</b>	Hospital	Hospital	SE
<b>Electrode Material</b>	Medical Grade 304VM (vacuum melt) Stainless Steel	Medical Grade 304VM (vacuum melt) Stainless Steel	SE
<b>Electrode Length</b>	7mm to 22 mm	7mm to 22 mm	SE
<b>Electrode Diameter</b>	0.4mm	0.4mm	SE
<b>Leadwire Material</b>	PVC-coated Electrical Wire	PVC-coated Electrical Wire	SE
<b>Leadwire Length</b>	1.0m to 3.0m	1.0m to 3.0m	SE
<b>Connector</b>	1.5mm DIN 42 802 pin touch proof connector	1.5mm DIN 42 802 pin touch proof connector	SE

<b>Characteristic</b>	<b>Subject Device: Guardian Needle™ Electrodes</b>	<b>Predicate Device: Subdermal Needle Electrodes</b>	<b>Substantial Equivalence (SE)<sup>1</sup></b>
<b>510(k) Number</b>	K200984	K022914 Subdermal Needle Electrodes	
<b>Use of Adhesive Tape</b>	YES	NO	SE; Material change has been validated by testing and does not adversely impact safety or effectiveness
<b>Needle Cover</b>	Thermoplastic polyester elastomer sheath	PVC tube	SE; Material and dimensional change has been validated by testing and does not adversely impact safety or effectiveness
<b>Operation of Applying the Device</b>	Applied subcutaneously to the skin	Applied subcutaneously to the skin	SE
<b>Prescription Use</b>	YES - Inserted by a licensed physician or practitioner	YES - Inserted by a licensed physician or practitioner	SE
<b>Targeted Procedures</b>	IONM, EMG, EP, EEG	IONM, EMG, EP, EEG	SE
<b>Compatibility with other devices</b>	Neurophysiology Monitors with a 1.5mm DIN 42 802 pin touch proof connector	Neurophysiology Monitors with a 1.5mm DIN 42 802 pin touch proof connector	SE
<b>Packaging</b>	24 electrodes sealed in a Tyvek pouch and placed inside a sealed, labeled box	24 electrodes sealed in a Tyvek pouch and placed inside a sealed, labeled box	SE
<b>Sterilization Method</b>	Supplied EtO Sterile	Supplied EtO Sterile	SE
<b>Single Patient Use</b>	YES – disposable	YES - disposable	SE
<b>Electrical Safety</b>	Connectors comply with IEC 60601-1 (1988) sub clause 56.3(c) per CFR 898.12	Connectors comply with IEC 60601-1 (1988) sub clause 56.3(c) per CFR 898.12	SE

<b>Characteristic</b>	<b>Subject Device: Guardian Needle™ Electrodes</b>	<b>Predicate Device: Subdermal Needle Electrodes</b>	<b>Substantial Equivalence (SE)<sup>1</sup></b>
<b>510(k) Number</b>	K200984	K022914 Subdermal Needle Electrodes	
<b>Mechanical Safety</b>	Leadwires are soldered to electrode using Tin/Silver solder and covered with heat shrink	Leadwires are soldered to electrode using Tin/Silver solder and covered with heat shrink	SE
<b>Duration of use</b>	≤24 hours	≤24 hours	SE
<b>Standards</b>	<ul style="list-style-type: none"> <li>• IEC 60601-1-1: 1988/a1: 1991/A2 :1995 § 56.3(c) per CFR 898.12</li> <li>• DIN 42802</li> <li>• ISO 14971</li> <li>• ISO 15223-1</li> </ul>	<ul style="list-style-type: none"> <li>• IEC 60601-1-1: 1988/a1: 1991/A2 :1995 § 56.3(c) per CFR 898.12</li> <li>• DIN 42802</li> <li>• ISO 14971</li> <li>• ISO 15223-1</li> </ul>	SE
<b>Biocompatible</b>	Yes	Yes	SE

<sup>1</sup>Only differences are indicated in this column.