



May 10, 2021

Anhui Bytech Medical Supplies Co., Ltd.  
% Ivy Wang  
Technical Manager  
Shanghai SUNGO Management Consulting Co., Ltd.  
13th F, 1500# Century Avenue  
Shanghai, Shanghai 200122  
China

Re: K201504

Trade/Device Name: Single-use medical poly (vinyl chloride) examination glove  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-Powdered Patient Examination Glove  
Regulatory Class: Class I, reserved  
Product Code: LYZ  
Dated: June 5, 2020  
Received: June 5, 2020

Dear Ivy Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmnmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray III, Ph.D.  
Acting Assistant Director  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K201504

Device Name

Single-use medical poly (vinyl chloride) examination glove

Indications for Use (Describe)

The Single-use medical poly (vinyl chloride) examination glove is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**K201504**

## 510(K) Summary

K201504

Date of preparation: 05 -10-2021

### A. Applicant

Name: ANHUI BYTECH MEDICAL SUPPLIES CO., LTD.

Address: LingBi Economic Development Zone (North), SuZhou City, Anhui Province

Name: Mr. Cheng Wang

Tel: +86-0557-6602888

Fax: +86-0557-6602888

Submission Correspondent:

Primary contact: Ms. Ivy Wang

Shanghai SUNGO Management Consulting Co., Ltd.

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-58817802

Email: [haiyu.wang@sungoglobal.com](mailto:haiyu.wang@sungoglobal.com)

Secondary contact: Mr. Raymond Luo

Room 1401, Dongfang Building, 1500# Century Ave., Shanghai 200122, China

Tel: +86-21-68828050

Email: [fda.sungo@gmail.com](mailto:fda.sungo@gmail.com)

### B. Device

Trade Name: Single-use medical poly (vinyl chloride) examination glove

Common Name: Vinyl Patient Examination Gloves (Powder Free)

Model(s): XS, S, M, L, XL

#### Regulatory Information

Classification Name: Vinyl Patient Examination Gloves (Powder Free)

Classification: Class I.

Product code: LYZ

Regulation Number: 21 CFR 880.6250

Review Panel: General Hospital

### C. Predicate device

<b>Sponsor</b>	HUIFU TRADING CO., LTD.
<b>Device Name</b>	Vinyl Examination Glove (Clear, Non-Colored)

## K201504

<b>510(k) Number</b>	K180849
<b>Product Code</b>	LYZ
<b>Regulation Number</b>	880.6250
<b>Regulation Class</b>	I

### D. Indications for use of the device:

The Single-use medical poly (vinyl chloride) examination glove is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.

### E. Device Description:

The proposed device is Powder Free Vinyl Patient Examination Gloves which is made of poly vinyl chloride. The proposed device is clear, non-colored. The design of proposed device is addressing the standards as ASTM D6124, ASTM D5151, and ASTM D5250. The proposed device is non-sterile.

### F. Comparison with predicate device

**Table 1 General Comparison**

ITEM	Proposed Device (K201504)	Predicate Device (K180849)	Remark
Product Code	LYZ	LYZ	SAME
Regulation No.	21 CFR 880.6250	21 CFR 880.6250	SAME
Class	I	I	SAME
Intended Use	The Single-use medical poly (vinyl chloride) examination glove is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.	The Vinyl Patient Examination Glove(Clear, Non-Colored) is a disposable device intended for medical purposes that is worn on the examiner's hands to prevent contamination between patient and examiner.	SAME
Powdered or Powered free	Powdered free	Powdered free	SAME
Design Feature	ambidextrous	ambidextrous	SAME
Surface Feature	Smooth	Smooth	SAME
Device Materials	Vinyl	Vinyl	SAME
Color	Clear	Clear	SAME

## K201504

Labeling Information	Single-use indication, powder free, device color, device name, glove size and quantity, Vinyl Examination Gloves, Non-Sterile	Single-use indication, powder free, device color, device name, glove size and quantity, Vinyl Examination Gloves, Non-Sterile	SAME
----------------------	---	---	------

**Table 2 Device Dimensions Comparison**

ITEM	Proposed Device (K201504)					Predicate Device (K180849)				Remark
Size	XS	S	M	L	XL	S	M	L	XL	similar
Length, mm	Minimum 230					Minimum 230	Minimum 235	Minimum 245	Minimum 245	similar
Width, mm	78±3	85±3	95±3	105±3	115±3	85±5	95±5	105±5	115±5	similar
Thickness, mm	Finger: minimum 0.05					Finger: minimum 0.08				Different
	Palm: minimum 0.08					Palm: minimum 0.08				Same

Discussion: The proposed device has slight difference on the dimensions of Length, width and finger thickness to the predicate device. However, the testing results meet the specification as stated in ASTM D5250. This difference will not affect the safety and effectiveness of the proposed device.

**Table 3 Performance Comparison**

ITEM		Proposed Device (K201504)	Predicate Device (K180849)	Remark
Colorant		Clear, Non-Colored	Clear, Non-Colored	Same
Physical Properties	Before Aging	Tensile Strength	11 MPa, min	Similar
		Ultimate Elongation	300 % min	
	After Aging	Tensile Strength	11 MPa, min	
		Ultimate Elongation	300 % min	
Freedom from Holes		Be free from holes when tested in accordance with ASTM D5151 AQL = 1.5	Be free from holes when tested in accordance with ASTM D5151 AQL =1.5	Same

## K201504

Powder Content	Meets Applicable Definition for Powder Free: $\leq 2$ mg per glove	0.5 mg	Similar
----------------	--	--------	---------

Discussion: The proposed device has slight difference on the performance of physical properties and powder content to the predicate device. However, the testing results meet the specification as stated in ASTM D5250 and ASTM D6124. This difference will not affect the effectiveness of the proposed device. The proposed device use the same AQL=1.5 as the predicate device, which is the requirement for European market and also meet the specification as stated in ASTM D5151. This will not affect the safety and effectiveness of the proposed device.

**Table 4 Safety Comparison**

ITEM		Predicate Device (K180849)	Proposed Device (K201504)	Remark
Biocompatibility	Irritation	Comply with ISO 10993-10	Under the conditions of the study not an irritant	Similar
	Sensitization		Under conditions of the study, not a sensitizer.	
	Cytotoxicity	Comply with ISO 10993-5	Under conditions of the study, did not show potential toxicity to L-929 cells.	
Label and Labeling		Meet FDA's Requirements	Meet FDA's Requirements	Same

### G. Summary of non-clinical performance testing

Test item		Test method	Pass criteria	Test results/ Verdict
Dimension (mm)	Overall Length	ASTM D5250	Min 230	>230/ Pass
	Width		XS: 76 $\pm$ 6 S: 85 $\pm$ 5 M: 95 $\pm$ 5 L: 105 $\pm$ 5 XL: 115 $\pm$ 5	XS: 79-81/ Pass S: 85-87 /Pass M: 95-97/ Pass L: 105-107/ Pass XL: 115-117/ Pass
	Finger thickness		0.05	0.09-0.12/ Pass
	Palm thickness		0.08	0.08-0.09/ Pass
Freedom from Holes		ASTM D5151	Meet the requirements of ASTM D5151 AQL 1.5	0/125, 1/ 125, 4/315, 5/500, 3/200 leaks / Pass
Powder Content		ASTM D6124	Meet the requirements of ASTM D6124 < 2.0mg	0.3 – 0.42/ Pass
Physical properties	Tensile Strength (MPa)	ASTM D412	Meet the requirements of ASTM D412 Min 11	14-19/ Pass
	Ultimate Elongation (%)		Meet the requirements of ASTM D412 Min 300%	390 – 430/ Pass
Biocompatibility	Irritation	ISO 10993-10	Non-irritating	Under the conditions of

## K201504

				the study, not an irritant/ Pass
	Sensitization	ISO 10993-10	Non-sensitizing	Under conditions of the study, not a sensitizer./ Pass
	Cytotoxicity	ISO 10993-5	Non-cytotoxic	Under conditions of the study, did not show potential toxicity to L-929 cells./ Pass

### Biocompatibility testing:

- **Skin sensitization & Irritation** -- ISO 10993-10:2010 Biological Evaluation of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.
- **Cytotoxicity** -- ISO 10993-5:2009 Biological Evaluation of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity

### Performance testing:

- **Residue Powder** -- ASTM D6124-06 (Reaffirmation 2017), Standard Test Method for Residual Powder on Medical Gloves
- **Freedom from holes** -- ASTM D5151-06 (Reapproved 2015), Standard Test Method for Detection of Holes in Medical Gloves.
- **Physical properties & Dimensions** -- ASTM D5250-06 (Reapproved 2015), Standard Specification for Poly (vinyl chloride) Gloves for Medical Application.
- ISO 2859-1:1999, "Sampling Procedures for Inspection by Attributes – Part I: Sampling Plans Indexed by Acceptable Quality Level (AQL) for Lot-by-Lot Inspection.

### H. Summary of Clinical Performance Test

No clinical study is included in this submission.

### I. Conclusion

The conclusion drawn from the nonclinical tests demonstrates that the subject device in 510(K) submission K201504, the Single-use medical poly (vinyl chloride) examination glove is as safe, as effective, and performs as well as or better than the legally marketed predicate device cleared under K180849.