



August 9, 2021

Edma Group, LLC  
% Boyle Wang  
Official Correspondent  
Shanghai Truthful Information Technology Co., Ltd.  
RM. 1801, No. 161, Lujiazui East Rd., Pudong  
Shanghai, 200120  
China

Re: K211540

Trade/Device Name: Edma Synthetic Nitrile Examination Gloves  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LZA  
Dated: May 12, 2021  
Received: May 19, 2021

Dear Boyle Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray III, Ph.D.  
DHT4B: Division of Infection Control  
and Plastic Surgery Devices  
OHT4: Office of Surgical  
and Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K211540

Device Name  
Edma Synthetic Nitrile Examination Gloves

Indications for Use (Describe)

A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hands or finger to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

---

---

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

---

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# 510(k) Summary

## (K211540)

This summary of 510(k) is being submitted in accordance with 21 CFR 807.92.

### **1.0 Submitter's Information**

Name: Edma Group, LLC.  
Address: 3634 E Piccadilly Rd, Phoenix, AZ 85018.  
Contact: Mr. Vio Cretu  
Date of Preparation: 07/03/2021

### **Designated Submission Correspondent**

Mr. Boyle Wang  
Shanghai Truthful Information Technology Co., Ltd.  
Room 1801, No. 161 Lujiazui East Rd., Pudong Shanghai, 200120 China  
Tel: +86-21-50313932  
Email: [Info@truthful.com.cn](mailto:Info@truthful.com.cn)

### **2.0 Device Information**

Trade name: Edma Synthetic Nitrile Examination Gloves  
Common name: Patient Examination Glove  
Classification name: Non-powdered patient examination glove  
Model(s): M, L, XL

### **3.0 Classification**

Production code: LZA  
Regulation number: 21CFR880.6250  
Classification: Class I  
Panel: General Hospital

### **4.0 Predicate Device Information**

Manufacturer: Careglove Global Sdn Bhd  
Device: Powder Free Nitrile Examination Gloves, Blue (colored)  
510(k) number: K172015

### **5.0 Indication for Use**

A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hands or finger to prevent contamination between patient and examiner.

### 6.0 Device Description

The subject device is powder free synthetic Nitrile patient examination gloves. The subject device is blue. The design of subject device is addressing the standards as ASTM D6124, ASTM D5151, and ASTM D6319. The subject device is non-sterile.

### 7.0 Technological Characteristic Comparison Table

Item	Subject device	Predicate device	Comparison			
510(k) number	K211540	K172015	--			
Product Code	LZA	LZA	Same			
Regulation No.	21CFR880.6250	21CFR880.6250	Same			
Class	I	I	Same			
Intended Use	A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hands or finger to prevent contamination between patient and examiner.	A patient examination glove is a disposable device intended for medical purposes that is worn upon the examiner's hand or finger to prevent contamination between patient and examiner.	Same			
Powdered or Powdered free	Powdered free	Powdered free	Same			
Main Material	Nitrile	Nitrile	Same			
Colorant	Blue	Blue	Same			
Design Feature	Ambidextrous	Ambidextrous	Same			
Single Use	Yes	Yes	Same			
Sterility status	Non-Sterile	Non-Sterile	Same			
Dimensions(mm)	Length: $\geq 230$ ; Width: M: $95 \pm 10$ ; L: $110 \pm 10$ ; XL: $120 \pm 10$ ;	Length: XS/S: $\geq 220$ ; M/L/XL: $\geq 230$ ; Width: XS: $70 \pm 10$ ; S: $80 \pm 10$ ; M: $95 \pm 10$ ; L: $110 \pm 10$ ; XL: $120 \pm 10$ ;	Similar			
Physical	Before	Tensile	14MPa,min	Tensile	14MPa,min	Same

Properties	Aging	Strength		Strength		
		Ultimate Elongation	500%min	Ultimate Elongation	500%min	Same
	After Aging	Tensile Strength	14MPa,min	Tensile Strength	14MPa,min	Same
		Ultimate Elongation	400%min	Ultimate Elongation	400%min	Same
Freedom from Holes	Be free from holes when tested in accordance with ASTM D5151 AQL=2.5		Meet AQL 1.5 Meet AQL 2.5		Same	
Powder Content	<0.07 mg per glove. Meet the requirements of ASTM D6124		Meet 2mg/glove max. Meet the requirements of ASTM D6124		Same	
Biocompatibility	Irritation: Under the conditions of the study, not an irritant or a sensitizer.		Under the conditions of this study the test material did not cause an irritant response.		Same	
	Sensitization: Under conditions of the study, not a sensitizer.		Under the conditions of this study, the test material did not produce a skin sensitization effect			
	Cytotoxicity: Under conditions of the study, did not show potential toxicity to L-929 cells.		/			
Labeling Information	Single use, powder free, device color, device name, glove size and quantity, Non-Sterile		Single use, powder free, device color, device name, glove size and quantity, Non-Sterile		Same	

## **8.0 Summary of Non-clinical Testing**

Non clinical tests were conducted to verify that the subject device met all design specifications as was Substantially Equivalent (SE) to the predicate device. The subject device were evaluated according to the following standards:

ASTM D6319-19, *Standard Specification for Nitrile Examination Gloves for Medical Application*.

ASTM D6124-06 (Reapproved 2017), *Standard Test Method for Residual Powder on Medical Gloves*

ASTMD5151-19, *Standard Test Method for Detection of Holes in Medical Gloves.*

**Biocompatibility Testing:**

The biocompatibility evaluation for the subject device were evaluated according to the following standard:

ISO 10993-5:2009 *Biological Evaluation Of Medical Devices - Part 5: Tests For In Vitro Cytotoxicity.*

ISO 10993-10:2010 *Biological Evaluation Of Medical Devices - Part 10: Tests For Irritation And Skin Sensitization.*

Biocompatibility testing including cytotoxicity test, sensitization test and irritation test according to ISO 10993-1 standards, have been conducted on the Edma Synthetic Nitrile Examination Gloves.

Under the parameters of the tests it is concluded that they are biocompatible, and that there are no new issues of biocompatibility regarding their use as intended.

Table 2: Performance Characteristic Comparison

Test Methodology	Purpose	Acceptance Criteria	Results
ASTM D6319	Physical Dimensions Test	Length(mm):≥230; Width(mm): M: 95±10; L: 110±10; XL: 120±10	Length: > 230 Width: M: 96-99 L: 106-109 XL: 114-118 <u>Pass</u>
		Thickness (mm): Finger: ≥0.08 Palm: ≥0.08	Finger: 0.08-0.13 Palm: 0.08-0.09 <u>Pass</u>
ASTM D5151	Watertightness Test for Detection of Holes	Meet the requirements of ASTM D5151 AQL 2.5	0/125 leaks <u>Pass</u>
ASTM D6124	Powder Content	Meet the requirements of ASTM D6124 < 2.0mg	0.07 mg <u>Pass</u>
ASTM D412	Physical properties	Before Aging   Tensile Strength   ≥14MPa	15-18.5 <u>Pass</u>

			Ultimate Elongation	≥500%	506-576 <u>Pass</u>
		After Aging	Tensile Strength	≥14MPa	14-17.6 <u>Pass</u>
			Ultimate Elongation	≥400%	400-522 <u>Pass</u>
ISO 10993-5	Cytotoxicity	Non-cytotoxic			Under conditions of the study, did not show potential toxicity to L-929 cells. <u>Pass</u>
ISO 10993-10	Irritation	Non-irritating			Under the conditions of the study, not an irritant. <u>Pass</u>
ISO 10993-10	Sensitization	Non-sensitizing			Under conditions of the study, not a sensitizer. <u>Pass</u>

### **9.0 Summary of Clinical Testing**

Clinical testing is not needed for this device.

### **10.0 Conclusion**

The conclusions drawn from the nonclinical tests demonstrate that the proposed device, K211540 is as safe, as effective, and performs as well as or better than the legally marketed predicated device.