



July 9, 2021

Magic Mobility
% Michelle Rubin-Onur
Senior Regulatory Specialist
AcKnowledge Regulatory Strategies, LLC
2251 San Diego Ave
Suite B-257
San Diego, California 92110

Re: K211574

Trade/Device Name: Extreme X8 Power Wheelchair
Regulation Number: 21 CFR 890.3860
Regulation Name: Powered Wheelchair
Regulatory Class: Class II
Product Code: ITI
Dated: May 20, 2021
Received: May 21, 2021

Dear Michelle Rubin-Onur:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531 - 542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Heather Dean, PhD
Assistant Director, Acute Injury Devices Team
DHT5B: Division of Neuromodulation
and Physical Medicine Devices
OHT5: Office of Neurological
and Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K211574

Device Name

Extreme X8 Power Wheelchair

Indications for Use (Describe)

The Extreme X8 Power Wheelchair is a battery-operated device with wheels that is intended for medical purposes to provide mobility to persons restricted to a sitting position who have the capability of operating a power wheelchair.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”

510(k) Summary

K211574

DATE PREPARED

July 9, 2021

MANUFACTURER AND 510(k) OWNER

Magic Mobility

3 International Court, Scoresby, VIC, 3179 Australia

Telephone: +61 3 8791 5600

Official Contact: Alex Suen, Quality Assurance and Regulatory Affairs Manager

REPRESENTATIVE/CONSULTANT

Michelle Rubin-Onur, Ph.D.

Matthieu Kirkland, M.S.

AcKnowledge Regulatory Strategies, LLC

Telephone: +1 (619) 458-9547

Email: mrubin@acknowledge-rs.com

DEVICE INFORMATION

Proprietary Name/Trade Name: Extreme X8 Power Wheelchair

Common Name: Wheelchair, Powered

Regulation Number: 21 CFR 890.3860

Class: Class II

Product Code: ITI

Premarket Review: Neuromodulation and Physical Medicine Devices (DHT5B)

Review Panel: Physical Medicine

PREDICATE DEVICE IDENTIFICATION

The Extreme X8 Power Wheelchair is substantially equivalent to the following predicates:

<i>510(k) Number</i>	<i>Predicate Device Name / Manufacturer</i>	<i>Predicate Relationship</i>
K090350	Extreme 4x4-X8 / Innovation In Motion	Predicate Device
K172384	Quickie® Q700-UP M / Sunrise Medical (US) LLC	Reference Device

DEVICE DESCRIPTION

The Extreme X8 Power Wheelchair is designed for everyday use for both indoor and outdoor environments including care facilities and private residences. The subject device is intended to provide mobility to persons who are restricted or limited to a sitting position.

The Extreme X8 Power Wheelchair is a battery powered, electric motor driven device that can be used on both indoor and outdoor surfaces (i.e., concrete, asphalt, indoor flooring such as carpet, gravel, grass, and bark/woodchips). The Extreme X8 Power Wheelchair offers two basic seating options: MPS and Rehab. The MPS option is more static and does not allow for additional aftermarket cushions. The Rehab option is more adjustable to adhere to recommendations from the user's physician or physical therapist.

The Extreme X8 Power Wheelchair includes the following accessories:

- Extra spreader bar
- Slide in table
- Lights
- Luggage rack
- Accessory charger
- Posture belt
- Roho cushion
- Jay cushion
- MPS push rail
- MPS peg push handle
- Scooter stopper
- Retractable docking pin
- Fold forward kit
- Steering lock

INTENDED USE

Magic Mobility power chairs are designed for the exclusive use of people (adults and children) who are unable to walk or have limited mobility and have the cognitive, physical and visual ability to control the vehicle safely. The Magic Mobility Extreme X8 is intended to be self-propelled on a range of surfaces. The All-Terrain tires can be used on both indoor and outdoor surfaces (i.e., concrete, asphalt, indoor flooring such as carpet, gravel, grass, and bark/woodchips).

INDICATIONS FOR USE

The Extreme X8 Power Wheelchair is a battery-operated device with wheels that is intended for medical purposes to provide mobility to persons restricted to a sitting position who have the capability of operating a power wheelchair.

COMPARISON OF TECHNOLOGICAL CHARACTERISTICS

Magic Mobility believes that the Extreme X8 Power Wheelchair is substantially equivalent to the predicate device based on the information summarized here:

The subject device has a similar design and dimensions and uses similar or identical materials as the device cleared in K090350. The subject device has the same intended use and similar technological characteristics (i.e., base technology and OEM joystick control) to the devices cleared in K093050 and K172384. The subject device has the same intended use environment, including off road capabilities, as the device cleared in K090350. The subject device uses the same software as the device cleared in K172384. The Extreme X8 Power Wheelchair has undergone testing to ensure that any differences in technological characteristics (i.e., battery, castor wheels, and no anti-pitch mechanism) do not affect safety and effectiveness when compared to the predicate device.

	Subject Device	Predicate Device	Reference Device	
	Magic Mobility Extreme X8 Power Wheelchair	Innovation In Motion Extreme 4x4-X8 K090350	Sunrise Medical (US) LLC Quickie® Q700-UP M K172384	Statement of Equivalence
Indications for Use	The Extreme X8 Power Wheelchair is a battery-operated device with wheels that is intended for medical purposes to provide mobility to persons restricted to a sitting position who have the capability of operating a power wheelchair.	The intended use of the X8 power wheelchair is to provide mobility to persons limited to a sitting position, who have the capability of operating a power wheelchair.	The Sunrise Medical Quickie® Q700-UP M power wheelchairs are battery operated devices, that are indicated for medical purposes to provide mobility and repositioning of the user, including a stand-up feature.	Identical to predicate device. No impact on safety or effectiveness.
Product Codes / Regulation Number	ITI /21 CFR 890.3860	ITI / 21 CFR 890.3860	IPL / 21 CFR 890.3900	Identical to predicate device. No impact on safety and effectiveness.
Regulation Description	Powered Wheelchair	Powered Wheelchair	Standup Wheelchair	Identical to predicate device. No impact on safety and effectiveness.
Technical Specifications				
General				
Maximum User Weight (lbs.)	400	400	265	Identical to predicate device. No impact on safety and effectiveness.
Storage Temperature (°C)	-40 to 70°C	-40 to 70°C	-40 to 70	Identical to predicate and reference devices. No impact on safety and effectiveness.
Location for Use	Indoors and outdoors including care facilities, residences, and soft/rough terrain.	Indoors and outdoors including care facilities, and residences	Indoors and outdoors including care facilities, and residences	Identical to predicate and reference devices. No impact on safety and effectiveness.
Frame Material	Steel	Steel	Steel and aluminum	Identical to predicate and reference devices. No impact on safety and effectiveness.
Biocompatibility	Uses materials common to many wheelchairs	Uses materials common to many wheelchairs	Uses materials common to many wheelchairs	Identical to predicate and reference devices. No impact on safety and effectiveness.
Ground Clearance	4"	4"	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Base				
Overall Dimensions (length by width; inches)	45.25"x37"	45.25"x37"	25"x36"	Identical to predicate device. No impact on safety and effectiveness.

	Subject Device	Predicate Device	Reference Device	
	Magic Mobility Extreme X8 Power Wheelchair	Innovation In Motion Extreme 4x4-X8 K090350	Sunrise Medical (US) LLC Quickie® Q700-UP M K172384	Statement of Equivalence
Rolling Base Weight (lbs.)	145	145	152	Identical to predicate device. No impact on safety and effectiveness.
Power Source	Batteries	Batteries	Batteries	Identical to predicate and reference devices. No impact on safety and effectiveness.
Battery Details	2- 73 AMP/hour AGM Group 24 12-volt Weight 55 lbs	2- 73 AMP/hour AGM Group 24 12 volt Weight 55 lbs	24V (2x12V) / 73 Ah/20h	Identical to predicate device. No impact on safety and effectiveness.
Battery Charger Input Voltage (volts)	110	110	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Motors	4-24V, 4 pole direct drive gear in line motors/gearbox. Freewheeling lever/motor lock releases & engages rear motor brakes	4-24V, 4 pole direct drive gear in line motors/gearbox. Freewheeling lever/motor lock releases & engages rear motor brakes	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Wheel Size (inches)	4"x14"	4"x14"	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Turing Radius (inches)	52"	52"	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Range (miles)	15.5	15.5	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Maximum Speed (mph)	6.2	6.2	6 (with an option of 8)	Identical to predicate device. No impact on safety and effectiveness.
Anti-pitch Mechanism for Climbing	None	None	Additional anti-pitch lock out	Identical to predicate device. No impact on safety and effectiveness.
Lift Range (inches)	0-12	0-12	0-12	Identical to predicate device. No impact on safety and effectiveness.
Tilt Range (degrees)	0-50	0-50	0-50	Identical to predicate device. No impact on safety and effectiveness.
Recline Range (degrees)	0-170	0-170	0-172	Identical to predicate device. No

	Subject Device	Predicate Device	Reference Device	
	Magic Mobility Extreme X8 Power Wheelchair	Innovation In Motion Extreme 4x4-X8 K090350	Sunrise Medical (US) LLC Quickie® Q700-UP M K172384	Statement of Equivalence
				impact on safety and effectiveness.
Suspension	None	None	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Braking System	Electromagnetic, regenerative brakes with a free-wheeling mode	Electromagnetic, regenerative brakes with a free-wheeling mode	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Minimum Braking Distance at Maximum Speed (meters)	1.7	1.7	Unknown	Identical to predicate device. No impact on safety and effectiveness.
User Controller	Joystick and hand control buttons	Joystick and hand control buttons	Joystick and hand control buttons	Identical to predicate device. No impact on safety and effectiveness.
Joystick Mount	Fixed mount, height adjustable, swing-away	Fixed mount, height adjustable, swing-away	Fixed mount, height adjustable, swing-away	Identical to predicate device. No impact on safety and effectiveness.
Software				
Software	R-Net from PGDT	Dynamic DX2-PMA70L electronic control and joystick	R-Net from PGDT	Identical to reference device. No impact on safety and effectiveness.
Seat/Armrest/Footrest				
Seat Height (minimum, inches)	19"	19"	16.2"	Identical to predicate device. No impact on safety and effectiveness.
Seat Width (inches)	16"-24"	16"-24"	16"-22"	Identical to predicate device. No impact on safety and effectiveness.
Armrest	Height adjustable, removable, flip up option	Height adjustable, removable, flip up option	Unknown	Identical to predicate device. No impact on safety and effectiveness.
Footrest	1 or 2 pieces, fixed or flip up, angle and height adjustable rigid footplates, and others	1 or 2 pieces, fixed or flip up, angle and height adjustable rigid footplates, and others	Unknown	Identical to predicate device. No impact on safety and effectiveness.

SUMMARY OF NON-CLINICAL TESTING

The following tests were performed to demonstrate substantial equivalence based on current industry and FDA recognized standards. The results of these tests indicate that the Extreme X8 Power Wheelchair is substantially equivalent to the predicate device.

PERFORMANCE

- Static stability (per ISO 7176-1)
- Dynamic stability (per ISO 7176-2)
- Effectiveness of brakes (per ISO 7176-3)
- Energy consumption (per ISO 7176-4)
- Dimensions, mass, and maneuvering space (per ISO 7176-5)
- Maximum speed, acceleration, and deceleration (per ISO 7176-6)
- Measurement of seat and wheel dimensions (per ISO 7176-7)
- Static, impact, and fatigue (per ISO 7176-8)
- Climatic test (per ISO 7176-9)
- Obstacle climbing (per ISO 7176-10)
- Test dummies (per ISO 7176-11)
- Power and control systems for power wheelchairs (per ISO 7176-14)
- Documentation and labeling (per ISO 7176-15)
- Resistance to ignition (per ISO 7176-16)
- Dynamic Test (per ISO 7176-19)
- Vocabulary (per ISO 7176-26)

EMC AND ELECTRICAL SAFETY

- EMC testing (per ISO 7176-21)
- Batteries and chargers per (per ISO 7176-25)

BIOCOMPATIBILITY

- Evaluation and testing within a risk management process (per ISO 10993-1)
- *In vitro* cytotoxicity (per ISO 10993-5)

The materials used in the subject device are identical in chemical composition, formulation processing, sterilization, and geometry as the materials found in the Extreme 4x4-X8 Power Wheelchair by Magic Mobility LLC (K090350) and the Quickie Q700-UP M by Sunrise Medical (US) LLC (K172384) device. Furthermore, the Extreme X8 Power Wheelchair has the same nature of tissue contact and contact duration as the predicate and reference devices. Therefore, based on previous use and the cytotoxicity testing conducted, the subject device is considered to have met the requirements of ISO 10993-1 and FDA's Guidance for Industry and Food and Drug Administration Staff – Use of International ISO 10993-1, "Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.

SOFTWARE

- Software life cycle process (per IEC 62304)

CONCLUSION

Based on the testing performed (including wheelchair dynamic testing and flammability) it can be concluded that the subject device does not raise new issues of safety or effectiveness compared to the predicate device. The similar indications for use, technological characteristics, and performance characteristics for the proposed Extreme X8 Power Wheelchair are assessed to be substantially equivalent to the predicate device.