



July 28,2021

Taiji Medical Supplies, Inc.
% Abdel Halim
President, Global and Regulatory Affairs
Global Quality and Regulatory Services
10 Scenic Way
Monroe, New Jersey 08831

Re: K211861

Trade/Device Name: TAIJI Professional Series Level 3 Surgical Face Mask
Regulation Number: 21 CFR 878.4040
Regulation Name: Surgical apparel
Regulatory Class: Class II
Product Code: FXX
Dated: July 21, 2021
Received: July 23, 2021

Dear Abdel Halim:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

For Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K211861

Device Name

TAIJI Professional Series Level 3 Surgical Face Mask

Indications for Use (Describe)

The TAIJI Professional Series Level 3 Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluid. This is a single use, disposable device, provided non-sterile.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



510(k) Summary

This 510(k) Summary is being submitted in accordance with the requirements of *Guidance for Industry and FDA Staff: Surgical Masks – Premarket Notification [510(k)] Submissions*, Issued March 5, 2004 (corrected July 14, 2004).

1. Device and Company Details:

510(k) Number: K211861

Submitter / Owner

Taiji Medical Supplies, Inc.
3211 Progress Drive
Lincolnton, NC 28092
T: (828) 310-2802
E: tiffany@taijimedical.com

Contact Person

Ms. Katherine L. Giannamore, Esq.
Shehadeh Giannamore, PLLC
396 Alhambra Circle, Suite 100A
Coral Gables, FL 33134
T: (305) 507-9843
E: katherine@sglawfl.com

Date of Preparation of this Summary: June 10, 2021

Device Trade or Proprietary Name: TAIJI Professional Series Level 3 Surgical Face Mask

Device Common/Usual Name: Surgical Mask

Regulation Name/Device Class/
Product Code: Surgical Apparel/Class II/FXX

Predicate Device Name/K Number: Premier Guard USA 3 Layer Ear Loop
ASTM Level 3 Surgical Face Mask/K202595

2. Device Description: The TAIJI Professional Series Level 3 Surgical Face Mask is in accordance with the FDA Guidance Document, Surgical Masks - Premarket Notification [510K] Submissions issued on March 5, 2004. TAIJI Professional Series Level 3 Surgical Face mask is a flat-pleated style mask with elastic ear loops to secure it over the users' mouth and face. The mask consists of three-layers. The inner facing layer is white and is manufactured from spunbond polypropylene (three layers of nonwoven polypropylene). The inner filter material is made of meltblown fiber. The outer facing layer is blue and is manufactured from spunbond polypropylene (three layers of nonwoven polypropylene).

The mask is a single use, disposable device, provided nonsterile.



Intended Use: The TAIJI Professional Series Level 3 Surgical Face Mask is a non-sterile, single-use, disposable mask, intended to be worn in healthcare and surgical procedures.

3. Technological Characteristics Comparison:

The following is a summary of the technological characteristics of the TAIJI Professional Series Level 3 Surgical Face Mask as compared to the predicate device.

Items	Subject Device TAIJI Professional Series Level 3 Surgical Face Mask	Predicate Device K202595	Comparison
Manufacturer	TAIJI Medical Supplies, Inc.	Premier Guard USA, LLC	N/A
510(k) Number	K211861	K202595	N/A
FDA Product Code	FXX	FXX	Same
Indications for Use	The TAIJI Professional Series Level 3 Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluid. This is a single use, disposable device, provided non-sterile	The Premier Guard USA 3 Layer Ear Loop ASTM Level 3 Surgical Face Masks are intended to be worn to protect both the patient and healthcare personnel from transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and body fluid. This is a single use, disposable device, provided non-sterile	Same

Materials	Inner and Outer Layers	Spun-bond polypropylene	Spun-bond polypropylene	Same
	Middle Layer	Melt blown polypropylene filter	Melt blown polypropylene filter	Same
	Ear Loops	Polyester with Spandex inner-core	Nylon and Spandex	Similar
	Nose Piece	Malleable polyethylene with aluminum wire	Polyethylene laminated soft annealed carbon steel wire	Similar
Dimensions		17.5cm length x 9.5cm height	17.5cm length x 9.5cm height	Same
Mask Style		Flat Pleated	Flat Pleated	Same
Design Features		Malleable nosepiece, flat-pleated elastic ear loops	Malleable nosepiece, flat-pleated elastic ear loops	Same
Sterility		Non-sterile	Non-sterile	Same
Use		Single Use, Disposable	Single Use, Disposable	Same
Color		Blue and White	Blue and White	Same
ASTM F2100 Level		Level 3	Level 3	Same
Biocompatibility				
Cytotoxicity Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity		Under conditions of the study, the device considered is cytotoxic (neat and 1:2 extracts were scored 3)	Under conditions of the study, the device is Non-cytotoxic	Fail - Acute Systemic Study was recommended by 3PR and agreed upon by FDA
Systemic Study / Systemic Injection Test / ISO GLP ISO 10993-11: Biological evaluation of medical devices - Part 11: Tests for systemic toxicity		Under conditions of the study, the device is non-cytotoxic	N/A	Pass No systemic toxicity: device is not toxic
Irritation		Under conditions of the study, the device is Non-irritating	Under conditions of the study, the device is Non-irritating	Same

Sensitization	Under conditions of the study, the device is Non-sensitizing	Under conditions of the study, the device is Non-sensitizing	Same
Fluid Resistance Performance (mmHg)	Pass at 160 mmHg	Pass at 160 mmHg	Same
Particulate Filtration Efficiency Performance (%)	Pass at $\geq 98\%$	Pass at $\geq 98\%$	Same
Bacterial Filtration Efficiency Performance (%)	Pass at $\geq 98\%$	Pass at $\geq 98\%$	Same
Differential Pressure (Delta-P) (mm H ₂ O/cm)	Pass at ≤ 6.0	Pass at < 6.0	Same
Flammability class Class 1	Pass Class 1	Pass Class 1	Same

4. Summary of Non-Clinical Testing

Per FDA document Guidance for Industry and FDA Staff: Surgical Masks - Premarket Notification [510(k)] Submissions, the below testing has been completed on the subject device:

Standard/Test Methodology	Purpose	Acceptance Criteria	Results
ASTM F1862: Standard Test Method for Resistance of Medical Face Masks to Penetration by Synthetic Blood	Barrier Testing	At least 29 out of 32 specimens show passing results at 160 mmHg	Samples met the predetermined acceptance criteria
ASTM F2101-19: Standard Specification for Performance of Materials Used in Medical Face Masks	Barrier Testing	≥98%	Samples met the predetermined acceptance criteria 96 Samples, 3 non-consecutive lots of 32 each
EN 14683:2019 Annex C and ASTM F2100-19: Standard Specification for Performance of Materials Used in Medical Face Masks	Physical Testing	<6.0 mm H ₂ O/cm ²	Samples met the predetermined acceptance criteria 96 Samples, 3 non-consecutive lots of 32 each
ASTM F2100-19: Standard Specification for Performance of Materials Used in Medical Face Masks	Barrier Testing	≥98%	Samples met the predetermined acceptance criteria 96 Samples, 3 non-consecutive

			lots of 32 each
21 CFR 1610	Safety Testing	Class 1, Does not ignite	Samples met the predetermined acceptance criteria 96 Samples, 3 non-consecutive lots of 32 each
ISO 10993-5: Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity	Safety Testing	Non-cytotoxic	Fail
ISO 10993-11: Biological evaluation of medical devices - Part 11: Tests for systemic toxicity	Safety Testing	Under conditions of the study, the device is non-cytotoxic	Pass No systemic toxicity: device is not toxic
ISO 10993-10: Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization	Safety Testing	Non-irritating	Pass
ISO 10993-10: Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization	Safety Testing	Non-sensitizing	Pass



5. Conclusion:

The conclusions drawn from the nonclinical tests demonstrate that the subject device K211861 TAIJI Medical Professional Series Level 3 Surgical Face Mask is as safe, as effective and performs as well or better than the legally marketed predicate device, K202595 Premier Guard USA 3 Layer Ear Loop ASTM Level 3 Surgical Face Mask.