



June 24, 2022

3M Company
Yumi Wackerfuss
Senior Regulatory Affairs Associate
2510 Conway Avenue, Bldg. 275-5W-06
Saint Paul, Minnesota 55144-1000

Re: K220942

Trade/Device Name: 3M Attest Steam Chemical Integrators (1243A, 1243B, 1243RE, 1243RES)
Regulation Number: 21 CFR 880.2800
Regulation Name: Sterilization Process Indicator
Regulatory Class: Class II
Product Code: JOJ
Dated: May 25, 2022
Received: May 27, 2022

Dear Yumi Wackerfuss:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part

801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Clarence W. Murray, III, PhD
Assistant Director
DHT4B: Division of Infection Control
and Plastic Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)
K220942

Device Name
3M™ Attest™ Steam Chemical Integrators (1243A, 1243B, 1243RE, 1243RES)

Indications for Use (Describe)

The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles.

The integrating indicator is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles:

Cycle Type	Temperature	Exposure Time
Gravity	250°F/121°C	30 minutes
Gravity	270°F/132°C	3 minutes
Gravity	270°F/132°C	4 minutes
Gravity	270°F/132°C	10 minutes
Gravity	270°F/132°C	15 minutes
Gravity	270°F/132°C	25 minutes
Gravity	275°F/135°C	3 minutes
Gravity	275°F/135°C	10 minutes
Dynamic Air Removal	250°F/121°C	30 minutes
Dynamic Air Removal	270°F/132°C	4 minutes
Dynamic Air Removal	270°F/132°C	5 minutes
Dynamic Air Removal	270°F/132°C	6 minutes
Dynamic Air Removal	270°F/132°C	7 minutes
Dynamic Air Removal	270°F/132°C	8 minutes
Dynamic Air Removal	270°F/132°C	9 minutes
Dynamic Air Removal	270°F/132°C	10 minutes
Dynamic Air Removal	273°F/134°C	3 minutes
Dynamic Air Removal	273°F/134°C	4 minutes
Dynamic Air Removal	275°F/135°C	3 minutes

Minimum Stated Values for 3M™ Attest™ Steam Chemical Integrators as determined in a resistometer:

250°F/121°C 270°F/132°C 273°F/134°C 275°F/135°C
16.5 Minutes 2.0 Minutes 1.4 Minutes 1.2 Minutes

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



**510(k) Summary
for
3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES**

Sponsor Information:

3M Company
3M Center, Bldg. 275-5W-06
St. Paul, MN 55144-1000

Contact: Yumi Wackerfuss
Senior Regulatory Affairs Associate
Phone Number: (651) 737-3556
Fax Number: (651) 737-5320

Date of Summary: March 29th, 2022



**3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES**

**510(k)
Summary**

1. Device Name and Classification:

Common or Usual Name Chemical Indicators
 Proprietary Name: 3M™ Attest™ Steam Chemical Integrators
 Classification Name: Physical/chemical sterilization process indicator
 Device Classification: Class II, 21 CFR § 880.2800(b)
 Product Code: JOJ

2. Predicate Device:

K193254, 3M™ Attest™ Steam Chemical Integrators

3. Description of Device:

3M™ Attest™ Steam Chemical Integrators are chemical indicators consisting of a paper wick and a steam and temperature sensitive chemical pellet contained in a paper/film/foil laminate. The chemical pellet melts and migrates as a dark color along the paper wick. The migration is visible through a window marked ACCEPT or window marked REJECT; the extent of migration depends on steam, time, and temperature.

4. Indications for Use

The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles. The integrating indicator is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles:

Cycle Type	Temperature	Exposure Time
Gravity	250°F/121°C	30 minutes
Gravity	270°F/132°C	3 minutes
Gravity	270°F/132°C	4 minutes
Gravity	270°F/132°C	10 minutes
Gravity	270°F/132°C	15 minutes
Gravity	270°F/132°C	25 minutes
Gravity	275°F/135°C	3 minutes
Gravity	275°F/135°C	10 minutes
Dynamic Air Removal	250°F/121°C	30 minutes
Dynamic Air Removal	270°F/132°C	4 minutes
Dynamic Air Removal	270°F/132°C	5 minutes
Dynamic Air Removal	270°F/132°C	6 minutes
Dynamic Air Removal	270°F/132°C	7 minutes
Dynamic Air Removal	270°F/132°C	8 minutes
Dynamic Air Removal	270°F/132°C	9 minutes
Dynamic Air Removal	270°F/132°C	10 minutes



**3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES**

**510(k)
Summary**

Cycle Type	Temperature	Exposure Time
Dynamic Air Removal	273°F/134°C	3 minutes
Dynamic Air Removal	273°F/134°C	4 minutes
Dynamic Air Removal	275°F/135°C	3 minutes

Minimum Stated Values for 3M™ Attest™ Steam Chemical Integrators as determined in a resistometer:

250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C
16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes



510(k) Summary
3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES

510(k) Summary

5. Summary of Technological Characteristics compared to Predicate Device

Element	Subject Devices: (this submission) 3M™ Attest™ Steam Chemical Integrators	Predicate Device (K193254: Cleared on Dec 23 rd , 2019) 3M™ Attest™ Steam Chemical Integrators	Comparison																																																																											
Device Models	1243A, 1243B, 1243RE and 1243RES	1243A and 1243B	1243A and 1243B are substantially equivalent as 1243RE and 1243RES (K191236).																																																																											
Indications for use	<p>The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles. The integrating indicator is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles:</p> <table border="1" data-bbox="373 776 968 1424"> <thead> <tr> <th>Cycle Type</th> <th>Temperature</th> <th>Exposure Time</th> </tr> </thead> <tbody> <tr><td>Gravity</td><td>250°F/121°C</td><td>30 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>3 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>4 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>10 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>15 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>25 minutes</td></tr> <tr><td>Gravity</td><td>275°F/135°C</td><td>3 minutes</td></tr> <tr><td>Gravity</td><td>275°F/135°C</td><td>10 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>250°F/121°C</td><td>30 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>4 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>5 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>6 minutes</td></tr> </tbody> </table>	Cycle Type	Temperature	Exposure Time	Gravity	250°F/121°C	30 minutes	Gravity	270°F/132°C	3 minutes	Gravity	270°F/132°C	4 minutes	Gravity	270°F/132°C	10 minutes	Gravity	270°F/132°C	15 minutes	Gravity	270°F/132°C	25 minutes	Gravity	275°F/135°C	3 minutes	Gravity	275°F/135°C	10 minutes	Dynamic Air Removal	250°F/121°C	30 minutes	Dynamic Air Removal	270°F/132°C	4 minutes	Dynamic Air Removal	270°F/132°C	5 minutes	Dynamic Air Removal	270°F/132°C	6 minutes	<p>The 3M™ Attest™ Steam Chemical Integrators are designed to respond to all critical parameters over a specified range of steam sterilization cycles. The integrating indicator is intended to be placed in each pack, pouch, container, tray or other containment device to function as an independent monitor of critical parameters for the following sterilization cycles:</p> <table border="1" data-bbox="1071 813 1665 1393"> <thead> <tr> <th>Cycle Type</th> <th>Temperature</th> <th>Exposure Time</th> </tr> </thead> <tbody> <tr><td>Gravity</td><td>250°F/121°C</td><td>30 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>3 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>4 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>10 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>15 minutes</td></tr> <tr><td>Gravity</td><td>270°F/132°C</td><td>25 minutes</td></tr> <tr><td>Gravity</td><td>275°F/135°C</td><td>3 minutes</td></tr> <tr><td>Gravity</td><td>275°F/135°C</td><td>10 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>250°F/121°C</td><td>30 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>4 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>10 minutes</td></tr> </tbody> </table>	Cycle Type	Temperature	Exposure Time	Gravity	250°F/121°C	30 minutes	Gravity	270°F/132°C	3 minutes	Gravity	270°F/132°C	4 minutes	Gravity	270°F/132°C	10 minutes	Gravity	270°F/132°C	15 minutes	Gravity	270°F/132°C	25 minutes	Gravity	275°F/135°C	3 minutes	Gravity	275°F/135°C	10 minutes	Dynamic Air Removal	250°F/121°C	30 minutes	Dynamic Air Removal	270°F/132°C	4 minutes	Dynamic Air Removal	270°F/132°C	10 minutes	<p>Similar. This submission proposes the additional exposure times for 5, 6, 7, 8 and 9 minutes for Dynamic Air Removal, 270°F/132°C.</p>
Cycle Type	Temperature	Exposure Time																																																																												
Gravity	250°F/121°C	30 minutes																																																																												
Gravity	270°F/132°C	3 minutes																																																																												
Gravity	270°F/132°C	4 minutes																																																																												
Gravity	270°F/132°C	10 minutes																																																																												
Gravity	270°F/132°C	15 minutes																																																																												
Gravity	270°F/132°C	25 minutes																																																																												
Gravity	275°F/135°C	3 minutes																																																																												
Gravity	275°F/135°C	10 minutes																																																																												
Dynamic Air Removal	250°F/121°C	30 minutes																																																																												
Dynamic Air Removal	270°F/132°C	4 minutes																																																																												
Dynamic Air Removal	270°F/132°C	5 minutes																																																																												
Dynamic Air Removal	270°F/132°C	6 minutes																																																																												
Cycle Type	Temperature	Exposure Time																																																																												
Gravity	250°F/121°C	30 minutes																																																																												
Gravity	270°F/132°C	3 minutes																																																																												
Gravity	270°F/132°C	4 minutes																																																																												
Gravity	270°F/132°C	10 minutes																																																																												
Gravity	270°F/132°C	15 minutes																																																																												
Gravity	270°F/132°C	25 minutes																																																																												
Gravity	275°F/135°C	3 minutes																																																																												
Gravity	275°F/135°C	10 minutes																																																																												
Dynamic Air Removal	250°F/121°C	30 minutes																																																																												
Dynamic Air Removal	270°F/132°C	4 minutes																																																																												
Dynamic Air Removal	270°F/132°C	10 minutes																																																																												



510(k) Summary
3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES

510(k) Summary

Element	Subject Devices: (this submission) 3M™ Attest™ Steam Chemical Integrators	Predicate Device (K193254: Cleared on Dec 23 rd , 2019) 3M™ Attest™ Steam Chemical Integrators	Comparison																														
	<table border="1"> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>7 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>8 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>9 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>270°F/132°C</td><td>10 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>273°F/134°C</td><td>3 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>273°F/134°C</td><td>4 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>275°F/135°C</td><td>3 minutes</td></tr> </table>	Dynamic Air Removal	270°F/132°C	7 minutes	Dynamic Air Removal	270°F/132°C	8 minutes	Dynamic Air Removal	270°F/132°C	9 minutes	Dynamic Air Removal	270°F/132°C	10 minutes	Dynamic Air Removal	273°F/134°C	3 minutes	Dynamic Air Removal	273°F/134°C	4 minutes	Dynamic Air Removal	275°F/135°C	3 minutes	<table border="1"> <tr><td>Dynamic Air Removal</td><td>273°F/134°C</td><td>3 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>273°F/134°C</td><td>4 minutes</td></tr> <tr><td>Dynamic Air Removal</td><td>275°F/135°C</td><td>3 minutes</td></tr> </table>	Dynamic Air Removal	273°F/134°C	3 minutes	Dynamic Air Removal	273°F/134°C	4 minutes	Dynamic Air Removal	275°F/135°C	3 minutes	
Dynamic Air Removal	270°F/132°C	7 minutes																															
Dynamic Air Removal	270°F/132°C	8 minutes																															
Dynamic Air Removal	270°F/132°C	9 minutes																															
Dynamic Air Removal	270°F/132°C	10 minutes																															
Dynamic Air Removal	273°F/134°C	3 minutes																															
Dynamic Air Removal	273°F/134°C	4 minutes																															
Dynamic Air Removal	275°F/135°C	3 minutes																															
Dynamic Air Removal	273°F/134°C	3 minutes																															
Dynamic Air Removal	273°F/134°C	4 minutes																															
Dynamic Air Removal	275°F/135°C	3 minutes																															
Indicator Agent	Proprietary formulation.	Proprietary formulation.	Identical																														
Sterilization method and cycles	Steam sterilization processes 250°F to 275°F (121°C to 135°C)	Steam sterilization processes 250°F to 275°F (121°C to 135°C)	Identical																														
Endpoint Specifications (Minimum Stated Values)	<p>The minimum stated values for the 3M™ Attest™ Steam Chemical Integrators as determined using a resistometer are provided in the table below.</p> <table border="1"> <tr> <td>250°F/121°C</td> <td>270°F/132°C</td> <td>273°F/134°C</td> <td>275°F/135°C</td> </tr> <tr> <td>16.5 Minutes</td> <td>2.0 Minutes</td> <td>1.4 Minutes</td> <td>1.2 Minutes</td> </tr> </table>	250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C	16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes	<p>The minimum stated values for the 3M™ Attest™ Steam Chemical Integrators as determined using a resistometer are provided in the table below.</p> <table border="1"> <tr> <td>250°F/121°C</td> <td>270°F/132°C</td> <td>273°F/134°C</td> <td>275°F/135°C</td> </tr> <tr> <td>16.5 Minutes</td> <td>2.0 Minutes</td> <td>1.4 Minutes</td> <td>1.2 Minutes</td> </tr> </table>	250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C	16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes	Identical														
250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C																														
16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes																														
250°F/121°C	270°F/132°C	273°F/134°C	275°F/135°C																														
16.5 Minutes	2.0 Minutes	1.4 Minutes	1.2 Minutes																														
Shelf life	Three (3) years	One (1) year for 1243A and 1243B	Extended via identical protocol. The																														



510(k) Summary
3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES

510(k) Summary

Element	Subject Devices: (this submission) 3M™ Attest™ Steam Chemical Integrators	Predicate Device (K193254: Cleared on Dec 23rd, 2019) 3M™ Attest™ Steam Chemical Integrators	Comparison
			1243RE and RES was cleared for 6 months shelf life originally.

6. Nonclinical Comparison to the Predicate Device

3M™ Attest™ Steam Chemical Integrators are identical to the previously cleared devices which are sold under the same tradename 3M™ Attest™ Steam Chemical Integrators (K193254) and was tested using the identified test methodology shown below.

The results of performance testing on 3M™ Attest™ Steam Chemical Integrators demonstrate the device performs and meet the acceptance criteria shown below.

Table 6.1: Summary of Nonclinical Testing

Test Name	Purpose	Acceptance Criteria		Results
Stated Value (SV) Testing	To identify the critical parameters required to achieve a stated inactivation, by referring to a stated test organism with stated D and z values. The integrator must turn to “ACCEPT” end point at the stated value time and must also remain “REJECT” when exposed to conditions of -1°C/-15% set point of the SV time. All testing is completed in a saturated steam resistometer. Integrator temperature coefficient and correlation coefficient are calculated to confirm alignment to biological indicator performance.	SV at 250°F/121°C	≥ 16.5 minutes	Pass
		SV at 270°F/132°C	≥ 2.0 minutes	
		SV at 273°F/134°C	≥ 1.4 minutes	
		SV at 275°F/135°C	≥ 1.2 minutes	
		Integrator temperature coefficient	10 - 27°C	
		Correlation coefficient	≥ 0.9	
Health Care Facility Simulated Use Testing	Confirm integrators provide acceptable performance in cleared customer use sterilization cycles.	Device reaches “ACCEPT” endpoint reaction when exposed to customer use cycles. Device does not reach endpoint (“REJECT”) when exposed to failing conditions in customer use cycles.		Pass
Dry Heat Testing	Verify device requires the presence of saturated steam to turn to reach endpoint.	Endpoint must not be met following dry heat exposure at 140°C for 30 min.		Pass
Side-by-Side Testing with Biological Indicator	Confirm integrators are parallel in performance to biological indicators (BI).	Chemical integrator parallels performance of BI and does not reach endpoint before BI is inactivated.		Pass
Endpoint Color Stability	Confirm endpoint color stability for samples exposed to passing and failing	Endpoint decision must remain unchanged after 6 months.		Pass



510(k) Summary
3M™ Attest™ Steam Chemical Integrators
1243A, 1243B, 1243RE and 1243RES

510(k)
Summary

Test Name	Purpose	Acceptance Criteria	Results
	conditions in a steam resistometer.		

7. Conclusion

Based on the intended uses, technological characteristics and non-clinical performance data, 3M™ Attest™ Steam Chemical Integrators are as safe, as effective, and performs as well as or better than the legally marketed predicate device, Attest™ Steam Chemical Integrators cleared under K193254, Class II (21 CFR 880.2800(b)), product code JOJ.